

Tooth Brushes



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**



Tooth Brushes



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Minty Finger Tooth Wipes (PacDent)

Textured teeth wipes that can be used anytime, anywhere to quickly clean the plaque off of teeth and freshen breath without the use of water or rinsing. Sugar-free and individually wrapped.



100 x individually wrapped finger tooth wipes

54-45491 **28.30**

Ultra Water Jet Flosser (Waterpik)



- Unique combination of water pressure and pulsation reduces gingivitis, prevents gum disease
- Contemporary compact design
- Includes six tips: two standard, one Pik-Pocket tip, one tongue cleaner, one orthodontic, and one toothbrush

Each

46-91281

49.95

Water Pik Flosser (Waterpik)

- Quick, easy, convenient.
- Gentle, yet effective.
- Tip action of 10,000 strokes per minute effectively removes plaque above and below the gum line.
- Single filament easily and comfortably reaches between the teeth.
- Ideal for use with braces, crowns, bridges and other dental work
- Contemporary ergonomic design.
- One-handed flossing.
- Portability on the road.
- Developed by a dentist.
- Flosser is battery operated w/replacement tip cartridge.



Flosser, Each

46-91271

26.95

Replacement Tip Cartridge, 30/Pkg

46-912712

6.95

ClearGrip Adult Toothbrush



- Two-tone soft end rounded bristles with two rows of blue brush use bristles to indicate replacement
- Oval head and narrow neck
- Clear plastic handles with elastomeric grips in assorted colors

Compact Head, 72/Box

54-91601

54.95

Full Head, 72/Box

54-916012

54.95

Full Head w/Angled Handle, 72/Box

54-9160123

54.95

RotaPower Toothbrush (Gum)



- Compact brush head
- Two circles of bristles that are of different lengths
- Battery operated

Benefits:

- Makes it easier to clean harder-to-reach areas of the mouth
- The different bristles clean both the tooth and along the gum line
- Easy to use and store, requires no charging

Power Toothbrush, 12/Pkg

54-11640

69.95

Bright Tooth Brush (TC Dental)



Beautifully designed with an assorted attractive colors. Soft and comfortable bristles provides extra cleaning. Handles are designed with good and firm grip. Each children's toothbrush imprinted with attractive cartoon. Each is individually weal wrapped in see through plastics.

For adult, 72 pcs/box

54-11116

26.99

Child Toothbrush

3-Row, 27 Tufts



Child Soft, 72/Box

54-91641

39.95

Interdental Flossing Brushes (PlasDent)



2 mm-3 mm Tapered Asst Lite Colors 100/Pkg

54-11114

31.15

2 mm-4 mm Tapered Asst Neon Colors 100/Pkg

54-11113

31.15

3 mm-7 mm Tapered Asst Neon Colors 100/Pkg

54-11112

31.15



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Disposable Pasted Toothbrush

- 4 colors, popular handle, U S formulated mint paste on soft bristles
- Often called a 'bathroom brush' before seeing dentist/hygienist
- Individually cello wrapped for hygiene & to avoid flat bristles



100/Box 54-91591 19.95

Disposable Toothbrush (ATW)



- 4 colors, popular handle, U S formulated soft bristles
- Often called a 'bathroom brush' before seeing dentist/hygienist
- Individually cello wrapped for hygiene & to avoid flat bristles

100/Box 54-91581 18.95

Fresh R Tongue Cleaner (Gum)



Brush and scraper combo for better cleaning and better control of bad breath.

Fresh R Tongue Cleaner, 6/Box 46-11244 8.69

Go Between (Gum)

- Made with birchwood to avoid splintering
- Anatomically designed to be gentle on gums
- Small cylindrical shape fits in spaces as small as 0.8mm

- Perfect for patients with dental restorations or trouble flossing
- Antibacterial coating keeps bristles cleaner for up to two weeks



Angled Cleaner, 36/Box 54-11210 26.95

Moderate Taper, 36/Box 54-11213 26.95

Tooth Brushes

Infant Toothbrush (PlasDent)



Clear Silicone Material Dishwasher Safe, each 54-11111 1.45

Interproximal Brushes (Prestige)

Individually Wrapped, 100/Box 54-91391 24.95

Interproximal Brushes (DA)



- Elastomeric Handle w/extended handle cap
- Brushes in individually sealed bags
- Soft bristles on coated thin wire
- Travels easily with patient

Brush Cylinder, Blue, 72/Box 54-91381 27.95

Brush Cone, Green, 72/Box 54-91382 27.95

Standard Denture Carebrush (PlasDent)



Firm, Blue Color, 12/box 41-97991 11.95

Firm, Magenta Color, 12/box 41-97992 11.95

Firm, Assorted Colors, 12/box 41-97993 11.95

Toothpaste Squeezers (PlasDent)



Asst 6 Bright Colors, 144 pcs/bag 54-11119 22.46

Tooth Brushes



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Nylon Clear Toothbrushes



Similar to Oral B

Soft-Youth , 3 Rows 27 Tufts, 72/Box	54-91611	39.95
Adult-Soft , 4 Rows 40 Tufts, 72/Box	54-91612	39.95
Adult-Soft , 4 Rows 48 Tufts, 72/Box	54-91613	39.95

Premium Adult Toothbrushes



Removes plaque from hard-to-reach places

- Color-coded round front tip
- Rubber grips prevent slipping
- Assorted colors
- Compares to Oral B Advantage

Clear , 72/Box	54-91621	49.95
Opaque , 72/Box	54-916212	49.95

Proxabrush Handles (Gum)

- Raised Thumb grip improves patient control
- Familiar EEZ-LOK design for easy brush changes
- Assorted colors



Single-Ended , 12/Pkg	54-11314	11.95
Dual-Ended , 12/Pkg	54-11315	11.95
Snap On Handles , each	54-11316	11.95

Toothbrush Covers (PlasDent)



Asst Colors 144 pcs/bag	54-11118	22.46
--------------------------------	----------	--------------

Smile Tooth Brush (TC Dental)



Beautifully designed with an assorted attractive colors. Soft and comfortable bristles provides extra cleaning.

Handles are designed with good and firm grip. Each children's toothbrush imprinted with attractive cartoon.

Each is individually weal wrapped in see through plastics.

For adult , 72 pcs/box	54-11117	26.99
-------------------------------	----------	--------------

Soft Tooth Brush (TC Dental)



Beautifully designed with an assorted attractive colors. Soft and comfortable bristles provides extra cleaning. Handles are designed with good and firm grip. Each children's toothbrush imprinted with attractive cartoon. Each is individually weal wrapped in see through plastics.

For children , 72 pcs/box	54-11115	26.99
----------------------------------	----------	--------------

Supreme Toothbrushes



This massager-style toothbrush offers:

- Tapered head with power tip bristles
- Narrow neck for better visibility and flexibility
- Elastomeric fingers which massage the gums while brushing
- Ribbed tongue cleaner on back of brush head
- Ergonomic handle for better gripping
- Imprinting available

Assorted Colors , 72/Box	54-91631	49.95
---------------------------------	----------	--------------



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Proxabrush SNAP-ONS Intro Kits (Gum)



Intro Kit

1 handle and 2 refills, Taper, each 54-11250 **19.95**

Intro Kit

1 handle and 2 refills, Cylinder, each 54-1125 **19.95**

Toothpaste (Oraline)

Formulas are equivalent to popular branded flavored toothpastes.



Fluoride Bubblegum Flavored

Toothpaste For Kids of All Ages, 0.85oz. 144/Pkg 46-46251 **67.95**

Fluoride Mint Flavored Toothpaste For

Adults, 0.85oz. 144/Pkg 46-46252 **67.95**

Orthodontic Toothbrushes

V-Cut, rounded edges, assorted colors, sturdy nylon bristles.



Soft, Opaque, 6 Dozen 54-98171 **29.95**

Angle Toothbrush (Butler)



Contemporary style, ribbed comfort, longer handle angled at 10 degrees.

Tooth Brushes

Compact w/34 Tufts, Soft, 12/Pkg	54-98211	11.95
Full Size w/42 Tufts, Soft, 12/Pkg	54-98212	11.95

Premium Denture Carebrush (PlasDent)



Angle Brush & Ergonomic Handle

Soft Muave/White Color, 12/box	41-98001	14.95
Soft Tangerine/White Color, 12/box	41-98002	14.95
Soft Blue/White Color, 12/box	41-98003	14.95
Soft Assorted Colors, 12/box	41-98004	14.95

Aquafresh Tooth Brushes



Soft, 12/Pkg	54-98191	10.95
Infant, Each	54-98192	6.49
Kids, X-Soft, Each	54-98193	2.95

Delicate Post-Surgical Brush (Butler)



Delicate, with ultra gentle 0.004 bristles Ideal for post-surgical cleaning or in the presence of aphthous ulcers, mouth irritations, extractions, implants and grafts.

Brush w/31 Tufts, 12/Pkg	54-98201	14.95
---------------------------------	----------	--------------



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

