


CAMMY™

USB INTRAORAL CAMERA



CAMMY™
USB Intraoral Camera
VERSION 2.0 (a)
Serial No: 7074-0514-0690-1041
DENTAMERICA®

CAMMY Software

- Directly connects to a computer via the USB port.
- Designed with a freeze button to capture images and a joystick to control digital zoom.
- Self-sufficient, CCD technology, digital image processing.
- Four white LED lamp illumination to provide maximum realistic color reproduction.
- Built-in lens protector designed to prevent dust and scratches to optical components.
- Mini-head design allows for easy access to the third molar.
- Lightweight (only 92 grams), ergonomically designed handpiece.
- User-friendly software. Easy to use and integrate into patient management software.
- Digital images can be easily stored, retrieved and transmitted electronically.
- Disposable sleeves included to prevent cross contamination.

ISO 9001
ISO 13485



DENTAMERICA®

