


CAMMY™

USB INTRAORAL CAMERA



CAMMY Software

- Directly connects to a computer via the USB port.
- Designed with a freeze button to capture images and a joystick to control digital zoom.
- Self-sufficient, CCD technology, digital image processing.
- Four white LED lamp illumination to provide maximum realistic color reproduction.
- Built-in lens protector designed to prevent dust and scratches to optical components.
- Mini-head design allows for easy access to the third molar.
- Lightweight (only 92 grams), ergonomically designed handpiece.
- User-friendly software. Easy to use and integrate into patient management software.
- Digital images can be easily stored, retrieved and transmitted electronically.
- Disposable sleeves included to prevent cross contamination.

ISO 9001
ISO 13485



DENTAMERICA®



USB INTRAORAL CAMERA

1/4" HIGH RESOLUTION COLOR CCD

LENS PROTECTOR TO PREVENT ACCIDENTAL LENS SCRATCH

4 LED WHITE LIGHT, COLOR TEMPERATURE APPROX. 5000°K

ERGONOMIC COMFORT GRIP OFFERS SMOOTH AND EASY HANDLING

MINI-HEAD DESIGN ALLOWS FOR EASY ACCESS TO THE THIRD MOLAR

ERGONOMIC COMFORT GRIP OFFERS SMOOTH AND EASY HANDLING

ZOOM TOGGLE

FREEZE BUTTON

ON/OFF BUTTON

LIGHTWEIGHT HANDPIECE

SPECIALY DESIGNED HANDPIECE CONTOURED TO INTEGRATE INTO ALL DENTAL UNITS

HIGHLY FLEXIBLE CABLE (NO FRAGILE FIBER OPTIC CABLE) SINGLE PLUG MAKES CAMERA PORTABLE AND MOBILE



Reorder NO. 336
ZOOM LENS

- Use the zoom lens to magnify and examine caries, root canal access etc.
- Attach the surface viewer ring to examine occlusal surfaces.

SPECIFICATIONS

INTRAORAL USB CAMERA

Interface	USB 2.0
Sensor	1/4" SONY CCD Sensor
Effective Pixels	1.3Mega Pixels
Max. Resolution	1280*960
Resolution Combination	1280*960,640*480,352*288,320*240, 176*144,160*120
Color Depth	24-16bit(16.7 million color)
Max. Frame Rate	30fps (640*480:30fps,1280*960:15fps)
White Balance	Auto
Auto Exposure (AE)	Auto

Motion JPEG Compression	Yes(Smooth Video)
Rating Voltage	5V
Rating Power Consumption	1050mW
Rating Input Current	210mA
USB Cable Length	260 cm
USB Type	USB_A
Weight(g)	92g
Dimensions(L)X(W)X(H)(mm)	182mm*24mm*20mm
Working Temperature	5°~40°C
Working Humidity	20%~80%



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455

DENTAMERICA®