

# CAMREX™ 192

WIRELESS INTRAORAL CAMERA  
WITH MULTIPLE IMAGE MEMORY



- Built-in wireless receiver with 4 separate frequencies for optimal signal reception.
- Video output connects directly to a TV or computer monitor through AV or VGA input.
- Self-sufficient, CCD technology, digital image processing.
- Four white LED lamps to provide maximum realistic color reproduction.
- Built-in lens protector designed to prevent dust and scratches to optical components.
- Mini-head design allows for easy access to the third molar.
- Disposable sleeves included to prevent cross contamination.
- Rotates between live and freeze frame images up to 4 channels.
- Simultaneously view all four stored images in quadrant format or each individual image full screen.
- Easy and convenient battery replacement.
- Large capacity lithium battery provides 3 hours of continuous use on a full charge.
- NTSC or PAL available.


ISO 9001  
ISO 13485



**DENTAMERICA®**

# CAMREX™ 192


## WIRELESS INTRAORAL CAMERA WITH MULTIPLE IMAGE MEMORY



Multiple unit operation




VGA video output



Computer

RCA video output



TV

### SPECIFICATIONS

#### Wireless Intraoral Camera

	NTSC	PAL
Pick Up Element	1/4" SONY Interline Super HAD CCD	
Total Number of Pixels	410K	470K
Resolution	480TV Line	
Electronic Shutter	Automatic, 1/60(1/50)-1/100,000sec.	
Gain Control	Automatic	
Video Out	1.0Vp-p Composite/75ohm	
Power Supply	DC 3.7V	
Dimensions	L227 x W27 x H32.9mm	
Weight	107.3g	

1/4" HIGH RESOLUTION COLOR CCD •

4 LED WHITE LIGHT, COLOR •  
TEMPERATURE APPROX. 5000°K

#### Wireless Receiver

TV System	NTSC or PAL
Video In	1.0Vp-p, 75ohm Impedance
Video Out	1.0Vp-p, 75ohm Impedance
VGA Out	1024*768
Digital Resolution	8 bit 256 Grad, 512 x 1024 pixels
Digital I/O	16 bits
Singal	52dB
Power Source	DC 9V
Freeze Mode	Frame
Antenna Dimensions	L151 x W13mm
Dimensions	L172 x 108.5 x 31mm
Weight	264.6g
RCA/VGA Video cable length	2m

LIGHTWEIGHT AND  
ERGONOMIC COMFORT  
OFFERS SMOOTH  
& EASY HANDLING

4 FREQUENCY  
SETUP

#### Handpiece Charging Station

Power Source	DC 9V
Dimensions	L119 x W100 x H74 mm
Weight	103.5g


• MINI-HEAD DESIGN ALLOWS  
FOR EASY ACCESS TO THE  
THIRD MOLAR

• LOW BATTERY INDICATOR  
• CHANNEL BUTTON  
• FREEZE BUTTON  
• ON/OFF BUTTON

• REPLACEABLE BATTERY

• MODERN HANDPIECE  
CHARGER

• LED CHARGING  
INDICATOR



## Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455

DENTAMERICA®