

CELUX™

DENTAL OPERATORY LIGHT BY DENTAMERICA®



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455

- Sharply defined rectangular light pattern is evenly distributed at the oral cavity without causing discomfort to the patient's eyes.
- 3800° to 4500° Kelvin temperature is ideal for shade matching and examination.
- Easy reach ON/OFF/INTENSITY switch for both doctor and assistant.
- Optimum size reflector with dichroic coating filters out infrared rays resulting in a cool and clear working condition.
- 12V 55W tungsten halogen H-3 bulb with twin filaments.
- Easy bulb replacement.
- Easy open safety shield allows aseptic cleaning.
- Two-way switch for low (12000 LUX) and high (24000 LUX) intensity.

