

LITEX™ 685

PLASMA ARC CURING LIGHT

- State-Of-The-Art, Plasma-Arc-based, Composite Curing Light and Whitening Unit
- Select 1-3 Second Curing Time or Optional Step-Cure Cycle (Low/High Intensity)
 - Easily Rotated, Magnetically Held, Fiber Optic Probe
 - Handpiece Actuated Fingertip Pressure Ring
 - High-Resolution, Flexible Liquid Light Guide
 - Optional Whitening Fiber Optic Probe Available
- Eye Protection Glasses and Disposable Light Shields included
- Available for 100/110 Volts 50/60 Hz or 220/240 Volts 50 Hz Operation



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455