



# Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455



## ROTEX™ 782

### COMPACT ELECTRIC HANDPIECE UNIT

- Ball bearing high torque electric motor.
- 1,000 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Standard package includes Variable Speed Foot Control, Handpiece and Handpiece Rest.
- Circuit overload protection.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.
- Unit is available for 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001  
ISO 13488

