


# SCALEX™ 800

BY DENTAMERICA®

## ULTRASONIC SCALER

- 25,000 scaling strokes per second, solid state power unit.
- Streamlined compact size and slim light weight handpiece.
- Convenient foot control.
- No power loss at maximum pressure.
- Unit is available in 100/110V 50/60Hz, or 220/240V 50Hz operation.
- Dimensions 22L x 19W x 7.5H cm, weight 3.6kg.



**Prestige Dental Products, Inc.**

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455