

SCALEX™ 850

SELF-CONTAINED ULTRASONIC SCALER



Self-contained solution tank

SCALEX™ 850

- Self-contained solution tank allows for convenient portability.
- One touch selection between medicament tank and external water source.
- 25k/30k scaling strokes per second, solid state power unit.
- Innovative water-cooled circuitry designed for greater stability during operation.
- Available selection of a wide variety of high quality inserts to correspond with specific clinical applications.
- Unit is available in 100/110V 50/60Hz or 220/240V 50Hz operation.

ISO 9001
ISO 13485



DENTAMERICA®

SCALEX™ 850

ULTRASONIC SCALER

※ULTRASONIC INSERTS

Reorder No.80005
Universal Insert 25K
(External Water)



Reorder No.83005
Universal Insert 30K
(External Water)



Reorder No.80006
Universal Insert 25K
(Internal Water)



Reorder No.83006
Universal Insert 30K
(Internal Water)




Reorder No.80007
Chisel Type Insert 25K
(External Water)



Reorder No.83007
Chisel Type Insert 30K
(External Water)



Reorder No.80008
Chisel Type Insert 25K
(Internal Water)



Reorder No.83008
Chisel Type Insert 30K
(Internal Water)



Reorder No.80011
Endodontic Insert 25K
(Canalsonic Tip)



Reorder No.88060

Foot control

Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.

SPECIFICATIONS:

- * FREQUENCY: 25,000 or 30,000 Hz/sec.
- * DIMENSIONS: D260xW140xH210 mm
- * VOLUME: 1,000 ml
- * WEIGHT: 3.3 kg
- * INPUT: 115V/60Hz or 230V/50Hz
- * OUTPUT: AC 24V 3A



Reorder No.85003

Self-contained solution tank

Easily change the solutions according to different applications including periodontics, endodontics, and prophylaxis.



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455

DISTRIBUTOR

DENTAMERICA®