

SCALEX™ 880 PLUS

PIEZO ULTRASONIC SCALER



- Auto-sourcing the frequency of any compatible tip to achieve smooth operation.
- Load-adjusting function provides stable and effortless scaling.
- Special TURBO function helps remove especially hard calculi by increasing power relative to the hardness.
- Streamlined compact size and slim lightweight handpiece.
- Innovative circuitry designed for greater stability during operation.
- Extremely gentle on the patient.
- Available selection of a wide variety of high quality tips to correspond with specific clinical applications.
- Unit is available for 100-240V 50/60Hz operation.
- EMS type tips compatible.

ISO 9001
ISO 13485



0434

GMP

DENTAMERICA®

SCALEX™ 880 PLUS

