

SCALEX™ 880 PLUS

PIEZO ULTRASONIC SCALER



- Auto-sourcing the frequency of any compatible tip to achieve smooth operation.
- Load-adjusting function provides stable and effortless scaling.
- Special TURBO function helps remove especially hard calculi by increasing power relative to the hardness.
- Streamlined compact size and slim lightweight handpiece.
- Innovative circuitry designed for greater stability during operation.
- Extremely gentle on the patient.
- Available selection of a wide variety of high quality tips to correspond with specific clinical applications.
- Unit is available for 100-240V 50/60Hz operation.
- EMS type tips compatible.

ISO 9001
ISO 13485



0434

GMP

DENTAMERICA®

SCALEX™ 880 PLUS

SCALEX™ 880 PLUS

PIEZO ULTRASONIC SCALER

※PIEZO SCALER TIPS



Reorder No. 88006
Perio



Reorder No. 88004
Universal B



Reorder No. 88005
Chisel



Reorder No. 88003
Universal A



Reorder No. 88002
Endo Chuck

SPECIFICATIONS:

- * INPUT: 100V~240V 50/60Hz 0.8A
- * OUTPUT: DC 24V 1.5A
- * FREQUENCY: 27,000 to 32,000 Hz/sec
- * DIMENSIONS: W140xD140xH50mm
- * WEIGHT: 1,000g



Reorder No. 88015
Handpiece

Autoclavable handpiece with piezo ceramic crystals and titanium tip socket structure.



Reorder No. 88060
Foot Control

Unique compact foot control designed to be activated at any angle and accessible from 360° degrees.



Reorder No. 88009
Torque-limiting
Tip Wrench

Special autoclavable wrench prevents damage to scaler tip when tightening.



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455

DENTAMERICA®