




Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455



TELEREX™

VIDEO X-RAY FILM VIEWER

- Special CCD & LED technology offers instant viewing of dental X-ray film on a TV or computer monitor.
- View old X-ray film or digitize and archive them to your computer.
- Optional USB Video Grabber enables the Telerex to connect to the computer via the USB port.
- Improves communication between the doctor and patient with enlarged viewing.
- Compact and stylish design with spring loaded film tray, only 9.4x9.2x7.2 cm.
- Special film adapter included for Pedo size X-ray film.
- 100-240V universal voltage, NTSC or PAL available.

ISO 9001
ISO 13488

