

FENDERMATE[®] Matrix

The World's Fastest Matrix.

*Greatly
Improved!*

**F
A
S
T


A
N
D

E
A
S
Y**

Insert the one piece matrix FenderMate as a wedge. FenderMate with its optimal matrix curvature is placed in only 5 seconds.




Fill cavity and get a restoration with a tight contact and a tight cervical margin. Finish the restoration.



Design and patent pending

Available in Assortment Kit or Refills in Right and Left (Regular or Narrow)

FENDERMATE[®] offers a unique solution by combining a wedge and a matrix in its design for fast and efficient restorations.



- Flexible wing separates the teeth and firmly seals the cervical margin.
- Optimal Matrix Curvature.
- Pre-shaped Contact.
- Ready for Restoration.
- No Ring Needed.



Product	Item#
Asst. Kit (72pcs)	6027-70
Left Regular (dark green - 18pcs)	6027-71
Left Narrow (light green - 18pcs)	6027-72
Right Regular (dark blue - 18pcs)	6027-73
Right Narrow (light blue - 18pcs)	6027-74

Food approved plastic (Polypropylene) and color.

Use Left and Right for a MOD restoration



We Care About How You Practice Dentistry



Prestige Dental Products, Inc

3031 E. Coronado St., #D Anaheim, CA 92806

Phone # (877) 772-3888, Fax # (310) 388-8455

JS DENTAL
MANUFACTURING, INC.