

Air Abrasion Systems



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Air Abrasion Systems



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Aquacut Quattro (Velopex)

The **AQUACUT's** unique fluid curtain suppresses the powder emissions common with air abrasion units from the past, while simultaneously increasing cutting efficiency by constantly removing debris created during the cutting process. The use of fluid reduces the need for high air pressure level, thereby increasing patient comfort. When your patient sits down in the dental chair, you are able to perform a fast and effective cleaning of the oral area. This will present you with a clearer indication of the patients dental condition. At the turn of a switch, the **AQUACUT** changed into a powerful, cutting tool that will in many instance cancel the need to use the dental turbine. It is a non-contact instrument, there is no vibration, heat generation or turbine noise. Most cases there is no need for an anesthetic



The Aquacut Quattro is a unique new patented generation of cutting tool, taking the original concept of air abrasion and with the addition of fluid, creating Aquabrasion, a far more efficient and gentle technology. Besides functioning as a cutting tool, the Aquacut also cleans, polished and dries.

Unit	11-11118
Powder, 4/pk 29, micron	11-11117
Powder, 4/pk 53, micron	11-11116
Powder, 4/pk, sodium/bicarbonate	11-11115
Replacement Handpiece Blue, 0.5mm	11-11113
Replacement Handpiece Gold, 0.8mm	11-11111
Replacement Handpiece Green, 0.4mm	11-11112
Disposable Tips Universal, 50/pkg	AB-10001

PrepStart (Danville)

Air abrasion cavity prep unit

The PrepStart is full featured



PrepStart w/ Foot Pedal Installation 11-109751

Kit: Handpiece & nozzle kit aluminum oxide 27 & 50 micron video of clinical use, 1 year warranty.

Bonus Restorative Kit: Starflow flowable composite, microprime, caries finder, 2 disposable mirrors & sure etch

Handpiece Kit 1	11-109752
Handpiece Kit 2	11-109753
.015 Tip (80°), 3/pk	11-109754
.019 Tip (80°), 3/pk	11-109755

PrepStart H2O Hydro Abrasion Cavity Prep (Danville)

Available as a completely new system or as an add-on accessory. The H2O may be added to an existing PrepStart in about 15 minutes with a screwdrivers.

The H2O is attractively styled to look integrated with the PrepStart



No electric or water connections are used. Water is added to a removable tank

Complete Unit	11-919131
Accessory Unit	11-919132
Handpiece & .019 Nozzle	11-919133
Handpiece, only	11-919134
.019 Tip (80°), 1/pk	11-919135



Best Quality

Dental Supplies



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

