

Cosmetic & Restorative



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Light Cure Composites

Herculite XRV Unidose (Kerr)

It is a micro-hybrid dental composite. Light cure composites work for both posterior and anterior high wear resistance.
- Long lasting luster and smoothness
- Ideal Vita shade match & shade selection.

15 years clinical efficacy.



A1 Enamel Unidose , 20 x 0.25gm	23-11572
A2 Enamel Unidose , 20 x 0.25gm	23-11571
A3 Enamel Unidose , 20 x 0.25gm	23-11581
A3.5 Enamel Unidose , 20 x 0.25gm	23-11580
B1 Enamel Unidose , 20 x 0.25gm	23-11570
B2 Enamel Unidose , 20 x 0.25gm	23-11579

Activa Bioactive Kids (Pulpdent)

ACTIVA KIDS is an esthetic, BioACTIVE resin restorative that stimulates apatite formation and the natural remineralization process with release and recharge of calcium, phosphate and fluoride. It is an opaque white shade ideally suited for pediatric dentistry. ACTIVA is the first bioactive restorative with an ionic resin matrix, a shock-absorbing resin component, and bioactive fillers that mimic the physical and chemical properties of natural teeth. It delivers all the benefits of glass ionomers and has the esthetics, strength, and durability of composites.



Starter Kit: 5ml/8gm syringe opaque white shade, ACTIVA-SPENSER, + 20 automix tips with bendable 20-gauge metal cannula	23-110501
Single Refill: 5ml/8gm syringe opaque white shade, + 20 automix tips with bendable 20-gauge metal cannula	23-110502
Value Refill: 2 x 5ml/8gm syringes opaque white shade, + 40 automix tips with bendable 20-gauge metal cannula	23-110503

Beautiful II Syringes (Shofu)

A universal nanohybrid composite incorporating S-PRG technology that contains a stable phase of glass ionomer for aesthetic and long-lasting restorations that release and recharge fluoride. Beautiful II responds to the concentration of fluoride in the mouth while maintain the strength and stability of a composite resin. High filler load and high flexural strength ensure long-term stability under rigorous intraoral conditions for anterior and posterior restorations, direct laminate veneers and core build-ups



A1 Syringe Refill , 4.5g	23-105101
A2 Syringe Refill , 4.5g	23-105102
A3 Syringe Refill , 4.5g	23-105103
A3.5 Syringe Refill , 4.5g	23-105104
B1 Syringe Refill , 4.5g	23-105105
B2 Syringe Refill , 4.5g	23-105106
C2 Syringe Refill , 4.5g	23-105107

Beautiful II Unidose (Shofu)

- Ideal light transmission and optical characteristics impact excellent shade reproduction with a chameleon effect by stimulating the internal structure of natural teeth
- Anti-plaque effect minimized plaque adhesion, inhibits bacterial colonization and plaque accumulation
- Fast & easy polishing creates an outstanding surface luster that remains over time



A1, 20 Tips Refill , .25g	23-105401
A2, 20 Tips Refill , .25g	23-105402
A3, 20 Tips Refill , .25g	23-105403
A3.5, 20 Tips Refill , .25g	23-105404
A4, 20 Tips Refill , .25g	23-105405
B1, 20 Tips Refill , .25g	23-105406
B2, 20 Tips Refill , .25g	23-105407
C2, 20 Tips Refill , .25g	23-105408
Incisal, 20 Tips Refill , .25g	23-105409
Bleach White, 20 Tips Refill , .25g	23-105410
A20, 20 Tips Refill , .25g	23-105411
A30, 20 Tips Refill , .25g	23-105412

Clearfil Majesty ES Flow (Kuraray)



Light cure universal flowable composite that contains 75 wt% filler and is indicated for the following restorative applications.

- Direct restorations for all cavity classes, cervical lesions, tooth wear, and tooth erosion
- As a cavity base/liner
- Corrections of tooth position and shape
- Intraoral repair of fractured restorations

A1, 20 Tips Refill , .25g	23-105401
A2, 20 Tips Refill , .25g	23-105402
A3, 20 Tips Refill , .25g	23-105403
A3.5, 20 Tips Refill , .25g	23-105404
A4, 20 Tips Refill , .25g	23-105405

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

B1, 20 Tips Refill, .25g 23-105406

Clearfil Majesty ES-2 Premium PLT (Kuraray)

A redefining composite technology. With the only VITA approved shade concept, Kurarays new composite material covers 15 VITA classical shades with only five enamel and dentin shades.



A1E, 10/pkg 23-11169
A2E, 10/pkg 23-11168
A3E, 10/pkg 23-11167

Denfil Composite Capsules (Vericom)

- Denfil is a light-cured restorative micro hybrid composite resin for use in both posterior and anterior restorations
- It is designed for minimizing polymerization shrinkage by inorganic fillers high loading, this Denfil has excellent marginal sealing and surface texture
- By adopting optimal polymerization initiator system, Denfil has high strength and abrasion resistance
- As Denfil contains radiopaque filler, any diagnosis using X-ray is convenient
- Also it has optimal viscosity for restoration, Denfil shows excellent workability



A1 Capsule Refill Shade, 20 x .25g 23-11206
A2 Capsule Refill Shade, 20 x .25g 23-11192
A3 Capsule Refill Shade, 20 x .25g 23-11204
A3.5 Capsule Refill Shade, 20 x .25g 23-11203
A4 Capsule Refill Shade, 20 x .25g 23-11202
B1 Capsule Refill Shade, 20 x .25g 23-11201
B2 Capsule Refill Shade, 20 x .25g 23-11200
C2 Capsule Refill Shade, 20 x .25g 23-11198
D3 Capsule Refill Shade, 20 x .25g 23-11197

Denfil Composite Syringes (Vericom)

- Denfil is a light-cured restorative micro hybrid composite resin for use in both posterior and anterior restorations
- It is designed for minimizing polymerization shrinkage by inorganic fillers high loading, this Denfil has excellent marginal sealing and surface texture
- By adopting optimal polymerization initiator system, Denfil has high strength and abrasion resistance
- As Denfil contains radiopaque filler, any diagnosis using X-ray is convenient



- Also it has optimal viscosity for restoration, Denfil shows excellent workability

A1 Syringe Refill, 4gm 23-11229
A2 Syringe Refill, 4gm 23-11230
A3 Syringe Refill, 4gm 23-11231
A3.5 Syringe Refill, 4gm 23-11233
A4 Syringe Refill, 4gm 23-11232
B1 Syringe Refill, 4gm 23-11236
B2 Syringe Refill, 4gm 23-11228
C2 Syringe Refill, 4gm 23-11226
Enamel Syringe Refill, 4gm 23-11224

Embrace WetBond Class V (Pulpdent)

- Bonds to the slightly moist tooth
- Self-etching to dentin
- Treats and desensitizes abfractions and exposed roots
- Tough. Resilient. Durable
- Excellent marginal adaptation, adhesion and retention
- Cures with all lights

Compressive Strength: 38,000 (1700) p.s.i. 260 (10) MPa
Diametral Tensile Strength: 8,650 (670) p.s.i. 60 (5) MPa
% Filler by Weight: 65%
% Solubility: 0.1%
Film Thickness: 18 microns



Class V Kit: 1.2ml syringe each A2, A3.5, B2, D2, Seal-n-Shine + accessories 23-88491
Class V Refill: 1.2ml syringe A2 + accessories 23-88492
Class V Refill: 1.2ml syringe A3.5 + accessories 23-88493
Class V Refill: 1.2ml syringe B2 + accessories 23-88494
Class V Refill: 1.2ml syringe D2 + accessories 23-88495

Estelite Flow Quick Syringe (Tokuyama)

ESTELITE FLOW QUICK and ESTELITE FLOW QUICK HIGH are universal soprano filled resin composites offering a choice of two optimal flowabilities, superior esthetics, ideal handling, and excellent gloss and polish ability, all in one system.



Building on the proven success of the spherical filler and the patented radical amplified photopolymerization initiator technology, Estelite Flow Quick and Estelite Flow Quick HF provide an extended working time and quicker curing time.

Low Flow Syringe (A1), 3.6gm + tips 23-11245

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Low Flow Syringe (A2), 3.6gm + tips	23-11258
Low Flow Syringe (A3), 3.6gm + tips	23-11259
Low Flow Syringe (A3.5), 3.6gm + tips	23-11260
Low Flow Syringe (A4), 3.6gm + tips	23-11261
Low Flow Syringe (B1), 3.6gm + tips	23-11262
Low Flow Syringe (B2), 3.6gm + tips	23-11267
Low Flow Syringe (B3), 3.6gm + tips	23-11266
Low Flow Syringe (B4), 3.6gm + tips	23-11265
Low Flow Syringe (C1), 3.6gm + tips	23-11264
Low Flow Syringe (C2), 3.6gm + tips	23-11257
Low Flow Syringe (C3), 3.6gm + tips	23-11256
Low Flow Syringe (OA1), 3.6gm + tips	23-11255
Low Flow Syringe (OA2), 3.6gm + tips	23-11246
Low Flow Syringe (OA3), 3.6gm + tips	23-11247
Low Flow Syringe (CE), 3.6gm + tips	23-11248
Low Flow Syringe (BW), 3.6gm + tips	23-11249
High Flow Syringe (A1), 3.6gm + tips	23-11250
High Flow Syringe (A2), 3.6gm + tips	23-11251
High Flow Syringe (A3), 3.6gm + tips	23-11252
High Flow Syringe (OPA2), 3.6gm + tips	23-11253
Syringe Tips, 50/box (Metal)	23-11254
Syringe Tips, 50/box (Plastic)	23-11263

Estelite Sigma Quick PLT (Tokuyama)



PLT (A1), 0.2gm x 20	23-11280
PLT (A2), 0.2gm x 20	23-11268
PLT (A3), 0.2gm x 20	23-11281
PLT (A3.5), 0.2gm x 20	23-11282
PLT (A4), 0.2gm x 20	23-11283
PLT (A5), 0.2gm x 20	23-11284
PLT (B1), 0.2gm x 20	23-11285
PLT (B2), 0.2gm x 20	23-11286
PLT (B3), 0.2gm x 20	23-11287
PLT (B4), 0.2gm x 20	23-11279
PLT (C1), 0.2gm x 20	23-11278

PLT (C2), 0.2gm x 20	23-11277
PLT (C3), 0.2gm x 20	23-11269
PLT (OA1), 0.2gm x 20	23-11270
PLT (OA2), 0.2gm x 20	23-11271
PLT (OA3), 0.2gm x 20	23-11272
PLT (OPA2), 0.2gm x 20	23-11276
PLT (CE), Clear enamel, 0.2gm x 20	23-11275
PLT (BW), Bleach white, 0.2gm x 20	23-11274
PLT (WE), White enamel, 0.2gm x 20	23-11273
Shade Guide	23-11288

Estelite Sigma Quick Syringe (Tokuyama)



Syringe (A1), 3.8gm	23-11308
Syringe (A2), 3.8gm	23-11297
Syringe (A3), 3.8gm	23-11295
Syringe (A3.5), 3.8gm	23-11296
Syringe (A4), 3.8gm	23-11294
Syringe (A5), 3.8gm	23-11293
Syringe (B1), 3.8gm	23-11292
Syringe (B2), 3.8gm	23-11291
Syringe (B3), 3.8gm	23-11290
Syringe (B4), 3.8gm	23-11307
Syringe (C1), 3.8gm	23-11298
Syringe (C2), 3.8gm	23-11299
Syringe (C3), 3.8gm	23-11306
Syringe (OA1), 3.8gm	23-11305
Syringe (OA2), 3.8gm	23-11304
Syringe (OA3), 3.8gm	23-11303
Syringe (OPA2), 3.8gm	23-11302
Syringe (CE), Clear enamel, 3.8gm	23-11301
Syringe (BW), Bleach white, 3.8gm	23-11300
Syringe (WE), White enamel, 3.8gm	23-11289

www.prestigedentalproducts.com

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Filtek Z250 Capsule (3M)

Low polymerization shrinkage for less stress on the tooth and less potential for patient sensitivity.

Speed and convenience: easy to use, fast placement technique, and easy to finish to a smooth surface

Excellent handling with exceptional resistance to slumping

Strong, long-lasting & natural looking results

Outstanding versatility for any area of the mouth

A1 Refill Capsule, 20 x 0.2gm	23-11325
A2 Refill Capsule, 20 x 0.2gm	23-11323
A3 Refill Capsule, 20 x 0.2gm	23-11324
A3.5 Refill Capsule, 20 x 0.2gm	23-11327
A4 Refill Capsule, 20 x 0.2gm	23-11326
B1 Refill Capsule, 20 x 0.2gm	23-11329
B2 Refill Capsule, 20 x 0.2gm	23-11328
C2 Refill Capsule, 20 x 0.2gm	23-11330



Filtek Z250 Syringe (3M)

Low polymerization shrinkage for less stress on the tooth and less potential for patient sensitivity.

Speed and convenience: easy to use, fast placement technique, and easy to finish to a smooth surface

Excellent handling with exceptional resistance to slumping

Strong, long-lasting & natural looking results

Outstanding versatility for any area of the mouth

A1 Refill Syringe, 4gm	23-11331
A2 Refill Syringe, 4gm	23-11333
A3 Refill Syringe, 4gm	23-11332
A3.5 Refill Syringe, 4gm	23-11334
B1 Refill Syringe, 4gm	23-11336



G-aenial Sculpt (GC America)

Universal composite with high density

Uniform dispersion *nano-filler* that creates a *self-shining* effect enabling a user to complete a restoration by doing simple rough finishing with a bur to bring the restoration to a brilliant luster and enables increased polish over time with an excellent shade match to enamel and to G-aenial Universal Flo when used as a system.

Intro Kit: 1 syringe of each A1, A2 & A3 23-108031 (2.0ml & 1 G-Premio bond 5ml kit)



Gradia Direct (GC America)

A visible light-cured micro-filled hybrid resin (MFR) composite restorative.

Recognizing the differences in radiopacity, polish ability, wear & durability, GC Gradia Direct is available in anterior & posterior versions

It provides the best option for optimal esthetics with its excellent shade matching ability

Unitip Refill, #BW, 10/pkg 23-108211



Herculite XRV Syringes (Kerr)

A micro hybrid dental composite developed using the Vita shade system.



Light-Cure composites for both posterior & anterior

A1 Enamel Syringe, 5gm	23-11578
A2 Enamel Syringe, 5gm	23-11577
A3 Enamel Syringe, 5gm	23-11576
A3.5 Enamel Syringe, 5gm	23-11575
B1 Enamel Syringe, 5gm	23-11574
B2 Enamel Syringe, 5gm	23-11573

Nano Hybrid Composite (Made in USA)

- Easy to control during placement
- Not sensitive to ambient light
- Adapts well to cavity preparations
- Does not stick to the placement instruments
- Can be contoured to the desired shape before curing
- Highly rated for blending with tooth structure, radiopacity, and ease of finishing
- Finished and polished to a high luster
- Fluoride release



A1 Syringe Refill, 4.5gm	23-11406
A2 Syringe Refill, 4.5gm	23-11405
A3 Syringe Refill, 4.5gm	23-11404
B1 Syringe Refill, 4.5gm	23-11402
A1 Compules Refill, 15 x 0.25g	23-11400
A2 Compules Refill, 15 x 0.25g	23-11399
A3 Compules Refill, 15 x 0.25g	23-11398
A3.5 Compules Refill, 15 x 0.25g	23-11397
B1 Compules Refill, 15 x 0.25g	23-11396
C1 Compules Refill, 15 x 0.25g	23-11408

Cosmetic & Restorative

Nexcomp (Meta Biomed)

- Excellent physical properties
- Less shrinkage
- Color match to the Vita shade guide
- Low 1.6% shrinkage for superior marginal integrity
- Filled with 40 nano sized filler for excellent polish ability
- Optimum handling - Excellent sculptability for restoring natural shape and function. Non slumping. Non-sticky

A1, 4g 1 syringe	23-11409
A2, 4g 1 syringe	23-11411
A3, 4g 1 syringe	23-11412

Prestige Nanohybrid Composite

- Prestige** is an ideal micro hybrid composite for all cavity classes. It is easy to polish and has good wear resistance. It is very color stable and matches Vita shades. Handling qualities that are superior and is not sticky, so application is fast & easy
- Universal for all cavity classes I-V
 - Non-sticky
 - High translucency
 - Excellent Polish ability
 - High color stability
 - Great aesthetics

A1, 4gm syringe	23-110461
A2, 4gm syringe	23-110462
A3, 4gm syringe	23-110463
A3.5, 4gm syringe	23-110464
B1, 4gm syringe	23-110465
B2, 4gm syringe	23-110466
OA1, 4gm syringe	23-110467
OA2, 4gm syringe	23-110468
OB1, 4gm syringe	23-110469

Spectrum TPH 3 Compules Composite Restorative



A1 Shade, compules refill 20/box	23-110661
A2 Shade, compules refill 20/box	23-110662
A3 Shade, compules refill 20/box	23-110663
A3.5 Shade, compules refill 20/box	23-110664



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Z100 Restorative Capsule (3M)

A1 Refill Capsule, 18 x 0.2gm	23-11582
A2 Refill Capsule, 18 x 0.2gm	23-11583
A3 Refill Capsule, 18 x 0.2gm	23-11584
A3.5 Refill Capsule, 18 x 0.2gm	23-11585



Z100 Restorative Syringe (3M)



Has a single filler which is **100% zirconia/silica**

A1 Refill Syringe, 4gm	23-11589
A2 Refill Syringe, 4gm	23-11590
A3 Refill Syringe, 4gm	23-11591
A3.5 Refill Syringe, 4gm	23-11592

Bonding Agents

G-Bond (GC America)

G-Bond is a new revolutionary **7th generation adhesive** that takes the guesswork out of bonding. The unique combination of phosphoric acid ester monomer and 4-MET adhesive technology creates superior etch and adhesion to enamel in addition to providing chemical and mechanical seal to dentin - referred to as the Nano Interaction Zone



Refill, 5ml bottle	23-21357
--------------------	----------

Optibond Solo Plus (Kerr)

5th generation bonding
Single liquid adhesive system

- One-step bonding
- Direct & indirect bonding
- Filled 15% with 0.4 micron glass
- Releases fluoride



Refill, 5ml bottle	23-11439
Unidose Refills 100-10ml, optibond solo plus unidose pkts.	23-11438

Clearfil SE Bond (Kuraray)

6th generation bonding. Etch, prime, and desensitize in 20 sec.

Light cure bonding system containing a water based primer

New, unique dispenser allows



Cosmetic & Restorative

single handed delivery of both SE Primer and SE Bond

Self-etching bonding systems virtually eliminate sensitivity.

No mixing, No rinsing, No sensitivity!!

Standard Kit: 1 primer 6ml, 1 bond 5ml 23-11173

Value Kit Contains: 3 - 5ml (bond), 3 - 6ml (primer) 23-11172

Clearfil SE Bond 2 (Kuraray)

- Universal clinical applications with new Clearfil DC activator
- New catalyst system for stronger bond strengths
- Exceptional, unparalleled, long-term, class-V, clinical research with Clearfil self-etch products
- Low film - Thickness with indirect restorations
- Ideal Bulk - Fill adhesive (Absorbs stress & maintains seal on Dentin)
- Ideal adhesive layer for direct restorations
- No-shaking, no-mixing, no-scrubbing, no multiple applications
- SE Primer 2 treats zirconia & alumina restorations
- No light-curing inside of indirect restorations



Kit: 1 primer (6ml), 1 bond (5ml), 1 dispensing dish, 1 light-blocking plate, 1 outer case, 1 applicator brush (fine 50pcs) 23-105851

G-aenial Bond (GC America)

A 7th generation one-step self-etch light-cured bonding agent

Highest bond strengths to enamel of all leading one-step and two-step self-etch bonding agents

Strengths to both enamel and dentin



5ml Refill, bottle 23-105631

Kit: 5ml with 50 tips, 20 trays 23-105632

Adper Single Bond (3M)

5th generation bonding, total-etch adhesive. Adper single bond plus adhesive, proven by 10 years of clinical performance.

- Excellent shear bond strength provides exceptional bonding performance
- Easy-to-use unit dose - simply press one chamber and use
- Nanofiller is stable and will not settle out of dispersion
- Unique, easy-dispensing bottle permits visual inspection of exactly how much adhesive remains in the bottle
- Cap snaps securely shut minimizing evaporation
- No shaking required - use straight out of the vial
- Ethanol/water-based adhesive maintains its shear bond strength over time to yield an outstanding shelf life.



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

6gm Refill, bottle 23-11112

Adper Prompt L-Pop Giant, 100/pkg 23-11111

BC Plus (Vericom)

BC-Plus is the 5th generation single component bonding agent designed to bond all classes of direct composite restorations to dentin, enamel, as well as procedures involving cast metals, composite, treated porcelain and set amalgam.



- BC-Plus is a light activated adhesive resin base on an ethanol formulation.
- BC-Plus allows for indirect restoration procedures, due to its minimal and uniform film thickness

5ml bonding agent, bottle 23-11117

Beautibond (Shofu)

Beautibond is ideal for highly aesthetic and minimally invasive restorations.

With excellent biocompatibility and bonding durability, Beautibond is an all-in-one adhesive that enables etching, priming and bonding in one simple step.



Kit: 50 - 0.1ml unit dose, 50 micro brush 23-16941 fine

Bond Force (Tokuyama)

7th generation, single-component bonding agent is self-etching and fluoride releasing. Bond force also contains a patented self-reinforcing (SR) monomer that creates a 3D link between the tooth and enamel and dentin. Bond force is indicated for bonding of light- or dual-cured composite material to cut or uncut enamel and dentin and to repair fractured porcelain or composite. It is available in both unit doses and bottles



Refill 5ml, bottle 23-11133

Standard Kit: 50 unit dose. 50 brushes 23-11135

Bond Force Pen (Tokuyama)

Bond Force Pen is easy to use and **can eliminate sensitivity issues** when used as recommended.

The bonding agent is available in a pen dispenser with a simple 2-click releaser



Bond Force Pen, Intro kit w/ tips 23-11596

Bond Force Pen, twin refill 2 x 2ml 23-11597

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Brush & Bond (Parkell)

Easy-to-use, high-strength bonding agent that protects and desensitizes



- Bonds direct and indirect composite restorations to all tooth surfaces
- Permanently seals the tubules to prevent sensitivity
- Protects the pulp against microleakage and bacterial invasion
- Creates an acid and wear-resistant shield for the tooth
- Doubles the surface hardness of the preparation
- Cures with all currently available curing lights
- Low film thickness won't affect crown fit

Kit w/ Standard Activator Brushes 23-16961

Include: 3ml dentin/enamel bonding liquid, mixing well plus, & 100 standard activator brushes

Bond Only, 3ml dentin/enamel bonding liquid 23-16962

Standard Activator Brushes, 100/pcs 23-16963

Clearfil Repair Kit (Kuraray)

A light-cure bonding system for repairing fractured porcelain, all ceramics or composite resin facing crowns & bridges. This system is also used for aesthetic repair intraorally.



Features & Benefits

- Simple, fast, reliable repair of chipped or fractured porcelain, composite resin, or metal oxide restoration
- Very high bond strength even to metal oxides, like zirconia
- No toxic hydrofluoric acid needed
- Long, stable shelf-life (two years from manufacturing)
- Re-glazing porcelain and composite restorations

Kit Contains: SE primer, SE bond, porcelain bond activator, clearfil ST opaquer, k-etchant gel & accessories 23-88461

Porcelain Bond Catalyst, liquid 23-88462

Porcelain Bond Universal, liquid 23-88463

Porcelain Bond Activator 23-88464

Clearfil Alloy, primer 23-88465

Clearfil S3 Bond Plus (Kuraray)

A light-cure, self-etching, one-step, fluoride releasing bonding agent. Based on the same progressive adhesive technologies from the gold standard Clearfil SE Bond, you can expect minimal postoperative sensitivity.



With its one-bottle liquid, Clearfil S3 Bond Plus dramatically simplifies bonding procedures. Along with Kuraray's MDP

technology, a new catalyst has been adopted; thereby, providing higher bond strength with improved working time.

Set Contains: 1 bond (4ml), 1 applicator brush (fine 50pcs), 1 dispensing dish, 1 light-blocking plate 32-50681

Value Pack Contains: 3 x 4ml bond 32-50682

Single Dose Pack Contains: 50 unit-dose 0.1ml/tip 32-50683

Single Dose Value Pack Contains: 100 unit-dose 0.1ml/tip 32-50684

Clearfil SE Protect 6th Generation Bonding (Kuraray)

Contains the same proprietary self-etching, light cure technology, plus two additional proprietary technologies; long-term fluoride release and the MDPB monomer, which has antibacterial cavity cleansing effect. Clearfil SE Protect, while having the same physical properties, features and as SE Bond, is actually more than an adhesive bonding agent... it is a revolutionary adhesive



Features & Benefits

- Antibacterial cavity cleansing properties
- Low post operative sensitivity
- High bond strength
- Fast & simple procedure (no rinsing)
- Fluoride releasing properties

Kit Contains: 1 clearfil SE protect primer 6ml, 1 clearfil SE protect bond 5ml, 2 brush holders (white & black), 1 mixing dish, 1 light-blocking plate, 1 outer case, 2 disposable brush tips (50pcs) 23-50691

Value Kit: 5ml (bond) x 3, 6ml (primer) x 3 23-50692

Bond Refill, 5ml 23-50693

Primer Refill, 6ml 23-50694

DenTASTIC All-Purpose Dental Adhesive System (Pulpdent)

This 4th generation adhesive system offers light cure and self cure options. The adhesive primer and resin bonding agent are separate components.



Patented chemistry features PMGDM technology with the magnesium salt of NTG-GMA, which provides higher bond strengths than the older PMDM formulas.

System Kit: 6ml primer part A, B, & C, 6 ml resin bonding agent part 1 & 2, etch-all, porcelain etch gel, silane, brush tips, mixing wells, applicator tips 23-88101



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Cosmetic & Restorative

Refill Kit: 6ml primer part A, B, & C, 6ml resin bonding agent part 1 & 2, etch-all	23-88102
6ml Primer Part A , dentin/enamel initiator	23-88103
6ml Primer Part B , universal base	23-88104
6ml Part C , metal initiator	23-88105
6ml Unfilled Resin Bonding Agent Part 1 , base, light activated	23-88106
6ml Unfilled Resin Bonding Agent Part 2 , catalyst	23-88107
Accessory Kit: 50 mixing wells, 20 etch gel applicator tips, 100 brush tips, 2 brush handles	23-88108

DenTASTIC UNO (Pulpdent)

A **one step wet bonding** for all direct bonding **light cure applications**

5th generation light cure; single component dental adhesive.

DenTASTIC UNO combines the patented PMGDM adhesive primer with hydrophilic resins and acetone solvent into an exceptionally strong single component adhesive.



Refill, 6ml bottle 23-11237

Standard Kit Includes: 2-6ml bottles Dentastic Uno, 1-5ml etch-rite etching gel, 20 applicator tips, 50 brush tips, 1 brush handle, and instructions for use 23-11238

DenTASTIC UnoDuo (Pulpdent)

5th generation dual cure, simplified dual cure adhesive for all direct and indirect clinical applications.

Pulpdent has developed DUO as the dual cure catalyst for UNO. UNO-DUO is the simplest dual cure adhesive available and is the top rated 5th generation adhesive for use with self-cure and dual cure core materials



Kit Contains: 6ml UNO, 3ml DUO, 5ml etch-rite, 20 applicator tips, 50 brush tips, brush handle 23-11239

E-BOND (Danville)

Light-curable, bonding agent for use on enamel, composite and etched silanated porcelain, unfilled resin, with ultra-thin film thickness



10ml, bottle 23-88421

Embrace Wetbond Rest. & PFM Repair (Pulpdent)

Solve Big Problems...Fast!!

Primes, protects, opaque's, seals, finishes & polishes; **no bonding agents required**

The latest generation of porcelain repair products uses first-coat in place of silane. First-coat is a solvent-free resin with strong bonds to ceramic and metal. Testing shows better results than silane.



Kit: 1.2ml syringe each first-coat, seal-n-shine, opaquer, porcelain etch gel, kool-dam + accessories 23-11243

IntegraBond (Premier)

This 5th generation bonding agent works with all composite materials and resin cements, light-cured, dual-cured and auto-cured. Easy-to-use, one-step, universal primer/bonding agent seals the dentin tubules to eliminate sensitivity while protecting the tooth.



Kit: 4ml integrabond, 4ml auto-cure activator, 2-1.2ml primer etch syringes, 50 microbrushes, 20 disposable etch tips, mix wells 23-16951

Univ. Bonding Agent, 7ml each 23-16952

Auto-Cure Activator, 7ml each 23-16953

Integrabond, 4ml each 23-16954

One-Step Bonding Adhesive (USA)

A **5th generation**, single component, multi-use adhesive, which is hydrophilic and suitable for use on moist dentin (wet bonding technique)



The adhesive incorporates aspects of a primer and bonding resin in a single bottle designed to bond composites and compomers to dents and enamel as well as treated amalgam, metal and ceramic surfaces

After etching, penetrations and saturation of the dentin and enamel are accomplished in 2 to 3 coats

One-Step Bonding Adhesive, 7ml 23-11437

One-Up Bond F (Tokuyama)

Light-cured bonding agent incorporates a color-chromatic agent that is visible during application and invisible when fully polymerized. Eliminating post operative sensitivity, the self-etching, one-step, agent releases fluoride and exhibits a bond strength of 20MPa to enamel and 18MPa to dentin.



Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

One-Up Bond F Liquid A, 5ml	23-11434
One-Up Bond F Liquid B, 5ml	23-11435
One-Up Bond F Applicator Brushes, 100/pkg	23-11598

One-Up Bond F Plus (Tokuyama)

A light-cured, self-etching bonding agent that permits a simplified procedure by eliminating the pretreatment of the tooth when performing direct composite restorations. One-Up Bond F Plus offers a quick, single application, with no etching, no rinsing and no drying. It requires 3 simple steps. (1) mix (2) apply (3) light-cure and each step can be visually confirmed by the built-in color coding, eliminating guesswork



One-Up Bond F Plus Kit: bonding agent A & B, 5ml bottle each, mixing well, applicator tips & handle 23-11433

Disposable Applicator Tip & Handle Set, 4 handles & 500 tips 23-11432

Disposable Brush Tips, 2 handles & 100 tips 23-11436

Optibond All-In-One (Kerr)

7th generation bonding, self etch essential ingredients needed etching, priming, bonding into a single solution.



Adhesive Kit: 1 bottle of adhesive (5ml), 3 unidose devices (0.18ml each), 50 kerr applicators, 25 disposable mixing wells, directions 23-99071

Optibond FL (Kerr)

4th generation bonding

- Light cure adhesive offering a simple two bottle process with only one light curing step
- Fluoride releasing



Prime, 8ml bottle 23-105711

Adhesive, 8ml bottle 23-105712

P & Bond (Meta Biomed)

Advantages

- Superior bonding strength
- **5 generations adhesive**
- One step priming & bonding
- Ethanol base adhesive
- Wet bonding type



Meta P & Bond, 5g 23-11363

Photo Bond Dual-Cure Bonding Agent (Kuraray)

Used for the total-etch technique. It was developed for use with light-cure restorative materials. It adheres to metal and amalgam, as well as to porcelain when mixed with porcelain bond activator

Features & Benefits

- High bond strength to enamel and dentin
- Dual-cure
- Simple and fast procedure
- Economical



Complete Kit: 1 catalyst 6ml, 1 universal 6ml, 1 k-etchant gel 6ml, 1 brush, 1 mixing dish, 1 sponge pledget 23-88451

Universal Liquid, 6ml 23-88452

Catalyst, 6ml 23-88453

Prelude One Bond (Danville)

7th generation bonding agent
A light-curable **single component self-etch adhesive** designed to bond composite resins to most dental surfaces



These surfaces include **dentin, enamel, zirconia, alumina, metals, porcelain, lithium disilicate, and composite restorations.**

As a universal adhesive, prelude one is truly compatible with light cured, dual cured, and auto cured composites (no activator required)

5ml, bottle 23-109141

Prelude SE (Danville)

Independent tests confirm that prelude has the very high 5 minute and 24 hour bond strengths. It also has the smallest standard deviation, proving Prelude to be very technique insensitive.



Prelude is extremely forgiving in the presence of moisture. Too wet, or too dry is not a concern. Its water/ethanol formulation tolerates an environment that would adversely affect other products. The remarkably low 5 micron film thickness means fit is not affected and it will not present a radiolucent layer. Bond strengths to cut and un-cut enamel are extremely high. It is ideal for Cerec** and other milled systems.

Self Etch Kit: Primer (5ml), link (5ml), mixing wells (10), microbrushes (40) 23-11510

Total Etch Kit: Adhesive (5ml), link (5ml), microprime (5ml), sureEtch (5gm), mixing wells (10), microbrushes (40) 23-11509

Primer/Adhesive, 2/pk (5gm) 23-11508

Prime, (5gm) 23-11511

Adhesive, (5gm) 23-11512

Link, (5gm) 23-11513

Cosmetic & Restorative

Proline Bonding Agent

5th generation

Light cure dentin/enamel one step bonding

Made in USA

7ml, bottle 23-151171



ProLink (Silmet)

6th generation bonding, single step to dentin, enamel & all popular light, dual or auto-cure composites

5ml, bottle 23-11516

10ml, bottle 23-11518



S-Bond (Danville)

Highly stable, **prehydrolyzed silane**. Provides full strength bond link to ceramics. Use on etched porcelain, and composite posts

10ml, bottle 23-88431



Scotchbond Multi Purpose (3M)

Dental Adhesive System
A light-cure system

Excellent bonds to a variety of surfaces
Fast & easy to apply

Adhesive, 8ml 23-108021

Primer, 8ml 23-108022



Touch & Bond (Parkell)

7th generation bonding
No-etch bonding agent
Requires a halogen curing Light
Creates strong bonds to dentin & prepared enamel. Quick & easy; only take **1 minute**. Can be used w/ light-cure or self-cure composites



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Kit: One bottle of touch & bond liquid (5.5ml), 175 activator pledgets, precision applicator microbrushes and organizer 23-110371

Transbond XT LC Brackets Adhesive (3M)

Light cure adhesive **bonds metal and ceramic brackets** to tooth surfaces.



The viscosity of Transbond XT adhesive was designed to prevent adhesive run-on and bracket skating that saves you money and reduces adhesive waste. 3M's adhesive filled capsules and dispensing gun also provide for easy application, increased control and convenience.

4g Syringe, only 44-108192

6ml Transbond XT Primer, only 44-108193

Kit: 2 syringes (4g each), 6ml transbond XT primer, 1 brush holder, 60 brush inserts 44-108191

Brush Holder & 60 Brush Inserts 44-108195

Transbond XT Primer (3M)

Light cure primer, bond metal and ceramics brackets to tooth surfaces

6ml, bottle 23-110531

Value Pack: 3 x 6ml bottle 23-110532



Bulk Fill Composites

BulkEZ (Danville)

Dual cure dental restorative composite designed for direct restorations.

When allowed to self cure, BulkEZ features excellent physical properties for posterior restorations and offers an unlimited depth of cure. Light curing is only needed to harden the air inhibited top layer after self-curing. With no layering required, time savings over light cured bulk composites are over 60%. Unlike other dual cure composites, BulkEZ is compatible with all bonding agents. The directed and "right" shrinkage of BulkEZ offers the best interproximal margins. No voids or gaps are formed at the bottom of the restoration as a result of curing, limiting the risk for post-operative sensitivity or recurrent decay, which is often seen at the interproximal gingival margin.

A1 Refill, 2/pkg 23-108531

A2 Refill, 2/pkg 23-108532

A3 Refill, 2/pkg 23-108533



Cosmetic & Restorative

Trial Kit Contains: 3 x 6-gram syringes (A1, A2, A3), 20 tips, 5 ZNano unidose A2 23-108534

Tips Refill 17 ga, 25/pkg 23-108535

Tips Refill 19 ga, 25/pkg 23-108536

Estelite Bulk Fill Flow Syringe (Tokuyama)

Resin-based dental restorative material; Estelite Bulk Fill Flow is a low viscosity, light-cured, radiopaque composite resin. This low stress flowable material can be placed 4mm increments. Estelite Bulk Fill Flow incorporates radical amplified photopolymerization initiator technology. RAP technology facilitates a shortened light curing time and ample working time



Syringe A1, 3gm 23-110861

Syringe A2, 3gm 23-110862

Syringe A3, 3gm 23-110863

Syringe B1, 3gm 23-110864

Syringe Universal, 3gm 23-110865

Composite Finishing

Luminescence Plus (Premier)

The first diamond polishing paste to contain potassium nitrate as a topical desensitizer. A single step, single-grip paste, offering superior intra-oral polishing on all restorative materials



Gel, 3gm 23-108131

Spirals Diamond Polishers

Flexible, diamond impregnated polishers for finishing and polishing composites, porcelain, glass ionomers, and ceramic materials designed to reach all tooth surfaces, including difficult to reach areas.



Multi-use, long-lasting, and autoclavable

RA-Type

All you need to quickly and efficiently produce the optimal polish, on any surface and on almost any material.

Coarse, 3/pk 32-110571

Medium, 3/pk 32-110572

Fine, 3/pk 32-110573

Kit, 6/pcs 32-110574

Call Us Today 833-733-6825



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Super Snap Discs (Shofu)

50/Box



Contouring Black - Coarse, #L506 32-109271

Finishing Violet/Medium, #L509 32-109272

Polishing Green/Fine, #L501 32-109273

Contouring Red/Superfine, #L502 32-109274

Contouring Black - Mini Disc, safe side for composite, #L507 32-109275

Finishing Dark Violet/Medium Regular Disc, safe side, #L508 32-109276

Polishing Green/Fine Mini Disc, safe side, #L519 32-109277

Polishing Red/Super Fine Mini Disc, safe side, #L521 32-109278

Polishing Dark Violet/Medium Mini Disc, double sided, #L522 32-109279

Finishing Dark Violet/Medium, double sided disc, #L528 32-109280

Super Snap Mandrels (Shofu)

The improved design of the mandrel for Super-Snap Disks and Buff Disks creates a secure seating of the disk onto the mandrel



- Stainless steel available in HP, CA, & FG Shanks

- Durable single-use plastic mandrels in CA shanks only

Shank CA, 6/pkg 32-111371

Flowable Composites

Revolution Formula 2 (Kerr)



A1 Refill Kit: 4-1gm syringes + 20 tips 23-11520

A2 Refill Kit: 4-1gm syringes + 20 tips 23-11522

A3 Refill Kit: 4-1gm syringes + 20 tips 23-11521



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Cosmetic & Restorative

B2 Refill Kit: 4-1gm syringes + 20 tips 23-11519

A3.5 Refill Kit: 4-1gm syringes + 20 tips 23-11523

StarFlow (Danville)

StarFlow's superior flow characteristics and esthetics make ideal results easy. Sub-micron particle size gives superior polish ability and resistance to staining. A leading clinical reviewer showed StarFlow to have the highest compressive strength of any flowable composite. True Vita shade match, radiopaque and non-sticky. It offers fluoride slow release and natural fluorescence.



Shade A1, 5gm syringe 23-11535

Shade A2, 5gm syringe 23-11536

Shade A3, 5gm syringe 23-11534

Shade A3.5, 5gm syringe 23-11533

Shade B1, 5gm syringe 23-11530

Intro Kit: (3-5gm) A2, A3.5, incisal 23-11532

Trial Kit: (4-1gm) A1, A2, A3.5, incisal 23-11531

Accolade SRO (Danville)

Thixotropic Flowable Composite. With a radiopacity of 320% of aluminum, (similar to titanium) radiographs were never easier to read! Recently confirmed by a leading clinical reviewer, Accolade SRO is among the most radiopaque flowable on the market, yet is value priced. As a 0.7 micron hybrid microfill, Accolade SRO also offers exceptional strength, esthetics and the polish ability of StarFlow and Accolade.



A1, 5gm 23-892412

A5, 5gm 23-892413

Extra Light, 5gm 23-892414

Aria Microfill Flowable Composite (Danville)

Aria's microfill formulation offers polish ability and flexibility no hybrid can match! Ideal for facials and class V. High translucency gives Aria "chameleon" qualities, necessitating only a limited number of shades. Aria seems to disappear into the surface making it an ideal companion to Accolade and StarFlow when less opacity is desired.



Shade A-2, 5gm syringe 23-11114

Order Now ➡➡

Denfil Flowable Composite (Vericom)

Light-Cured Radiopaque Flowable Composite Resin, has a low shrinkage and good mechanical properties

- Denfil is designed to provide optimum flow ability and handling can restore narrow deep cavity as well as cervical area

- Denfil containing micro-hybrid fillers provides excellent color adaptability, nature-like esthetics and high-gloss polish ability

- Vericom has supplied two other distinct grades besides universal viscosity type for the limited clinical application

Indications

- Restorations of class III, Class V, smaller Class IV

- Base/liner in the Class I, Class II restorations

- Repair of resin, porcelain, and acrylic temporary materials

- Pit and fissure sealant

- Undercut blackout

A1 Sys Kit: 4 x 2gm, each 23-21218

A2 Sys Kit: 4 x 2gm, each 23-21219

A3 Sys Kit: 4 x 2gm, each 23-21220

A3.5 Sys Kit: 4 x 2gm, each 23-21221

B1 Sys Kit: 4 x 2gm, each 23-21222

B2 Sys Kit: 4 x 2gm, each 23-21223



Flowable Composite (Made in USA)

Highly polish able, 67% filled, esthetic composite used for all types of cavity preparations. The composite has a nominal particle size of 0.7 micron. Ideal choice for anterior and posterior applications as well as use as a pit and fissure

sealant. Also ideal for repairing small defects and filling in cervical gaps caused by toothbrush erosion and receding gums. Also used to eliminate undercuts prior to taking impressions for composite inlays and to line cavities and fill in difficult access areas prior to placing condensable composites. Easy to polish to high luster. Has high compressive strength as well as low shrinkage and resist wear and staining. Radiopaque for easy identification in radiographs

Benefits

- The micro-hybrid formula provides strength and polish ability

- Fast, deep light activation in only 30sec. for lighter shades & 40sec. for darker shades

- Fluoride-releasing

- Low shrinkage

- High compressive strength 290 MPa

- Weight/volume: 67%/54%

- Diametral tensile strength: 47 MPa

- Flexural strength: 125 MPa

A1 Refill, 2g syringe 23-11345

A2 Refill, 2g syringe 23-11344

A3 Refill, 2g syringe 23-11343



Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

A3.5 Refill, 2g syringe	23-11342
B1 Refill, 2g syringe	23-11341
B2 Refill, 2g syringe	23-11340

G-Aenial Universal Flo (GC America)

- Flowable composite
- Handles like a low-flow flowable and performs like a restorative
- Highly thixotropic
- Recommended for Class I, II, III, IV and V restorations
- Bis-GMA free



A2, 2ml syringe	23-109811
------------------------	-----------

N'Durance Dimer Flow (Septodont)

This light-cured, low-shrinkage, flowable restorative material is indicated for use in both anterior and posterior restorations

- Non-sticky and easy to handle
- Low viscosity and non-slumping to provide excellent handling characteristics
- Outstanding wear resistance and polish ability
- Extremely hydrophobic compatible with all bonding agents
- Highly radiopaque



A1, Refill, 2 x 2g syringes	23-105821
A2, Refill, 2 x 2g syringes	23-105822
A3, Refill, 2 x 2g syringes	23-105823
A3.5, Refill, 2 x 2g syringes	23-105824
B1, Refill, 2 x 2g syringes	23-105825
C1, Refill, 2 x 2g syringes	23-105826
Bleach White, Refill, 2 x 2g syringes	23-105828

Nexcomp Flow (Meta Biomed)

Light curing time
20s ISO4049-Resin
based dental
restorative
material.



Indications

- restoration of minimally invasive cavity (Class V & Class III)
- Pits and fissure sealant
- Base and liner for direct restoration
- Repair acrylic temporary materials
- Repair of esthetic indirect restorations

Key Features

- Excellent biocompatibility with dentin and pulp
- Light curing typer cavity liner, cure extremely hard
- Chemically bonds to adhesives and composites
- Contains calciumhydroxyphosphate
- Releases calcium, hydroxyl, phosphate and fluoride ions
- Fluoride releasing
- Insoluble in water and oral fluids
- Radiopacity

- Flow-on demand handling and easy dispensing

A1, 2g x 2 syringes, 10 disposable tips	23-11417
A2, 2g x 2 syringes, 10 disposable tips	23-11418
A3, 2g x 2 syringes, 10 disposable tips	23-11419
B1, 2g x 2 syringes, 10 disposable tips	23-11422

OptiFlow II (Pacdent)

Fluoride-releasing, resin based, low viscosity flowable light-cure composites.



Natural translucency, high radiopacity and great controlled flowability

Also available with high viscosity

A2 4-Pack: 4 x 1.5gm syringes & 20 blue flo-tips	23-11446
A3.5 4-Pack: 4 x 1.5gm syringes & 20 blue flo-tips	23-11459
B1 4-Pack: 4 x 1.5gm syringes & 20 blue flo-tips	23-11457
A1 4-Pack: 4 x 1.5gm syringes & 20 blue flo-tips	23-11447
A3 4-Pack: 4 x 1.5gm syringes & 20 blue flo-tips	23-11460

Palfique Estelite LV Flow Syringe (Tokuyama)

An esthetic flowable resin composite, which offers 3 different viscosities; high flow, medium flow and low flow, to suit different clinical cases and the preferences of the clinicians. Just as other Estelite brand of esthetic resin composites, every inorganic filler contained in Palfique Estelite LV is a submicron spherical filler, providing the same benefits as Estelite Sigma, such as excellent esthetics, great self-shade matching, superior polish ability and gloss retention



Low Flow Syringe (A1), 3gm + tips	23-11462
Low Flow Syringe (A2), 3gm + tips	23-11469
Low Flow Syringe (A3), 3gm + tips	23-11470
High Flow Syringe (A1), 3gm + tips	23-11471
High Flow Syringe (A2), 3gm + tips	23-11472
High Flow Syringe (A3), 3gm + tips	23-11467
Medium Flow Syringe (A1), 3gm + tips	23-11466
Medium Flow Syringe (A2), 3gm + tips	23-11465
Medium Flow Syringe (A3), 3gm + tips	23-11464

Cosmetic & Restorative

Medium Flow Syringe (A3.5), 3gm + tips 23-11463

Black Syringe Tips, 50/box 23-11468

Vertise Flow (Kerr)

Light cured self-adhering, resin based flowable composite dental restorative designed for direct placement



A2 Refill, 2 x 2g syringes, 20 bendable brushes and dispensing tips 23-109282

A1 Refill, 2 x 2g syringes, 20 bendable brushes and dispensing tips 23-109281

Glass Ionomer Restoratives

Fuji IX GP Capsule (GC America)

- Posterior glass ionomer restorative
- Ideal alternative to amalgam and composite
- For geriatric & pediatric applications
- Uses as an IRM, for long-term temporary restorations, as a core build-up material, or restorative base



A2 Regular, 50 capsules/box 23-98431

Fuji IX GP Fast (GC America)

- Extremely hard
- Very wear resistant



A3, 50 cap/box 23-105121

GI Restorative Radiopaque Glass Ionomer Restorative Cement (Prime Dental)

- A self-cured glass ionomer restorative
- Its high surface hardness provides for a durable restoration
- Radiopaques. Fluoride release
- Superior tear & bond strengths



Powder and liquid, 18g powder bottle, 17.5ml liquid bottle and accessories 23-11357




Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Condensable Composite

Filtek P60 Condensable - Packable Composite (3M)

- Fastest placement technique, 20 seconds to cure an increment of 2.5mm
- Posterior restorative
- Excellent packability, resists slumping
- Ideal flow and viscosity for all Class I & II restorations
- Doesn't stick



A3, 4gm syringe refill 23-105491

Etching Agents

EtchVance Gel Jumbo Size 70g Syringe

37% phosphoric acid etching gel used for etching enamel and dentin prior to bonding procedures



Stays where placed and will not run

Does not dry out and has long shelf life when properly stored

70g Etching, blue gel 23-11322

Best Etch 38% (Vista)

38% phosphoric solution for use on dentin and enamel.

Material will stay and not flow onto surround tissue.

Leaves no residue, rinses away clean



Bulk Etch Refill, 1-30ml syringe & 1 syringe docking port 23-11128

Bulk Etch Kit: 1-30ml filled syringe, 1 syringe docking port, 20 stat-flo tips, 20 of 1.2cc empty syringes 23-11131

Standard Kit: 4 of 1.2cc pre-filled syringes w/25 of 22ga pre-bent applicator tips 23-11130

Value Pack Kit: 12-3cc pre-filled syringes w/75 of 22ga pre-bent applicator tips 23-11132

Denfil Etching Gels (Vericom)

The strong point of this etchant is **thickened with natural polymer materials** so that it does **not leave any white residues on the surface** of etched tooth like silica thickened etchant; the white residues **reduce the bonding strength**. The numbers mean the concentration of phosphoric acid



Cosmetic & Restorative

12 sys x 3ml, each w/ 72 pre-bent tips 23-11208

Kit: 4 sys x 1.2ml each w/ 72 pre-bent tips 23-11209

Etch-Rite Etching Gel (Pulpdent)

Etch-rite is a superior quality dental etching gel specially formulated for maximum working and handling characteristics

- Thixotropic
- Stays where placed - will not run onto exposed tissue
- Washes off easily without leaving any residue
- 15-second etch
- Syringe dispensing
- Precise placement through small-gauge pre-bent applicator tips

Kit: 4-1.2 ml syringes gel & 8 bent needles 23-11316

Jumbo Refill: 2 x 25 ml syringes 23-11318

Etching Gel Bulk Pack: 24 x 1.2 ml syringes 23-11317



Etchening Gel (Prime Dental)

Prime-Dental Etchant Gel/Liquid with 37% phosphoric acid is used for etching enamel and dentin prior to bonding procedures

Prime-Dental Etchant Gel is thixotropic, stays where placed and will not run

Does not dry out and has long shelf life when properly stored

Etching Gel, 50ml syringe 37% blue 23-11320

Etching Gel, 4 x 1.2ml syringe & 8 dispensing tips, 37% blue 23-11319

Etching Gel, 20 x 1.2ml syringe & 40 dispensing tips, 37% blue 23-11321



Porcelain Etch (Vista)



Kit: 2-1.2ml pre-filled etch syringes, 2-1.2ml pre-filled silane syringes with 25-25ga micro flowthru brush tips 23-11487

Refill, 4-1.2ml pre-filled silane syringes 23-11488

Porcelain Etch Gel (Pulpdent)

A one-minute application etches porcelain crowns, veneers, inlays one fractured porcelain surfaces prior to bonding. This superior quality porcelain etch gel will not stain porcelain, composites or resin cement.



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Kit Contains: 4-1.2ml syringes etch gel; 8 dispenser tips 23-11490

Porcelain Prep Kit (Pulpdent)

Contains Porcelain Etch Gel, Silane, Kool-Dam and Dry-Rite drying agent

Protect soft tissue with Kool-Dam light cure liquid dam. Apply pulpdent porcelain etch gel to ceramic surface for one minute. The one-minute application will not stain the ceramic. For best results, dry the ceramic surface completely with Dry-Rite before applying silane.

Kit Contains: 1x1.2ml syringes porcelain etch gel, 1x1.2ml syringes kool-dam, 1x1.2ml syringe dry-rite, 1x1.2ml syringe silane, 4 dropper tips & 8 applicator tips 23-11492



Silane Bond Enhancer (Pulpdent)

For bonding composites and resin cements to porcelain

A single component material that increases the bond strength of organic resin, such as composites and cement resin to porcelain.

Apply silane to porcelain after etching with porcelain Etch Gel and drying dry-rite

Conveniently dispenses one drop at a time from the syringe

Kit Contains: 4 x 1.2ml syringes + 8 dropper tips 23-11527



Composite Sealant

Eco-S Pit & Fissure Sealant (Vericom)

Light-cured dental material designed for sealing the enamel pit and fissure of teeth for the prevention of caries.

Has good properties like shear bond strength, micro-leakage and wear resistance.

Has an appropriate viscosity, to enable the sealant to penetrate the depths of pits

Provides good efficacy as well as work ability

4 sys x 1.2ml sealant with 20 tips 46-92762



Embrace Wetbond Pit & Fissure Sealant (Pulpdent)

Key Features

- Contains no Bisphenol A, No Bis-GMA, No Bis-DMA
- Wet Bonding - Moisture tolerant resin chemistry
- Tooth Integrating - Seals against microleakage
- Margin-Free - No chipping, pitting, or staining
- No Bonding Agents Required - Fewer steps - Saves



Cosmetic & Restorative

time & money

The BEST sealant EVER, bonds to the moist tooth

Kit: 4 x 1.2ml syringes, natural shade + 46-11114
20 applicator tips

Bulk Pack: 20 x 1.2ml syringes, natural shade + 100 applicator tips 46-11116

Kit: 4 x 1.2ml syringes, off-white shade + 20 applicator tips 46-11117

Seal-N-Shine (Pulpdent)

Finishes and polishes restorations - seals composites and margins - save times

- Tough, resilient, long-lasting resin provides smooth, durable, protective finish and seals against microleakage
- Eliminates final finishing and polishing steps
- Does not alter occlusion or anatomy
- Bonds to the moist tooth
- Cures with all lights
- Heals white line - Prevents marginal staining
- Cures clear. No yellow tint
- Available in bottles or syringes



Bottle Kit: 6ml bottle clear shade, brush handle, 100 brush tips 23-11525

Syringe Kit: 2-1.2ml syringes, 40 flocked tips 23-11526

Seal-Rite Pit & Fissure Sealant (Pulpdent)

Seal-Rite is 34.4% filled and flows nicely from cusp to cusp. Its pearlescent color is easily seen in the mouth. **Seal-Rite Low Viscosity** is 7.7% filled for those practitioners who prefer a more flowing sealant.

Off-white color contrasts with tooth enamel.

- Fluoride release - fluoride uptake by adjacent enamel
- Visible light cure
- Syringe dispensing



Kit: 4 x 1.2ml syringes seal-rite + 8 applicator tips 46-11136

Procedure Kit: 1-2ml syringe each: 2 x seal-rite, etch-rite, dry-rite + 12 applicator tips 46-11138

Composite Applicators

Pre-Bent Needle Tips

Disposable pre-bent dispensing tips for use with etchants, resins and flowable composites



18G - Pink

20G - Black

20G - Yellow

22G - Gray

25G - Blue



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

18ga., Pink, 100/pkg 23-99181

20ga., Black, 100/pkg 23-99182

20ga., Yellow, 100/pkg 23-99183

22ga., Grey, 100/pkg 23-99184

25ga., Blue, 100/pkg 23-99185

Bendable Irrigation Needles (DZ)

- 1" multi-purpose bendable irrigation needle tip
- Excellent flexibility for easy access to any canal
- Use for root canal irrigation and periodontal pockets
- Universal luer-lock design
- Color coded
- Non-sterile



18G, Pink, 1.2 x 25mm, 100/pkg 29-11175

19G, Cream, 1.1 x 25mm, 100/pkg 29-11176

20G, Yellow, 0.9 x 25mm, 100/pkg 29-11177

22G, Black, 0.7 x 25mm, 100/pkg 29-11178

25G, Orange, 0.5 x 25mm, 100/pkg 29-11179

27G, Medium Gray, 0.4 x 25mm, 100/pkg 29-11180

Bendable Needle Tip (Vista)

- Endo/perio
- Blunt-end bendable dead-soft needles resist crimping
- Designed for use in evacuation and irrigation
- Use for dispensing of flowable composites, etchants, resins, and dentin primers
- Extended 3/4" length to reach pockets, fistulas and surgical sites
- 6 gauge sizes
- Universal luer-style design
- Color coated for easy identification



18 Gauge (Pink), 100/pkg 23-11118

25 Gauge (Purple), 100/pkg 23-11127

18 Gauge (Pink), 100/pkg 23-11126

19 Gauge (Brown), 100/pkg 23-11125

19 Gauge (Brown), 100/pkg 23-11124

20 Gauge (Yellow), 100/pkg 23-11123

20 Gauge (Yellow), 100/pkg 23-11122

22 Gauge (Black), 100/pkg 23-11121

22 Gauge (Black), 100/pkg 23-11120

25 Gauge (Purple), 100/pkg 23-11119

Dental Supplies

Cosmetic & Restorative

Brush Applicators



Color Coded



Regular Size, Blue Color, 100/tube, 4 tubes/box	23-105881
Fine Size, Green Color, 100/tube, 4 tubes/box	23-105882
Ultrafine Size, Purple Color, 100/tube, 4 tubes/box	23-105883
Cylinder Size, White Color, 100/tube, 4 tubes/box	23-105884

Composite Brush Handle (Plasdent)



Blue, each	23-11569
White, each	23-11568
Green, each	23-11567
Red, each	23-11566
Yellow, each	23-11565
Brush Tips, 100pcs/bag	23-11564
Micro Tips, 100pcs/bag	23-11562
Kit: 50 brush tips, 50 micro tips & 5 handles	23-11563

Composite Brush Tips

Tapered DuPont nylon bristles ensure even and smooth application of liquid, etchants, resins and other monomers. **Fits all** leading brands of brush holders snugly. Regular black bristle



Black Brush Tips, Regular, 100/pkg	23-11178
White Brush Tips, Fine, 100/pkg	23-11177
Brush Tips Kit: 500 black regular brush & 2 handles	23-11176



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Brush Tips Kit: 500 white fine brush & 2 handles 23-11179

Composite Brush Tips & Handle (3M)



Brush Tip, 60/pkg	23-111191
Applicator Handle, each	23-111192
Standard Kit: 2 handles & 120 tips	23-111193

Composite Brushes & Handles (Vista)



Brushes Fine White, 100/pkg	23-92551
Brushes Regular Black, 100/pkg	23-92552
Handle, 2/pkg	23-92553

Delivery Tips (Pac-Dent)



Pre-bent Needle Tip w/ Brush-padded end, 100 tips/pkg	23-25611
Pre-bent Needle Tip w/ Micro Brush, 100 tips/pkg	23-25612
Pre-bent Needle Tip w/ Yellow I/O Connection, Green, 100 tips/pkg	23-25613
Pre-bent Needle Tip w/ Teal I/O Connection, Green, 100 tips/pkg	23-25614
I/O Tips w/ Cap, 100 tips/pkg	23-25615
I/O Tips w/ Snap-Off Cap, 100 tips/pkg	23-25616
Curved Black w/ Brush Tips, 100 tips/pkg	23-25617
Curved Black w/ No Brush, 100 tips/pkg	23-25618

www.prestigedentalproducts.com

Cosmetic & Restorative

Disposable Etch Tips

- Universal luer-lock design
- Curved needle tip for use of placing viscous to thin materials
- Color code



Blue Color, 100/pkg

23-11313

Disposable Flow Tips

- Universal luer-lock design
- Curved needle tip for use of placing viscous to thin materials
- Color code



Pre-Bent, Black Color, 100/pkg

23-11240

Flowthru Dispensing Tips (Microbrush)

Disposable plastic applicators are designed for use in place of standard luer lock needle tips. The plastic tip can bend up to 90° without pinching and greatly reduces needle phobia. A threaded locking system ensures that the applicators are securely fastened to the syringe when dispensing material. Applicators have a protective ring along the base to prevent waste material from reaching the threads, reducing the risk of cross contamination



18ga, 100/pkg

23-11352

22ga, 100/pkg

23-11354

25ga, 100/pkg

23-11355

Flowthru Microbrush-Style w/ Foam Tips (Microbrush)

Same feature as the standard-style with the addition of non-linting, non-absorbent fibers that form a perfect sphere of solution (as small as 1/8 of a drop) without spilling or dripping



18ga, 100/pkg

23-11348

20ga, 100/pkg

23-11349

22ga, 100/pkg

23-11350

25ga, 100/pkg

23-11351

Maxbrush Applicators (Plasdent)

Assorted, 4" Long, 100/pkg

23-98081

Assorted, 4" Long, 400/pkg

23-98082



Maxext Micro Applicators - Extended

Black, 100/pkg

23-98091


Black, 400/pkg

23-98092



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Maxmicro Micro Applicators (Plasdent)



Regular Size 2.0mm, 100 purple, 100 peach, 100 blue & 100 green color, 400/pkg

23-97981

Fine Size 1.5mm, 200 pink & 200 yellow color, 400/pkg

23-97982

Super Fine Size 1.0mm, white color, 400/pkg

23-97983

Ultra Fine Size 0.50mm, black color, 400/pkg

23-97984

Super Fine Size 1.0mm, white color, 100/pkg

23-97985

Ultra Fine Size 0.50mm, black color, 100/pkg

23-97986

Micro Applicators (Plasdent)



Blue, 144/box

23-108231

Yellow, 144/box

23-108232

Green, 144/box

23-108233

Red, 144/box

23-108234

Purple, 144/box

23-108235

Black, 144/box

23-108236

MicroBrush (MicroBrush)

Disposable applicators the working end of the microbrush is comprised of non-linting, non absorbent fibers arranged in a special spherical shape, designed to hold (in suspension) solutions as small as 1/8 of a drop



This applicator eliminates dripping, spilling, and waste. Can be bent, and will hold to any shape.

Regular, Assortment Package, 400/pkg

23-11378

Fine, 400/pkg

23-12379

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Super Fine, 400/pkg 23-12377

MicroBrush Plus (MicroBrush)



In addition to original line of Microbrushes. Dispenser series utilizing a patented bending portion, which makes bending these applicators easier than ever before.

W/ Dispenser Regular, 400/pkg 23-11384

W/ Dispenser Fine, 400/pkg 23-11387

Regular Refill (400) 23-11381

Fine Refill (400) 23-11382

Super Fine Refill (400) 23-11383

W/ Dispenser White S. Fine, 400/pkg 23-11386

MicroBrush X (MicroBrush)

A unique new extended and extra-thin applicator. These patented applicators are designed for an even distribution of material to inner surfaces for application of material in confined spaces. Applicators bend firmly up to 90 degrees. Extra short fibers allow for scrubbing



Kit: Dispenser, w/ 100 applicators 23-11379

Refill (100) 23-11380

MicroStix Adhesive Tip Applicator (MicroBrush)



Original, Blue, 64/pkg 23-98422

Mini Micro Applicators

Regular, 400/pkg 23-98791

Fine, 400/pkg 23-98792

Super Fine, 400/pkg 23-98793

Prebent Irrigation Needles (Dentazon)

- Universal luer-lock design
- Curved needle tip for use of placing viscous to thin materials
- Color coded
- Non-sterile



18G: Pink, 1.2 x 25mm, 100/pkg 23-11506

20G: Yellow, 1.1 x 25mm, 100/pkg 23-11505

20G: Black (for curing light composite resin), 0.9 x 25mm, 100/pkg 23-11504

22G: Dark Gray, 0.7 x 25mm, 100/pkg 23-11503

25G: Blue, 0.5 x 25mm, 100/pkg 23-11507

Spira-Flo Brush Tips (Vista)

Features

- Curved cannula for optimal entry
- Internal spiral channel
- Adjustable brush tip fiber length
- Nylon brush fiber for greater controlled delivery



100/pkg 23-26231

500/pkg 23-26232

Stat-Flo Syringe Applicator Brush Tips (Vista)

Nylon bristle end allows for effective brushing of solutions.



19ga., Yellow, 100/pkg 23-94241

Ultrabrush 1.0 (MicroBrush)

These applicators combines the precision of a smaller bristle brush, satisfying the needs of the Ultrabrush Tube Series and Ultrabrush 2.0 applicator design



Kit: Blue (100) 23-11554

Refill: Orange/blue (200) 23-11553

Ultrabrush 2.0 (MicroBrush)

Applicator features chemical resistant fibers to hold up to ever the harshest of materials



Kit: (100) 23-11555

Refill: (200) 23-11556

No Stick Lube

No Stick Lube (Danville)

Glide over composite without sticking and without reducing bond or composite strength



10ml, bottle 23-11431

Cosmetic & Restorative

Shade Guides

Rite-Lite LED Shade Matching Device (AdDent)

L.E.D. shade matching light designed to eliminate costly remakes, the Rite-Lite shade matching light incorporates LED technology to provide simulated north sky daylight with a standardized color temperature of 5500K. The small, handheld, portable unit comes with a neutral color tabs for placement behind the teeth to eliminate environmental color distractions.



Rite-Lite Device 23-11524

Sunlight LED Illuminated Magnifier Lighting Unit (Meta)

It illuminates and magnifies with the most natural color used for measuring accurate shading level for porcelain works.

Features

1. Most natural LED lighting source
2. Surface emitting light; illuminates equally that prevent shadows.
3. Detachable 2x & 3x magnifier with preventable shape distortion.
4. Three step dimmer makes it easy to control the brightness
5. Use USB for easy charging and usage
6. Most comfortable observation rate for both eyes; 16:9 rate rectangle.



Kit: Light unit, 2x magnifier lens, 3x magnifier lens, USB charging cable and manual 23-98811

Vita Shade Guide (VITA)



One universal shade determination system for all fields of dental prosthesis includes: A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4

Kit: 16 shades 23-88201

Prestige Dental Products



**Excellent Qualities
at Affordable Prices**

877-772-3888

www.prestigedentalproducts.com



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Composite Caddies

Ultraviolet Covers (Plasdent)

For protection of light-curing materials from ultraviolet light



(7 1/2" L x 4 1/2" W x 1 1/2" H) 23-11559

(3 1/4" L x 3" W x 1" H) 23-11558

Composite Dispensing Gun

Cavifil Composite Gun Dispenser (Ivoclar Vivadent)



Composite Unidos gun for application of dental restorative materials

Made in Germany

Each 23-110761

Comax Composite Dispensing Gun (AdDent)

- Composite dispenser for all composite compules ergonomic
- Handle and trigger high mechanical advantage mechanism.
- Easy one hand unloading high-tech plastic & stainless steel materials
- Completely autoclavable

2 year warranty

Each 23-11175



Composite Dispensing Gun

- Unique dispenser for the delivery of all composite compules
- Ergonomically designed dispensing gun provides greater leverage for dispensing both high and low viscosity materials
- Compule dispenser gun
- Multi-use compatible with most unidose tips

Each 23-11190



**Order
Online!**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Cosmetic & Restorative

Composite Dispensing Gun (DA)

For use with all Unidose compules



Each

23-91961

Composite Dispensing Gun (Pac-Dent)

- Ergonomically designed dispensing gun
- Compatible with most unidose tips
- ideal for all viscosity materials
- Slot is designed for fast loading



Metal tip

1 Plastic Composite, dispensing gun 23-11180

1 Composite, dispensing gun w/ metal tip 23-11181

Composite Gun (Plasdent)

- Ergonomic design
- Accept most carpules
- Autoclavable
- Fast loading
- For all viscosity material



Blue Ocean color 23-88571

Green Grass color 23-88572

Orange Yellow color 23-88573

Light Gray color 23-88574

Purple Gala color 23-88575

Black Knight color 23-88576

Composite Warmer

Calset Composite Warmer (AdDent)



110007-02



110007

Warming composites couldn't be easier with AdDents new Calset Tri Tray. Not only will it warm your compule in a dispenser gun, it warms 2 placement instruments that will enable you to sculpt the composite easier and with more precision. It also holds, and warms 3 extra compules. Using the AdDent CoMax dispenser, the compule can be picked up directly from the tray without even touching the compule, however other composite dispensers can be used for warming process.



110007-12



110040

Features

- Heats composite material to 98°F
- Maintains constant temperature
- Increase degree of cure
- Warms finishing instruments

Benefits

- Shortens curing time by over 80%
- Improves composite flow by 68%
- Reduces microleakage
- Improves physical properties
- Sculpting and getting into the marginal ridges made easier

Advantages

- Calset Tri-Tray warms composite
- Warms 2 finishing instruments, making it easier to sculpt the composite
- Holds 3 spare compules
- Warmed composite provide quicker cures
- Less microleakage
- Stronger restoration

Calset composite heater, 120 volts AC power supply, with standard composite compule tray 23-11144

Standard composite compule tray 23-11143

Anesthetic tray (holds 3 carpules) 23-11141

7 hole syringed composite tray 23-11140

Porcelain veneer tray 23-11139

Calset Tri-Tray composite heater, 120 volts AC power supply, for use with CoMax, Caulk, Kerr (original model), SDI and compo jet. also holds 2 finishing instruments 23-11138

Calset composite heater, 120 volts, AC power supply, with dispenser tray for use with centrix, 3M, venus, kerr (newer model) guns 23-11137

Calset heater for syringed composites, 120 volts AC power supply, with 7 hole syringed composite tray 23-11136

Follow Us & Like Us



Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Composite Accessories

Bondwe-II Mixing Wells



2-Well, White, 200pcs/box	23-16601
4-Well, White, 200pcs/box	23-16602

Bondwell Mixing Wells (Plasdent)



4-Well, Black, 200pcs/box	23-11389
4-Well, Asst in White, Blue, Yellow and Red, 200pcs/box	23-11390
4-Well White, 200pcs/box	23-11391
4-Well White, 500pcs/box	23-11392
2-Well, Asst in White, Blue, Yellow and Red, 200pcs/box	23-11393
2-Well White, 200pcs/box	23-11394
2-Well White, 500pcs/box	23-11395

Brush Applicators



Color Coded



Regular Size, Blue Color, 100/tube, 4 tubes/box	23-105881
Fine Size, Green Color, 100/tube, 4 tubes/box	23-105882
Ultrafine Size, Purple Color, 100/tube, 4 tubes/box	23-105883
Cylinder Size, White Color, 100/tube, 4 tubes/box	23-105884

Composite Brushes & Handles (Vista)



Brushes Fine White, 100/pkg	23-92551
Brushes Regular Black, 100/pkg	23-92552
Handle, 2/pkg	23-92553


Composite Syringe Organizer (Plasdent)



Large, (8" W x 5 1/8" H x 4 1/4" D)	26-11240
Small, (4" W x 5 1/8" H x 4 1/4" D)	26-11247

Cosmetic Aligner (IndiGreen)

After tooth preparation, the anterior midline and the horizontal plane are often difficult to determine, yet they are critical to communicate to the lab for accurate fabrications of esthetic restorations.



The Cosmetic Aligner is easier, faster and more economical than sophisticated articulation equipment. Together with any quick-set bite material, the compact one-piece device registers the anterior midline between central incisors and the horizontal plane. By only recoding centrals, it reduces material waste and cumbersome, inaccurate cross-stick bite impression methods

50/box	23-87971
--------	----------

Mixing Sticks (Dentsply)

Mini Mixing Sticks: Small design makes it easy to handle and it allows for more control when applying adhesive

Standard Mixing Sticks: Perfect for mixing banding and bonding cements. Large size facilitates easy control



Mini Blue Sticks, 2.2" long, 100/pkg	45-12965
Standard White, 4.4" long, 100/pkg	45-12964
Standard White, 4.4" long, 1,000/pkg	45-12963

Cosmetic & Restorative



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Spatulas (Plasdent)



Mixtix Mixing, White Color, 100/pkg	21-11210
Nylon Spatula Firm, neutral each	21-11230
Alginate Spatula Soft, white each	21-11251

Therma-Flo Step Down Tips (Vista)

Featuring a self-tapping design that will adapt to most standard compules our patent pending step down tips combine the ultimate in convenience with precise placement of composite materials

Features

- Narrow profile for optimal visualization of work area
- Extended length for deeper access
- Narrow profile allows for optimal visualization
- Designed for easy extrusion of highly-filled composites
- Available in 2X and 4X reduction for precision placement
- Orifice reduction for greater precision during micro-restorative procedures



Developed for use with the Therma-Flo composite warming kit, but can be used with any standard compules or compule guns

Step Down Tips 2X, 40 cup/pkg 29-26211

Step Down Tips 4X, 40 cup/pkg 29-26212

Transbond XT Primer (3M)

Light cure primer, bond metal and ceramics brackets to tooth surfaces



6ml, bottle 23-110531

Value Pack: 3 x 6ml bottle 23-110532

Universal Wells (Plasdent)



Blue Ocean, each	23-98761
Green Grass, each	23-98762
Orange Yellow, each	23-98763
Light Gray, each	23-98764
Purple Gala, each	23-98765
Black Knight, each	23-98766

Tubes and Plugs

Applicator Tips with Needle



100 Tips + 100 Plugs 23-107801

ProDose Disposable Tubes & Plugs (Pac-Dent)

- Available in low viscosity, 20ga needle tube, and 19ga needle tube
- Can be used with syringes on the market
- Disposable so they eliminate cross-contamination and cleanup
- Easy to load
- Dispense cement or composite
- Restore with minimal voids and porosities
- Photo-blocking brown tint blocks light, extending the life of light-activated materials
- Material mixture accessible, eliminating pad-to-mouth motion



Needle Tubes & Plugs, Teal, Low Viscosity, 19ga., 100/pk 23-11560

Needle Tubes & Plugs, Photo-Blocking Brown, Low Viscosity, 20ga., 100/pk 23-11561

Separate Tubes & Plugs, Photo-Blocking Brown, Low Viscosity, 100/pk 23-21562

Separate Tubes & Plugs, Clear, Medium, Viscosity, 100/pk 23-21562

Black, High Viscosity, 100/pk 23-21564

Brown, High Viscosity, 100/pk 23-21565

Tubes & Plugs - Separate (LA)

Clear tube & Gray rubber plug

Clear Tubes, 100/pkg & Gray Plugs, 100/pkg 23-11551

Clear Tubes, 500/pkg & Gray Plugs, 500/pkg 23-11552





Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

