

Impression Material



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Vinyl Polysiloxane Impression Material

BiteVance Bite Registration Material (ATW)

- Optimal hardness, 80+ durometer
- Thixotropic for easy application
- Non-slumping consistency
- Fluffy, mousse-like
- Excellent dimensional stability
- Improves fit & alignment of crowns & bridges
- Used for bite registration, quick die for temporary inlays & onlays, copying transfers, and primary impression with bite tray
- Won't compress or flex
- 30 seconds working time, and 50 seconds setting time



Fast Set 2 x 50ml Cart., + 6 tips 16-11129

Fast Set 10 x 50ml Cart., no tips 16-12129

Fast Set 60 x 50ml Cart., no tips 16-84651

Exaflex Putty (GC America)

For preliminary impressions and bite registration

Made in Japan



Std Pkg: One 500gm jar ea. of base (blue) & catalyst (yellow), measurer & instructions 37-11230

Clinic Pkg: Five jars (2500gm) ea. of base (blue) & catalyst (yellow), scoops & instructions 37-11231

Examix NDS (GC America)

An advanced addition reaction (platinum catalyzed) VPS impression material that demonstrates outstanding physical properties plus a unique combination of innovative, clinical advantages for impression taking. It produces crisp, extraordinarily smooth and highly precise impressions for models, crown & bridges, dentures and inlays and onlays.



- Unsurpassed dimensional stability
- Patented hydrogen scavenger formula
- Excellent elastic recovery and tear strength
- Thixotropic ideal consistency and flow under pressure makes Examix NDS easy to handle and extremely accurate
- Non-slumping
- Virtually tasteless and odorless

Injection 2 x 48ml Cart., + 6 tips 37-11241

Heavy Body 2 x 48ml Cart., + 6 tips 37-11239

Injection 80 x 48ml Cart. 37-11236

Heavy Body 80 x 48ml Cart. 37-11235

Injection 32 x 48ml Cart. 37-11232

Heavy Body 80 x 48ml Cart. 37-11240

Express Putty (3M)

An excellent all-purpose impression putty that mixes easily without streaking, manipulates nicely, and provides a firm set. Refill includes 1 jar base and 1 jar catalyst-305ml each



Standard Contains: 2 jars (610ml total) 37-11243

ImpreVance VPS Impression Material (ATW)

- High viscosity
- High hydrophilic
- Thixotropic snap-set
- Color coded
- 1:1 mixing ratio, easy to remove tray from the mouth
- Choose your preferred technique-laminar Thixotropic
- No transformation by dimensional stability
- Excellent elastic recovery
- High dimensional stability
- **Regular Set Working Time:** 1' 30" & **Setting Time:** 4' 00"
- **Fast Set Working Time:** 1' 00" & **Setting Time:** 2' 30"



Heavy Body Reg. Set 2 x 50ml Cart., + 6 tips 37-26441

Heavy Body Reg. Set 10 x 50ml Cart. 37-26442

Heavy Body Reg. Set 60 x 50ml Cart. 37-26443

Heavy Body Fast Set 2 x 50ml Cart., + 6 tips 37-26444

Heavy Body Fast Set 10 x 50ml Cart. 37-26445

Heavy Body Fast Set 60 x 50ml Cart. 37-26446

Light Body Reg. Set 2 x 50ml Cart., + 6 tips 37-26450

Light Body Reg. Set 10 x 50ml Cart. 37-26451

Light Body Reg. Set 60 x 50ml Cart. 37-26452

Light Body Fast Set 2 x 50ml Cart., + 6 tips 37-26453

Light Body Fast Set 10 x 50ml Cart. 37-26454

Light Body Fast Set 60 x 50ml Cart. 37-26455

Intraoral Yellow Tips, 100/pkg 37-26459

Intraoral Pink Tips, 100/pkg 37-26460

Intraoral Teal Tips, 100/pkg 37-26461

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Affinis (Coltene/Whaledent)

A perfect AFFINIS impression is the key for accurately fitting, laboratory fabricated restorations. It provides the ideal basis for successful teamwork between the dentist and the dental technician.



Light Body , single pack	37-11541
System 75 Heavy Body , single pack	37-11542
Rigid Body Micro , single pack	37-11543
Regular Body , single pack	37-11544
Light Body Micro , single pack	37-11545
Fast Putty , single pack	37-11546
System 75 Monobody , each	37-11547
Soft Putty , single pack, each	37-11548
Light Body 20 x 25ml Microsystem , bulk pack, each	37-11549
Light Body 20 x 25ml , bulk pack, each	37-11550
Heavy Body 20 x 75ml System 75 , bulk pack, each	37-11551
Monobody 20 x 75ml System 75 , bulk pack, each	37-11552
System 75 , single pack, heavy body, fast set pack	37-11553
System 75 , bulk pack, heavy body, fast set pack	37-11554
Microsystem , single pack, light body, fast set pack	37-11555
Light Body 50ml , single pack fast set	37-11556
Regular Body 50ml , single pack fast set	37-11557
System 360 , single pack, heavy body pack	37-11558
Heavy Body 50ml , single pack	37-11559
System 50 , single pack fast set	37-11560

Affinis Precious (Coltene/Whaledent)

Benefits

- 77% lower than the leading competitive brands
- Allows the impression material to easily capture all critical details around the margins of the preparation
- Produces accurate void free impressions
- Produces clear, easy to read margins
- Leaves you in control of the working time while providing an oral set time similar to fast set material



- Saving time and reducing the potential deformations that can occur when the tray is moved

Light Body 2 x 50ml , 12 mix tips pkg	37-11537
Regular Body 2 x 50ml , 12 mix tips pkg	37-11538
Light Body 4 x 50ml , 20 mix tips + 20 oral tips pkg	37-11539
Regular Body 4 x 50ml , 20 mix tips + 20 oral tips pkg	37-11540

Affinis System 360 Putty (Coltene/Whaledent)

- Excellent positioning pressure
- Instantly tactile, kneadable consistency
- Does not stick to gloves
- Economical 380ml cartridge



Putty Refill Contains: 2 base/catalyst 380ml each, 1 fixation ring	37-11534
Yellow Mixing Tips , 40/pkg	37-11535
Fixation Ring , each	37-11536

AlphaBite Bite Registration (Danville)

Highly stable and accurate bite registration material that is formulated to the specifications of *Blu Mousse*. It has the same 'standby' properties, yet is more value priced. It is also easy to carve without breaking. Alpha Bite is made from quality and silicones.

Work time is 45 sec., set is 75sec.



Fast Set , 2 - 50ml cart.	16-110421
----------------------------------	-----------

Cinch (Parkell)

- Exceptionally accurate, hydrophilic and strong
- Use as a final impression material or an alginate substitute
- Choices of viscosities, final hardness and setting times



Heavy Body Heavy , 90sec. set, 60sec. working time, 4 - 50ml cart. (200ml ttl.) and mixing tips	37-110361
Light , 90sec. set, 60sec. working time, 4 - 50ml cart. (200ml ttl.) and mixing tips	37-110362
Blue-Green Base Mixing Tip , 30 tips (for mucosoft, mucohard, sharp, medium & heavy, cinch putty & heavy)	37-110363

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Colorise Accessories (Zhermack)



Set Temperature Control Unit (110V), used to adjust setting and working times of Colorise impression material. Cool the material to extend the working time. Heat the material to speed up the set time.	37-11175
Universal Tray Adhesive, 1 x 10ml bottle	37-11166
Large Mixing Tips, 48x green mixing tips	37-11167
Small Mixing Tips, 48x yellow mixing tips	37-11169
Intra Oral Tips, 48x yellow tips	37-11171
Dynamic Mixers, 50x	37-11172
Bayonet Rings, 2x	37-11173
Dynamic/Static Mixing Tips, 50x	37-11175
Dynamic/Static Mixing Tip Locker, 2x	37-11165

Colorise Chromatic Color Change Impression Material (Zhermack)



- Utilizes A.M.D.A. - Advanced Moisture Displacement Action to insure accurate detail reproduction
- High flowability and tear strength
- **Fast set times (4 minutes)**
- Vanilla flavor

Intro Kit: 2x heavy body 50ml cart./1x light body 50ml cart./1x regular body 50ml cart./1x clear tray upper/1x clear tray lower/mixing tips/intra-oral tips	37-11177
Dual Arch Kit: 2x colorise rigid 50ml cart./1x colorise light body 50ml cart./6x green mixing tips/6x yellow mixing tips/6x yellow intra-oral tips, hi-tray light dual arch = 1x anterior	37-11190
Putty Blue 1x 300ml base/ 1x 300ml Catalyst/ 2x Spoons	37-11178
Maxi Putty 2x 380ml cart./ 15x Dynamic Mixers/ 2x Bayonet Rings	37-11179

Heavy Body Blue to White to Blue (blue at room temperature) 2x 50ml cart./ 6x **Green Mixing Tips** 37-11180

Heavy Body Maxi Blue to White to Blue (blue at room temperature) 2x 380ml cart./ 15x **Dynamic Mixers/ 2x Bayonet Rings** 37-11182

Rigid Heavy Body Blue to White to Blue (blue at room temperature) (for dual arch technique) 2x 50ml cart./ 6x **Green Mixing Tips** 37-11183

Monophase Blue to White to Blue (blue at room temperature) 2x 50ml cart./ 6x **Pink Mixing Tips** 37-11185

Monophase Maxi Blue to White to Blue (blue at room temperature) 2x 380ml cart./ 15x **Dynamic Mixers/ 2x Bayonet Rings** 37-11187

Regular Body Wash Material Changing Color Form Green to Yellow, 2x 50ml cart./ 12x **Pink Mixing Tips** 37-11188

Light Body Wash Material Changing Color Form Green to Yellow, 2x 50ml cart./ 12x **Yellow Mixing Tips** 37-11189

Extra Light Body Wash Material Changing Color Form Green to Yellow, 2x 50ml cart./ 12x **Yellow Mixing Tips** 37-11176

Defend VPS Impression Material

High performance cartridge system available in regular and fast set. All viscosities are available.

Excellent tear strength to maintain impression integrity upon removal and precise impressions no matter what technique you use.

Exceeds ISO 4380 standard which will help eliminate remakes.



Heavy Body Fast Set, 2 x 50ml cart. + 6 mixing tips	37-105901
Heavy Body Regular Set, 2 x 50ml cart.	37-105902
Light Body Fast Set, 2 x 50ml cart.	37-105903
Light Body Regular Set, 2 x 50ml cart.	37-105904
Heavy Body Fast Set, 4 x 50ml cart.	37-1059051
Heavy Body Regular Set, 4 x 50ml cart.	37-1059052
Light Body Fast Set, 4 x 50ml cart.	37-1059053
Light Body Regular Set, 4 x 50ml cart.	37-1059054

Elite Glass (Zhermack)

Creation of transparent matrices for intra and extra-oral polymerization of both light-cured & self-cured composites

Applied in a wide variety of clinical



Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Situations: presence of pigmentations, fixing diastema, correcting rotations, fractures, dental anomalies, processing of temporary crown & bridge, etc.

- Thixotropic, does not slump
- Clear VPS material
- Rigid after setting (70 shore-A)

Standard Pack Clear, 2x 50ml cart./ 6x 37-11195
Green Mixing Tips

Elite HD + (Zhermack)

The presence of the nanostructure increases the fluidity of the material, particularly when pressure is applied

- High flowability due to nano-technology
- Super ability to displace moisture resulting in excellent detail reproduction
- Excellent tear strength
- High deformation strength
- Excellent elastic recovery
- High dimensional stability



3 year guarantee

Fast Set - 4 minutes - Regular Set - 5 minutes 30 seconds

Tray Material - Heavy Body - Fast set 37-11203
orange 2x 50ml cart./ 6x green mixing tips

Wash Material - Light Body - Fast set 37-11201
green 2x 50ml cart./ 12x yellow mixing tips

Wash Material - Super Light Body- Fast 37-11199
set violet 2x 50ml cart./ 12x yellow mixing tips

Wash Material - Light Body - Regular 37-11207
set blue 2x 50ml cart./ 12x yellow mixing tips

Wash Material - Light Body - Regular 37-11208
set blue 1x 90ml tube base/ 1x 90ml tube catalyst

Wash Material - Regular Body - Regular 37-11218
set blue 2x 50ml cart./ 12x yellow mixing tips

Wash Material - Regular Body - Regular 37-11215
set blue 1x 90ml tube base/ 1x 90ml tube catalyst

Monophase - Medium Body - Regular 37-11213
set yellow 2x 50ml cart./ 12x yellow mixing tips

Monophase - Medium Body - Regular 37-11211
set maxi yellow 2x 380ml cart./ 15x dynamic mixers/ 2x bayonet rings

Monophase - Medium Body - Regular 37-11210
set maxi - economy pack yellow 6x 380ml cart.

Tray Material - Heavy Body - Fast set 37-11209
maxi orange 2x 380ml cart./ 15x dynamic mixers/ 2x bayonet rings

Tray Material - Heavy Body - Fast set 37-11197
maxi - economy pack orange 6x 380ml cart.

Elite HD + Accessories (Zhermack)



D2 Dispenser (1:1) 37-11229

Replacement Plunger for Universal Dispenser, 50ml 37-11220

D3 Multi-Dispenser (1:1, 2:1, 10:1, 4:1) 37-11221
1x applicator/ 1x 4:1 plunger/ 1x 1:1/ 2:1 plunger/ 1x 10:1 plunger/ 1x 1:1/ 2:1 breech/ 1x 4:1/ 10:1 breech

Universal Tray Adhesive, 1x 10ml bottle 37-11222

Large Mixing Tips, 48x green mixing 37-11223
tips

Small Mixing Tips, 48x yellow mixing 37-11224
tips

Intra Oral Tips, 48x yellow tips 37-11225

Dynamic Mixers, 50x dynamic mixers 37-11226

Bayonet Rings, 2x bayonet rings 37-11227

Dynamic/Static Mixing Tips, 50x 37-11228

Dynamic/Static Mixing Tip Locker, 2x 37-11219

Elite HD + Soft Putty (Zhermack)

The presence of the nanostructure increases the fluidity of the material, particularly when pressure is applied

- High flowability due to nano-technology
- Super ability to displace moisture resulting in excellent detail reproduction
- Excellent tear strength
- High deformation strength
- Excellent elastic recovery
- High dimensional stability



3 year guarantee

Fast Set - 4 minutes - Regular Set - 5 minutes 30 seconds

Putty Soft - Regular Set orange 1x 37-110631
250ml tub base/ 1x 250ml tub catalyst/ 2x spoons

Putty Soft - Regular Set Economy Pack 37-110632
orange 2x 450ml tub base/ 2x 450ml tub catalyst/ 2x spoons

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Putty Soft - Fast Set yellow 1x 250ml tub 37-110633
base/ 1x 250ml tub catalyst/ 2x spoons

Putty Soft - Fast Set Economy Pack 37-110634
yellow 2x 450ml tub base/ 2x 450ml tub
catalyst/ 2x spoons

Extrude VPS Impression Material (Kerr)

Extrude is the proven standard against which other impression materials are judged

With five viscosities, including a putty, you have all the options you need for accurate results, regardless of your preferred technique



Heavy (Extra), 2 cart. w/o tips 37-99041

Heavy (Extra), 24 cart. w/o tips 37-99044

Wash, Blue, 2 cart. 50ml each, w/o tips 37-99045

Wash, Blue, 24 cart. w/o tips 37-99048

Medium, Green, 2 cart. 50ml each w/o tips 37-99049

Medium, Green, 24 cart. 50ml each w/o tips 37-99052

First Quarter VPS Impression Material (Danville)

1 ¾ Minute Set Time! First Quarter VPS is a super fast impression material great for smaller, triple tray type impressions

Works in a 1/4 of the time of traditional 6-minute materials, making single unit impression fast, easy and trouble-free



Light Body, 4 x 50ml cart. + 12 tips 37-11272

Heavy Body, 4 x 50ml cart. + 12 tips 37-11260

Light Body, 2 x 50ml cart. + 6 tips 37-11261

Heavy Body, 2 x 50ml cart. + 6 tips 37-11262

Genie (Sultan)

Ultra hydrophilic

Patient-preferred taste

Low extrusion force

Consistent quality



Light Body Fast Set, 2 x 50ml cart. 37-11303

Light Body Fast Set Bulk, pkg 60 x 50ml 37-11304
cart.

Light Body Regular Set, 2 x 50ml cart. 37-11305

Light Body Regular Set Bulk, pkg 37-11306
60 x 50ml cart.

Regular Body Fast Set, 2 x 50ml cart. 37-11302

Regular Body Fast Set Bulk, pkg 37-11301
60 x 50ml cart.

Regular Body Regular Set Bulk, pkg 37-11299
60 x 50ml cart.

Heavy Body Fast Set, 2 x 50ml cart. 37-11298

Heavy Body Fast Set Bulk, pkg 37-11297
60 x 50ml cart.

Heavy Body Regular Set, 2 x 50ml cart. 37-11296

Heavy Body Regular Set Bulk, pkg 37-11295
60 x 50ml cart.

Regular Body Regular Set, 2 x 50ml 37-11307
cart.

iMix Dynamic Cartridge VPS (Pac-Dent)

- Easy to use
- Easy dispensing
- Highly accurate impressions
- Exceptional coverage for precision
- Offers ample working time



Monophase Standard Set, 1 x 380ml 37-87491
cart.

Heavy Body Standard Set, 1 x 380ml 37-87492
cart.

Imprint 3 Garant (3M)

Improved technology features high hydrophilicity, tensile strength, elongation, toughness, and memory. For crown & bridge, implant, and orthodontic impression in all types of trays



Heavy Body Refill Contains: 4 x 50ml 36-35751
garant cart. & 10 green mixing tips

ModulMix (Zhermack)

Automatic A-Silicone mixer with a 5:1 mixing ratio. ModulMix, available in 6 colors, will fit perfectly with all types of surgery décor

ModulMix automatically mixes impression material for impressions prepared in 5:1 mixing ratio cartridges (362-380-390ml)



Impression Material

With respect to manual mixers, the Zhermack system automatic mixing is easy, constant and offers higher quality

Mixing phase is notably reduced to give more product working time



Two expulsion button located on each side of the machine to facilitate the filling of the impression tray with two speeds (70ml per minutes and 155ml per minute)

Motorized forward and reverse motion of the finger-safe pistons to expedite cartridge substitution (max 3 seconds)

- Cartridge level indicator
- Unit can be wall-mounted

2 year warranty

Ultimate Impression Package Includes: 37-11382
C309035 modulmix - 1x side panel (gray), 1x base support, 1x digital timer, 1x AC adapter, C300100 freealgin alginate replacement, maxi 2x 380ml cart., 15x dynamic mixers, 2x bayonet rings

Cover Light , blue (2x)	37-11381
Cover , orange (2x)	37-11380
Cover , blue (2x)	37-11379
Cover , green (2x)	37-11378
Cover , white (2x)	37-11377
Single Wall Mounting Unit	37-11376
Double Wall Mounting kit	37-11375
Charger , 115 volt	37-11374

Mojo II Syringes (Danville)

Disposable, single-use syringes that don't waste impression material



They have locking and storing capabilities that other material syringes don't have.

The locking positions are: fill, store, and activate

The 18g tip syringes work with 1:1 impression cartridges only
100/pkg 36-106401

PuttyVance Impression Material (ATW)

Used in putty wash technique

Non slumping, air-free filling to reduce air bubbles at margin

Setting Time: 2':30"
Should be mixed 1:1 ratio



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Regular Set, 2 Jars: 300ml base, 300ml catalyst & 2 spoons 37-11395

Fast Set, 2 Jars: 300ml base, 300ml catalyst & 2 spoons 37-11325

Star VPS Impression Material (Danville)

Extra smooth texture & ideal thixotropic flow properties result in excellent surface detail reproduction



Ultra Light , 4x50ml cart.	37-11399
Ultra Light Econo Pack , 24x50ml cart.	37-11410
Light Body , 2x50ml cart., + 6 tips	37-11411
Light Body , 4x50ml cart.	37-11412
Light Body Econo Pack , 24x50ml cart.	37-11413
Light Body Super Econo Pack , 96x50ml cart.	37-11414
Heavy Flex , 2x50ml cart., + 6 tips	37-11415
Heavy Flex , 4x50ml cart.	37-11416
Heavy Flex Econo Pack , 24x50ml cart.	37-11417
Heavy Flex Super Econo Pack , 96x50ml cart.	37-11409
Ultra Heavy , 4x50ml cart.	37-11408
Ultra Heavy Econo Pack , 24x50ml cart.	37-11400
Heavy Stiff Bite , 2x50ml cart. + 6 tips	37-11401
Heavy Stiff Bite , 4x50ml cart.	37-11402
Heavy Stiff Bite Econo Pack , 24x50ml cart.	37-11403
Heavy Stiff Bite Super Econo Pack , 96x50ml cart.	37-11404
Clear Bite , 2x50ml cart., + 6 tips	37-11405
Clear Bite , 4x50ml cart.	37-11406
Clear Bite Econo Pack , 24x50ml cart.	37-11407
Putty , 2 of 400gm tubs base & catalyst	37-11418

Vonflex Putty (Vericom)

Advantages

- Easy removal dosage
- Ideal consistency
- Easy handling of mixing
- Precise reproduction of the finest details
- Exceptional dimensional stability
- Optimal elastic recovery

Indications



Impression Material

- For individual tray
- For taking of crown & bridge

2 Jars of Base, 280ml, catalyst 280ml & 37-11425
2 measuring scoops

Vonflex VPS Impression Material (Vericom)

Indications

- For crown & bridge impression
- For inlay and onlay impression
- For all types of pickup impression

Advantages

- Homogeneous and bubble free material quality
- Easy & clean handling
- Easy & safe removal from mouth
- Ideal final hardness
- Two types; injection & tube



Heavy Body Regular Intro Kit: Set 2 x 23-98131
50ml cart., + 6 tips

Heavy Body Regular Value Kit: Set 12 x 23-98132
50ml cart.

Light Body Regular Intro Kit: Set 2 x 23-98133
50ml cart., + 6 tips

Light Body Regular Value Kit: Set 12 x 23-98134
50ml cart.

Polyether Impression Material

Impregum F Double Pack Tube (3M)

Handmix formulation, medium body polyether material with high-precision impressions.

Working Time (including mixing): 2:45min

Setting Time (from start of mix): 6:00min



Double Pack: 2 tubes base, 2 tubes catalyst & mixing pad 37-11340

Impregum Penta Medium Body (3M)

Impregum Penta material displaces moisture when it really counts - right from the first contact with tissue - for precise impressions.

Impregum Penta material produce finely detailed impressions for a better fitting restoration

Working Time (including mixing): 2:45min

Setting Time (from start of mix): 6:00min



Double Pack: 2-X 300ml foil bags, base plate & 2 x 60ml catalyst paste (red) 37-11343

Penta Syringe Tip, 12/pkg 37-11342

Red Mixing Tips, 50/pkg 9 (HB) 37-11341

Elastomer Syringe 37-12345



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Impregum Penta Soft Medium (3M)

Indications

- Crowns & bridges
- Inlays & onlays
- Functional impressions
- Implant abutments
- Quick-setting materials especially suited for single-unit & 2-unit impressions.



Double Pack: 2 x 300ml base plate, 2 x 60ml catalyst 37-11345

Mixing Tips, 50/pkg (HB) 37-11346

Cure Box (Meta)

- Cures light activated composted
- Includes four 9 watts bulb
- Heavy duty and durable for long lasting use
- Equipped w/ 5 and 10 minute timer
- LxWxH: 11" x 7.5" x 4.25"



Each

30-98841

Cure Box (Vericom)

Easy handling, powerful intensity & selectable curing time, light curing device for individual tray resin

- Short-time polymerization/ easy handling
- Very effective for polymerization
- Making precise individual trays (5min)
- 220V/60Hz, 110V/50Hz
- 4 lamps x 9w



Cure Box

30-50154

Impregum Penta Elastomer Syringe (3M)

Single Syringe, reusable 37-119651

Penta Elastomer Syringe, full set 37-109652

Syringe Tips, 12/pkg 37-109651



www.prestigedentalproducts.com

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Rubber Base Impression Material

Coe-Flex (GC America)

Impression material for inlays, crowns, bridges, dentures, rebasing, & relining



Regular Standard package	37-11163
Injection Standard package	37-11162
Heavy Body , standard package	37-11161
Fast Set , standard package	37-11164

Indurent Gel C-Silicone Impression Material (Zhermack)

High Quality Putty 85 Shore

The first silicone specifically researched for laboratory use, simplifying the dental technicians to a point where it is irreplaceable. Thanks to versatility, this is the most utilized laboratory silicone in the world.



Catalyst for Putty Pink , 1 x 60ml tube	37-11348
--	----------

Permlastic (Kerr)

Polysulfide tear resistant material for denture impressions.



Light Body , standard package in tubes	37-25161
Regular Body , standard package in tubes	37-25162
Heavy Body , standard package in tubes	37-25163

Alginate

Alghamix II (Zhermack)

- 2 bowl rotation speeds with ergonomic speed selector
- Available bowls: small 350cc for alginate, large 600cc for plaster or stone
- Bayonet bowl-fastener



Automated Mixer 1x alghamix II mixer/ 1x alghamix flexible mixing bowl (grey)/ 1x mixing spatula	37-11118
--	----------

Mixing Bowl - Light blue 1x mixing bowl (350CC)	37-11117
Mixing Bowl - Orange 1x mixing bowl (350CC)	37-11116
Mixing Bowl - Navy blue 1x mixing bowl (350CC)	37-11115
Mixing Bowl - Green 1x mixing bowl (350CC)	37-11114
Mixing Bowl - White 1x mixing bowl (350CC)	37-11113
Mixing bowl for stones/ Plasters 1x mixing bowl (600cc)	37-11112
Disc For Mixing Silicones , 1x pads (20 sheets per pad)	37-11111

Alginate Flavoring



Mint Flavor , 2oz	37-11128
Grape Flavor , 2oz	37-11120
Punch Flavor , 2oz	37-11121
Cherry Flavor , 2oz	37-11122
Orange Flavor , 2oz	37-11123
Blueberry Flavor , 2oz	37-11124
Raspberry Flavor , 2oz	37-11125
Bubblegum Flavor , 2oz	37-11126
Pina Colada Flavor , 2oz	37-11127
Strawberry Flavor , 2oz	37-11119

Alginate Flavoring Holder (Dentsply)

Manufactured from heavy duty acrylic. Will hold up to six bottles of alginate flavoring. Uniquely designed handgrip is perfect for carrying alginate flavoring between chairs or satellite offices.



Smoke Acrylic	37-11129
----------------------	----------



Find us on:
facebook®

Impression Material

Alginate Flavorings (Dentsply)

Creates great tasting alginate at an economical low prices. Highly concentrated for maximum flavor. Sugar free. 2oz bottle



Banana	37-11149
Blueberry	37-11139
Bubble Gum	37-11137
Cherry	37-11136
Chocolate	37-11135
Cinnamon	37-11134
Cola	37-11133
Crème Caramel	37-11132
Fruit Punch	37-11131
Grape	37-11138
Jaimaica Java Coffee	37-11140
Lemon	37-11148
Orange	37-11147
Peppermint	37-11146
Pina Colada	37-11145
Pineapple	37-11144
Raspberry	37-11143
Root Beer	37-11141
Strawberry	37-11142
Watermelon	37-11130

Alginate Mixing Bowls (Dentsply)

Made from extremely smooth and pliable plastic, for easy mixing and cleaning. PVD is not used in the manufacturing process, eliminating any danger of migrating plasticizers. Glossy inside surface, making alginate removal easy



Blue Dark	45-11221
Blue Ocean Splash	45-11208
Blue Pearl	45-11207
Blue Tropical	45-11209



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Green Kiwi Splash	45-11210
Green Pearl	45-11211
Green Tropical	45-11212
Orange Tangerine Splash	45-11216
Orange Tropical	45-11215
Patriotic Splash	45-11214
Pearlescent	45-11213
Pink Tropical	45-11218
Red	45-11217
Splash	45-11219
Yellow Tropical	45-11220

Bio Jel Alginate

Dustless alginate type I/ type II applications? Removable prosthetic? Fixed and removable prosthetic characteristics? Type I and type II? Apricot scent with a hint of mint advantages? Soft surface and high detail reproduction (0.050mm)? Excellent resistance (160.6PSI)? 3 layers packaging to guarantee reliable material.



Regular Set, 1lb	37-110561
Fast Set, 1lb	37-110562

Blendex Automatic Alginate Mixer (Dentamerica)

- Automatic trouble-free alginate mixing at the touch of a button
- Time saving process takes only seconds to complete
- Consistent, smooth, and bubble-free mixing each and every time
- Features automatic timer with audible tone
- Three customized memory presets available
- LED displays shows mixing time and memory presets
- Automatic door safety lock ensures safe operation

Specifications: 100/110V or 220/240V AC

Dimension: W205 x D245 x H295mm

Weight: 18kg

Accessories: Main unit, power cable, mixing bowl, spatula, mixing stick, water bottle, magnetic plastic basket

Unit	37-11157
-------------	----------



Cavex Colorchange Alginate

Cavex Colorchange, visible setting process, dimensionally stable for 5 days

- 9 days dimensional stability
- Fast set, dust-free & easy mixing
- Smooth surface, optimum gypsum compatibility
- Rapid, complete hardening reaction
- Highly elastic, ideal for orthodontic
- Good tear and deformation resistance



Impression Material

- Detailed reproduction at least 25µm

Fast Set, 1.1 lb bag (500g) 37-108431

Disposable Mixing Bowls



White Bowls , 36/pkg	37-99201
Neon Blue Bowls , 36/pkg	37-99202
Neon Yellow Bowls , 36/pkg	37-99203
Neon Green Bowls , 36/pkg	37-99204
Neon Pink Bowls , 36/pkg	37-99205
Neon Purple Bowls , 36/pkg	37-99206
Assorted Bowls , 36/pkg	37-99207
White Spatulas , 36/pkg	37-99208
Neon Blue Spatulas , 36/pkg	37-99209
Neon Yellow Spatulas , 36/pkg	37-99210
Neon Green Spatulas , 36/pkg	37-99211
Neon Pink Spatulas , 36/pkg	37-99212
Neon Purple Spatulas , 36/pkg	37-99213
Assorted Spatulas , 36/pkg	37-99214

Elite Mix (Zhermack)

Vacuum-sealed mixer for mixing; prosthetic plaster models, coatings with the lost-wax casting technique and silicones for duplication
2 year warranty



1x elite mix/ 1x power cord/ 1x 500ml lid/ 1x 500ml bowl	41-112201
Stand , 1x elite mix stand/ 1x wrench	41-112191
Lid (blade) 1x 200ml	41-112181
Lid (blade) 1x 500ml	41-112171
Lid (blade) 1x 1,000ml	41-112161
Bowl , 1x 200ml	41-112151
Bowl , 1x 500ml	41-112141
Bowl , 1x 1,000ml	41-112131



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Freealgin Alginate Replacement (Zhermack)

A medium viscosity addition silicone, used in substitution of traditional alginates for preliminary impression taking



One of its advantages is its high dimensional stability, which ensures flexibility in the impression casting times and great accuracy. Greater elasticity, a better recovery from deformation and a higher tear resistance also complement its characteristics.

For automatic mixing, the packaging is also available in a 5:1 ratio, which is quick and convenient.

Medium Body, 2x 50ml cart./ 4x lilac mixing tips 37-11284

Medium Body - Economy Pack, 24x 50ml cart. 37-11283

Medium Body - Maxi, 2x 380ml cart./ 15x dynamic mixers/ 2x bayonet rings 37-11282

Medium Body - Maxi - Economy Pack, 6x 380ml cart./ 2x bayonet rings 37-11281

Dispenser 1:1, each 37-11280

Replacement Plunger for Universal Dispenser, 50ml 37-11279

Multi-dispenser (1:1, 2:1, 10:1, 4:1) 1x applicator/ 1x 4:1 plunger/ 1x 1:1/ 2:1 plunger/ 1x 10:1 plunger/ 1x 1:1/ 2:1 breech/ 1x 4:1/ 10:1 breech 37-11278

Dynamic Mixers, 50x 37-11276

Bayonet Rings, 2x 37-11275

Universal Tray Adhesive, 1x 10ml bottle 37-11274

Hydrogum (Zhermack)

In addition to the alginates with specific characteristics, the wide Zhermack range also allots an important place to the universal alginates, to meet the many needs of daily clinical practice.



Like Hydrogum and Hydrogum Soft, that are particularly elastic and resistant. Or like the chromatic alginates Tropicalgin and Phase Plus, designed to highlight in an easy and immediate way the processing steps through their color change.

Refill Green, 1x 500g (1.1 lb) bag 37-11318




Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Hydrogum 5 (Zhermack)

Hydrogum 5 and Hydrocolor 5 are the high dimensional stability solutions studied by Zhermack, which will allow to have more time for the casting of the model after the impression taking, thus overcoming working time that is often too tight, and simplifying the daily work of the dentist. They enable the impression to be cast up to 5 days after the impression is taken, thus obtaining precise and reliable models




Hydrogum 5 is an extra-fast setting alginate, with a homogeneous consistency and an accurate reproduction of details. It ensure less stress for the patient, thanks to the time in mouth of only 45 seconds.

Refill Lilac, Mangustan Fruit Flavor, 1x 37-11322
453g (1 lb) bag

Refill Lilac, Mangustan Fruit Flavor, 37-11323
12x 453g (12 lb) bag

Hydrogum Soft (Zhermack)

In addition to the alginates with specific characteristics, the wide Zhermack range also allots an important place to the universal alginates, to meet the many needs of daily clinical practice.



Like Hydrogum and Hydrogum Soft, that are particularly elastic and resistant. Or like the chromatic alginates Tropicalgin and Phase Plus, designed to highlight in an easy and immediate way the processing steps through their color change.

Refill Light Blue, 1x 453g (1 lb) bag 37-11324

Identic Alginate (Dux)

- Dust-free alginates
- A full strength alginate that is dust-free
- Antimicrobial, easy to mix
- Pleasant tasting and produces an excellent stone surface
- Pink with cinnamon scent



Fast Set, 1 lb bag 37-109851

Fast Set, 8 x 1 lb bag 37-109852

Kromatica Alginate Color Change

Alginate mint, color changing alginate, 1 lb aluminum foil pouch. With the chromatic phase indicator you visually control the different working phases

Stable hydrocolloid class A alginate. Keeps its shape up to 100 hours



Thixotropic. Fast set, in-mouth setting time is about 30 seconds

Fast Set, 1 lb mint bag 37-110521

Regular Set, 1 lb mint bag 37-110522

Kromopan Color Changing Alginate Dust Free

Chromatic alginate for high precision and elastic impressions

Changes color from purple to white and will hold an accurate impression for up to 100 hours



Dust-free

Made in Italy

Fast Set, 1 lb/bag 37-108121

Regular Set, 1 lb/bag 37-108122

Max Print Alginate Dustless (MDC)

Applications

- Removable prosthetic
- Fixed and removable prosthetics

Advantages

- Less time in mouth (type I)
- Soft surface and high detail reproduction (0.050mm)
- 3 layer packaging to guarantee reliable material



Fast Set, Cherry Scent, 1 lb/bag 37-11354

Regular Set, Cherry Scent, 1 lb/bag 37-11355

Max Print Chromatic Alginate

Applications

- Removable prosthetics
- Fixed and removable prosthetics

Characteristics

- Chromatic alginate - three phases
- Violet - mixing phase
- Pink - insert into mouth
- White - remove from mouth
- Fast gel time type I

Advantages

- 130 hours without alterations and without losing and details
- 3 layer packaging to ensure the materials stability
- Less salivation from the patient
- Good resistance



Fast Set, Mint Scent, 1 lb/bag 37-11351

Max Print Cyan Alginate (MDC)

Applications

- Removable prosthetics
- Fixed and removable prosthetics
- Accurate impressions

Characteristics

- Extra fine powder
- Dustless

Advantages

- Made with hydrochloric, an antiseptic which eliminates cross contamination
- Smooth texture and high detail reproduction
- Excellent resistance (160.6 PSI) does not tear upon removal



Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Fast Set, Cherry Scent, 20 lb/box	37-11434
Regular Set, Cherry Scent, 20 lb/box	37-11435
Fast Set, Cherry Scent, 1 lb bag	37-11352

Max Print Dustless Orthodontic Alginate

Applications

- Orthodontic use

Characteristics

- Helps avoid tearing when taking dental impressions with braces
- Three layer packaging to guarantee material reliability



1 lb bag	37-11449
20 lb master box	37-11450

Max Print Pediatric Alginate

Applications

- Pediatric use

Advantages

- Less salivation from patient
- Less time in mouth
- Soft surface and high detail reproduction (0.050mm)
- 3 layer packaging to guarantee material reliability



Chocolate Scent, Extra Fast Set, 1 lb/bag	37-11445
Apple Scent, Extra Fast Set, 1 lb/bag	37-11446
Chocolate Scent, Extra Fast Set, 20 lb master box	37-11447
Apple Scent, Extra Fast Set, 20 lb master box	37-11448

Mixing Bowls

- Flexible plastic material.
- Three different sizes to fit different jobs
- Washable for repeated use



Mixing Bowls, Set of 3: Large, medium, & small	37-11373
---	----------

Mixing Bowls

Dark Green mixing bowls



Large 600c, each	37-108801
Medium 350cc, each	37-108802

Small 150cc, each	37-108803
--------------------------	-----------

Noginate Alginate Substitute Material (Danville)



2 x 50ml cart. + 6 short green mixing tips	37-110511
---	-----------

Orthoprint (Zhermack)

Orthoprint is a specific alginate for orthodontic use that combines precision and practicality



Highly elastic and suitable for impression taking even in the presence of brackets. It is characterized by extra-fast working and setting times to allow for the quick removal of the impression from the oral cavity, while recording even the smallest details

Its vanilla scent, with emetic reflex reduction, helps patients to overcome any fears.

Refill Yellow, Vanilla Flavor, 1.1 lb bag	37-11390
One Case Yellow, Vanilla Flavor, 12 x 1.1 lb bag	37-11391

Phase Plus (Zhermack)

In addition to the alginates with specific characteristics, the wide Zhermack range also allots an important place to the universal alginates, to meet the many needs of daily clinical practice.



Like Hydrogum and Hydrogum Soft, that are particularly elastic and resistant. Or like the chromatic alginates Tropicalgin and Phase Plus, designed to highlight in an easy and immediate way the processing steps through their color change.

Refill White, 1 lb bag	37-11394
-------------------------------	----------

Spatula for Alginate & Plaster

Spatula A



Spatula B

Spatula A	40-98911
Spatula B	40-98912

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Spatulas for Plaster & Alginate

Fig. 4



Fig. 5



Fig. 6



Fig. 4 40-13784

Fig. 5 40-13785

Fig. 6 40-13786

Spectrum Flowabowl Mixing Bowls (Plasdent)

Large



Medium

Spectrum Blue, large 600cc 37-11369

Spectrum Yellow, large 600cc 37-11370

Spectrum Green, large 600cc 37-11371

Spectrum Ruby Red, large 600cc 37-11356

Spectrum Purple, large 600cc 37-11368

Spectrum Pink, large 600cc 37-11367

Spectrum Orange, large 600cc 37-11366

Spectrum Blue, medium 350cc 37-11365

Spectrum Yellow, medium 350cc 37-11364

Spectrum Green, medium 350cc 37-11363

Spectrum Ruby Red, medium 350cc 37-11362


Spectrum Purple, medium 350cc 37-11361

Spectrum Pink, medium 350cc 37-11360

Spectrum Orange, medium 350cc 37-11359

Tropicalgin (Zhermack)

In addition to the alginates with specific characteristics, the wide Zhermack range also allots an important place to the universal alginates, to meet the many needs of daily clinical practice.



Like Hydrogum and Hydrogum Soft, that are particularly elastic and resistant. Or like the chromatic alginates Tropicalgin and Phase Plus, designed to highlight in an easy and immediate way the processing steps through their color change.

Refill Yellow, 1 lb bag

37-11424

Vacuum Mixer

Vacuum mixer is designed for mixing all types of gypsum, plaster and investment material

Features

- Powerful vacuum - The vacuum pumps quickly, creates high vacuum that eliminates bubbles in the mix
- Hands free mixing - Pre-mix materials (power & liquid) are easily mix with the alternating paddle direction
- Memory Functions - 11 pre-program speed and mix time for specific materials
- Mixing bowls - 500ml polycarbonate bowl that's easy to clean and that last a long time

Specifications

Dimensions: 350 x 160 x 220mm (H x W x D)

Weight: 18 lbs

Voltage: 100-240 voltage, 50/60 Hz

Mixing Speed: 100rpm-500rpm

Vacuum: 27.5 in Hg (930 mbar)

User Programs: 11 (P:00 - P:10)

Mixing Time: 5 minutes

Fuse Type: 4 Amp

Operational Altitude: Sea level - 2,000 meters (6,560 feet)

Operational Temperature: 5°C - 40°C

Each, contents 1 - bowl 500ml 1 - power cord 1 - user manual 1 - extra vacuum filter 1 - allen wrench base cover column/base assembly hardware



Mixing Tips

Affinis Mixing Tips (Coltene/Whaledent)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg 37-99179

Yellow, 40/pkg 37-99178

Pink, 40/pkg 37-99177

Affinity Multi-Prep Mixing Tips (Clinician's Choice)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg 37-109991

Impression Material

Yellow, 40/pkg 37-109992
Pink, 40/pkg 37-109993

Affinity VPS Mixing Tips (Clinician's Choice)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip
Pink: 5.4mm tip
Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110011
Yellow, 40/pkg 39-110012
Pink, 40/pkg 39-110013

Algin X Ultra Mixing Tips (Dentsply Caulk)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip



Teal/Green, 40/pkg 39-110021

Alginate Alternative Mixing Tips (H Shein)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Pink: 5.4mm tip



Pink, 40/pkg 39-110041

AlgiNot Mixing Tips (Kerr Restoratives)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip



Teal/Green, 40/pkg 39-110051

Contact Us



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Aquasil Mixing Tips (Dentsply Caulk)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip
Pink: 5.4mm tip
Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110031
Yellow, 40/pkg 39-110032
Pink, 40/pkg 39-110033

Chinch Mixing Tips (Parkell)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Pink: 5.4mm tip
Yellow: 4.2mm tip



Pink, 48/pkg 39-110061
Yellow, 48/pkg 39-110062

Colorise Mixing Tips (Zhermack)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip
Pink: 5.4mm tip
Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110071
Yellow, 40/pkg 39-110072
Pink, 40/pkg 39-110073

Delivery Tips (PacDent)

125 -

Brush-padded end

1169 -

Adjustable brush tip fiber length, tightly fitted fibers for low viscosity materials



1169N - 1169 with no fibers

126 - Snap-off tip provides reseal ability

126-2 - Included cap provides reseal ability

1357 - Brush-padded end, stand-out fibers

128 - Pre-bent needle tip, intra-oral connection w/ HP yellow mixers

Impression Material

Pre-bent Needle Tip w/ Brush-padded End, 100 tips/pkg	23-25611
Pre-bent Needle Tip w/ Micro Brush, 100 tips/pkg	23-25612
Pre-bent Needle Tip w/ Yellow I/O Connection, Green, 100 tips/pkg	23-25613
Pre-bent Needle Tip w/ Teal I/O Connection, Green, 100 tips/pkg	23-25614
I/O Tips w/ Cap, 100 tips/pkg	23-25615
I/O Tips w/ Snap-off Cap, 100 tips/pkg	23-25616
Curved Black w/ Brush Tips, 100 tips/pkg	23-25617
Curved Black Tips w/ No Brush, 100 tips/pkg	23-25618

Dimension Garant L Mixing Tips (3M)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit. **Yellow:** 4.2mm tip



Yellow, 40/pkg 39-110081

Dynamic Green Mixing Tips (Defend)

Ultra high quality mixing tips for use with automatic impression material mixing systems. Compatible with

- Splash fat pak/ discus dental
- Precision fat pak/ discus dental
- Genie/ sultan
- Exajet/ GC America
- Flexitime magnum 360/ heraeus Kulzer
- System 380/ coltene whaledent
- P2/ heraeus Kulzer
- VP mix 360/ henry shein



50/pkg 37-105381

Dynamic Green Mixing Tips (PacDent)

Ultra high quality mixing tips for use with automatic impression material mixing systems. Compatible with

- Splash fat pak/ discus dental
- Precision fat pak/ discus dental
- Genie/ sultan
- Exajet/ GC America
- Flexitime magnum 360/ heraeus Kulzer
- System 380/ coltene whaledent
- P2/ heraeus Kulzer
- VP mix 360/ henry shein



50/pkg 37-35502

Dental Supplies



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Dynamic Mixing Tips



Red, for 3M Penta 5:1, 50/pkg	37-99171
Blue, for Dentsply Duomix, 50/pkg	37-99172
Yellow, for Heraeus Kulzer, GC, Huge, Muller-Omicron, Zhermack, 50/pkg	37-99173
Grey, for Zhermack Modulmix 5:1, 50/pkg	37-99174
Blue, for Dentsply, (old style machine), 50/pkg	37-99175
Green, for Kettenbach, (old style machine) 50/pkg	37-99176

Dynamic Mixing Tips for Flexitime



50/pkg 37-111001
Bulk Pkg, 5 x 50/pkg 37-111002

Dynamic Red Mixing Tips (PacDent)



50/pkg 37-11193

Dynamic Teal Mixing Tips (PacDent)

Ultra high quality mixing tips for use with automatic impression material mixing systems. Compatible with

- Elite HD+ maxi/ zhermack
- Impregum 1 & 2/ 3M
- Virtual 380ml/ ivoclar
- Flexitime/ Heraeus
- Aquasil auto mix/ dentsply
- Colorise/ zhermack
- Correct plus mighty mix/ pentron
- Defend auto mix/ mydent
- Take 1 advanced, alginot FS/ kerr



Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

50/pkg 35-98591

Elite HD+ Mixing Tips (Zhermack)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110091

Yellow, 40/pkg 39-110092

Pink, 40/pkg 39-110093

Exa'lence Mixing Tips (GC America)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110101

Yellow, 40/pkg 39-110102

Pink, 40/pkg 39-110103

Exa'lence Mixing Tips (GC America)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110111

Yellow, 40/pkg 39-110112

Pink, 40/pkg 39-110113

Examix NDS Mixing Tips (GC America)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg 39-110121

Yellow, 40/pkg

39-110122

Pink, 40/pkg

39-110123

Express Mixing Tips (3M)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Yellow: 4.2mm tip



Yellow, 40/pkg

39-110131

Extrude Mixing Tips (Kerr Restoratives)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg

39-110141

Yellow, 40/pkg

39-110142

Pink, 40/pkg

39-110143

First Half VPS Mixing Tips (Danville)

The new mixing technology provides less material waste and consistently high mixing quality. The changeover to the new T-mix tips is seamless: 100% compatible with the existing MIXPAC systems is guaranteed - as before, all original accessories fit.

Teal/green: 6.5mm tip

Pink: 5.4mm tip

Yellow: 4.2mm tip



Teal/Green, 40/pkg

39-110151

Yellow, 40/pkg

39-110152

Pink, 40/pkg

39-110153

FujiCem Automix - Mixing Tips (GC America)

Use with GC FujiCem resin-modified glass ionomer cement in paste-pak dispensing system



20/pkg

37-98581

www.prestigedentalproducts.com

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Hub & Intra Oral Tip (Plasdent)



Tip For GC Syringe Clear, 100pcs/bag	37-11314
Hub For Free-Flo Syringe	37-11311
Tip For Free-Flo Syringe Clear, 100pcs/bag	37-11312
Hub For GC Syringe	37-11313

Impregum Penta - Mixing Tips



50/pkg 37-35061

Impregum Penta Elastomer Syringe (3M)



Single Syringe, reusable	37-119651
Penta Elastomer Syringe, full set	37-109652
Syringe Tips, 12/pkg	37-109651

Intra Oral XX-Fine Mixing Tips



50/pkg 37-105391

Intraoral Tips



Yellow, 100pcs/bag	37-11350
Black, 100pcs/bag	37-11490
Clear, 100pcs/bag	37-11491
Clear, 100pcs/bag	37-11492
Clear, 100pcs/bag	37-11493
Ribbon Tips, Clear, 100pcs/bag	37-11494

Mixing Tips for Cement

Works better, lasts longer, costs less

These mixing tips works for these core materials

- Adherence DC - septodont
- Beauticem SA - shofu
- Bifix QM - voco
- Breeze - pentron
- Cement it - pentron
- Clearfil esthetic cement ex - kuraray America
- G-cem linkace - GC America
- Geristore syringeable - denmat
- Natural elegance resi cement - h schein
- Ultra-bond plus - denmat



Brown Mixing Tips, 25/pkg 37-110171

Intra-Oral Tips, 25/pkg 37-110172

Mixing Tips for Core Materials

Works better, lasts longer, costs less

These mixing tips works for these core materials

- Clearfil DC core plus - kuraray America
- Compcore AF syringemix - premier dental
- Core material automix - h schein
- Core paste XP - denmat
- formcore DC - J morita
- Gradia core - GC America
- Multicore flow - ivoclar vivadent
- TI core auto E - essential dental systems
- TI core flow+ - essential dental systems



Brown Mixing Tips, 25/pkg 37-110161

Intra-Oral Tips, 25/pkg 37-110162

Mixing C&B Mixing Tips - Short

- Works better
- Lasts longer
- Cost less
- 100% compatible with current mixing system
- Improved mixing performance
- Latest technology in dental mixing tips
- Ratio 10:1



25/pkg

10-106691

Call Us Today 833-733-6825

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Mixpac Colibri Mix Tips (Pacdent)

Endodontical applications, core build-ups and precision impressions demand high standards of hygiene and accuracy. Colibri fulfills both of it. Mixers of the double syringe and cartridge systems have been selected in cooperation with dentists and combined with a needle made of medical stainless steel



It is not required anymore to attach an intraoral tip. Beside the existing advantages of the S- and L- system mixers, the Colibri offers a high level of flexibility and ductility but also a constant inner diameter. This allows not only a maximum material flow rate but also an easy accessibility even when working within a restricted intraoral situation and when an ergonomic application is needed.

Advantages

- High process safety
- All parts can be steam sterilized in the autoclave
- No risk of needle chipping off
- Time saving
- All parts can be steam sterilized in the autoclave

Brown 20ga, 40/pkg	37-107821
Brown 17ga, 40/pkg	37-107822
Light Blue 20ga, 40/pkg	37-107823
Light Yellow 17ga, 40/pkg	37-107824

Mixpac T-Mixer Mixing Tips

Innovative mixer design allows more compact shape reduced material waste volume

Advanced technology leads to more efficient material processing more homogeneous mixing results enhances application safety.



Fully compatible to MIXPAC S-System no change in working environment

Green 6.5mm, 40/pkg	37-106112
Pink 5mm, 40/pkg	37-106113
Yellow 4.2mm, 40/pkg	37-106111
Green/ Teal Intra-Oral Tips, Fit Green/ Teal Tips, 100/pkg	37-11612
Yellow Intra-Oral Tips, Fit Yellow Tips, 100/pkg	37-11613

Mojo II Single Use Disposable Syringe (Danville)



The material wasted inside a standard mixing tip (1.6cc) is enough to do three crowns with Danville's Mojo syringe (DMS) the small yet sensibly angled tip of our syringe gives you easy access even to the distals of second molars and allows the material to be applied right at the sulcus circumferentially, thus minimizing air bubbles and retakes. Avoid messy cleanups and save \$2 per impression or \$1,000 per year.

Regular 18g Tips, 100/pkg	37-13941
Endo 22g, 100/pkg	37-13942

Mixing Bowls

Alginate Mixing Bowls (Dentsply)

Made from extremely smooth and pliable plastic, for easy mixing and cleaning. PVD is not used in the manufacturing process, eliminating any danger of migrating plasticizers. Glossy inside surface, making alginate removal easy



Blue Dark	45-11221
Blue Ocean Splash	45-11208
Blue Pearl	45-11207
Blue Tropical	45-11209
Green Kiwi Splash	45-11210
Green Pearl	45-11211
Green Tropical	45-11212
Orange Tangerine Splash	45-11216
Orange Tropical	45-11215
Patriotic Splash	45-11214
Pearlescent	45-11213
Pink Tropical	45-11218
Red	45-11217
Splash	45-11219
Yellow Tropical	45-11220

Disposable Mixing Bowls



White Bowls, 36/pkg	37-99201
----------------------------	----------

Impression Material

Neon Blue Bowls, 36/pkg	37-99202
Neon Yellow Bowls, 36/pkg	37-99203
Neon Green Bowls, 36/pkg	37-99204
Neon Pink Bowls, 36/pkg	37-99205
Neon Purple Bowls, 36/pkg	37-99206
Assorted Bowls, 36/pkg	37-99207
White Spatulas, 36/pkg	37-99208
Neon Blue Spatulas, 36/pkg	37-99209
Neon Yellow Spatulas, 36/pkg	37-99210
Neon Green Spatulas, 36/pkg	37-99211
Neon Pink Spatulas, 36/pkg	37-99212
Neon Purple Spatulas, 36/pkg	37-99213
Assorted Spatulas, 36/pkg	37-99214

Mixing Bowls

Dark Green mixing bowls



Large 600cc, each	37-108801
Medium 350cc, each	37-108802
Small 150cc, each	37-108803

Mixing Bowls

- Flexible plastic material.
- Three different sizes to fit different jobs
- Washable for repeated use



Mixing Bowls, Set of 3: Large, medium, & small 37-11373

Plastic Spatulas (Dentsply)

Designed to fit all mixing bowls. Won't scratch the surface of bowls, and clean up easily. Perfect combination of flexibility and stiffness, for good control. Easy grip handle reduces slipping



Blue Autoclavable	45-13216
Ocean	45-13209
Blue Pearl	45-13208
Blue Tropical	45-13207
Kiwi	45-13212
Green Tropical	45-13210



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Orange Tropical	45-13213
Pearl	45-13211
Pink Tropical	45-13214
Red	45-13215
Yellow Tropical	45-13206
Splash	45-13207

Silicone Mini Bowls (Plasdent)



Medium 30cc, 2pcs/box	37-11397
Small 8cc, 3pcs/box	37-11396

Spatula

Sturdy, durable, and comfortable to hold on. Washable for repeated use



Spatula, 12/box	37-11398
-----------------	----------

Spectrum Flowabowl Mixing Bowls (Plasdent)



Spectrum Blue, large 600cc	37-11369
Spectrum Yellow, large 600cc	37-11370
Spectrum Green, large 600cc	37-11371
Spectrum Ruby Red, large 600cc	37-11356
Spectrum Purple, large 600cc	37-11368
Spectrum Pink, large 600cc	37-11367
Spectrum Orange, large 600cc	37-11366
Spectrum Blue, medium 350cc	37-11365
Spectrum Yellow, medium 350cc	37-11364
Spectrum Green, medium 350cc	37-11363
Spectrum Ruby Red, medium 350cc	37-11362

Impression Material

Spectrum Purple , medium 350cc	37-11361
Spectrum Pink , medium 350cc	37-11360
Spectrum Orange , medium 350cc	37-11359

Impression Accessories

Acumix (Plasdent)

1:1/2:1 , 50ml	37-11440
10:1/4:1 , 50ml	37-11439



Alginates Accessories (Zhermack)

Flexible Mixing Bowl - for alginate - rubber	37-11150
Alginates Measuring Set - Regular type 1x water measure/ 1x powder scoop	37-11151
Alginates Measuring Set - For hydrogum 5 1x water measure/ 1x powder scoop	37-11152
Mixing Spatulas - plastic, multicolor, autoclavable 6x spatulas	37-11153

C-Silicone Impression Material Accessories (Zhermack)

Mixing Spatula for silicone	37-11160
Mixing Pad , 20x sheets	37-11159
Putty Cut - Cutter for silicones	37-11158

Cartridge Dispensing Gun (Defend)

- Can be easy to automatic mixing and deliver the impression materials
- Used to extrude dental temporary crown material or dental impression materials and bite registration materials from cartridges.
- It consists of a release lever and a plunger, which are coordinated simultaneously to load and remove the cartridge
- Our dispensing gun can automatic mix and deliver any-brand dental silicone impression materials



1:1/2:1 , each	37-11437
-----------------------	----------




Like us on Facebook



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Disposable Impression Syringes (Plasdent)



Bendable Clear Tip , 50set/box	37-11192
Bendable 1 1/2" Longer Tip , 50set/box	37-11191

Freealgin Accessories (Zhermack)

D2 Dispenser (1:1)	37-11291
Replacement Plunger For Universal Dispenser , 50ml	37-11290
D3 Multi-Dispenser (1:1, 2:1, 10:1, 4:1) 1x applicator/ 1x 4:1 plunger/ 1x 1:1/ 2:1 plunger/ 1x 10:1 plunger/ 1x 1:1/2:1 breech/ 1x 4:1/ 10:1 breech	37-11289
Dynamic Mixers , 50x	37-11287
Bayonet Rings , 2x	37-11286
Universal Tray Adhesive , 1x 10ml bottle	37-11285

Hub & Intra Oral Tip (Plasdent)



Tip For GC Syringe Clear , 100pcs/bag	37-11314
Hub For Free-Flo Syringe	37-11311
Tip For Free-Flo Syringe Clear , 100pcs/bag	37-11312
Hub For GC Syringe	37-11313

Impression Cartridge Transfer Connector (Plasdent)



5pcs/bag	37-98532
----------	----------

Impression Material



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Impression Cartridges Connector (PacDent)


HP cartridges connector to transfer unused left-over materials to charge another cartridge



5/pkg

37-98531

Impression Guns Holder (PlasDent)



Upright Triple, 26-11392
(11"W x 7 1/4"H x 6 1/2" D)

Upright Dual, 26-11393
(7 1/4"W x 7 1/4"H x 6 1/2" D)

Straight Triple, 26-11394
(10 7/8"W x 5"H x 6" D)

Straight Dual, 26-11391
(7"W x 5"H x 6" D)

Impression Injector

Clear & Adjustable 37-11347
Tip, 50sets/box



Impression Organizer



(12 1/4"W x 9 1/2"H x 6" D)

26-11395



Impression Protector (Dentsply)

This elastomer strip that fits over the brackets covering up areas the alginate normally adheres to. The impression can be removed easily with minimal damage. Each strip is made to fit one full arch, but can be cut for patients who are not fully bonded. They are individually wrapped to prevent cross-contamination



50 sets of two

45-12452

Intra Oral XX-Fine Mixing Tips

50/pkg 37-105391



Mixpac Dispenser Gun (PacDent)

Elegant, simplified workflow self explain handling, outstanding ergonomics.



Reduced dispensing force fatigue-free usage and a precise work, even with high viscosity materials.

1:1/2:1 cartridges 37-105721

4:1/10:1 cartridges like crown & bridge 37-105722

Mojo II Syringes (Danville)

Disposable, single-use syringes that don't waste impression material. They have locking and storing capabilities that other material syringes don't have. The locking positions are: fill, store, and activate. The 18g tip syringes work with 1:1 impression cartridges only.




100/pkg

36-106401

Tac Tray Adhesive (Bosworth)

3.5oz, spray can

37-108461



Impression Material

Tray Adhesive

A paint on tray adhesive for use with all alginates, it helps to keep the alginate firm in the tray during removal from the mouth

The tray adhesive comes with convenient brush inside the cap.

Amber Tray, 2oz. bottle

37-11419

Blue Tray, 2oz. Bottle

37-11420




Universal Tray Adhesive (Zhermack)



1x 10ml bottle

37-20003

www.PrestigeDentalProducts.com



**Check our
website for
more items
and
variation**



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Find Us On Social Media

twitter 



Find us on
Facebook



Pinterest

Prestige Dental Products



**Excellent Qualities
at Affordable Prices**

877-772-3888

www.prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

