

Waste Management



Prestige Dental Products, Inc

Prestigedentalproducts.com

Tel: 877-772-3888 Fax: 714-632-8821

Waste Management

Amalgon Pre Paid Postage Amalgam Recycling (WCM)

Amalgon® is the most convenient way to properly manage amalgam/amalgam contaminated items. Amalgon® follows the ADA's establishes BMP (Best Management Practice) for the management of amalgam contaminated items. **Amalgon® is a mail back (pre paid postage) recycling & disposal product;** after you fill the container with amalgam contaminated waste, simply attach the pre paid shipping label and give the Amalgon® container to your UPS driver. The Amalgon® container is shipped to our authorized recycling center where the contents are promptly recycled.



1 Gallon Capacity 12-11124

Chemgon Fixer and Developer Disposal System

When equal amounts of developer and fixer are drained into the 5 gallon Chemgon container, the liquid mixes with Chemgon's special patented dry ingredient causing the contents to solidify. When full, simply throw the Chemgon container away with your normal trash.



5 Gallon (each) 59-11138

Eco-Star Amalgam Separator (Tech-West)

The Eco-Star Amalgam Separator systems are on the leading edge of technology in design and performance

- Designed to serve up to 10 chairs
- Easy filtration of amalgam, mercury, alloy and other waste particles
- Long filter life




Amalgam Separator 12-111841
Cartridge Filter 12-111842



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Foilgon Pre Paid Postage Lead Foil Recycling (WCM)

Includes storage container, prepaid UPS shipping label and recycling. Foilgon Lead Foil Recycling is a convenient way to properly manage lead foil



7lb. Capacity (each) 59-11208

Hg5 Amalgam Separation System (Solmetex)

Hg5 systems provide is outstanding performance, capturing amalgam and sediment from dental waste streams.



- Professionally engineered, robust, simple design
- Easy to operate and maintenance free
- Flexible for wet or dry vacuum systems
- System is clear, making it possible to visualize functionality
- Simple collection container change process include no tools, no mess; packaging for delivery to certified recycler, recycle certificate issued to practice, no middleman required for disposal, no additional charges for shipping to recycling facility
- No contracts w/ hidden fees
- No electricity required
- No decanting, pumps, timer

Warranty: 2 years

1 to 4 operatories 12-25412
1 to 10 operatories 12-25413
11 to 20 operatories 12-25414
Replacement Collection Container/ 12-25415
Recycle Kit: Complete kit
Recycle Kit Only: Mail in program, 12-25416
 shipping to the recycler, recycling and documentation

Call Us Today 833-733-6825

Waste Management



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Isolyser SMS Self Disposal (WCM)

Using Isolyser/
SMS is a time &
money saving
way of
managing sharps
waste. It
encapsulates
sharps waster
material after
our proprietary
catalyst packets
are added to our
patented odor-free liquid NMA base ingredient.



3 Liters

10 Liters

After Curing, the waste is solidified, rendered unrecognizable and unusable; making it safe to dispose of in your regular trash in most states

3 Liters	13-11162
10 Liters	13-11160-1

Medentotainer Waste Amalgam Solution (Medentex)

Specifically
designed for
the proper
hygienic
disposal and
recycling of
dental waxes,
Including

amalgam, lead foils, and
extracted teeth containing
amalgam. Medentex features a
full service disposal, recycling,
re-shipment, and creation of a
recycling certificate with our
medentotainer. Customers can
be assured that all of the waste
medentex receives will not end up in a landfill. Medentex recycles
100% of all amalgam waste at its own recycling facility



X-large 6.5, gallon container capacity	12-98691
Large 4.0, gallon container capacity	12-98692
Small 1.5, gallon container capacity	12-98693
Tooth box 0.35, gallon container capacity	12-98694
Cap guard 1.5, gallon container capacity	12-98695

Medentotainer Waste Amalgam Solution (Medentex)

For x-ray lead foils from dental practice

Lead Foils



1.5, gallon container capacity	59-98701
---------------------------------------	----------

Monoject Chimney Top Sharps Containers (Covidien)

- Polished viewing on 4 sides
- Translucent **RED**
- Clear lid w/ extra wide butterfly closure
- Can be autoclaved or incinerated



Small, 4 qt. size 6.75" W x 10.56" L x 7.08" H, each	13-89431
Medium, 8 qt. size 6.75" W x 10.56" L x 10.89" H, each	13-89432
Large, 14 qt. size 6.75" W x 10.56" L x 17.77" H, each	13-89433

Sharps Containers (Crosstex)

Dual access, clear top provides
visualization of collector
fill-levels



1.5qt., red each	13-88281
-------------------------	----------

Waste Bags



10AG Bio Hazard Waste Bag Red, (250pcs/box)	39-11371
16AG Bio Hazard Waste Bag Red, (250pcs/box)	39-11373
Bio Hazard Waste Bag 9" x 10", (200pcs/box)	39-11372
Bio Hazard Bag 6" x 6", (200pcs/box)	39-11374
Utility Waste Bag 9" x 10", (200pcs/box)	39-11375
Utility Waste Bag 6" x 6", (200pcs/box)	39-11376



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

