

Prestige Dental Products

Call: (877) 772-3888

Fax: (714) 632-8821


Telescopic Loupes

Extra wide field of view. Large depth of focus. High light efficiency. Advanced mounting system with 3 adjustments for individualizing posture, working angles and to achieve perfect eye alignment. Lightweight frames with soft silicone nose pad. Adjustable head strap. Removable side shields.




Magnification	Working Distance	Field of View
250	17"	4.5"
300	17"	3.5"
350	17"	2.5"


Centering Control



Rotate Center to adjust Pupil Distance



Adjust Angle Lens



Adjust Lens & Glasses Distance

