

MICROMAX™

Cordless Prophylaxis Handpiece




- Cordless convenience offers complete portability.
- Heavy-duty high torque brushless motor.
- Forward/reverse functions.
- Reciprocating setting available without use of additional attachment.
- Five rotation speeds available for different applications.
- Ergonomic design provides comfort and easy operation.
- High capacity lithium-ion battery provides stable and prolonged usage.
- Offers quick and easy coupling.
- Low battery warning light.
- Four different attachments available to choose from.

DENTAMERICA®

MICROMAX™

MICROMAX™

CORDLESS PROPHYLAXIS HANDPIECE



Optional Attachments



Latch Attachment



Screw Attachment



Snap-on Attachment



Straight Attachment



SPECIFICATIONS:

INPUT	AC 100V~240V 50/60Hz
FUSE	500mA 250V
BATTERY TYPE	3.7V 1620mAh
DIMENSIONS (HANDPIECE)	L192 xW26 xH29 mm
DIMENSIONS (CHARGING STATION)	L105xW85xH59mm
WEIGHT (HANDPIECE)	150g
WEIGHT (CHARGING STATION)	125g

Order No.

- 177401 Latch type
- 177402 Screw type
- 177403 Snap-on type
- 177408 Straight type

DENTAMERICA®

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.asia

DISTRIBUTOR

Prestige Dental Products

3031 E. Coronado St. #D
Anaheim, CA 92806
714-632-8817
www.prestigedental.com