

# Autoclaves & Sterilizers



**Prestige Dental Products, Inc**

**PrestigeDental.com**  
**877-772-3888**

# Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

**As an authorized Dealer of many brands in North America**, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

**Prestige Dental Products** offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

**Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and**

## Customer Information

### Order

**PHONE** To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

**FAX** You may FAX your orders to 714-632-8821.

**INTERNET** [www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



*Note:* New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



### Customize order for special customer

Email us to [mass@prestigedentalproducts.com](mailto:mass@prestigedentalproducts.com) a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

## Payment

**CREDIT CARD** We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

### Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

### MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

## Shipping

### GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

**EXPEDITED SHIPPING** UPS Next Day and Second Day Service for which additional shipping fees will be charged.

**SHIPPING HAZARDOUS ITEMS** All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

## State and federal laws

### DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog  
Are Subject To Change Without Prior Notice.**



# Autoclaves & Sterilizers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Autoclave, Sterilizer Unit (Tuttnauer 2340M)

### Tuttnauer 2340M Manual Autoclave, Sterilizer



Manual autoclave is designed to complement any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2340M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.

#### Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications
- **1 year warranty**

**Tuttnauer 2340M** 30-11145 **Call for Price**

## Statim 5000 Automatic Autoclave (SciCan)



Fully automatic. Simply select one of four sterilization programs and with a single touch of the keypad, begin the cycle. An advanced microprocessor controls all functions and maintains optimal sterilization conditions, from start to finish.

SciCan's Statim cassette sterilizers have a biological effectiveness that has been proven in tests for hollow and solid instruments at many internationally recognized institutes in Canada, the United States and Europe.

Cycle times do not include drying

#### Specifications:

- Unit Size: L 21.75" x W 16.25" x H 7.5"
- Cassette Internal Dimensions: L 15" x W 7" x H 3"
- Reservoir Capacity: 4.0 litres (distilled water) Approximately 40 cycles
- Weight without water: 33 kg / 72 lbs
- Power Consumption: 110-120 V, 60Hz, 1300 W

## Optional Stat printer:

- Type: Thermal Printer
- Print: 20 characters/line
- Speed: 1 line/sec
- Capacity: 80 cycles/roll

**Unit** 30-26011 **Call for Price**

## Dry Heat Sterilizer (CPAC)



### The Only F.D.A. registered dry heat sterilizer of its kind.

It takes about 15-30 minutes from normal room temperature it reaches desired sterilizing temperatures 160° to 180°C (320° to 360°F).

In addition to that--allow one(1) hour to sterilize at 320f.

Avoid sterilizing items likely to melt during heat process (like rubber, plastic, etc).

<b>Unit with 2 Draws, 110V</b>	30-92391	<b>649.95</b>
<b>Unit with 2 Draws, 220V</b>	30-92392	<b>695.95</b>
<b>Unit with 3 Draws, 110V</b>	30-92393	<b>795.95</b>
<b>Unit with 3 Draws, 220V</b>	30-92394	<b>895.95</b>

## Prestige Cassette Sterilizer



Stylish aesthetics, together with a host of enhanced features, build on the successful attributes of the original Series 2100 range, to ensure the 2100 'Classic' perfectly complements today's modern clinics and the exacting needs of the health-care professional.

Prestige Medical 210025 Dental Cassette Sterilizer, as an ISO 9001-2000 / EN13485 accredited company ensures all products including the 2100 Classic have been designed to



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

products including the 2100 Classic have been designed to the latest international standards for safety and performance and are fully CE marked to the Medical Devices Directive 93/42/EEC.

#### Specifications:

- Chamber Size: 9" (H) X 8" (W) 0" Deep
- Automatic (electronic) controlled
- Steam Sterilizer Power Supply: 110V
- 16 minutes cycle
- Table Style

**Unit** 30-11306 **1199.00**

#### Rapid COX Dry Heat Sterilizer (CPAC)



#### 6 Minutes cycle Low maintenance No supplies needed

- You can adjust the time to: 6 Minutes cycle time for un-wrapped Instrument, or 8 Minutes cycle time for Hand pieces, or 12 Minutes cycle time for wrapped Instrument
- Contaminated Implants or those needing alterations can be re-sterilized fast
- Avoids added risk and time of making new appointment and re-sedating patient if alteration is needed
- Quick cycle times which reduce needs for increased instrument inventory
- Excludes corrosion, rust and dulling associated with steam sterilizers
- No foul odors, chemical residues or possible health hazards attributed to chemical sterilizers
- Scheduling of your patients at your convenience instead of waiting for sterilization of long cycle sterilizers

**Unit, 110V** 30-92381 **3995.95**

**Unit, 220V** 30-92382 **4095.95**

## Autoclaves & Sterilizers

#### Dry Heat Sterilizer (DA)

- Cost efficient
- Sterilize trays with lid on
- Automatic timer
- Even heat, rapid heat-up
- Simple to operate
- Burnout-proof heating element
- Precision thermostat
- No pressure, no moisture



- Chamber sealing door
- Stainless steel, rust resistant
- 1-year warranty
- 2-Tray Sterilizer Dimensions, Outside: 15 1/2" W x 10 1/2" H x 9 1/2" D, Inside: 12 1/2" W x 5 1/2" H x 7 1/2" D, Tray: 10" x 6 1/2"

**2 Trays Unit Model 200, 450 Watts** 30-92401 **699.95**

**3 Trays Unit Model 300, 650 Watts** 30-92402 **899.95**



## BEST QUALITY DENTAL SUPPLIES



# Autoclaves & Sterilizers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Manual Kwiklave Sterilizer (Tuttnauer 2540MK)



The Tuttnauer 2540MK Series is also available in a quick cycle model, the Kwiklave (MK). This sterilizer offers all the proven, dependable features of our manual autoclave, with the added benefit of reducing the total heat up time. This autoclave is the perfect choice for quick instrument turn-around or a high volume office. Speed and reliability are accomplished without sacrificing load size. 1 year parts and labor warranty. Model 2540MK requires a voltage between 220V and 235V for proper operation.

### Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certification
- **1 year parts and labor warranty**

**Tuttnauer 2540MK** 30-11274 **3749.95**

## Large Capacity Fully Automatic Autoclave/Sterilizer (Tuttnauer 3870EA)



Our largest automatic table top autoclave. Features a closed

door active drying system. Perfect for multi-doctor practices, clinics or surgery centers. 15"x 30" chamber accommodates longer instruments or a large number of small instruments. Capacity for 15 full size IMS cassettes. 2 year parts and labor warranty.

### Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- Hepa filter 0.2 µm air filter (replaceable)
- Rapid cycles optimize the availability of ready-to-use sterile instruments
- Easy-to-use, single touch control system increases productivity
- 3 Pre-set sterilization programs
- Program and customize cycle parameters to suit your needs and maximize flexibility
- Integrated printer for cycle documentation (optional)
- Interdependent temperature and pressure monitoring
- Cycle information recovery in the case of power failure or cycle interruption
- Insufficient water detector with automatic shut-off

**Tuttnauer 3870EA** 30-11256 **Call for Price**

## Large Capacity Manual Autoclave/Sterilizer (Tuttnauer 3870M)



Our largest manual table top autoclave. Perfect for multi-doctor practices, clinics or surgery centers. 15" x 30" chamber accommodates longer instruments or a large number of small instruments. Capacity for 15 full size IMS cassettes. 1 year parts and labor warranty.

### Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications
- **1 year parts and labor warranty**

**Tuttnauer 3870M** 30-11257 **Call for Price**

# Autoclaves & Sterilizers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Tuttnauer Autoclave/Sterilizer (Tuttnauer 2540M)



The 2540M Tuttnauer manual autoclave is designed to complement any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2540M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.

### Features:

- Automatic shut off at the end of both the sterilization and dry cycles
- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- Dual safety thermostat to protect against overheating
- International certifications

**•1 year parts and labor warranty**

**Tuttnauer 2540M** 30-114111 **Call for Price**

## Valueklave 1730 Sterilizer (Tuttnauer)



- The Valueklave 1730 Autoclave will satisfy all of your basic sterilization needs at an economical price
- The compact size (**7" chamber**) makes it a perfect primary unit for small offices or a second autoclave for larger practices
- The fast cycle time** and heated drying delivers a quick turnaround for wrapped or unwrapped instruments
- The simple controls and easy closing door make this unit a breeze to operate
- High speed** and reliability are built into every Valueklave 1730

**Valueklave 1730** 30-114141 **1949.95**

## Delta Q Sterilizer (Pelton & Crane)



- Closed door filtered air drying ensures a higher quality of sterilization
- Sterilization cycle time of 15 minutes for unwrapped instruments
- Automatic chamber preheat
- Electronic temperature and time control monitor sterilization conditions in chamber
- One custom and four automatic preprogrammed sterilization settings
- Approved for hand piece sterilization
- Digital function display allows easy and accurate monitoring
- Designed and built to ASME specifications

Specifications	Delta Q 8	Delta Q 10
<b>Chamber Size</b>	8 1/4" x 14" (21 x 36cm.)	9 7/8" x 17 1/2" (25.1 x 44.5cm.)
<b>Chamber Capacity</b>	739cu.in. (12,110cu.cm.) Holds 2 European trays	1340cu.in. (22,019cu.cm.) Holds 3 large or 9 small Hu-Friedy cassettes or 3 European trays
<b>Start-Up Time (to 134°C)</b>	10-13min. Warm	12min. Warm
<b>Sterilization Time</b>	16min. warm (unwrapped)	16min. warm (unwrapped)
<b>Drying Cycle Time</b>	30min. preprogrammed, variable 0-99min.	30min. preprogrammed, variable 0-99min.
<b>Outside Dimensions (with stand)</b>	17 1/2" W x 12 1/4" H x 19 1/4" D (44 x 31 x 50cm.)	19 1/2" W x 14" H x 23" D (49 x 36 x 60cm.)
<b>Amp Circuit</b>	Requires dedicated 20-amp circuit	

**Delta Q 8, 8 in. 115V, each** 30-13201 **4694.95**

**Delta Q 10, 10 in. 115V, each** 30-13203 **4794.95**

**Delta Q 10, 10 in. 230V, each** 30-13204 **4794.95**



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Autoclaves & Sterilizers

### Elara 11 Class B Automatic Autoclave Sterilizer (Tuttnauer)

The Elara 11 Sterilizer is designed for sterilization of wrapped, porous and hollow loads. The chamber volume is large, external dimensions are small, making it especially suitable for smaller clinics requiring a high turnover. The control system is smart and intuitive, easy to use and has a large memory for consistent and reliable performance.

#### Efficient:

- Short cycle and high volume pre and post vacuum pump for perfect class B sterilization results
- 316L grade stainless steel chamber and door with electro-polish finish
- A built-in steam generator provides steam and facilitates a standby steam mode
- The door closes smoothly and locks automatically

#### Features:

- Large 7.5 Gallon chamber within a small foot print. A perfect fit for standard 24 inch shelves or counters
- Provided with 5 stainless steel trays and can also hold cassettes
- Water filling from the front economizes on office space by eliminating the need to keep the top of the autoclave accessible

#### Safety Features:

- A cycle will not start unless the door is properly locked
- An automated double locking safety device prevents the door from opening while the chamber is pressurized
- An automatic safety shut-off prevents overheating of the chamber
- The autoclave comes with safety pressure-release valves for the chamber and generator



**Elara 11 Class B**

30-11190

**Call for Price**

### EZ Automatic Autoclave/Sterilizer (Tuttnauer)

#### Fully Automatic Autoclaves EZ9TM &

#### EZ10

With the simplicity of one touch design, all your sterilization and drying needs are fulfilled. The Tuttnauer fully automatic autoclaves fill, sterilize, exhaust and dry at the touch of a button. FORGET IT ONCE YOU SET IT.



#### Features:

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- High efficiency air pump circulates hot air after the sterilization cycle for fast drying
- Hepa filter 0.2µm air filter (replaceable)
- Pre-programmed cycles for: unwrapped instruments; wrapped/packs; glassware and additional drying cycle
- All program parameters can be customized and stored
- Automatic shut off at the end of both the sterilization and dry cycles
- Audible and visual cycle interruption alert
- Double safety locking device prevents door from opening while chamber is pressurized
- Low water sensor prevents activation of sterilization cycle when there is insufficient water in the chamber or reservoir.
- Constant monitoring of temperature and pressure for sterilization
- Power outage recovery system
- Door design prevents steam from coming in contact with control panel
- Optional internal printer documents: date, time, temperature and pressure
- 510k clearance for dental hand pieces
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- International certifications

**• 2 year parts and labor warranty**

#### Fully Automatic Kwiklave™ EZ 10K

The Tuttnauer EZ series is also available in a quick cycle model, the Kwiklave (EZ10k). This sterilizer offers all the proven, dependable features of our fully automatic autoclaves, with the added benefit of reducing the total heat up time. This autoclave is the perfect choice for quick instrument turnaround or a high volume office. Speed and reliability are accomplished without sacrificing load size. 2 year parts and labor warranty. Model EZ10k requires a voltage between 220V and 235V for proper operation.

**Tuttnauer EZ10**

30-11205

**Call for Price**

**Tuttnauer EZ9**

30-11204

**Call for Price**

**Tuttnauer EZ10K**

30-11206

**Call for Price**

## Autoclaves & Sterilizers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888



[www.prestigedental.com](http://www.prestigedental.com)

# Order Now



877-772-3888



2340 M  
No Chemical

We Carry Brands Like:

**Tuttnauer**  
Your Sterilization & Infection Control Partners

## Buy any EZ Autoclave...



**Closed Door System Fully Automatic Autoclaves**

- Hepa Filtered Air Drying • Cassette Capacity 3 large + 3 half

**EZ Models come with a 2 Year Parts and Labor Warranty**

**cpac**  
Equipment Inc.

And More!!!



Dry Heat Sterilizer

Best  
Prices!!!





**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

#### ☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

#### ☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

