

# Gypsums



**Prestige Dental Products, Inc**

**PrestigeDental.com**  
**877-772-3888**

## Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

**As an authorized Dealer of many brands in North America**, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

**Prestige Dental Products** offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

**Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and**

### Customer Information

#### Order

**PHONE** To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

**FAX** You may FAX your orders to 714-632-8821.

**INTERNET** [www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



*Note:* New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



#### Customize order for special customer

Email us to [mass@prestigedentalproducts.com](mailto:mass@prestigedentalproducts.com) a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

#### Payment

**CREDIT CARD** We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

#### Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

#### MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

#### Shipping

##### GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

**EXPEDITED SHIPPING** UPS Next Day and Second Day Service for which additional shipping fees will be charged.

**SHIPPING HAZARDOUS ITEMS** All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

#### State and federal laws

##### DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog  
Are Subject To Change Without Prior Notice.**



# Gypsums



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Magnum Die Pink Stone (MDC)

### Magnum Nic Stone®

#### Class IV

Manufactured under ISO and ADA specifications  
Also ISO 9001 - 2008 standar.



#### CHARACTERISTICS

- Class IV super stone
- Regular setting
- Extra fine powder
- Optimal hardness

#### ADVANTAGE

- Smooth surface
- Creamy smooth rock
- High detail
- Low setting expansion
- Dimensionally stable

#### PHYSICAL PROPERTIES

Setting expansion 0.15%  
Pressure resistance 9500 psi  
Setting time 15 ± 2 min.  
Water/Powder ratio 20 cc / 100 gr



#### Presentation

Box 22 lbs / 100gr

#### Code

22DP02

Die Pink Stone, 22 lbs.

34-88691 24.95

## Magnum Orthodontic Stone (MDC)

### Magnum® Ortho stone

#### Class II

Snow white  
Specially designed for orthodontic models  
Super smooth and creamy



#### CHARACTERISTICS

- Class III super stone
- Regular setting
- Extra fine powder
- Optimal hardness

#### ADVANTAGE

- Smooth surface
- Creamy smooth gypsum
- High detail
- Low setting expansion
- Dimensionally stable

#### PHYSICAL PROPERTIES

Setting expansion 0.15%  
Pressure resistance 8000 psi  
Setting time 10 ± 2 min.  
Water/Powder ratio 32 cc / 100 gr



#### Presentation

Box 22 lbs / 10 kg

#### Code

22OP02

Magnum Orthodontic Stone, 22 lbs.

34-11120 22.95

## Magnum Plaster White Stone (MDC)

### Magnum Plaster

#### Class II



#### CHARACTERISTICS

- Class II plaster
- Regular setting
- Extra fine powder
- Optimal hardness

#### ADVANTAGE

- Smooth surface
- Creamy gypsum
- High detail
- Low setting expansion
- Dimensionally stable

#### PHYSICAL PROPERTIES

Setting expansion 0.20%  
Pressure resistance 5000 psi  
Setting time 6 ± 2 min.  
Water/Powder ratio 24 cc / 100 gr



#### Presentation

Box 22 lbs / 10 kg


#### Code

20PW02

Creamy White Stone, 22 lbs.

34-11121 22.95

## Model Carriers



### Carrier Set:

1 tray w/2 foams, 300 set/case

34-98041 194.95

### Foams for model carriers,

Charcoal, 1000/Pkg

34-98042 99.95

## Magnum Yellow & Blue Stones (MDC)

### Magnum®

#### Stone class III

Manufactured under ISO and ADA specifications  
Also ISO 9001 - 2008 standard.



#### CHARACTERISTICS

- Class III super stone
- Regular setting
- Extra fine powder
- Optimal hardness

#### ADVANTAGE

- Smooth surface
- Creamy smooth gypsum
- High detail
- Low setting expansion
- Dimensionally stable

#### PHYSICAL PROPERTIES

Setting expansion 0.15%  
Compression resistance 390 kg/cm2  
Setting time 10 ± 2 min.  
Water/Powder ratio 24 cc / 100 gr

#### Presentation

Box 22 lbs / 10 kg

#### Yellow

26Y02

#### Blue

22B02

Yellow, 22 lbs.

34-11122 21.95

Blue, 22 lbs

34-11123 21.95

## Disposable Model Base Formers (Dentsply)



Saves clean up time. Assures creation of consistent models.  
Durable plastic will not deform during model creation.  
Flexible plastic for easy separation.

### 60/Pkg,

20 each small, medium and large

34-11111 29.95

### 120/Pkg,

40 each small, medium and large

34-11112 54.95





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Green Rubber Base Formers for Plastic Articulators (Meta Dental)



<b>Q-Left</b>	34-11117	<b>5.95</b>
<b>Q-Right</b>	34-11116	<b>5.95</b>
<b>Small</b>	34-11115	<b>5.95</b>
<b>Medium</b>	34-11113	<b>5.95</b>
<b>Large</b>	34-11114	<b>5.95</b>
<b>X-Large</b>	34-11118	<b>5.95</b>

### Extend-a-Pour (Dux)



<b>12 fluid oz, each</b>	41-90521	<b>13.95</b>
--------------------------	----------	--------------

### Green Silicon Rubber Base Formers (Meta Dental)




Base formers made of silicon rubber (Left, Right, Small, Medium, Large). Simply pour stone, level off, then place the teeth cast on the form Easy to use and reusable

[TOP: Medium full arch; Large full arch]  
 [BOTTOM: Small full arch; Right & Left Quadrant]

<b>Set of 5/Box</b>	34-11119	<b>19.95</b>
---------------------	----------	--------------

## Gypsums

### Model Formers (PlasDent)



<b>Disposable Large/Deep, Clear,</b> 60 pcs/box	34-11132	<b>18.95</b>
<b>Disposable Large, Clear,</b> 60pcs/box	34-11131	<b>15.95</b>
<b>Disposable Medium Clear,</b> 60pcs/box	34-11130	<b>13.95</b>
<b>Disposable Asst Sizes 20 of</b> <b>each size, 60pcs/box</b>	34-11129	<b>15.95</b>
<b>Disposable Asst Sizes 40 of</b> <b>each size, 120pcs/box</b>	34-11128	<b>29.95</b>
<b>Large Blue,</b> 3pcs upper & 3pcs lower/box	34-11127	<b>9.95</b>
<b>Medium Green,</b> 3pcs upper & 3 pcs lower/box	34-11126	<b>8.95</b>
<b>Small Yellow,</b> 3pcs upper & 3pcs lower/box	34-11125	<b>7.95</b>
<b>Asst Sizes, 6pcs/box</b>	34-11124	<b>8.95</b>

### Gleco Plaster Trap (IndiGreen)



Trap System Invented by a plumber, the Gleco Trap brings revolution to the age-old chore of sink trap retrieval and plaster trap maintenance. No more leaky bags or bulky buckets!

Install the pre-assembled Gleco Trap once, then replacing the bottle takes



## Gypsums

less than one minute start to finish without tools. Simply turn side valve to drain off excess water, unscrew full bottle, and screw in a clean one.

- Side drain siphons out all water prior to removal, eliminating spillage
- Translucent bottle shows waste level.
- Adapts easily to existing plumbing — the Gleco Trap is the only plaster trap to meet US Uniform Plumbing Code!

System includes one trap assembly, two replacement bottles, all necessary 1 1/2" pipe components and simple instructions for a five-minute installation

<b>64 oz</b>	34-87931	<b>94.95</b>
<b>19 oz</b>	34-87932	<b>89.95</b>

### Gleco Trap HV System (IndiGreen)



Item #10012 64 oz. Gleco Trap System    Item #11010 3.5-Gallon Gleco Trap HV Kit  
Item #10010 19 oz. Gleco Trap System    Item #11011 5-Gallon Gleco Trap HV Kit

From the inventor of the original Gleco Trap, the Gleco Trap HV Sanitary High-Volume Kit brings the same clean, dry and simple maintenance to full-size plaster traps. Never again will you soak your cabinets and yourself with foul-smelling, contaminated wastewater when you change your plaster trap container.

The Gleco Trap HV features a specially-engineered air-tight lid with threaded inlet and outlet extensions that reserve an air pocket inside. When the trap is full, simply remove the air-release plug to drain all inlet and outlet plumbing into the bucket before disconnecting. Its patented design prevents the mess and contamination that occurs with the removal of traditional plaster traps.

Instead of leak-prone flexible hoses found on other traps, the Gleco Trap HV uses standard 1 1/2" threaded inlet and outlet ports to ensure a positively leak-proof assembly. Adapts easily to existing plumbing.

Each kit includes a heavy-duty container, all 1 1/2" PVC plumbing components to fit most sink configurations and instructions.

<b>1 gallon, Refill (4)</b>	34-87945	<b>36.95</b>
<b>32 oz. Refill (12)</b>	34-87946	<b>36.95</b>
<b>64 oz. Refill (6)</b>	34-87947	<b>36.95</b>
<b>GT-19 w/19oz. Bottle</b>	34-87948	<b>92.95</b>
<b>GT-2 w/64oz. Bottle</b>	34-87950	<b>96.95</b>
<b>Replacement Gasket Set</b>	34-87950	<b>7.49</b>



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Trap-Eze Disposable Plaster Traps (Buffalo)

#### Economical trap systems designed to slow water discharge flow

- Requires minimum 9" clearance between bottom of sink tail pipe and top of trap bucket for easy installation without tools
- All trap buckets feature high-strength plastic construction and permanently-sealed lids
- 3 1/2 gallon size (12 1/2" diameter x 11" H) ideal for intermediate applications
- 5 gallon size (12 1/2" diameter x 14 1/2" H) best for larger offices and labs with higher volume collection needs (space permitting)
- Two flexible hoses (1 1/2" ID x 1 3/4" OD x 20" long) conveniently connect trap to sink tail pipe and sewer tail pipe
- Flexible adaptor fitting for sewer tail pipe included for improved drain seal
- Complete kit includes sealed trap bucket, 2 hoses, and sewer tail pipe adaptor fitting
- Refill includes trap and two plastic caps for sealing full trap
- Fast and easy installation and trap replacement
- Made in the USA



#### 3 1/2 Gal. Trap-Eze Complete Kit

(11" Height Sealed Bucket, Hoses) 34-11133 **46.95**

#### 3 1/2 Gal. Trap-Eze Refill

(11" Height Sealed Bucket) 34-11134 **20.95**

#### 5 Gal Trap-Eze Complete Kit

(14.5" Height Sealed Bucket, Hoses) 34-11135 **50.95**

#### 5 Gal. Trap-Eze Refill

(14.5" Height Sealed Bucket) 34-11136 **22.95**

#### 3 1/2 Gal. Trap-Eze SS Self-Sanitizing Trap Complete Kit New

34-11137 **54.95**

#### 3 1/2 Gal. Trap-Eze SS Self Sanitizing Trap Complete Kit-New

34-11138 **28.95**

#### 5 Gal Trap-Eze SS Self-Sanitizing Trap Complete Kit-New

34-11139 **57.95**

#### 5 Gal Trap-Eze SS Self Sanitizing Trap Refill- New

34-11140 **30.95**

#### Hoses only for Trap-Eze, Set of 2

(18" Overall Length Each) Flexible 34-11141 **35.95**

#### Drain Fitting for Trap-Eze Exhaust Hose

34-11142 **4.95**

#### Trap-Eze BT Bottle Trap Complete Kit ( 1 1/2" Plumbing + One 64 oz. Bottle

34-11143 **73.95**


#### Trap-Eze BT Bottle Trap 64 oz. Refills w/Lids & Bibs, Pkg of 2

34-11144 **19.95**



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Gypsums



**IndiGreen**  
INNOVATIONS

Brighter Ideas for Better Dentistry



**Buffalo**  
Buffalo Dental Manufacturing Co., Inc.



Call Us To Place Your  
Order Or Visit Us  
Online  
[www.prestigedental.com](http://www.prestigedental.com)

# Order Now



Call us for Quarterly Promotion

[Prestigedentalproducts.com](http://Prestigedentalproducts.com)



**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	





## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

