

Handpieces



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**



Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Fiber Optic High Speed Handpieces

Chuck Style Fiber Optic Handpiece (Dentex)

Chuck Style Fiber
Optic SINGLE
PORT WATER
SPRAY Handpiece
(Dentex) , Fixed
Back End, 4-Hole and 5-Hole Connector



Standard Head Fiber Optic	35-111771	162.00
Mini Head Fiber Optic	35-111781	162.00

Chuck Style Fiber Optic Handpiece with Triple Port Water Spray (Dentex)

Chuck Style
Handpiece 4-
Hole and 5-Hole
Connector Fiber
Optic (Dentex),
Triple Port Water Spray, Fixed Back End



Standard Head Fiber Optic	35-111801	243.00
Mini Head Fiber Optic	35-11179	243.00

Gemini Fiber Optic Handpiece for Kavo (Johnson-Promident)



Handpiece for KaVo Super Torque 625/630/640 (large head)	35-25851	349.95
Handpiece for KaVo Bella Torque 632/642/643 (small head)	35-25852	349.95

Integrity Leverage Fiber Optic Handpiece with Lever Cap (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best high-speed handpiece performance and value at a fraction of the cost of the top leading brands.



- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip

Fiber Optic High-Speed Handpiece	35-117811	289.95
---	-----------	---------------

Integrity Leverage Fiber Optic Handpiece with Push Button Cap (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best high-speed handpiece performance and value at a fraction of the cost of the top leading brands.



- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip Grip design

Integrity High-speed Handpiece w/ Push Button Cap	35-117821	239.95
---	-----------	---------------

Integrity Leverage Fiber Optic Highspeed with Lever Cap (Osseo Scientific)



Integrity Leverage Fiber Optic High-speed w/ Lever Cap	35-117831	299.49
--	-----------	---------------

Push Button Fiber Optic Handpiece (Dentex)

Push Button Fiber
Optic SINGLE
PORT WATER
SPRAY 4-Hole
and 5-Hole Connector Handpiece
(Dentex) , Fixed Back End



Standard Head Fiber Optic	35-119391	208.49
Mini Head Fiber Optic	35-119381	208.49

Push Button Fiber Optic Handpiece Interchangeable with Kavo Coupler (Dentex)

Push Button Fiber
Optic Handpiece Inter-
changeable with Kavo
Coupler (Dentex),
Swivel Handpiece



Standard Head Fiber Optic	35-119401	324.00
Mini Head Fiber Optic	35-119411	324.00



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

KAVO Interchangeable Fiber Optic Handpiece (Sable)



**Sable ACCESS S4500KL KAVO In-
terchangeable: Mini Head (No Coupler)** 35-88651 **494.95**

**Sable ROTAMAX S6500KL KAVO
Interchangeable:
Torque Head (No Coupler)** 35-88652 **494.95**

**6 Hole KAVO Multiflex LUX Type
Coupler** 35-88653 **199.95**

**6 Hole KAVO Multiflex LUX Type
Coupler: LED Bulb Long Life** 35-88654 **229.95**

**Sable Light Control Set Complete
Includes:** Sable Light Control Board, for 1-2
Handpieces with adjustable time delay &
voltage feature Sable Power Transformer
120v, 60Hz, 24VDC 1 Amp 6 Pin, 7 Foot,
Straight GRAY Tubing 35-88655 **279.95**

**Sable Light Control Set Complete (5
Hole) Includes:** Sable Light Control Board,
for 1-2 Handpieces with adjustable time delay
& voltage feature Sable Power Transformer
120v, 60Hz, 24VDC 1 Amp 5 Pin, 7 Foot,
Straight Gray Tubing with Bulb (1000815) 35-88656 **294.95**

Push Button Fiber Optic Handpiece Interchangeable with Mid- west/Bien Air Coupler (Dentex)

Push Button Fiber Optic
Handpiece Interchangeable
with Midwest/Bien Air
Coupler (Dentex), Swivel
Handpiece, Triple Port
Water Spray



Standard Head Fiber Optic 35-119421 **324.00**

Mini Head Fiber Optic 35-119431 **324.00**

Handpieces

Push Button Fiber Optic Handpiece Interchangeable with NSK Swivel Coupler (Dentex)



Push Button Fiber Optic Handpiece Interchangeable with
NSK Swivel Coupler (Dentex), Triple Port Water Spray,
Swivel Handpiece

Standard Head Fiber Optic 35-119441 **324.00**

Mini Head Fiber Optic 35-119451 **324.00**

Push Button Fiber Optic Handpiece Interchangeable with Star Coupler (Dentex)



Push Button Fiber Optic Handpiece Interchangeable with
Star Coupler (Dentex), Triple Port Water Spray, Swivel
Handpiece

Standard Head Fiber Optic 35-119461 **324.00**

Mini Head Fiber Optic 35-119471 **324.00**

Push Button Fiber Optic Handpiece Interchangeable with W & H Coupler (Dentex)



Push Button Fiber Optic Handpiece Interchangeable
with W&H Coupler model RA-24 & RA-25 (Dentex),
Triple Port Water Spray, Swivel Handpiece

Standard Head Fiber Optic 35-119481 **324.00**

Mini Head Fiber Optic 35-119491 **324.00**

Push Button Fiber Optic Handpiece with Triple Port Wa- ter Spray (Dentex)



Handpiece 4-Hole and 5-Hole Connector, Push Button,
Triple Port Water Spray, Fixed Back End FIBER OPTIC

Standard Head Fiber Optic 35-119501 **285.00**

Mini Head Fiber Optic 35-119511 **285.00**

Handpieces

Quest Optic Handpieces (Lares)

American-Manufactured Highspeeds

- Quest highspeeds feature low-maintenance, long-life, ceramic bearings
- The Quest stainless steel exterior and subtle grip design allows for fast, easy cleaning and autoclaving/chemiclaving
- Quest highspeeds are priced for value-conscious dentists unwilling to compromise on quality and performance
- The Quest 100 Series features a compact head for easy posterior access
- These handpieces are fast cutting, and have exceptional power and visibility, easy posterior access
- The Quest 300 Series offers a full sized head for powerful torque and smooth, quiet operation
- Quest 300 models feature heavy-duty construction and multi-port water spray

Optic Ceramic Bearing:

Push Button, Small Head, 5 Holes, model 120 35-119751 **399.95**

Optic Ceramic Bearing: Push Button, Full Sized Head, 5 Holes, model 320 35-119761 **399.95**

TwinPower Turbine High Speed F.O. Handpiece (J-Morita)

J. Morita's TwinPower Turbine now features ceramic ball bearings. Forty percent lighter and 3 times harder than conventional bearings, they offer extended turbine life, reduced operation noise, and less vibration. For maximum torque control, TwinPower's unique double-impeller design provides up to 22 W of power. Other features include mini head option, zero suck back in the air line, rapid braking within 2 seconds, glass rod optics, and direct connection to popular coupling systems.



The TwinPower operates by air from the drive air nozzles powering the primary impeller. The exhaust air is then directed through fixed fins, which power the secondary impeller. The result is that the handpiece is more powerful and provides constant torque and controlled speed, even under load.

Fits Morita, Fiber Optic Standard

Head, Includes: Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool 35-121421 **864.95**

Fits Sirona, Fiber Optic Standard

Head, Includes: Lubricant Spray Nozzle, Chuck Grip Tester, Spray Hole Wire, Screwdriver, Cap Tool 35-121431 **864.95**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Handpieces Parts & Accessories

2-3 Holes Tubing Adapters (Parts)

HP Tubing Adaptor, 4-Hole HP To 2/3-Hole Tubing



ing: Adaptors

4-Hole HP to 2/3-Hole HP tubing, Provides Drive Air & Coolant Water, Does not provide Coolant Air.

HP Tubing Adaptor, 4-Hole HP To 3-Hole Tubing: Adapts 4-Hole HP to 3-Hole HP tubing only, Provides Drive air, water & air coolant

2/3-Hole HP To 4-Hole Tubing: Adapts 2 or 3-Hole HP to 4-Hole HP tubing

HP Tubing Adaptor, 4-Hole HP To 2/3-Hole Tubing 35-48311 **74.95**

HP Tubing Adaptor, 4-Hole HP To 3-Hole Tubing 35-48312 **74.95**

2/3-Hole HP To 4-Hole Tubing 35-48313 **74.95**

3 Holes Tubing Adapters (Parts)

3 Hole Borden Adaptor, Plastic: Lightweight plastic with brass fittings

3 Hole Borden Adaptor, Metal: Quality plated brass

3 Hole Borden Nut & Plastic Connector: Plated brass nut and plastic connector, Brass fittings

3 Hole Borden Nut & Metal Connector: Plated brass nut and connector

3 Hole Borden Adaptor, Plastic 35-44271 **6.95**

3 Hole Borden Adaptor, Metal 35-44272 **7.95**

3 Hole Borden Nut & Plastic Connector 35-44273 **11.49**

3 Hole Borden Nut & Metal Connector 35-44274 **12.49**

4 Hole Tubing Adapters (Parts)



4 Hole Midwest Adaptor, Plastic 35-44811 **6.95**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

4 Hole Midwest Adaptor, Metal	35-44812	7.95
4 Hole Midwest Nut & Plastic Connector	35-44813	11.49
4 Hole Midwest Nut & Metal Connector	35-44814	12.49

4 Holes Tubing Connectors (Parts)

- Handy for extending HP or Foot Control Tubing
- Plastic body with brass barbs



Each	35-4841	7.95
-------------	---------	-------------

Bearing Puller



TST710



TST713

Bearing Puller	35-111531	112.49
Heavy Duty	35-111521	239.99

Bien Air Cap Wrench (HPP)



Part #: HBA130



Part #: HBA230



Part #: HBP130



Part #: HBP230

Outer Cap Wrench for Black Pearl & Bora	35-97781	24.95
Inner Cap Wrench for Black Pearl & Bora	35-97782	49.95
Outer Cap Wrench for Prestige	35-97783	24.95
Inner Cap Wrench for Prestige	35-97784	49.95

Bur Changer (HB)



Tradition



Pana Air/Naminex/Glow



Quiet Air



Pana Air/Naminex/Glow Mini

Fits: Midwest Tradition, each	35-92231	11.95
Fits: Midwest Quiet Air, each	35-92232	11.95
Fits: Pana Air/Naminex/Glow, each	35-92233	11.95
Fits: Pana Air/Naminex/Glow Mini, each	35-92234	11.95

Handpieces

Cap Wrenches



WR100



WR100-1

Universal Back Cap Wrench	35-111481	0.99
Back Cap Wrench	35-111471	1.29

Handpiece Gaskets (HPP)



Part #: HHX4HG
4 HOLE GASKET



Part #: HHX5HG
5 HOLE GASKET



Part #: SST004
WATER LINE GASKET

4 Hole Gasket, each	35-117241	3.19
5 Hole Fiber Optic Gasket, each	35-117231	3.19
2 Hole Gasket, each	35-11722	3.19
4 Hole Gasket, 5/Pkg	35-11721	14.95
5 Hole Fiber Optic Gasket, 5/Pkg	35-11720	14.95
2 Hole Gasket, 5/Pkg	35-11719	14.95

Handpiece Pull Tester (HPP)

A pull tester is an absolute must when repairing/rebuilding turbines. Our pull testers are simple, but very effective.



Part #: THR008

A precision pin is mounted to a spring loaded shaft that exerts a pull of 8 lbs which more than exceeds the industry standard of 6 lbs.



Part #: THR808

The standard pull tester is to be used on a turbine that is not in motion.

The dynamic pull tester can be used like the standard, however the test pin is attached that allows you to perform the pull test on a handpiece while it's running.

This is the most accurate way to test the holding power of a spindle.

Standard	35-97681	39.95
Dynamic	35-97682	64.95

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Impact Air Cap Wrench (HPP)



Each 35-97751 10.95

Kavo Back Cap Flat Wrench



Each 35-10018 11.29

Kavo Cap Wrench (HPP)



Part #: HKC130

Part #: HKC230

Part #: HKC230

Cap Wrench 35-97721 9.95

Heavy Duty Cap Wrench 35-97722 13.95

Kavo 635, 647, 649 Cap Wrench 35-97723 13.95

Kavo Slinger and Spray Cover Tools (HPP)



Part #: HK5222

Part #: HK7222

Part #: HK5240

Kavo 635 Slinger Install Tool 35-97701 19.95

Kavo 647-649 Slinger Install Tool 35-97701 19.95

Kavo Dual-Sided Spray Cover Tool 35-97703 39.95

Kavo Spare Parts



TST609

TST657

TST658

TST659

Kavo Cap Wrench 35-117931 10.19

Kavo 625/630 Head Expander 35-117921 39.79

Kavo 634/642 Head Expander 35-11791 39.79

Kavo 640B Head Expander 35-11790 39.79

Lares Cap Wrench (HPP)



Part #: HL5135



Part #: HL5130

Heavy Duty Dual Cap Wrench 35-97741 10.95

Lares 557 Cap Wrench 35-97742 10.95

Midwest Cap Wrench (HPP)



Part #: HMQ130



Part #: HMQ230



Part #: HMT230



Part #: HMT330



Part #: 753333

Quiet Air Standard Cap Wrench 35-97711 9.95

Quiet Air Standard Cap H/D Cap Wrench 35-97712 9.95

Tradition Standard H/D Cap Wrench 35-97713 9.95

Multi-Tool Cap Wrench works on: XGT, Tradition (L&PB), Quiet Air (L&PB), Stylus PB 35-97714 9.95

Tradition-L, Quiet Air Lever Cap Wrench (OEM) 35-97715 13.95

Midwest Wrenches



WR-A141



WR 140



WR141

Midwest Quiet Air Back Cap Wrench 35-118781 11.95

Midwest Tradition Bur Wrench 35-118791 11.95

Midwest Quiet Air Bur Wrench 35-118771 11.95

NSK Chuck-Type Wrench



Each 35-119011 9.99



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Star Cap Wrench (HPP)



Each 35-97761 **9.95**

Swivel Handpiece Connector 360 Degree (Parts)

- Add a swivel to 4-hole H.P. tubing
- Threads in place between tubing nut and handpiece
- Allows for 360 degree rotation



Each 35-44821 **234.95**

Turbine Rebuild Kits

Midwest Type Quiet Air, All Models	35-121311	20.99
Midwest Type Tradition, All exc. XGT	35-121301	20.99
Star 430 Type, All Models	35-121291	29.95
Kavo 640/642 Type	35-121281	24.99
Kavo 635B Type	35-121271	43.95
Kavo 647/649B Type	35-121261	49.95
NSK Phatelus/Viper Type	35-121251	22.95

W&H Adec Cap Wrench (HPP)



Part #: HAS136



Part #: HAS137



Part #: HAS138

Synea TA-96 Cap Wrench	35-97771	29.95
Synea TA-97 Cap Wrench	35-97772	29.95
Synea TA-98 Cap Wrench	35-97773	29.95

Handpieces

Implant Handpieces

SD Torque Surgeons Tool

- SD-Torque is the surgeons last tool to use
- Surgeons are able to reach deep molar areas without difficulty compared to other torques.



Each 36-17111 **794.95**

Lab Hanpieces

Lab Highspeed Handpiece (Nakamura)

Made in Japan



Highspeed Laboratory Handpiece 320,000rpm 41-11415 **259.95**

Highspeed Laboratory Hanpiece Cartridge for TS-300Z 41-11414 **129.95**

Low Speed Hanpieces

Air Torx E-Type Air Motor (J-Morita)

Air Torx offers efficient, powerful, and constant torque with operational speeds up to 20,000 rpm. Its innovative fluid dynamics generate about



twice as much torque compared to a conventional low speed motor in the 5,000–10,000 rpm range. Air Torx is comfortable to operate - it is lightweight, perfectly balanced, and delivers smooth power control. Versatile, it can be used for a wide range of tasks such as grinding dentures, adjusting, finishing, and polishing restorations, and for tooth polishing with a prophyl angle. Air Torx is designed for enhanced durability and offers an extended working life. Other features include: forward/reverse drive with continuous speed control, double-lock connection, autoclave safe design, and compatibility with standard ISO-type attachments.

- Powerful torque even at a low speed.
- Long working life
- Perfectly balanced, lightweight design
- Forward/reverse drive with continuous speed control
- Safe and easy double-lock connection
- Autoclavable for complete infection control
- Compatible with all standard ISO-type attachments


Unit Includes: Lubricant Spray Nozzle 35-98621 **849.95**

Handpieces

Airmotor E-Type (Nakamura)

Max Air Pressure 50-60psi

Made in Japan



Midwest Type, E-Type, 20,000 rpm 4 Holes	35-111391	209.95
Midwest Type, E-Type, 5,000 rpm 4 Holes	35-111381	269.95
Borden Type, E-Type, 20,000 rpm 2 Holes	35-111371	234.95
Borden Type, E-Type, 5,000 rpm 2 Holes	35-111361	269.95

Contra Angle Handpiece NRS-E (NSK)

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 degrees C and highly durable.



Product Features:

- 10:1 Reduction Head
- For CA burs(Ø2.35)
- Clean Head System
- Push Button Chuck

NRS-E Contra Angle Handpiece	35-11186	159.95
-------------------------------------	----------	---------------

Contra Angle Sheath for Midwest Airmotor Type (Nakamura)

CONTRA ANGLE SHEATH FOR MIDWEST MOTOR



Midwest
Shorty Motor
Rhino XP Motor

Spare Head for SHS-MWC



SHS-MWC Contra Angle Sheath	35-120651	99.95
EH-50PK Spare Head	35-120531	44.95
EH-50PS Spare Head	35-120541	44.95
EH-20L Spare Head	35-120631	42.95



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

EH-20ML Spare Head	35-120621	45.95
EH-20PS Spare Head	35-120561	39.95
EH-20PK Spare Head	35-120551	39.95
EH-35BL Spare Head	35-120571	94.95
EH-30BL Spare Head	35-120581	129.95
EH-30BLP Spare Head	35-120591	129.95
EH-30TL Spare Head	35-120601	59.95
EH-20FG Spare Head	35-120611	49.95

Contra Angle with Push Button Head Handpiece FBP-E (NSK)

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common handpiece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 deg. and highly durable.



Features:

- For CA burs(Ø2.35)
- Clean Head System
- Ultra Push Chuck
- Connect directly to all brands of ISO E-type air or electric micromotors 1:1 transmission gear

1:1 Direct Drive: For CA burs(Ø2.35) - Clean Head System -Push Button Chuck - Ball Bearing -Max. Speed 40,000 min	35-111881	234.95
--	-----------	---------------

E-Type Contra Angle 4:1 (Nakamura)



EG-30BL Midwest Type	35-11233	148.79
EG-30TL Midwest Type	35-112341	148.79
EG-20L Midwest Type	35-112351	120.83
EG-20ML Midwest Type	35-112361	125.37
EG-20FG Midwest Type	35-112371	130.82
EG-20AP Midwest Type	35-11228	125.37




Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Handpieces

EF-20APS Midwest Type	35-11227	172.61
EG-20PS Midwest Type	35-11218	119.00
EG-20PK Midwest Type	35-11219	119.00
EG-50PS Midwest Type	35-11220	123.55
EG-50PK Midwest Type	35-11221	123.55
EGN-30BL NSK or STAR Type	35-11223	148.79
EGN-30TL NSK or STAR Type	35-11222	148.79
EGN-20L NSK or STAR Type	35-11224	120.83
EGN-20ML NSK or STAR Type	35-11225	125.37
EGN-20FG NSK or STAR Type	35-11226	130.82
EGN-20AP NSK or STAR Type	35-112381	125.37
EGN-20APS NSK or STAR Type	35-11232	172.61
EGN-20PS NSK or STAR Type	35-11231	119.00
EGN-20PK NSK or STAR Type	35-11230	119.00
EGN-50PS NSK or STAR Type	35-11229	123.55
EGN-50PK NSK or STAR Type	35-11217	123.55

E-Type Contra Angle Handpieces (Nakamura)




Made in Japan

EC-30BL Push-Up Latch, Ball-Bearing Head for CA bur, 30,000 rpm	35-11207	119.95
EC-30BLP Push Button, Ball-Bearing Head for CA bur, 30,000 rpm	35-11206	194.95
EC-20L Standard Latch Head for CA bur, 20,000 rpm	35-11208	99.95

ESG-30AR 4:1 Speed Reduction Straight Handpiece	35-112109	234.95
ES-30A Straight Handpiece 40,000 rpm	35-11210	94.95
ME-80B 2 Holes 5,000 rpm Airmotor	35-11211	269.95
ME-80M 4 Holes 5,000 rpm Airmotor	35-112131	269.95
ME-20B 2 Holes 20,000 rpm Airmotor	35-11204	234.95
EG-20L 4:1 Speed Reduction Standard Latch Head for CA bur	35-11203	234.95
EG-20PS 4:1 Speed Reduction Regular Prophecy Screw-On Head	35-11193	119.95
ME-20M 4 Holes 20,000 rpm Airmotor	35-112151	209.95

E-Type Durable Handpiece (Dentex)

Durable Stainless Steel



E-Type Low-Speed Handpiece System Durable Stainless Steel (Dentex)

20,000 RPM Motor	35-112481	252.00
1:1 Contra Angle	35-112461	86.99
4:1 Reduction Contra Angle	35-112451	129.99
Manual Latch Head	35-11242	36.59
Push Latch Head	35-11241	70.49
Prophy Head Screw Type	35-11240	34.99
Quarter Turn Endo Head	35-11239	70.49
5,000 RPM Air Motor	35-112471	309.99
1:1 Nose Cone	35-112441	140.99
4:1 Reduction Nose Cone	35-11243	184.49

Low Speed Handpieces (Nakamura)

Slim design, light weight, smooth run, minimum heat generation, achieve less fatigue for its user and guarantee trouble free operation.



50,000rpm, each	35-12851	349.95
20,000rpm, each	35-12852	319.95

Call us for Quarterly Promotion


Prestigedentalproducts.com

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

E-Type Handpiece Sets (Nakamura)



ME-20MS
includes
ES-30A Straight Handpiece
EC-20L Standard Latch Head
Contra Angle Handpiece
ME-20M 20,000rpm Airmotor

ME-30MS
includes
ES-30A Straight Handpiece
EC-30BL Push-Up Latch, Ball-Bearing
Head CA Handpiece
ME-20M 20,000rpm Airmotor

ME-80MS
includes
ES-30A Straight Handpiece
EC-20L Standard Latch Head
Contra Angle Handpiece
ME-80M 5,000rpm Airmotor

Made in Japan

Set:

ES-30A 1:1 Straight Nosecone HP, EC-20L
1:1 Standard Latch Head Contra Angle HP
& ME-20M 20,000 rpm Air motor 35-111421 **399.95**

Set: ES-30A 1:1 Straight Nosecone HP,
EC-30BL 1:1 Push-Up Latch, Ball-Bearing
Head CA HP & ME-20M 20,000 rpm Air
motor 35-111411 **419.95**

Set: ES-30A 1:1 Straight Nosecone HP,
EC-20L 1:1 Standard Latch Head Contra
Angle HP & ME-80M 5000 rpm Air motor 35-111401 **449.95**

E-Type Handpiece with Latch Head BB-E (NSK)

E-type External Water Spray – for Varying
Treatment Needs
NSK provides multiple types of contra-angle
handpieces for your varying treatment needs.
All NSK handpieces work perfectly with a
common handpiece shank, offering conven-
ience and flexibility to every clinician's di-
versified needs. Heads and shank are all
autoclavable at up to 135 deg. and highly
durable.



Features:


- For CA burs(Ø2.35)
- Clean Head System
- Swing latch type
- Connect directly to all brands of ISO E-type air or electric micro motors 1:1 transmission gear
- 40,000 rpm Ball bearing head for smoother and vibration free operation

NSK BB-E 1:1

w/Latch Head 40,000 rpm 35-112491 **174.95**

E-Type Ultra Light Handpiece (Dentex)

Ultra Light Aluminum



E-Type Low-Speed Handpiece System Ultra Light Aluminum

20,000 RPM Motor	35-112601	211.49
1:1 Contra Angle	35-112571	70.49
4:1 Reduction Contra Angle	35-112581	112.49
Manual Latch Head	35-11255	36.59
Push Latch Head	35-11254	70.49
Propy Head Screw Type	35-11253	35.19
Quarter Turn Endo Head	35-11251	70.49
5,000 RPM Air Motor	35-11256	269.99
1:1 Nose Cone	35-11252	126.99
4:1 Reduction Nose Cone	35-11259	155.99

Friction Grip Push Button Type Handpiece NFB-E (NSK)

E-type External Water Spray – for Varying
Treatment Needs

NSK provides multiple types of contra-angle
handpieces for your varying treatment needs.
All NSK handpieces work perfectly with a
common handpiece shank, offering conven-
ience and flexibility to every clinician's di-
versified needs. Heads and shank are all autoclav-
able at up to 135 deg. and highly durable.



Features:

- Standard Head
- Swing Latch type
- For friction grip burs (Ø1.6)
- Max.40,000 min-1
- Connect directly to all brands of ISO E-type air or electric micromotors 1:1 transmission gear


**NFB-E Friction Grip Push Button type
Handpiece**

35-112641 **214.95**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Hygienist Airmotor Handpiece (Nakamura)



Midwest Type 5,000 rpm 4 Holes
for Propy Angle use 35-117691 **229.95**

Regular Screw-In Propy Angle HP 35-117692 **14.95**

Regular Screw-On Propy Angle HP 35-117693 **14.95**

Sealed Screw-In Propy Angle HP 35-117694 **29.95**

Sealed Screw-On Propy Angle HP 35-117695 **29.95**

Integrity Low Speed Handpiece Set (Osseo Scientific)

Includes:

- 20,000 RPM E-Type Motor (NTE-20)
- 1:4 Attachment (NTE-14)
- Standard Latch Head (NTL-01)
- 1:1 Straight Nose Cone (NTS-11)



Low Speed Handpiece Combination Set
includes: 20,000 RPM E-Type Motor (NTE-20)
1:4 Attachment (NTE-14) Standard Latch Head
(NTL-01) 1:1 Straight Nose Cone (NTS-11) 35-117881 **427.95**

Low Speed Handpiece Combination Set
includes: 20,000 RPM E-Type Motor (NTE-20)
1:4 Straight Nose Cone (NTS-14) 1:1 Straight
Nose Cone (NTS-11) 35-117871 **427.95**

Low Speed Handpiece Combination Set
includes: 5,000 RPM E-Type Motor (NTE-5)
1:1 Attachment (NTE-11) Standard Latch Head
(NTL-01) 1:1 Straight Nose Cone (NTS-11) 35-11786 **427.95**

Low Speed Handpiece Combination Set
include: 20,000 RPM E-Type Motor (NTE-20)
1:4 Attachment (NTE-14) 1:1 Attachment (NTE-
11) 2 Standard Latch Heads (2NTL-01) 35-11785 **427.95**

Lynx Low Speed Handpiece (MTI)

Air Motor, Low Speed Handpiece System

- Length: 405"
- Weight: 2.7oz
- Speed: 2,000 to 20,000
- Operating Air Pressure: 35 to 50 psi.
- Air Consumption: 1.4 cfm.
- Connections: 4-, or 2-hole

2 Hole 35-118121 **264.95**

4 Hole 35-118111 **264.95**

Handpieces

Micromite Low Speed Handpiece (Micromite-Japan)

Autoclavable, High
torque

5,000RPM 4 Holes 35-118551 **379.95**

Micromite Type Low Speed Handpieces (Nakamura)

U-type Nose,
4-holes



Made in Japan

Midwest type Straight type 4 holes
5,000 rpm 35-118581 **314.95**

Borden type Straight type 2 holes
5,000 rpm 35-118571 **314.95**

Midwest type Straight type 4 holes
20,000 rpm 35-118561 **314.95**

Borden type Straight type 4 holes
20,000 rpm 35-218561 **314.95**

Prophy Angle E-Type 1:1 (Nakamura)



Midwest Type, EC-20PS 35-21546 **84.95**

Midwest Type, EC-20PK 35-21547 **84.95**

Midwest Type, EC-50PS 35-21548 **84.95**

Midwest Type, EC-50PK 35-21549 **84.95**

Push Button Head Contra Angle Handpieces (Nakamura)

30,000rpm, Ball-Bearing
Push-Button Type Head, For
CA latch bur



Made In Japan

E-type (ISO Type)
Connection Contra Angle 35-82361 **194.95**

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Star Titan-type Connection Contra Angle 35-82362 **269.95**

Midwest-type Connection Contra Angle 35-82363 **234.95**

Rotary Latch Files Head Handpiece NAC-E316R (NSK)

E-Type External Water Spray

The speed of all Endodontic low speed handpieces is determined by the speed of the micromotor to which they are connected. Micromotors can vary in speed ranges. Micromotor speed should be confirmed before selecting a reduction handpiece.



Features:

- For engine files
- 16:1 Reduction
- 360°rotation
- For Ni-Ti files Standard Head (NAC-Y) Latch Head, (E16R) Contra Angle
- Spray Nozzle (Optional)

NAC-E16R

Rotary Latch Files Head Handpiece 35-120141 **349.95**

Slow Speed E Type (Johnson-Promident)



Classic Slow Speed Set

(E-Motor + Nose + Contra Angle 30K) 35-27241 **438.95**

Classic "E" Type Motor (25,000 RPM) 35-27242 **264.95**

Classic "E" Type Motor (5,000 RPM) 35-27243 **318.95**

Star Titan Type Low Speed Handpieces (Nakamura)

Superior performance, Very low noise level, High durability, Various system choice, Easy & comfortable use, Autoclavable



Head Attachment

For Friction Grip bur 20,000rpm 35-92341 **64.95**

Head Attachment

For CA latch bur 20,000rpm 35-92342 **52.95**

Head Attachment For CA latch bur

Ball-Bearing 30,000rpm 35-92343 **89.95**

Head Attachment For CA latch bur

Ball-Bearing Push-Button 30,000rpm 35-92344 **119.95**

Head Attachment For Propy Snap-on

35-92345 **39.95**

Head Attachment

For Propy Screw-in 1-64 35-92346 **39.95**

1:1 Contra Angle Sheath Low Speed Attachment

35-92347 **139.95**

4:1 Speed Reduction Torque Multiplier Low Speed Attachment

35-92348 **149.95**

Straight Nosecone

Ball-Bearing, U-type Nose 40,000rpm Low Speed Attachment 35-92349 **264.95**

1:1 Contra Angle Sheath 12-teeth Low Speed Attachment

35-92350 **154.95**

20,000rpm Low Speed Air Motor

4-hole Connector 35-92351 **349.95**

5,000rpm Low Speed Air Motor

4-hole Connector 35-92352 **379.95**

Straight Durable Handpiece (Dentex)

Compact easy-to-use low-speed handpiece for use with doriot style contra angles and straight handpiece/lab burs. 20,000 rpm

- Straight nose swivels 360°, reducing tubing twist and hand fatigue
- Speedadjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

Durable Stainless Steel 35-120881 **340.49**

Straight Durable Handpiece M5 (Dentex)

Compact, easy-to-use low-speed handpiece for use with doriot-style contra angles, disposable and metal prophy angles and straight handpiece/lab burs



- Variable speed up to 5,000 rpm
- Straight nose cone swivels 360°, reducing tubing twist and hand fatigue



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888


- Speed-adjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

Durable Stainless Steel 35-120891 **392.99**

Straight Handpieces E-Type (Nakamura)



ES-30A



ES-35L

1:1, 40,000 rpm,
 3 Ball Bearings U Type Nose 35-120911 **169.95**

1:1, 35,000 rpm,
 2 Ball Bearings For Laboratory 35-120901 **239.95**

Straight Handpieces for Midwest Type Airmotor (Nakamura)

Made in Japan



MWS-30A
40,000rpm

40,000 rpm 3 Ball Bearings 35-120921 **194.95**

Straight Nose Cone Handpiece EX5 (NSK)

EX5 E-type External Water Spray –
 for Varying Treatment Needs

NSK provides multiple types of contra-angle handpieces for your varying treatment needs. All NSK handpieces work perfectly with a common hand-piece shank, offering convenience and flexibility to every clinician's diversified needs. Heads and shank are all autoclavable at up to 135 deg. and highly durable.



Features:

- For Prophylaxis Application
- For screw in-cups and brushes
- Connect directly to all brands of ISO E-type air or electric micro-motors 4:1 reducer

NSK EX5
 4:1 reduction Straight Handpiece 35-120931 **219.95**

Straight Nose Cone Handpiece EX6 (NSK)

NSK EX-6/ Propy Straight Hand-piece, EX-6: For HP burs(dia. 2.35), Non-Optic, Clean Head System, Spray Nozzle(Optional). Disposable Propy Heads can be attached to EX-6



EX6 1:1 Straight Nose Cone Handpiece 35-120941 **129.95**

Handpieces

Straight Ultra Light Handpiece (Dentex)

Compact easy-to-use low-speed handpiece for use with doriot style contra angles and straight handpiece/lab burs. 20,000 rpm



- Straight nose swivels 360°, reducing tubing twist and hand fatigue
- Speed adjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

Ultra Light Aluminum 35-120951 **292.49**

Straight Ultra Light Handpiece M5 (Dentex)

Compact, easy-to-use low-speed handpiece for use with doriot-style contra angles, disposable and metal proph angles and straight handpiece/lab burs.



- Variable speed up to 5,000 rpm
- Straight nose cone swivels 360°, reducing tubing twist and hand fatigue
- Speed-adjusting ring controls speed and direction of motor
- No lubrication required
- Fully autoclavable

Ultra Light Aluminum 35-120961 **342.00**

Non-Fiber Optic High Speed Handpieces

Highspeed Handpiece (Nakamura)

Manual Chuck Type



TC-40PM (4 holes)
TC-40PB (2 holes)
 Standard Size Head, 400,000rpm



TC-35YM (4 holes)
TC-35YB (2 holes)
 Large Torque Head, 350,000rpm



Quick-Joint Model
TC-80QM (4 holes)
TC-80QB (2 holes)
 Standard Size Head, 380,000rpm



TCQM2950
 360° Quick Connect Swivel for TCP-70QM, TC-80QM

Made in Japan

Standard Size Head, 400,000rpm,
 4 holes 35-97541 **169.95**

Standard Size Head, 400,000rpm,
 2 holes 35-97542 **169.95**

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Large Torque Head, 350,000rpm, 4 holes	35-97543	179.95
Large Torque Head, 350,000rpm, 2 holes	35-97544	179.95
Standard Size Head with swivel, 380,000rpm, 4 holes	35-97545	219.95
Standard Size Head with swivel, 380,000rpm, 2 holes	35-97546	219.95
360 Quick Connect Swivel for TC-80QM & TC-80QB	35-97547	74.95

Highspeed Handpiece- Pana-MAX Push Button


Made in Japan



Each	35-16531	249.95
-------------	----------	---------------

Highspeed Push Button Handpiece (Nakamura)

Push-Button Chuck Type



TCP-450M (4 holes)
TCP-450B (2 holes)
Standard Size Head, 400,000rpm

TCP-350M (4 holes)
TCP-350B (2 holes)
Large Torque Head, 350,000rpm

Quick-Joint Model
TCP-70QM (4 holes)
TCP-70QB (2 holes)
Standard Size Head, 380,000rpm

Standard Size Head, 400,000rpm, 4 holes	35-11764	179.95
Standard Size Head, 400,000rpm, 2 holes	35-11763	179.95
Large Torque Head, 350,000rpm, 4 holes	35-117651	194.95
Large Torque Head, 350,000rpm, 2 holes	35-117661	194.95
Standard Size Head with swivel, 380,000rpm, 4 holes	35-117662	229.95
Standard Size Head with swivel, 380,000rpm, 2 holes	35-117663	229.95

Integrity Leverage Classic Highspeed with Push Button Cap (Osseo Scientific)



Integrity Leverage Classic Highspeed with Push Button Cap	35-117801	179.49
--	-----------	---------------

Integrity Leverage Handpiece with Lever Cap (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best highspeed handpiece performance and value at a fraction of the cost of the top leading brands.



- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip Grip design

Handpiece with Lever Cap	35-117841	279.95
---------------------------------	-----------	---------------

Integrity Push Button Handpiece (Osseo Scientific)

Like top brand handpieces, the Integrity™ is an exceptionally designed handpiece that offers smooth, quite controlled cutting action with high torque through the entire operating speed range, receive the best highspeed handpiece performance and value at a fraction of the cost of the top leading brands.



- Light Weight
- Whisper Quiet
- Small Head Size
- Push Button Bur Changing
- Ergonomic, Comfort For Smooth Control Attractive surface with No-Slip Grip design

Highspeed Handpiece Push Button	35-117891	129.95
--	-----------	---------------

Prophy Angle U-Type Right Angle HP (Nakamura)

RA-10PS



RA-10PS	35-119321	13.96
RA-10PK	35-119331	13.95
RA-50PS	35-119341	19.95
RA-50PK	35-119351	19.95



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Handpieces

Quest Push Button for Swivel Handpieces (Lares)



100 Series
 Power and Durability in a Compact Head

300 Series
 Power for Your Toughest Cutting Requirements



- Small, compact head for excellent visibility and accessibility
- High speed for fast cutting
- Reduced contra angle for easier posterior access
- Easy, secure push button bur changing
- Chair-side rotor and chuck replacement
- Optional bright twin beam fiber optics
- Optional 360° swivel and quick disconnect
- Autoclavable and chemiclavable
- Low maintenance Ceramic Bearings



- Full-size head for powerful torque
- Heavy duty construction
- Smooth, quiet operation
- Optional 360° swivel and quick disconnect
- Optional bright twin beam fiber optics
- Easy, secure push button bur changing
- Chair-side rotor and chuck replacement
- Autoclavable and chemiclavable
- Low maintenance Ceramic Bearings
- Multiple Port Water Spray

Quest 150 LO	Height	Diameter	Weight
PB SW MX	11.9 mm	10.2 mm	50 g
11320	Speed	Water Spray	Torque
	450,000 rpm	Single Port	15.3 watts

Quest 350 LO	Height	Diameter	Weight
PB SW MX	14.65 mm	13 mm	56 g
11322	Speed	Water Spray	Torque
	350,000 rpm	Multi-Port	20.8 watts

American Manufactured highsPEEDS

- Quest highsPEEDS feature low-maintenance, long-life, ceramic bearings. The Quest stainless steel exterior and subtle grip design allows for fast, easy cleaning and autoclaving/chemiclaving
- Quest highsPEEDS are priced for value-conscious dentists unwilling to compromise on quality and performance
- The Quest 100 Series features a compact head for easy posterior access
- These handpieces are fast cutting, and have exceptional power and visibility, easy posterior access
- The Quest 300 Series offers a full sized head for powerful torque and smooth, quiet operation
- Quest 300 models feature heavy-duty construction and multi-port water spray

130 Quest Non Fiber Optic Small Head
 HP 360 Swivel, Push Button, Multiflex
 Compatible, Ceramic Bearings, w/o Swivel 35-119801 **479.95**

330 Quest Non Fiber Optic Standard Head
 HP 360 Swivel, Push Button, Multiflex Compatible, Ceramic Bearings, w/o Swivel 35-119791 **479.95**

4 Holes Coupler,
 Multiflex Compatible Coupler 35-119781 **194.95**

5 Holes Coupler, Multiflex Coupler 35-119771 **295.95**

Highspeed Basic Non-Push Button Handpiece (Alfa Air)

Handpiece Non-Push 35-117561 **99.99**

Turbine 35-117551 **29.99**

Turbine Deluxe, Made in Japan 35-11754 **59.99**

Highspeed Basic Push Button Handpiece (Alfa Air)

Handpiece Push Button, Made in Japan 35-117581 **169.99**

Turbine Deluxe, Made in Japan 35-117571 **69.99**

Quick Joint Type High Speed Handpieces (Nakamura)



Made in Japan

Standard Chuck 380,000rpm ISO Con 35-82351 **259.95**

Push-Button (Auto)
 Chuck 380,000rpm ISO Con 35-82352 **279.95**

360° Quick Connect Swivel for TC-80QM, TCP-70QM 35-82353 **64.95**

E-Type Straight Handpiece (Nakamura)



ES-30A 35-112501 **152.28**

Handpieces

Micromotor Handpiece and E-Type Micromotor (Nakamura)

2 LB-25HM



LB-20EM
E-type Micromotor
Speed Range : 3,000 - 20,000rpm



LB-25HM
Micromotor Handpiece
Speed Range : 3,000 - 20,000rpm



LB-30HM	35-118591	447.87
LB-35EM	35-118601	290.71
LB-25HM	35-118611	284.35
LB-20EM	35-118621	194.71

Micromotor Handpiece Set for Laboratory (Nakamura)



LX-30FS

Made in Japan



LXC-25FS

LX-30FS Set Includes: LX-35V, LB-30HM, AND FC-500	35-118641	1062.90
LX-35FS	35-118651	1050.18
LXC-20FS	35-118661	580.51
LXC-25FS Set Include: LXC-12V, LB-25H, FCC-150 and cradle	35-118631	557.79

Portable Micromotor System (Nakamura)

BD-550 (Battery)



BD-20FS



Set Includes: BD-550, LB-20EM, ES-30A, CHR-80A, Cradle	35-119171	619.20
Set Includes: BD-550, LB-25HM, CHR-80A, Cradle	35-119181	604.00
Battery-Pack Controller	35-119191	245.60



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Table-Top Type Controller (Nakamura)

1 LXC-12V



LB-35V	35-121141	466.95
LXC-12V	35-121151	199.86

U-Type Contra Angle Sectional (Nakamura)

DU-30BL



DU-30TL



DU-20L



DU-20ML



DU-20FG



DU-30BL	35-121451	81.76
DU-30TL	35-121461	81.76
DU-20L	35-121471	54.51
DU-20ML	35-121481	59.96
DU-20FG	35-121491	65.41

U-Type Contra Angle Standard (Nakamura)

DC-10L



DC-10APS



DC-10L	35-121501	28.95
DC-10FG	35-121511	64.51
DC-10APS	35-12152	99.94

Chuck Style Handpiece (Dentex)

Chuck Style Fiber
Optic SINGLE
PORT WATER
SPRAY Hand-
piece (Dentex),
Fixed Back End



Standard Head Non-Optic	35-111821	134.95
Mini head Non-Optic	35-111811	134.95



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Chuck Style Triple Port Water Spray Handpiece (Dentex)

Chuck Style Handpiece 4-Hole and 5-Hole Connector Non-Optic (Dentex), Triple Port Water Spray, Fixed Back End



Standard Head	35-111831	171.00
Mini Head	35-111841	171.00

Push Button Handpiece (Dentex)

Handpiece 4-Hole and 5-Hole Connector Non-Optic (Dentex), Push Button, Triple Port Water Spray, Fixed Back End, Highspeed



Standard Head Non-Optic	35-119521	213.00
Mini Head Non-Optic	35-119531	213.00

Push Button Handpiece (Sable)

SABLE HIGH QUALITY
 PUSH BUTTON,
 QUICK CONNECT,
 360° SWIVEL,
 HANDPIECES.

ROTAMAX™ II
 Standard Head Handpiece
 w/ 4 Hole Coupler
 2000006

ROTAMAX™ II
 Standard Head
 Handpiece
 2000004

ACCESS™ II
 Mini Head Handpiece
 w/ 4 Hole Coupler
 2000005

ACCESS™ II
 Mini Head
 Handpiece
 2000003

4 Hole / Quick Connect / 360° Swivel
 2000101

2 Hole / Quick Connect / 360° Swivel
 2000102

1 YEAR WARRANTY
 Sable handpiece turbines are
 manufactured to meet and exceed
 ISO 7785-1 handpiece standards.

ACCESS II: PB Mini Head, Handpiece Only (No Coupler)	35-88661	298.95
ACCESS II: PB Mini Head 360/Quick Connect Coupler Included	35-88662	329.95
ROTAMAX II: PB Standard Head, Handpiece Only (No Coupler)	35-88663	298.95

Handpieces

ROTAMAX II: PB Standard Head, 360 Quick Connect Coupler Included	35-88664	329.95
---	----------	---------------

4 Hole 360 Quick Connect Coupler To Fit Access II & Rotamax II	35-88665	59.95
---	----------	--------------

2 Hole 360 Quick Connect Coupler To Fit Access II & Rotamax II	35-88666	59.95
---	----------	--------------

Push Button Handpiece Interchangeable with Kavo Coupler (Dentex)

Push Button Handpiece Interchangeable with Kavo Coupler (Dentex), Swivel Handpiece



Standard Head Non-Optic	35-119541	259.95
Mini Head Non-Optic	35-119551	259.95

Push Button Handpiece Interchangeable with Midwest/ Bien Air Coupler (Dentex)

Push Button Handpiece Interchangeable with Midwest/Bien Air Coupler (Dentex), Swivel Handpiece Triple Port Water Spray



Standard Head Non-Optic	35-119561	252.00
Mini Head Non-Optic	35-119571	252.00

Push Button Handpiece Interchangeable with NSK Swivel Coupler (Dentex)

Push Button Handpiece Interchangeable with NSK Swivel Coupler (Dentex) Swivel Handpiece, Triple Port Water Spray



Standard Head Non-Optic	35-119581	252.00
Mini Head Non-Optic	35-119591	252.00

Push Button Handpiece Interchangeable with Star Coupler (Dentex)

Push Button Highspeed Handpiece, Triple Port Water Spray, Swivel Handpiece Interchangeable with Star Coupler (Dentex)



Standard Head Non-Optic	35-119601	252.00
Mini Head Non-Optic	35-119611	252.00

Call us for Quarterly Promotion

Prestigedentalproducts.com

Handpieces

Push Button Handpiece Interchangeable with W&H Coupler (Dentex)

Push Button Handpiece Interchangeable with W&H Coupler model RA-24 & RA-25 (Dentex) Triple Port Water Spray, Swivel Handpiece



Standard Head Non-Optic 35-119621 **252.00**

Mini Head Non-Optic 35-119631 **252.00**

Push Button Handpiece M1 (Dentex)

Push Button SINGLE PORT WATER SPRAY 4-Hole and 5-Hole Connector Handpiece (Dentex) , Fixed Back End



Standard Head Non-Optic 35-119641 **149.95**

Mini Head Non-Optic 35-119651 **149.95**

Quest Handpieces (Lares)

American Manufactured Highspeeds

- Quest highsPEEDS feature low-maintenance, long-life, ceramic bearings. The Quest stainless steel exterior and subtle grip design allows for fast, easy cleaning and autoclaving/chemiclaving
- Quest highsPEEDS are priced for value-conscious dentists unwilling to compromise on quality and performance
- The Quest 100 Series features a compact head for easy posterior access
- These handpieces are fast cutting, and have exceptional power and visibility, easy posterior access
- The Quest 300 Series offers a full sized head for powerful torque and smooth, quiet operation
- Quest 300 models feature heavy-duty construction and multi-port water spray




Push Button Ceramic Bearing: Small

Head Handpiece, 4 Head, model 110 35-119701 **369.95**

Push Button Ceramic Bearing: Full

Sized Head Handpiece, 4 Head, model 310 35-119691 **369.95**

Sheath of E-Type Contra Angle (Nakamura)



SHS-EC (1:1) Midwest Type 35-120381 **58.14**




Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

SHS-ECG (4:1) Midwest Type	35-120401	89.94
SHS-EGG (16:1) Midwest Type	35-120391	131.73
SHS-ECN (1:1) NSK or STAR Type	35-120411	58.14
SHS-EGN (4:1) NSK or STAR Type	35-120421	89.94
SHS-EGGN (16:1) NSK or STAR Type	35-120431	131.73

Spare Head with Elbow for Star-Titan Type Contra Angle (Nakamura)



MH-30BL Midwest Type	35-120821	67.23
MH-30TL Midwest Type	35-120831	67.23
MH-20L Midwest Type	35-120841	39.97
MH-20ML Midwest Type	35-120851	44.51
MH-20FG Midwest Type	35-120861	49.97
MH-20AP Midwest Type	35-120771	44.51
MH-20PS Midwest Type	35-120761	38.16
MH-20PK Midwest Type	35-12067	38.16
MH-50PS Midwest Type	35-12068	42.70
MH-50PK Midwest Type	35-120691	42.70
MH-35BL Midwest Type	35-120701	84.49
MHN-30BL NSK or STAR Type	35-120711	67.23
MHN-30TL NSK or STAR Type	35-120721	67.23
MHN-20L NSK or STAR Type	35-120731	39.97
MHN-20ML NSK or STAR Type	35-12074	44.51
MHN-20FG NSK or STAR Type	35-12075	49.97
MHN-20AP NSK or STAR Type	35-120871	44.51
MHN-20PS NSK or STAR Type	35-120811	38.16
MHN-20PK NSK or STAR Type	35-120801	38.16
MHN-50PS NSK or STAR Type	35-120791	42.70
MHN-50PK NSK or STAR Type	35-120781	42.70
MHN-35BL NSK or STAR Type	35-12066	84.49



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Cartridge Turbine for Handpiece (Nakamura)

CARTRIDGE FOR AIRTURBINE HANDPIECE



Made in Japan

Standard Chuck , for TC-40PM / 40PB & TC-80QM / 80QB	35-11176	74.95
Large Size, Standard Chuck , for TC-35YM / 35YB	35-111751	89.95
Push Button Chuck , for TC-450M / 450B & TCP-70QM / 70QB	35-111741	89.95
Push Button Chuck , for TC-350M / 350B	35-111731	94.95
Cartridge for TS-300ZT	35-111732	119.95

Prophy Handpieces

Hygiene Handpiece (MTI)

- Excellent for hygienic procedures, pins, post excavations and prophys
- This high torque, low speed handpiece is perfectly balanced to reduce hand and arm fatigue
- Fits all straight burs, doriot attachments and disposable prophy heads
- Fixed Forward 5,000 rpm



Each	35-89791	239.95
-------------	----------	---------------

Hygiene Prophy AirLite Handpiece with Friction Grip (MTI)

- Fixed forward 3,000 rpm
- Contoured for comfort
- Super quiet vibration free
- No more torn gloves with our Friction Grip System
- Your choice of attractive colors
- Now with a **2-Year Warranty**




Blue Color	35-89801	239.95
Purple Color	35-89802	239.95
Pink Color	35-89803	239.95

Handpieces

Hygienist Airmotor Handpiece (Nakamura)

AIRMOTOR HANDPIECE FOR HYGIENIST



Midwest Type 5,000 rpm 4 Holes for Prophy Angle use	35-117691	229.95
Regular Screw-In Prophy Angle HP	35-117692	14.95
Regular Screw-On Prophy Angle HP	35-117693	14.95
Sealed Screw-In Prophy Angle HP	35-117694	29.95
Sealed Screw-On Prophy Angle HP	35-117695	29.95

Hygienist Durable Prophy Handpiece (Dentex)

Durable Stainless Steel M5P Prophy Handpiece (Dentex)

- High Torque Motor
- Designed for use with disposable and metal autoclavable prophy angles
- Quiet, compact and light weight-minimizes hand and wrist fatigue
- Nose Cone swivels 360°
- No Lubrication required
- Fully Autoclavable



0-5,000 Durable Stainless Steel M5P Hygienist Prophy Handpiece (Dentex)	35-117671	319.49
--	-----------	---------------

Hygienist Handpiece (Johnson-Promident)

Compare to Midwest RDH

HYGIENE HANDPIECE

Compare to Midwest RDH



- Same design, performance, and warranty as the original for just half the price
- Lightweight ergonomic design for better positioning and reduced fatigue
- Smooth 360 degree swivel for better comfort and access
- Removable autoclavable color bands for handpiece identification

Each	35-117958	389.99
-------------	-----------	---------------

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Hygienist Propy Handpiece (Dentex)

Ultra Light Aluminum M5P Hygienist Propy Handpiece (Dentex)

- High Torque Motor
- Designed for use with disposable and metal autoclavable proph angles
- Quiet, compact and light weight-minimizes hand and wrist fatigue
- Nose Cone swivels 360°
- No Lubrication required
- Fully Autoclavable



0-5,000 Ultra Light Aluminum M5P Hygienist Propy Handpiece (Dentex) 35-117681 **269.99**

iBrite Deplaquer Polisher (Pac-Dent)

A convenient cordless stain remover that cleans better than just brushing and flossing. Its unique fusion of professional technology and ease of use allows for the ultimate removal of tobacco, coffee, tea, wine, and other types of stains even years after they form. The results will be visible immediately after use, and you will be able to smile more confidently. The kit comes with several different types and sizes of cleaning heads to ensure the cleaning cups reach all tooth crevices and gaps, and effectively remove the stains accumulated over time. The included polishing paste contains fluoride for cavity prevention.



Advantages:

- Cleans better than brushing and flossing alone
- Effective removal of tobacco, coffee, tea, wine, and other types of stains
- Immediate results
- Included paste contains fluoride for cavity prevention

1x iBrite Deplaquer/Polisher Unit: 1x Handpiece, 1x Remi desensitizing paste pen, 4x AA batteries, 1x Shade guide, 6x Propy pastes, 5x Bids with adhesive stickers, 1x Instruction, 6x Propy angles (2-Pedo cups, 3-ProAngle Ergo Torque cups, 1-Brush cup) 35-55410 **199.45**

Integrity Hygiene Handpiece (Osseo Scientific)

Lightest Swivel Hygiene Handpiece!
Only 2.3oz (60g)

- Friction grip chucking, just push in or pull out the propy head
- Smooth, vibration-free operation up to 4,000 RPM
- Ergonomic design provides perfect balance



- 360° swivel for superb maneuverability
- Compatible with all standard doriot-style disposable metal angles

Brand	Weight (oz)	Length (in)
Prophyhygiene™	2.3	3.9
NSK Hygiene Pro	2.5	4.0
Kavo PROPHYwiz	2.7	4.3
Midwest RDH	3.0	4.4

Integrity Prophyhygiene 35-117791 **344.95**

Micromax Cordless Prophylaxis Handpiece (Dentamerica)

- Cordless** convenience offers complete portability
- Heavy-duty high torque motor**
- Forward/reverse** functions
- Reciprocating setting available without use of additional attachment
- Nine rotation speeds** available for different applications
- Ergonomic design provides comfort and easy operation
- High capacity lithium-ion battery** provides stable and prolonged usage
- Offers quick and easy coupling
- Low battery warning light



Specifications:

- 100-240V AC 50/60Hz
- Dimensions: L30 x W27 x H197 mm
- Weight: 150g
- Accessories: Charging base, Handpiece, Power Cable

Unit with Latch Attachment	35-118541	674.95
Unit with Screw Attachment	35-118531	674.95
Unit with Snap-On Attachment	35-118521	674.95
Unit with Straight Attachment	35-118511	674.95
Latch Attachment	35-118501	184.95
Screw Attachment	35-118491	184.95
Snap-On Attachment	35-118481	184.95
Straight Attachment	35-118471	184.95

ProMate E Portable Hygiene Handpiece (Pac-Dent)

A completely portable and re-chargeable handpiece.

This portable hygienist's handpiece is the ultimate in convenience, comfort, and control.





Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

No more cords to mess with, lightweight, fully autoclavable, rechargeable, comfortable grip and contour, as well as minimal vibration and chatter. Smoothly control speed with a simple twist to your desired setting anywhere between 2,000 and 20,000 rpm.

Advantages:

- Cordless, portable, and rechargeable
- Textured silicon grip reduces vibration and chatter for enhanced comfort & non-slip control
- Fully autoclavable and lightweight
- Consistent, powerful torque delivery
- Variable speed between 2,000 and 20,000 rpm
- Interchangeable autoclavable grips allow color-coding, personalization & easy maintenance
- Works with any standard disposable prophylaxis angle on the market

1 X Portable Hygiene Handpiece
 with 3 X silicon grips 35-55400 **484.45**

Dusty Rose silicon grips for hygiene handpiece, 10/pk 35-55401 **20.95**

Green silicon grips for portable handpiece, 10/pk 35-55402 **20.95**

Blue silicon grips for portable handpiece, 10/pk 35-55403 **20.95**

Lavender silicon grips for portable handpiece, 10/pk 35-55404 **20.95**

ProMate™ Hygiene Prophylaxis Handpiece (Pac-Dent)

Reduce fatigue and daily micro-traumas on the hand and arm muscles



Reduce fatigue and daily micro-traumas on the hand and arm muscles with our new ProMate™ hygiene handpiece. It has a lightweight, ergonomic design that minimizes vibration and chatter as well as a textured silicon grip and contoured shape that provides superior comfort and non-slip control. Maintain the ProMate™ easily by autoclaving the entire handpiece with the grip, or remove and replace the grip for color coding and personalization. Use with your favorite standard disposable prophylaxis angles today and see the difference! The grips are available in blue, green, lavender and dusty rose.

Advantages:

- Very economical, best buy for price & quality
- Soft, non-slip grips help reduce fatigue and daily micro-traumas
- Color-coded grips are ideal for personalization
- Easy maintenance with the entire handpiece being autoclavable
- Ergonomic lightweight design ensures superior comfort
- Compared to the best brands on the market
- Works with any standard disposable prophylaxis angle

1 X Hygiene Handpiece with 3 X silicon grips 35-55391 **398.95**

Handpieces

Green silicon grips for hygiene handpiece, 10/pk 35-55392 **20.95**

Blue silicon grips for hygiene handpiece, 10/pk 35-55393 **20.95**

Lavender silicon grips for hygiene handpiece, 10/pk 35-55394 **20.95**

Dusty Rose silicon grips for hygiene handpiece, 10/pk 35-55395 **20.95**

Surgical Handpieces

Impact Air 45 Surgical Handpiece (Palisades)

- Impact Air 45 is designed for use in periodontal and endodontic procedures
- Unique 45 Degree angled head, allowing maximum access and visibility
- Available with or without fiber optics
- Highspeed Air Powered Handpiece
- Fully Autoclavable



4-Hole Surgical Handpiece 35-117771 **719.95**

Fiber Optic 5-Hole Surgical Handpiece 35-117761 **869.99**

2-Hole Surgical Handpiece 35-11775 **789.95**

Model 9 Push-Button Conversion Kit:
 (includes turbine, head cap and cap wrench) 35-11774 **299.95**

Model 9 Push-Button Turbine (replacement) 35-11773 **266.99**

Model 9 Push-Button Chuck 35-11772 **99.95**

Model 9 Push-Button Head Cap with Cap Wrench 35-11771 **48.99**

IMPACT AIR 45 Autoclavable Lubricant 35-11770 **78.99**

Surgical Push Button Handpiece



4-Hole 35-16511 **199.95**

Surgical Push Button Handpiece Pana MAX



45 Degree Angle

Made in Japan

Swivel adapter included 35-16521 **394.95**

Handpieces

Surgical Push Button Handpiece 45 Degree

6 Month
Warranty



4-Hole Surgical Handpiece with Swivel 35-121131 **299.95**

Twin Power Surgical Handpiece (J-Morita)

- Safe, rear-facing exhaust vents
- High torque: up to 20 watts of power
- Quiet & long-lasting ceramic ball bearings



Handpiece 4H PAR-4HEX-O-45
(Morita Couplings) 35-98821 **979.95**

Handpiece 4H PAR-4HEX-O-KV-45
(Kavo MULTiflex LUX Coupling) 35-98822 **964.95**

Turbine 4HEX-45 Capsule Cartridge
(Ceramic Bearing) 35-98823 **379.95**

Electric Motors

Rotex 782 Electric Micro Motor (Dentamerica)



- Ball bearing high torque electric motor
- 1,000 to 35,000 R.P.M. variable speed
- Double insulation heavy duty construction
- Automatic cooling fan for continuous high speed performance
- Forward/Reverse switch Standard package includes Foot Switch, Handpiece, and Handpiece Rest
- Circuit overload protection
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics, and steel
- Unit is available for 100/110V 50/60 Hz, or 220/240V 50 Hz operation

•1 year warranty

Unit with Electric Handpiece and Variable Speed Foot Control 35-120211 **454.95**

Unit with E-Type Straight, Contra Angle Attachment and Variable Speed Foot Control 35-120201 **519.95**

Handpiece only with cable for Rotex 782 35-120191 **294.95**

Handpiece only with cable for Rotex 782E 35-12018 **199.95**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

M35 Compact Electric Lab Handpiece System (Buffalo)

New High Torque Variable Speed System Features a Compact and Versatile Design



Console: 5 in. (128 mm) Width x 6.5 in. (166 mm) Depth x 3.25 in. (82 mm) Height x 2.9 lb. (1.3 kg) Weight

- High Torque Cutting Power from 0 – 35,000 RPM New high torque handpiece and control system provide excellent cutting power and smooth control in a compact system
- Improved Versatility System is designed for 120V AC operation but can be easily switched to 220V AC via setting on the underside of the Console
- Improved Console and Controls Improved controls provide fully variable speed via console or foot control, simplify operation and protect the controls and handpiece
- Excellent High and Ultra-Low Speed Operation System excels at providing a full range of usable operating speeds from ultra-low to 35,000 RPM
- Heavy Duty Construction Components and system designed to provide reliable and trouble-free operation. **1 year warranty**
- Sizes:
Handpiece: 6.25 in. (159 mm) Long x 1.1 in. (28.5 mm) Back Diameter .79 in. (20 mm) Front Diameter x 8 oz. (.23 kg) Weight

120V/220V AC
(Console+MG Gray HP+FC) 35-118131 **419.95**

Magic Wand Lab Handpiece System (Buffalo)

Features/Benefits:

- Compact electric handpiece system features fully variable cutting speeds from 0 to 35,000 RPM
- Premium black electric handpiece provides smooth, quiet and comfortable operation with a history of long term durability
- Reliable, no-frills, made-in-USA Console
- Fully variable speed conveniently controlled by a heavy-duty, all-metal Variable Speed AC Foot Control
- The convenience and benefits of an electric handpiece without the added extras of higher priced systems
- Sizes



Handpiece: 6 in. (153 mm) Long x 1 1/8 in. (29 mm) Major Diameter, 5/8 in. (16 mm) Nose Diameter x 9.5 oz. (270 gm) Weight

Console: 4.25 in. (108 mm) Width x 5.25 in. (133 mm) Depth x 3 1/8 in. (80 mm) Height, 2 lb. 7 oz. (1.1 kg)

120V AC (Console + Black HP + FC) 35-118151 **449.95**

220V AC (Console + Black HP + FC) 35-118141 **449.95**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

V35 Electric Lab Handpiece System 110V (Buffalo)

Features/Benefits:

- Reliable Gray Electric Handpiece, New design and high-torque motor provide smooth operation and powerful performance with less operator fatigue. Larger bearings, extra-strong chucking mechanism, and proven design promote increased reliability
- Variable Speed Foot Control, Variable speed foot control provides fully variable speed operation from 0 – 35,000 rpm and allows the system to excel at both higher speed and ultra-low speed (0 – 500 rpm) operation
- Dependable V35 Control Console, Four multi-speed options on the control console (14k, 22k, 27k 35k rpm) simplify console electronics, improve system reliability, and allow ideal speed selection for most cutting needs. Simplified design creates higher power transfer through the console and more torque to the handpiece. Console **Made in USA!**



120V AC

Console + Gray HP + Foot Control 35-12154 **449.95**

220V AC

Console + Gray HP + Foot Control 35-12153 **449.95**

X35 Premium Electric Lab Handpiece System (Buffalo)

•Upgraded Handpiece:

The new black signature handpiece features an improved heavy-duty spindle, upgraded bearings, and more durable square carbon brushes, providing smoother operation and increased reliability

•**Improved Console:** The upgraded control console features dual power outlets allowing fast and easy switching between two separate handpieces as desired (one handpiece included with system), saving time-wasting tool changes. Added overload protection helps protect the handpiece from damage during extended heavy grinding

•**Superior Control:** The system provides exceptional control at ultra-low Rotational speeds, high torque cutting power, and fully variable speed from 0-35,000 rpm. The X35 comes complete with a fully variable speed foot control for added versatility

•**Heavy Duty Construction:** Components and system designed to provide reliable and trouble-free operation. **1 year warranty**

•Sizes

Handpiece: 6 in. (153 mm) Long x 1 1/8 in. (29 mm) Major Diameter, 5/8 in. (16 mm) Nose Diameter x 9.5 oz. (270 gm) Weight

Console: 7.25 in. (184 mm) Width x 7 in. (177 mm) Depth x 3.5 in. (89 mm) Height x 6.3 lb. (2.9 kg) Weight

120V AC

(Console + Black HP + Foot Control) 35-12156 **639.95**

220V AC

(Console + Black HP + Foot Control) 35-12155 **639.95**



Handpieces

X40 Premium Electric Lab Handpiece System (Buffalo)

New 40,000 RPM High-Torque System Updates Our Classic X35 and Provides Additional Speed and Torque



- New 40,000 RPM Handpiece 40k blue electric handpiece features an improved heavy-duty spindle, upgraded bearings, and more durable square carbon brushes, providing higher speeds, smoother operation, and increased reliability
- Improved Console** The upgraded control console features dual power outlets allowing fast and easy switching between two separate handpieces as desired (one handpiece included with system), saving time-wasting tool changes. Higher voltage output and additional overload control help protect the handpiece from damage during extended heavy grinding
- Superior Control:** The system provides exceptional control at ultra-low rotational speeds, high torque cutting power, and fully variable speed from 0-40,000 rpm. The X40 comes complete with a fully variable speed foot control for added versatility
- Heavy Duty Construction:** Components and system designed to provide reliable and trouble-free operation. **1 year warranty**
- Sizes**

Handpiece: 6.25 in. (159 mm) Long x 1 1/8 in. (29 mm) Major Diameter x 5/8 in. (16 mm) Front Diameter x 8.6 oz. (244 gm) Weight

Console: 7.25 in. (184 mm) Width x 7 in. (177 mm) Depth x 3.5 in. (89 mm) Height x 6.3 lb. (2.9 kg) Weight

120V AC

(Console + 40k rpm HP + Foot Control) 35-12158 **669.95**

220V AC

(Console + 40k rpm HP + Foot Control) 35-12157 **669.95**

X40 40,000 rpm Blue Electric Handpiece Only- New!

35-63021 **369.95**

X50 Brushless Electric Lab Handpiece System (Buffalo)

- Ultra-High Torque Cutting Power Brushless high-torque motor and control system provide brute-strength 7.8 N. cm. torque and precise control from 1,000 to 50,000 rpm with no carbon brush maintenance



- New Whisper-Quiet Motor New brushless motor features ultra-quiet and extremely smooth operation. Improved bearings and balance provide highly-concentric, vibration-free operation in a sealed, maintenance-free design. Convenient twist lock/unlock chuck provides safe and secure holding of HP-sized instruments
- New Advanced Console Design with Microcomputer Control New intelligent microprocessor control provides fully variable speed via console or foot control. Console provides precise speed control, digital rpm display, and controls for hand/foot and forward/reverse operation.

Handpieces

- Foot Control "Smart Cruise" Feature New smart cruise control circuitry automatically senses when the Foot Control speed is held constant and auto-activates precise rpm control with an audible two beep tone. Feature immediately deactivates with next touch of Foot Control and ready for resetting as desired
- Heavy Duty Construction Components and system designed to provide safe, reliable and trouble-free operation. High impact housing is tough and easy to clean. **1 year warranty**
- Improved Safety Intelligent safety circuitry is designed to protect the computer controls and the handpiece. The safety-lock on the console speed control dial helps prevent unintentional higher speed operation

•Sizes

Handpiece: 6.25 in. (159 mm) Long x 1.125 in. (29 mm) Back Diameter, .75 in. (19 mm) Front Diameter x 8.2 oz. (232 gm) Weight

Console: 7 in. (178mm) Width x 9.5 in. (241 mm) Depth x 4.25 in.

120V AC

50,000 RPM Electric Handpiece System 35-12160 **899.95**

220V AC

50,000 RPM Electric Handpiece System 35-12159 **929.95**

207S Series (Saeshin)



207S SET B 106HP 35K Variable 35-120351 **295.00**

207S SET B 106HP 45K Variable 35-120341 **295.00**

207S SET B 107HP 35K Variable 35-120361 **295.00**

207S SET B OPTIMUS HP 35K Variable 35-120371 **295.00**

207S SET B OPTIMUS HP 45K Variable 35-120331 **295.00**

107S Control box only 35-120321 **90.00**

E-File Portable (Saeshin)



RPM Portable E-File Set in Gray Willow Set 35-120241 **299.00**

Ergo Slim Handpiece Only in Gray Willow, Each 35-120231 **199.00**

RPM Ergo Slim E-Type Motor in Gray, Each 35-120251 **149.99**

RPM On/Off Pedal for Table, Each 35-120221 **10.00**



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Handy I Portable Micromotor (Saeshin)

Handy I Portable Micromotor Set, 20K Set 35-11739 **339.99**

Handy I Control Box, 20K Each 35-11742 **125.00**

Handy E-Type Handpiece, 20K RPM Each 35-11741 **140.00**

1:1 Straight Nose Cone, 35K Each 35-117441 **105.00**

Handy I Control Box, 20K with Charger Each 35-11743 **145.00**

Handy I Charger for Handy I Controller, Each 35-117451 **22.00**

Charger for Handy I Controller, 230V Each 35-11740 **22.00**

Pouch for Handy I Controller Each 35-11738 **15.00**

Micro Motor 210 (Saeshin)

- Variable Speed dial system
- ON/Off Variable Speed foot control
- Right & Left turning ability
- Max, 45,000 rpm
- Torque 280 gf.cm
- Input 110 / 230 Volts (50 / 60Hz)
- Output DC 0-32W
- Weight: Control box: 2.7 Lb Handpiece: 7.1 oz. (Strong 102L)



Ram 210 35K Control Box, Tech2000 HP ON/OFF 35-119901 **185.00**

Ram 210 35K Control Box, Tech2000 HP Variable 35-119891 **210.00**

Ram 210 35K Control Box, Tech2000 HP None 35-119911 **175.00**

Ram 210 35/45 Control Box Only 35-119921 **75.00**

Ram 210 45K Control Box, 45000 HP ON/OFF 35-119931 **185.00**

Ram 210 45K Control Box, 45000 HP Variable 35-119881 **210.00**

Ram 210 45K Control Box, 45000 HP None 35-119871 **175.00**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Handpieces

Micro-Cad Pocket Micromotor (Saeshin)

Saeshin Product

MICRO-CAD
POCKET MICROMO-
TOR - 12V - 20,000
RPM (RAM)



Complete set includes controller, handpiece, rubber stand and charger

Controller:

- Forward/Reverse
- Variable speed dial
- Runs 5 - 7 hours on a single charge
- Charging time is 5 - 7 hours
- Battery recharge - 500 times +
- 110 volts

Handpiece:

- Length - 5" x Height - 1 1/3" x Width - 2 1/3"
- 3,000 - 20,000 rpms

Micro-Cad E-Type With Straight Nose Cone Complete Set, 110V	35-118331	464.95
Micro-Cad E-Type With Straight Nose Cone Complete Set, 220V	35-118341	464.95
Micro-Cad Complete With Lab Style Handpiece Set 110V	35-118321	399.95
Micro-Cad Complete With Lab Style Handpiece Set, 220V	35-118311	399.95
Micro-Cad Controller Only Each	35-118301	132.95
Micro-Cad E-type Motor Only Each	35-118291	209.95
One-Piece Micro-Cad Lab Style Handpiece Only, Each	35-118371	279.95
1:1 Straight Nose Cone 35K, Each	35-118351	104.95
Micro-Cad Charger (110V Only), Each	35-118361	21.95
110V-230V Micro-Cad Universal Charger, Each	35-118281	24.95

Microlab Brush Motor (Saeshin)

- Feedback control system
- Variable Speed dial system
- Right & Left turning ability
- Hand/Foot speed control system
- Automatic overload protection system
- 45,000 rpm high-powered carbon brush motor
- ON/Off and variable speed foot control
- Torque 450 gf.cm
- Input 110 / 230 Volts (50 / 60Hz)



- Output DC 0-32W
- Weight Control box: 3.82 Lb Handpiece: 8.0 oz. (Strong 106)

MicroLab 350 Control Box, Tech2000 HP ON/OFF	35-118411	259.00
MicroLab 350 Control Box, Tech2000 HP Variable	35-118401	285.00
MicroLab 350 Control Box, Tech2000 HP None	35-118391	249.00
MicroLab 350/450 Control Box Only	35-118421	120.00
MicroLab 450 Control Box, 45,000 RPM HP ON/OFF	35-118431	259.00
MicroLab 450 Control Box, 45,000 RPM HP Variable	35-11844	285.00
MicroLab 450 Control Box, 45,000 RPM HP None	35-118451	249.00
MicroLab 350/450 CB, E-Type HP, 1:1 Straight Nose Cone ON/OFF	35-118461	305.00
MicroLab 350/450 CB, E-Type HP, 1:1 SNC, 1:1 CA ON/OFF	35-118381	360.00

Microtorque Motor (Saeshin)



Microtorque Control Box, Tech2000 Handpiece Variable	35-118731	495.00
Microtorque Control Box, Slim Nose Handpiece Variable	35-118741	495.00
Microtorque Control Box, 45,000 RPM Handpiece Variable	35-118721	495.00
Microtorque Control Box, Tech2000 Handpiece ON/OFF	35-118711	435.00
Microtorque Control Box, Slim Nose Handpiece ON/OFF	35-118701	435.00
Microtorque Control Box, 45,000 RPM Handpiece ON/OFF	35-118691	465.00
35,000 rpm Microtorque Control Box Only (110V)	35-118681	175.00
35,000 rpm Microtorque Control Box Only (230V)	35-118761	175.00
45,000 rpm Microtorque Control Box Only (110V)	35-118751	175.00
45,000 rpm Microtorque Control Box Only (230V)	35-118671	175.00

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Motor 210 (Saeshin)

RPM max. 45.000 rpm
Torque 280 gf.cm
Input 110/ 230 Volts (50/60Hz)
Output DC 0-32V
Weight: Control box: 2.7 lb
Handpiece: 7.1oz (Strong 102L)
•Right & Left turning ability
•On/Off and variable speed foot control



Saeshin 210 35K Control Box,
Tech2000 HP ON/OFF 35-118811 **185.00**

Saeshin 210 35K Control Box,
Tech2000 HP Variable 35-118821 **210.00**

Saeshin 210 35K Control Box,
Tech2000 HP None 35-118831 **175.00**

Mystisa Pocket Controller Motor (Saeshin)

Includes:

- MYSTISA POCKET CONTROL BOX
- MYSTISA HANDPIECE
- Length - 5" x Height - 1 1/3" x Width - 2 1/3"
- A/C charger (110v or 230v)
- Leather pouch with belt clip
- Handpiece rubber stand
- Carry Bag
- Under Nail Carbide Bur
- Large Barrel Carbide Bur
- 10 Sanding Bands
- Sanding Band Mandre



Mystisa Pocket Micromotor Set Complete Set 35-118851 **399.95**

Mystisa Handpiece Only Each 35-118861 **224.95**

Mystisa Pocket Controller Only ,
Gray Each 35-118871 **134.95**

110V-230V Universal Charger
For Mystisa Pocket Cont Each 35-118881 **24.95**

Gray Leather Pouch
For Mystisa Controller Each 35-118841 **15.95**

Mystisa Pocket Controller, Thumb Lab HP 35-120000 **399.95**

Mystisa Pocket Controller, Thumb E-Type HP 35-120002 **449.95**

Mystisa Pocket Drill (Saeshin)

- Hands Free Portability
- Speed range from 3,000 to 30,000 Rpm
- Runs on AC or battery power
- Runs 5 to 7 hours on a single charge
- Forward/Reverse direction
- 24V/30K Rpm



Mystisa Table Micromotor Set Complete Set 35-118891 **429.95**

Mystisa Handpiece Only Each 35-118911 **225.00**

Mystisa Table Control Box Only,
Blue Each 35-118921 **140.00**

Mystisa Table Control Box Only,
Pink Each 35-118931 **140.00**

Mystisa Table Control Box Only,
White Each 35-118941 **140.00**

Charger For Mystisa Table Control Box
Each 35-118901 **29.95**

Nextgen Brushless Motor (Saeshin)

Nextgen Brushless Motor (Ram)

- Ergonomic design, Excellent durability, Minimal noise and vibration
- Right & Left turning ability
- Hand/Foot speed control system
- Built-in micro processor
- Self diagnosis, Error display, Work memory storage and handpiece auto recognition function
- Variable speed foot control (Optional)
- Max, 40,000 rpm
- Input 110 / 230 Volts (50 / 60Hz)
- Torque 4 N.cm
- Output 150W
- Weight Control box: 4.2Lb Handpiece: 8.0 oz. (Forte 300II)



NextGen Micromotor Set with
Brushless Handpiece Variable 35-119001 **695.00**

NextGen Micromotor Set with
E-type Handpiece Variable 35-118981 **765.00**

NextGen Control Box Only Each 35-118971 **365.00**

NextGen Variable Foot Pedal Only Each 35-118991 **85.00**

NextGen Brushless Handpiece Only Each 35-118961 **340.00**

NextGen E-Type Handpiece Only Each 35-118951 **255.00**

Optimus 35/Optimus 45 HP (Saeshin)



206 Control Box,
OPTIMUS 35 HP ON/OFF 35-119071 **420.00**

206 Control Box,
OPTIMUS 35 HP Variable 35-119081 **410.00**

206 Control Box, OPTIMUS 35 HP None 35-119091 **395.00**

206 Control Box,
OPTIMUS 45 HP ON/OFF 35-119101 **420.00**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Handpieces

206 Control Box, OPTIMUS 45 HP Variable	35-119111	410.00
206 Control Box, OPTIMUS 45 HP None	35-119061	395.00

OZ Plus Brushless Motor (Saeshin)

RPM max.

50.000 rpm

Torque 7.8

N.cm

Input 110/ 230

Volts (50/60Hz)

Output 230W

Weight: Control box: 5.5

lb Handpiece: 8.9oz (Forte100)

•Right & Left turning ability

•Hand/Foot speed control system

•Self-diagnosis, Error display, Work memory

•Storage function

•Variable speed foot control (optional)

OzPlus Set With

Sensor 50K Set B Variable " with Pedal " 35-119141 **795.00**

Oz Plus Set With Sensor 50K Set C None "
without Pedal " 35-119131 **765.00**

Oz Plus Controller Only (Sensor 50K)
Each 35-119151 **410.00**

Oz Plus 50K Handpiece w/Sensor Each 35-119161 **340.00**

Variable Pedal Only For Oz Plus Each 35-119121 **50.00**



Power Brush Type Motor (Saeshin)

•Variable speed dial system

•Right & Left turning ability

•Hand/Foot speed control system

•Automatic overload protection system

•Reset apparatus ability

•ON/Off and variable speed foot control

•Max, 45,000 rpm

•Torque 450 gf.cm

•Input 110 / 230 Volts (50 / 60Hz)

•Output DC 0-32W

•Weight Control box: 4.8Lb Handpiece:
6.7 oz. (strong 110LN)

Rampower 35 Control Box,
Tech2000 HP ON/OFF 35-120021 **259.00**

Rampower 35 Control Box,
Tech2000 HP Variable 35-120011 **285.00**

Rampower 35 Control Box,
Tech2000 HP None 35-120001 **249.00**

Rampower 35/45 Control Box Only 35-119991 **135.00**

Rampower 45 Control Box,
45,000 RPM HP ON/OFF 35-119981 **259.00**

Rampower 45 Control Box,
45,000 RPM HP Variable 35-119971 **285.00**



Ram Power 45 Control Box,
45,000 RPM HP None 35-119961 **249.00**

Ram Power 35 CB with
Optimus 35K HP ON/OFF 35-119951 **375.00**

Ram Power 35 CB with
Optimus 35K HP Variable 35-120061 **385.00**

Ram Power 35 CB with
Optimus 35K HP None 35-120031 **299.00**

Ram Power 45 CB with
Optimus 45K HP ON/OFF 35-120041 **375.00**

Ram Power 45 CB with
Optimus 45K HP Variable 35-120051 **385.00**

Ram Power 45 CB with
Optimus 45K HP None 35-119941 **299.00**

RamPro500N Series (Saeshin)

•Computer controlled microprocessor

•Soft touch keypad control digital
speed display

•Hand/foot mode and forward/reverse
direction

•0 - 35,000 or 45,000 RPMs

•Auto cruise control on/off button for
continuous power in foot pedal mode

•Four memory speed reset buttons

•Dual handpiece output jacks

•Safety overload protection

•Forward/Reverse direction

•Hand/Foot selection



RamPro500N Control Box,
Tech2000 HP ON/OFF 35-120081 **425.00**

RamPro500N Control Box,
Tech2000 HP Variable 35-120091 **450.00**

RamPro500N Control Box,
Tech2000 HP None 35-120101 **399.00**

RamPro500N Control Box Only 35-120111 **185.00**

RamPro500N Control Box,
45,000 RPM HP ON/OFF 35-120121 **435.00**

RamPro500N Control Box,
45,000 RPM HP Variable 35-120131 **460.00**

RamPro500N Control Box,
45,000 RPM HP None 35-120071 **409.00**

Strong 204 Brush Motor (Saeshin)

Control Box:

•35,000 rms

•Torque: 280 gf.cm

•Input 110 - 230 volts

•Compact design

Tech2000 Handpiece:

•35,000 RPM

•Twist type chuck release



Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

- Improved dust seal protection
- Includes 3/32" chuck
- Multi-chuck feature (easy to change chuck)
- Optional chuck sizes available: 1/8", 1/16", or 3MM

On/Off Pedal

- Maintains constant speed

204 35K Control Box,
Tech2000 HP ON/OFF 35-120971 **179.95**

204 35K Control Box,
Tech2000 HP None 35-120981 **169.95**

204 45K Control Box,
45000 RPM HP ON/OFF 35-120991 **179.95**

204 45K Control Box,
45000 RPM HP None 35-121001 **172.95**

Strong 207A Brush Motor (Saeshin)

Control Box

- 35,000 rpms
- Torque: 300 gf.cm
- DC 0-32V
- Input 110 - 230 volts
- Overload protection

Tech2000 Handpiece

- 35,000 RPM
- Twist type chuck release
- Improved dust seal protection
- Includes 3/32" chuck
- Multi-chuck feature (easy to change chuck)
- Optional chuck sizes available: 1/8", 1/16", or 3MM

On/Off Pedal

- Maintains constant speed

207A 35K Control Box,
Tech2000 HP ON/OFF 35-121021 **189.95**

207A 35K Control Box,
Tech2000 HP Variable 35-121031 **204.95**

207A 35K Control Box,
Tech2000 HP None 35-121041 **179.95**

207A 35/45 Control Box Only 35-121051 **89.95**

207A 45K Control Box,
45000 RPM HP ON/OFF 35-121061 **189.95**

207A 45K Control Box,
45000 RPM HP Variable 35-121071 **204.95**

207A 45K Control Box,
45000 RPM HP None 35-121011 **184.95**

Strong 90 Carbon Brush Motor (Saeshin)

Basic model, Mini-Compact design,
Less vibration standard carbon brush
motor

- Variable Speed dial system
- Right & Left turning ability
- ON/Off Foot control
- Max, 35,000 rpm
- Torque 280 gf.cm



- Input 110 / 230 Volts (50 / 60Hz)
- Output DC 0-32W
- Weight: Control box: 2.9 Lb Handpiece: 7.1 oz. (Strong 102)

SS90 Control Box,
Tech2000 HP ON/OFF 35-121081 **175.00**

SS90 Control Box,
5000 RPM HP ON/OFF 35-121101 **175.00**

SS90 Control Box Only 35-121111 **80.00**

SS90 On/Off Foot Pedal Only ON/OFF 35-121091 **21.00**

Thumb Portable Micromotor (Saeshin)

Saeshin Product

Thumb Portable Micromotor 24V
30,000 Rpm (Ram)



Set Consists of:

- Thumb Control Box, Charger,
- Thumb Pouch, Either an E-Type Handpiece with a 1:1 Straight Nose Cone or Lab Style Handpiece
- Recharger Battery takes 5 Hours to Fully Charge, Up to 8 Hours Battery Life, Variable Speed Control
- Light Weight Handpiece, Chuck Locking, Easy Bit Change, High Torque, Quiet Operations, No Vibration

Thumb Lab Type Portable Set 35-121201 **399.95**

Thumb E-Type Portable Set With
Straight Nose Cone Set 35-121191 **459.95**

Thumb Lab Style Motor Handpiece Only 35-121181 **279.95**

Thumb E-Type Motor Handpiece Only 35-121171 **259.95**

Thumb Control Box Only Each 35-121161 **194.95**

Thumb Charger/Adaptor Each 35-121211 **24.95**

Thumb Pouch Each 35-121231 **24.95**

Thumb Battery Each 35-121221 **74.95**

Xi Brush Type Motor (Saeshin)

- Compact / Eegonomic design
- Vertical / Horizontal placement
- Right & Left turning ability
- Hand/Foot speed control system
- Circuit safety system
- Variable Speed dial system (209A)
- ON/Off and variable speed foot control
- Max, 35,000 rpm
- Torque 300 gf.cm
- Input 110 / 230 Volts (50 / 60Hz)
- Output DC 0-32W
- Weight: Control box: 2.6 Lb Handpiece: 8.2 oz. (Strong 110LN)



XI Control Box, Tech2000 HP ON/OFF 35-12165 **185.00**

XI Control Box, Tech2000 HP Variable 35-12164 **199.00**

XI Control Box, Tech2000 HP None 35-12163 **175.00**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Control Box Only Gray 32V 35/45K	35-12162	90.00
XI Control Box , 45000 RPM HP ON/OFF	35-12161	185.00
XI Control Box , 45000 RPM HP Variable	35-12168	199.00
XI Control Box , 45000 RPM HP None	35-12167	175.00
XI Control Box , 45000 RPM HP, 2 Collets None	35-12166	199.00

Handpiece Bulbs

LED Bulbs (Chapman-Huffman)



LED Bulb 16-120-00
 Fits following Swivel Couplers Chapman-Huffman 20-079-00KAVO MultiflexMK-DENT QC5016, QC5016KKAVO.



LED Bulb 16-126-00
 Fits KAVO Motors



LED Bulb 16-130-00
 Fits SIRONA Quick Connectors



LED Bulb 16-134-00
 Fits W&H Turbines



LED Bulb 16-132-00
 Fits SIRONA Turbines
 and SIRONA Motors



LED Bulb 16-122-00
 Fits BIEN-AIR Motors



LED Bulb 16-124-00
 Fits BIEN-AIR Quick Connectors



LED Bulb 16-128-00
 Fits NSK Quick Connector



LED Bulb 16-136-00
 Fits STAR Quick Connector

Fits following: Swivel Couplers
 Chapman-Huffman 20-079-00
 KAVO Multiflex MK-DENT
 QC5016, QC5016K KAVO cross
 ref. # 1007.5372

	35-118061	Call for Price
Fits KAVO Motors	35-118051	Call for Price
Fits SIRONA Quick Connectors	35-118041	Call for Price
Fits SIRONA Turbines and SI- RONA Motors	35-118031	Call for Price
Fits W&H Turbines	35-118021	Call for Price
Fits BIEN-AIR Motors	35-118011	Call for Price
Fits BIEN-AIR Quick Connec- tors	35-118001	Call for Price
Fits NSK Quick Connector	35-117991	Call for Price
Fits STAR Quick Connector	35-118071	Call for Price

Handpieces

Handpieces Turbines & Couplers Swivels

Athena Champion Push Button Turbine (HPP)



Cap Wrench	35-112711	99.95
Push Button cap	35-112721	14.95
Complete Turbine Assembly	35-112731	9.95

Athena Champion Standard Turbine (HPP)



Cap Wrench	35-11266	6.95
Cap	35-112671	6.95
Bur Wrench	35-112691	9.95
Complete Turbine Assembly	35-112681	49.95

Bien Air Black Pearl Friction Grip Turbine (HPP)



Complete Turbine Assembly	35-112771	89.95
----------------------------------	-----------	--------------

Bien Air Black Pearl Lab Turbine (HPP)



Bearing Kit	35-112821	24.95
Complete Turbine Assembly	35-112831	99.95

Bien Air Black Pearl Push Button Turbine (HPP)



Complete Turbine Assembly	35-113011	89.95
----------------------------------	-----------	--------------

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Bien Air Bora S36-S36L Push Button Turbine (HPP)



Bearing Kit	35-112801	24.95
Complete Turbine Assembly	35-112811	94.95

Bien Air Ondine Friction Grip Turbine (HPP)



Bearing Kit	35-112951	28.95
Complete Turbine Assembly	35-112971	89.95

Bien Air Ondine Push Button Turbine (HPP)



Bearing Kit	35-112891	29.95
Complete Turbine Assembly	35-112881	109.95

Bien Air Prestige Push Button Turbine (HPP)

Bearing Kit	35-97572	29.95
Complete Turbine Assembly	35-97571	109.95

Canister Turbine (HPP)



Push Button Turbine	35-11308	88.00
Push Button Turbine, Angular H	35-11306	94.85

Push Button Turbine, Angular M	35-11307	99.85
Push Button Turbine, Phenolic	35-11305	98.85
Push Button Turbine, Angular Phenolic H	35-11303	99.85
Push Button Turbine, Angular Phenolic M	35-11304	99.85
Push Button Turbine, Silver / Maintenance free	35-113101	99.85
Push Button Turbine, Maintenance Free	35-113091	109.85
Push Button End Cap	35-113111	96.90
Cap Wrench	35-11302	5.85

Canister Turbine Mini (HPP)



Mini Screw Type Turbine H	35-11320	35.99
Mini Screw Type Turbine M	35-11319	35.99
Mini Screw Type Turbine, Angular H	35-11317	39.99
Mini Screw Type Turbine, Angular M	35-11318	47.99
Mini Screw Type Turbine, Phenolic H	35-11316	41.85
Mini Screw Type Turbine, Phenolic M	35-11314	41.85
Mini Screw Type Turbine, Angular Phenolic H	35-11315	49.85
Mini Screw Type Turbine, Angular Phenolic M	35-11313	49.85
Mini Screw Type Turbine, Silver / Maintenance Free	35-113241	49.85
Mini Screw Type Turbine, Maintenance free	35-11321	64.85
Mini Bur Wrench	35-113221	8.00
Mini Screw Type End Cap	35-113231	5.85
Cap Wrench	35-11312	6.00

Canister Turbine Standard (HPP)




Screw Type Turbine H	35-113471	39.99
Screw Type Turbine M	35-113451	39.99



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Screw Type Turbine, Angular H	35-11344	44.99
Screw Type Turbine, Angular M	35-11343	48.99
Screw Type Turbine, Phenolic H	35-11341	43.99
Screw Type Turbine, Phenolic M	35-11342	22.99
Screw Type Turbine, Angular Phenolic H	35-11340	49.99
Screw Type Turbine, Angular Phenolic M	35-11338	49.99
Screw Type Turbine, Silver / Maintenance Free	35-11339	49.99
Screw Type Turbine, Maintenance Free	35-11336	63.99
Bur Wrench	35-11337	7.90
Screw Type End Cap	35-113461	4.99
Cap Wrench	35-11335	5.99

Cartridge Turbine for Handpiece (Nakamura)



Made in Japan

Standard Chuck, For TC-40PM / 40PB & TC-80QM / 80QB	35-11176	74.95
Large Size, Standard Chuck, For TC-35YM / 35YB	35-111751	89.95
Push Button Chuck, For TC-450M / 450B & TCP-70QM / 70QB	35-111741	89.95
Push Button Chuck, For TC-350M / 350B	35-111731	94.95
Cartridge For TS-300ZT	35-111732	119.95

Impact Air Push Button Turbine (HPP)



Complete Turbine Assembly	35-97581	99.95
Complete Turbine Assembly with Cap	35-97582	109.95
Bearing Kit	35-97583	29.95

Kavo 6000 Series Push Button Turbine (HPP)



Complete Turbine Assembly	35-97611	109.95
---------------------------	----------	---------------

Handpieces

Kavo 625/630/640 Standard Head Friction Grip Turbine (HPP)



Bearing Kit	35-114021	34.95
Complete Turbine Assembly	35-114031	89.95

Kavo 625/630/640 Standard Head Push Button Turbine (HPP)



Bearing Kit	35-97592	34.95
Complete Turbine Assembly	35-97591	99.95

Kavo 632/633 Miniature Head Friction Grip Turbine (HPP)

Bearing Kit	35-97602	34.95
Complete Turbine Assembly	35-97601	89.95

Kavo 632/642/643 Miniature Head Push Button Turbine (HPP)



Bearing Kit	35-113651	34.95
Complete Turbine Assembly	35-113661	104.95

Kavo 634 Miniature Head Friction Grip Turbine (HPP)



Bearing Kit	35-113811	34.95
Complete Turbine Assembly	35-113801	94.95

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Kavo 634 Miniature Head Push Button Turbine (HPP)



Bearing Kit	35-11388	34.95
Complete Turbine Assembly	35-113901	104.99

Kavo 635 Mira Lux 3 Push Button Turbine (HPP)



Push Button Cap	35-11391	34.95
Bearing Kit	35-113941	44.95
Complete Turbine Assembly	35-113931	104.95
Cap Wrench	35-11392	10.95

Kavo 647/649 Lux 3 Standard Head Push Button Turbine (HPP)



Push Button Cap	35-113713	28.95
Bearing Kit	35-113711	38.95
Complete Turbine Assembly	35-113721	109.95
Cap Wrench	35-113712	10.95

KAVO 650 Type Ceramic Turbine (Sable)

12 Month Warranty

Made in Canada



Each	35-54251	224.95
------	----------	--------

Kavo 6500 Series Push Button Turbine (HPP)

Complete Turbine Assembly	35-97621	109.95
---------------------------	----------	--------

Kavo 659 Compact Push Button Turbine (HPP)



Complete Turbine Assembly	35-114141	104.95
---------------------------	-----------	--------

Kavo Style Coupler (Dentex)

Non-Optic



Standard Head Non-Optic	35-117941	150.00
Mini Head Non-Optic	35-117951	150.00

Kavo Style Fiber Optic Coupler (Dentex)



Standard Head Fiber Optic	35-117961	195.00
Mini Head Fiber Optic	35-117971	195.00

Lares 557 Small Head Push Button Turbine (HPP)



Push Button Cap	35-11416	31.95
Bearing Kit	35-114171	33.95
Complete Turbine Assembly	35-114181	109.95
Back Cap Wrench	35-11415	10.95

Lares 757 Large Head Push Button Turbine (HPP)



Complete Turbine Assembly	35-114211	109.95
Bearing Kit	35-114201	33.95
Back Cap	35-11419	54.95
Cap Wrench	35-12419	10.95



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Handpieces

Midwest 8000/8000i Standard Turbine (HPP)



Complete Turbine Assembly	35-114341	79.95
Bearing Kit	35-11432	24.95
Bur Wrench	35-114331	10.95
Cap Wrench	35-11422	8.95

Midwest Quiet Air Lever & Push Button Turbine (HPP)



Complete Turbine Assembly	35-114451	99.95
Complete Turbine Assembly with P.B. Cap	35-11442	109.95
Cap Wrench Lever Style	35-11443	11.95
Push Button Cap	35-114441	24.95
Cap Wrench	35-11436	9.95
Bearing Kit	35-11437	24.95

Midwest Quiet Air Standard Turbine (HPP)



Complete Turbine Assembly with Cap	35-114701	79.95
Bearing Kit	35-11459	19.95
Bur Wrench	35-11468	10.95
Cap Wrench	35-11458	10.95

Midwest Stylus Push Button Turbine (HPP)



Complete Turbine Assembly	35-114801	109.95
Push Button Cap	35-11477	35.95

Cap Wrench	35-11478	10.95
Bearing Kit	35-114791	53.95

Midwest Tradition Lever Turbine (HPP)



Complete Turbine Assembly	35-114911	99.95
Bearing Kit	35-11482	24.95
Cap Wrench Lever Style	35-11481	14.95

Midwest Tradition Standard Turbine (HPP)



Complete Turbine Assembly with Cap	35-115161	79.95
Bearing Kit	35-11514	24.95
Bur Wrench	35-115151	10.95
Cap Wrench	35-11504	9.95

Midwest Tradition XGT Push Button Turbine (HPP)



Complete Turbine Assembly	35-115281	104.95
Complete Turbine Assembly with Cap	35-11525	114.95
Bearing Kit	35-11526	26.95
Push Button Cap	35-115271	24.99
Cap Wrench XGT Lever	35-11517	9.95

Ney Hurricane Friction Grip Turbine (HPP)

Complete Turbine Assembly	35-97631	69.95
Bearing Kit	35-97632	39.95

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

NSK 75NL Push Button Turbine (HPP)

Complete Turbine Assembly 35-97671 **109.95**

NSK Mini Phatelus Push Button Turbine (HPP)



Complete Turbine Assembly 35-116011 **99.95**

Bearing Kit 35-116021 **28.95**

NSK Pana Air Mini Canister (HPP)



Complete Turbine 35-115411 **39.95**

Bearing Kit 35-11537 **19.95**

Bur Wrench 35-11539 **10.95**

Back Cap 35-115401 **8.95**

Back Cap Wrench 35-11529 **9.95**

NSK Pana Air Mini Head Push Button Canister (HPP)



Complete Turbine 35-113341 **79.95**

Bearing Kit 35-113321 **19.95**

Back Cap 35-113331 **13.95**

Cap Wrench 35-11325 **9.95**

NSK Pana Air Standard Head Canister (HPP)



Complete Turbine Assembly 35-115511 **39.95**

Bearing Kit	35-11548	19.95
Bur Wrench	35-11549	10.95
Back Cap	35-115501	8.95
Cap Wrench	35-11543	9.95

NSK Pana Air Standard Head Push Button Canister (HPP)



Complete Turbine Assembly 35-97641 **79.95**

Bearing Kit 35-97642 **22.95**

Cap Wrench 35-97643 **9.95**

Cap 35-97644 **14.95**

NSK Push and Swivel Push Button Turbine (HPP)



Complete Turbine Assembly 35-11559 **74.95**

Bearing Kit 35-115601 **19.95**

NSK Push and Swivel Standard Turbine (HPP)

Complete Turbine Assembly 35-97651 **54.95**

Bearing Kit 35-97652 **24.95**

NSK Style Coupler (Dentex)

Non-Optic



Standard Head Non-Optic 35-119031 **76.49**

Mini Head Non-Optic 35-119021 **76.49**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Handpieces

NSK Style Fiber Optic Coupler (Dentex)



Standard Head Fiber Optic	35-119041	96.00
Mini Head Fiber Optic	35-119051	96.00

NSK Super Grade Push Button Canister (HPP)



Complete Turbine Assembly	35-115711	74.95
Bearing Kit	35-11562	19.95

NSK Super Grade Standard Canister (HPP)



Complete Turbine Assembly	35-115811	44.99
Bearing Kit	35-115821	19.95
Bur Wrench	35-11572	10.95

NSK Torque Phatelus Push Button Turbine (HPP)




Complete Turbine Assembly	35-97661	94.95
Bearing Kit	35-97662	28.95

NSK Viper Phatelus Push Button Turbine (HPP)



Complete Turbine Assembly	35-115921	99.95
Bearing Kit	35-11590	28.95

Quest Highspeed Handpiece Couplers (Lares)



6 Pin Fiber Optic MX Coupler,

Multiflex Compatible Coupler (with 0.7 Amp Bulb)

35-119711 **259.95**

5 Holes MX Coupler, Multiflex Coupler with Fiber Optic Bundle

35-119721 **259.95**


4 Holes Coupler, Multiflex Compatible Coupler

35-119731 **194.95**

MX Coupler Fiber Optic (1.5 Amp Bulb) (Lares Power Supply Compatible)

35-119741 **259.95**

Star Concentrix Standard Turbine (HPP)



Complete Turbine Assembly	35-116211	39.95
Bur Wrench (Hex)	35-116212	10.95
Cap	35-116213	9.95

Star Futura I and II Standard Turbine (HPP)



Complete Turbine Assembly (Square)	35-116321	39.95
Complete Turbine Assembly (Hex)	35-116322	39.95
Bur Wrench (Square)	35-116323	10.95
Bur Wrench (Hex)	35-116324	10.95

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Star Ultra Standard Turbine (HPP)



Complete Turbine Assembly (Square)	35-116381	44.95
Complete Turbine Assembly (Hex)	35-116382	44.95
Bur Wrench (Square)	35-116383	10.95
Bur Wrench (Hex)	35-116384	10.95

Star Vista 430 Solara Push Button Turbine (HPP)



Complete Turbine Assembly (Greased Bearing)	35-116641	109.95
Complete Turbine Assembly (Lube-Free)	35-116642	119.95
Bearing Kit (Angular Greased Bearings)	35-116643	32.95
Bearing Kit (Lube-Free)	35-116644	54.95
Push Button Cap	35-116645	24.95

Star Vista 430 Standard Turbine (HPP)



Complete Turbine Assembly (Square)	35-116551	84.95
Complete Turbine Assembly (Hex)	35-116552	84.95
Bearing Kit (Angular Greased Bearings)	35-116553	32.95
Bearing Kit (Lube-Free)	35-116554	54.95
Bur Wrench (Square)	35-116555	9.95
Bur Wrench (Hex)	35-116556	54.95

W&H Adec 890 Series Large Push Button Turbine (HPP)



Complete Turbine Assembly	35-117001	99.95
Complete Turbine Assembly (Phenolic Bearing)	35-117002	104.95
Bearing Kit	35-117003	29.95

W&H Adec 890 Series Small Push Button Turbine (HPP)



Complete Turbine Assembly	35-117071	99.95
Complete Turbine Assembly (Phenolic Bearing)	35-117072	104.95
Bearing Kit	35-117073	11.95

W&H Adec Push Button Turbine (HPP)



Complete Turbine Assembly	35-116721	99.95
Complete Turbine Assembly (Phenolic Bearing)	35-116722	104.95
Bearing Kit	35-116723	29.95

W&H Adec Synea Large Push Button Turbine (HPP)



Complete Turbine Assembly	35-116861	99.95
Complete Turbine Assembly (Phenolic Bearing)	35-116862	104.95
Bearing Kit	35-116863	29.95



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

W&H Adec Synea Small Push Button Turbine (HPP)



Complete Turbine Assembly	35-116931	99.95
Complete Turbine Assembly (Phenolic Bearing)	35-116932	104.95
Bearing Kit	35-116933	29.95

Handpieces Lubricants & Cleaners

EasyOil Lubrication Adapter (Sable)

The best lubrication adapter. Ever.



The universal solution for a trouble-free handpiece care. EasyOil is the universal support that can conveniently gather up to five different adapters to meet all your dental handpiece maintenance needs. An effective handpiece care results in fewer failures and great savings on repair costs. EasyOil allows you to take care of your handpieces in an easy, fast and comfortable way because all you need is already there, readily at hand. With just a few, simple gestures EasyOil will help you make the life of your handpieces longer

Patented-100% **Made in Italy**

Kavo Easy Adapter	35-98771	49.95
W&H Easy Adapter	35-98772	49.95
Bien Air Easy Adapter	35-98773	49.95
Sirona Easy Adapter	35-98774	49.95
NSK Easy Adapter	35-98775	49.95
Midwest Easy Adapter	35-98776	49.95
Handpiece Maintenance Multi-adapter	35-98777	219.95



Like us on Facebook

Call us for Quarterly Promotion

Handpieces

Handpiece Cleaning & Lubrication System (TPC)

Handpiece Cleaning and Lubrication System

Product Features

- Economical, efficient and time-saving
- A fail-safe maintenance solution for your office
- Filtered air used to rotate handpieces for safety
- Care for most handpieces and all handpieces brands
- Disposable mist filtration system
- Solution level indicator
- Detachable door for easy cleaning of unit
- Easy replenishment of solution tank
- Air Pressure gauge
- Oil flow control valve

- **Saves Time**
- **Simple to Operate**
- **Optimizes Performance**



H6000
(3 high speed positions)



H6005
• 2 high speed positions
• 1 low speed E-type motor



H6025
• 1 high speed positions
• 2 low speed E-type motor



Content: 3 high speed handpiece positions, Mist filters (pack of 12), Lubrication fluid, Funnel & Air filter

35-12580	794.95
----------	---------------

Content: 2 high speed, 1 low speed E-type Motor positions, Mist filters (pack of 12), Lubrication fluid, Funnel & Air filter

35-12581	794.95
----------	---------------

Content: 1 high speed, 2 low speed E-type Motor positions, Mist filters (pack of 12), Lubrication fluid, Funnel & Air filter

35-12582	794.95
----------	---------------

Lubrication Fluid

35-12583	39.95
----------	--------------

Mist Filter 12/Pkg

35-12584	18.95
----------	--------------

Oil absorbent door strips 20/Pkg

35-12585	15.95
----------	--------------

Kavo Type HP adapter

35-12586	94.95
----------	--------------

NSK Type HP adapter

35-12587	94.95
----------	--------------

Bien Air Type HP adapter

35-12588	94.95
----------	--------------

W&H Type HP adapter

35-12589	94.95
----------	--------------

Handpiece Lubricant (Kavo)



500ml Bottle	35-16581	39.95
Nozzle, each	35-16582	9.95

Prestigedentalproducts.com

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Handpiece Lubricant Oil (HPP)

General purpose lubricating oil, for use in highspeed and slowspeed applications.



Each 35-97731 14.95

Handpiece Lubrication and Cleaning System (ATW)

Features:

- High security, Easy to operate the dental handpiece lubrication system
- 3 Dental Handpieces can be Oiled and Cleaned One Time
- Remained oil can be cleared by pressing the air key after handpiece lubrication
- The mist filters attached behind the door can keep mist leakage to minimum
- Mode Time of cleaning and oiling Time of cleaning
- Short Mode Almost 35 Seconds Almost 60 Seconds
- Long Mode Almost 45 Seconds Almost 60 Seconds
- Long Mode Almost 50 Seconds Almost 90 Seconds
- The use of dental handpiece can improve the maintenance drive 40 percent of life



Unit for: 2 High Speed Handpieces & 1 Low Speed Handpiece, comes with Mist Filter 12/Box 35-117271 799.00

Lubricant Fluid, 1Q 35-117261 39.99

Mist Filter, 12/Pkg 35-11725 18.99

Handpiece Lubrication Oil Spray (Yoshida)

Made in Japan



Each 35-97551 59.95

Handpiece Purge Station with Hepa Filter (Parts)

- Purges excess lubricant and water from handpiece
- Patented Bio-Cleanse HEPA exhaust filter
- Hepa Filter helps protect operator from air-borne pathogens
- Provides for quick and easy handpiece maintenance
- Exhaust reservoir safely collects excess lubricant



- Purges handpiece coolant air and water lines
- Easy view gauge with adjustable air regulator
- Air supply line with 1/4" male quick-disconnect included
- Compact size of a foot control - requires little counter space
- Handpiece adapter required from below

Each 35-117371 164.95

Midwest Adapter 35-117361 29.95

Star Adapter 35-11735 29.95

Kavo Adapter 35-11734 37.95

NSK, Kinetics, Healthco, Viper, Atech Adapter 35-11733 42.95

Borden Adapter 35-11732 37.95

Midwest XGT Adapter 35-11731 42.95

Lares Adapter 35-11730 37.95

W&H Adapter 35-11729 32.95

Single Hepa Filter 35-11728 7.99

Lubrex Handpieces Cleaning and Lubricating System (Dentamerica)

- Automatic cleaning, purging and lubricating system
- Ergonomic front-loading, easy access design
- Can be used for all high speed, low speed, air scaler, and air motor handpieces
- Up to 4 handpieces can be processed simultaneously
- Quick release air connector
- Automatic fan and filter for vapor removal
- Large reservoir capacity with visible oil level indicator
- Specifications : 100~240V 50/60 Hz AC
- Dimensions : W330 x H430 x D430 mm
- Weight : 13 kg
- Accessories : Main unit, Lubricant, Quick Connector, Hose, Power Cord



Unit 35-118081 989.95

Lubricant & Clean Spray

Effective dental handpieces lubrication 2 to 3 seconds twice a day.

Spray Can 18oz, 550ml 35-118091 19.99



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Lubricant Spray for Kavo (DA)

For Kavo attachments



12oz Spray Bottle 35-99031 **12.95**

Lubrinsa Automatic Handpiece Lubrication System (J-Morita)

Lubrinsa cleans and purges handpieces, and a variety of other equipment, via a rapid air blow system. The speed of air flow has been substantially increased resulting in more power. It delivers superior performance in the same amount of time – just 25 seconds per handpiece. Lubrinsa also offers a long blow out mode, if needed, to help ensure proper maintenance and maximize service life. Rated as highly versatile in 2009 for its setting options, four ports, and compatibility with multiple handpiece and lubricant brands, Morita continues to improve Lubrinsa's versatility. In addition to 4-hole and 6-pin, it is now compatible with a standard 5-hole connection. Lubrinsa is also currently the ONLY system on the market that accepts multiple brands of spray lubricant including the use of two different sprays in the same cycle. Powerful, versatile, and easy to operate, Lubrinsa simplifies the maintenance process and keeps handpieces running smoothly.



Flexible

Only unit that accepts multiple brands of spray
 Allows two different sprays per cycle
 Accommodates wide variety of attachments & handpiece types
 Short and long mode options


Saves time & Money

Uses only 1/3 the lubricant (compared to manual lubrication)
 Services up to 4 handpieces in only 100 seconds
 Protects your equipment, helps keep it running smoothly


Unique

Built-in chuck maintenance system
 Air suction system eliminates mist for clean operation


Lubrinsa is compatible with the following connections:



6-Pin Connection



4-Hole Connection



5-Hole Connection

Unit Includes: AR Spray (420ml), AR
 Spray Can Stand, Oil Absorption Pad (5)
 Oil Absorption Case, Tray, Tray Filter, Air
 Tube, Power Cord

35-118101 **1999.95**

Lubrinsa 4H (Morita Type) Coupling 35-118102 **137.95**

Lubrinsa TR (E-Type Joint) Coupling 35-118103 **137.95**

Handpieces

Lubrinsa Alpha Coupling	35-118104	137.95
Lubrinsa Root ZX II Coupling	35-118105	189.95
Lubrinsa Tri Auto ZX Coupling	35-118106	189.95
Lubrinsa W&H Coupling	35-118107	189.95
Lubrinsa Kavo Coupling	35-118108	189.95
Lubrinsa Sirona Coupling	35-118109	189.95

Once-A-Day Handpiece Spray (Palmero)

TopCat Once-A-Day 1 Second Spray,
 For use in air turbine hand pieces,
 proph angles and contra angles. Can
 also be used as a lubricant/cleaner for
 hand pieces. Spray daily. No corrosive
 (contains hexane). 8.8 oz. aerosol can
 with extension tube.



8.8 oz can 35-99021 **10.95**

Pana Spray Handpiece Lubricants and Conditioners (NSK)

Handpiece Lubricants and Conditioners
 Spray (NSK America)



Spray 16oz 35-94231 **39.95**

Spray Lube Surgical Instrument Lubricant (Miltex)

Oil 8oz Anti-Corrosive



8oz oil/bottle 35-99011 **10.95**

Handpieces



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Waterline Treatment

BluTab Waterline Maintenance Tablets (Confirm Monitoring)

- Specially designed to be continuously present in your waterlines
- Eliminates the need for weekly shock treatments
- Odorless, tasteless, non-toxic, and safe for equipment
- No need to purge lines at the end of the day



50 tabs/Box for 2 liters of water

35-94152

36.95

Citrisil Dental Waterline Cleaner (Sterisil)

- Non-toxic, non-corrosive, no taste, no odor
- Just drop tablets into water reservoir and fill with water
- Long lasting
- Eliminates the need to empty bottles and purge and air-dry



20 Maintenance Tablets & 1 Shock Tablet: White For 0.7-1 liter bottles Mfg. Part #: C20-W

35-94156

18.95

20 Maintenance Tablets & 1 Shock Tablet: Blue For 0.7-1 liter bottles Mfg. Part #: C20-B

35-94157

18.95

Mint-A-Kleen Waterline Disinfectant (Anodia Systems)

- EPA/FDA Registered.
- Fast and ready to use—no mixing or diluting necessary.
- 0.12% Chlorhexidine Gluconate—active ingredient directly contacts waterline walls and bottles for thorough treatment, cleaning, and protection of waterlines and equipment.
- Formulated without harsh chemicals.
- Compatible with all dental units.
- Proven mint-blue, liquid solution—does not corrode, clog, or leave residue in waterlines, bottles, or equipment.

Mint-A-Kleen's weekly treatment prevents bad taste in the water, unpleasant odors, slimy and clogged waterlines, bottles, handpieces and syringes. **Available for sale in ALL 50 States.**



Case of 10 x 16 oz Bottle

39-11261

115.99

1 x 16oz bottle

39-11262

14.99



Order Online at
PrestigeDental.com





Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Handpieces

Handpiece Bulbs

LED Bulbs (Chapman-Huffman)



LED Bulb 16-120-00

Fits following Swivel Couplers Chapman-Huffman 20-079-00 KAVO Multiflex MK-DENT QC5016, QC5016K KAVO



LED Bulb 16-126-00

Fits KAVO Motors



LED Bulb 16-130-00

Fits SIRONA Quick Connectors



LED Bulb 16-134-00

Fits W&H Turbines



LED Bulb 16-132-00

Fits SIRONA Turbines and SIRONA Motors



LED Bulb 16-122-00

Fits BIEN-AIR Motors



LED Bulb 16-124-00

Fits BIEN-AIR Quick Connectors



LED Bulb 16-128-00

Fits NSK Quick Connector



LED Bulb 16-136-00

Fits STAR Quick Connector

Fits following Swivel Couplers Chapman-Huffman 20-079-00 KAVO Multiflex MK-DENT QC5016, QC5016K KAVO cross ref. # 1007.5372

Fits KAVO Motors

35-118061 **70.80**

Fits SIRONA Quick Connectors

35-118051 **70.80**

Fits SIRONA Turbines and SIRONA Motors

35-118041 **70.80**

Fits W & H Turbines

35-118031 **70.80**

Fits BIEN-AIR Motors

35-118021 **70.80**

Fits BIEN-AIR Quick Connectors

35-118011 **70.80**

Fits NSK Quick Connector

35-118001 **70.80**

Fits STAR Quick Connector

35-117991 **70.80**

35-118071 **70.80**

Variety of Handpieces to choose from:





Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

