

Laboratory



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**





Prestige Dental Products, Inc.
Tel: 877-772-3888

Lab Plasters

Magnum Plaster White Stone (MDC)



Creamy White Stone, 22LB 34-11121 **22.95**

Lab Stones

Magnum Orthodontic Stone (MDC)



Die Pink Stone, 22LB 34-88691 **24.95**

Magnum Super Stone Buff (MDC)



Blue, 22lbs 34-11123 **21.95**

Yellow, 22lbs 34-11122 **21.95**

Model Formers

Disposable Model Base Formers (Dentsply-Glenroe)



120/Pkg 34-11112 **54.95**

60/Pkg 34-11111 **29.95**

Green Rubber Base Formers for Plastic Articulators (Meta Dental)



Green Rubber Base Former Q-Left 34-11117 **5.95**

Green Rubber Base Former Q-Right 34-11116 **5.95**

Green Rubber Base Former (Small) 34-11115 **5.95**

Green Rubber Base Former (Medium) 34-11113 **5.95**

Green Rubber Base Former (Large) 34-11114 **5.95**

Green Rubber Base Former (X-Large) 34-11118 **5.95**

Green Silicon Rubber Base Formers (Meta Dental)



Green Silicon Base Formers (Set of 5/Box 34-11119 **19.95**

Laboratory

Model Formers (Plasdent)



Asst Sizes 34-11124 **8.95**

Medium Green 34-11126 **8.95**

3pcs Upper and 3 pcs Lower per Box

Small Yellow 34-11125 **7.95**

3pcs Upper and 3 pcs Lower per Box

Large Blue 34-11127 **9.95**

3pcs Upper and 3 pcs Lower per Box

Disposable Asst Sizes 120 Pcs/Box 34-11128 **29.95**

Disposable Asst Sizes 60 Pcs/Box 34-11129 **15.95**

Disposable Medium 60 Pcs/Box 34-11130 **13.95**

Disposable Universal 60 Pcs/Box 34-11131 **15.95**

Disposable Universal Deep 60 Pcs/Box 34-11132 **18.95**

Orthodontic Plaster

Magnum Die Pink Stone (MDC)



22Lb 34-11120 **22.95**

Plaster Traps

Gleco Plaster Trap (IndiGreen)



Trap System Invented by a plumber, the Gleco Trap brings revolution to the age-old chore of sink trap retrieval and plaster trap maintenance. No more leaky bags or bulky buckets!

64 oz. Gleco 34-87931 **94.95**

19 oz. Gleco 34-87932 **89.95**

It's easy to join our mailing list!

Just send your email address by text message:

Text
PRESTIGEDENTAL
to **22828** to get started.



Message and data rates may apply.

Call us for Quarterly Promotion

Prestigedentalproducts.com

Laboratory

Trap-Eze™ Disposable Plaster Traps (Buffalo)



Requires minimum 9" clearance between bottom of sink tail pipe and top of trap bucket for easy installation without tools
All trap buckets feature high-strength plastic construction and permanently-sealed lids
3 1/2 gallon size (12 1/2" diameter x 11" H) ideal for intermediate applications
5 gallon size (12 1/2" diameter x 14 1/2" H) best for larger offices and labs with higher volume collection needs (space permitting)
Two flexible hoses (1 1/2" ID x 1 3/4" OD x 20" long) conveniently connect trap to sink tail pipe and sewer tail pipe
Flexible adaptor fitting for sewer tail pipe included for improved drain seal
Complete kit includes sealed trap bucket, 2 hoses, and sewer tail pipe adaptor fitting

Trap-Eze BT Bottle

Trap 64 Oz. Refills w/ Lids & Bibs, Pkg. of 2 34-11144 **19.95**

Trap-Eze BT Bottle Trap Complete Kit (1 1/2"

Plumbing + One 64 oz. Bottle 34-11143 **73.95**

Drain Fitting for Trap-Eze Exhaust Hose 34-11142 **4.95**

Hoses Only for Trap-Eze, Set of 2
(18" Overall Length Each) Flexible 34-11141 **35.95**

5 Gal.
Trap-Eze SS Self-Sanitizing Trap Refill - New! 34-11140 **30.95**

5 Gal.
Trap-Eze SS Self-Sanitizing Trap Complete Kit- New! 34-11139 **57.95**

5 Gal. Trap-Eze Refill
(14.5" Height Sealed Bucket) 34-11136 **22.95**

5 Gal. Trap-Eze Complete Kit
(14.5" Height Sealed Bucket, Hoses) 34-11135 **50.95**

3 1/2 Gal.
Trap-Eze SS Self-Sanitizing Trap Refill - New! 34-11138 **28.95**

3 1/2 Gal.
Trap-Eze SS Self-Sanitizing Trap Complete Kit - New! 34-11137 **54.95**

3 1/2 Gal. Trap-Eze Refill
(11" Height Sealed Bucket) 34-11134 **20.95**

3 1/2 Gal. Trap-Eze Complete Kit
(11" Height Sealed Bucket, Hoses) 34-11133 **46.95**



Prestige Dental Products, Inc.
Tel: 877-772-3888

Aluminum Oxide (Deldent)



In 2.2lb. Plastic Bottles 50 Micron for general etching and blasting.

34-11216.....**28.95**

90 Micron for heavier work.

Use with the Miniblaster™, Minietcher™,
and all other clinical sandblasting units.

41-20000**28.95**

Aluminum Oxide

99.6% pure
Non-contaminating
Uniform particles
Desiccant-packaged
for moisture prevention
25lb.

50 Micron White
41-20001**85.95**



Aluminum Oxide (Danville)

1lb.

27 Micron 41-20002 **20.95**

50 Micron White 41-20003 **13.95**

90 Micron Tan 41-20004 **13.95**





Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

