

# Matrix & Wedges



**Prestige Dental Products, Inc**

**PrestigeDental.com**  
**877-772-3888**

# Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

**As an authorized Dealer of many brands in North America**, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

**Prestige Dental Products** offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

**Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and**

## Customer Information

### Order

**PHONE** To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

**FAX** You may FAX your orders to 714-632-8821.

**INTERNET** [www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.

*Note:* New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



### Customize order for special customer

Email us to [mass@prestigedentalproducts.com](mailto:mass@prestigedentalproducts.com) a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

### Payment

**CREDIT CARD** We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

#### **Bank to Bank payment through High Secure Intuit Payment-Network**

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

#### **MONTHLY BILLING "Net 30 days payment"**

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

### Shipping

#### **GROUND DELIVERY**

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

**EXPEDITED SHIPPING** UPS Next Day and Second Day Service for which additional shipping fees will be charged.

**SHIPPING HAZARDOUS ITEMS** All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

### State and federal laws

#### **DEA CERTIFICATE & LICENSE**

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog  
Are Subject To Change Without Prior Notice.**





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Matrix & Wedges

### Automatrix (Dentsply)



AutoMatrix is a Retainerless Metal Matrix which is time-saving, gives you a clear field of vision and makes it easier on your patient. AutoMatrix Retainerless Metal Matrix is indicated where placement of a matrix band around a prepared tooth will facilitate contouring of a direct placement of restorative material (class II preparations).

Med/Regular Refills, 1/4" x .002" 72/Pkg	42-11113	<b>54.99</b>
Narrow/Regular Refill, 3/16" x .002" 72/Pkg	42-11112	<b>54.99</b>
Med/Thin Refills, 1/4" x .0015", 72/Pkg	42-11111	<b>54.99</b>

### Tofflemire Matrix Band (Arel)



Band #13: Pedo - .0015" thick, 144/Pkg	42-11184	<b>6.95</b>
Band #3 - .0015" thick, 144/Pkg	42-11181	<b>6.95</b>
Band #2 - .0015" thick, 144/Pkg	42-11182	<b>6.95</b>
Band #1 - .0015" thick, 144/Pkg	42-11183	<b>6.95</b>

### Plastic Wedges (Plasdent)



Green, 12mm, 100 pcs/bag	42-11178	<b>7.02</b>
Yellow, 14mm, 100 pcs/bag	42-11177	<b>7.02</b>
Blue, 10mm, 100 pcs/bag	42-11179	<b>7.02</b>
Purple, 16mm, 100 pcs/bag	42-11176	<b>7.02</b>
Asst 4 Sizes, 100 pcs/bag	42-11180	<b>7.02</b>

### Sectional Matrix

#### Contact Matrix Rings (Danville)



**MegaRings** make perfect, tight contacts even easier. They apply the greatest force and resist deformation to keep their strength. Application is easy as ever, with the same forceps. Design allows easy 'stacking'. Converging tines prevent 'pop-offs' and the widest tines in the market to span the widest preps. MegaRings are made of a unique, ultra 'springy' alloy and also features a graduated, higher profile towards the back of the ring. This is engineered to spread the force of bending over the entire arch of the ring. **Contact Rings** hold the matrices in place and cause the teeth to spread apart, creating tight, reliable contacts. The wide blades on the contact rings span wide preps where competitive brands fall in. Rings are also easily stacked for multiple proximals. Converging tines prevent pop-offs.

<b>Trial Kit Contact Rings,</b> 1 Inward & 1 OutwardThin Flex Matrices (50 Small & 50 Large)	42-11142	<b>61.95</b>
<b>Mega Ring Trial Kit</b>	42-11148	<b>79.95</b>
<b>Mega Ring Clinical Kit</b> PlusRing Pliers Plus & Mega Grip Forceps	42-11149	<b>259.95</b>
<b>Clinical Kit Contact Rings,</b> 3 Inward & 3 OutwardThin Flex Matrices (100 Small & 100 Large)	42-11141	<b>124.95</b>
<b>Ring Refill 1ea</b>	42-11137	<b>21.95</b>
<b>Thin Flex Matrices/Large 100/pk</b>	42-11136	<b>37.95</b>
<b>Ultra Thin Flex Matrices/Small 100/pk</b>	42-11135	<b>47.69</b>
<b>Ultra Thin Flex Matrices/Large 100/pk</b>	42-11134	<b>47.69</b>
<b>Stiff Flex Matrices/Large 100/pk</b>	42-11133	<b>37.95</b>
<b>Thin Flex Matrices/Small 100/pk</b>	42-11138	<b>37.95</b>
<b>Ring Refill 2 Outward (gold)</b>	42-11139	<b>21.95</b>
<b>Ring Refill 2 Inward</b>	42-11140	<b>21.95</b>
<b>Mega Ring Refill 2 Inward</b>	42-11132	<b>37.95</b>
<b>Subgingival Matrices 50/pk</b>	42-11147	<b>34.95</b>
<b>Mega Ring Clinical Kit3 Each Mega Rings, 100 Each S/L Matrices</b>	42-11146	<b>149.95</b>
<b>Mega Ring Refill 1 Inward &amp; 1 Outward Ring</b>	42-11145	<b>37.95</b>
<b>Mega Ring Refill 2 Outward</b>	42-11144	<b>37.95</b>
<b>Stiff Flex Matrices/Small 100/pk</b>	42-11143	<b>37.95</b>



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Matrix & Wedges


### FenderMate (Directa)



Right Regular	Dark Blue	18pcs/Box	42-11164	<b>55.27</b>
Right Narrow	Light Blue	18pcs/Box	42-11165	<b>55.27</b>
Left Narrow	Light Green	18pcs/Box	42-11168	<b>55.27</b>
Assortment Kit	72pcs/Box (4 x 18Pcs each size)		42-11167	<b>164.95</b>
Left Regular	Dark Green	18pcs/Box	42-11166	<b>55.27</b>


### CoForm Directa (Directa)

Transparent Corner Matrices for class IV fillings, restoration of incisal edges, restoration of corners, support and form for core build-ups, support and form when making pontics in bonded composite bridges.



64pc Kit (11,12,21,22 M & D)		42-11114	<b>65.25</b>
Size: 41-42, Mesial: , 8/Pkg		42-11115	<b>9.53</b>
Size: 31-32, Mesial: , 8/Pkg		42-11116	<b>9.53</b>
Size: 31-32, Distal: , 8/Pkg		42-11117	<b>9.53</b>
Size: 23, Distal: , 8/Pkg		42-11118	<b>9.53</b>
Size: 23, Mesial: , 8/Pkg		42-11119	<b>9.53</b>
Size: 22, Distal: , 8/Pkg		42-11120	<b>9.53</b>
Size: 22, Mesial: , 8/Pkg		42-11121	<b>9.53</b>
Size: 21, Distal: , 8/Pkg		42-11122	<b>9.53</b>
Size: 21, Mesial: , 8/Pkg		42-11123	<b>9.53</b>
Size: 13, Distal: , 8/Pkg		42-11124	<b>9.53</b>
Size: 13, Mesial: , 8/Pkg		42-11125	<b>9.53</b>
Size: 12, Distal: , 8/Pkg		42-11126	<b>9.53</b>
Size: 12, Mesial: , 8/Pkg		42-11131	<b>9.53</b>
Size: 11, Distal: , 8/Pkg		42-11130	<b>9.53</b>
Size: 11, Mesial: , 8/Pkg		42-11129	<b>9.53</b>
128pc Kit (1pkg of each size)		42-11128	<b>123.30</b>
Size: 41-42, Distal: , 8/Pkg		42-11127	<b>9.53</b>

### FaciForm Directa (Directa)



Matrices for forming composite facings, pontics, cervical fillings, incisal contouring.

Size: 11, 8/Pkg	42-11156	<b>9.53</b>
56pc Kit (1pkg of each size)	42-11157	<b>65.25</b>
Size: 12, 8/Pkg	42-11158	<b>9.53</b>
Size: 13, 8/Pkg	42-11159	<b>9.53</b>
Size: 21, 8/Pkg	42-11160	<b>9.53</b>
Size: 22, 8/Pkg	42-11161	<b>9.53</b>
Size: 23, 8/Pkg	42-11162	<b>9.53</b>
Size: 32-42, 8/Pkg	42-11163	<b>9.53</b>

### Wedges

#### Contact Wedge (Danville)



ContactWEDGE Refill 85 Small	42-11150	<b>39.95</b>
ContactWEDGE Refill 85 Large	42-11152	<b>39.95</b>
ContactWEDGE Refill 85 Medium	42-11151	<b>39.95</b>
Mega Grip Forceps	42-11153	<b>53.95</b>
Elastic Wedge System Contact WEDGE Kit Small, Medium, Large 85ea.	42-11154	<b>96.95</b>

#### Wykle Wedges (Wykle)

The patented anatomical design is contoured for optimum performance. Crafted from Lindenwood, the Wykle Wedge prevents injury to gingival tissue and remains where placed.



Wykle Wedges Refill Blue 15mm, 500/Pkg	42-11190	<b>23.49</b>
Wykle Wedges Refill Yellow 13mm, 500/Pkg	42-11191	<b>23.49</b>
Wykle Wedges Refill Natural 11mm, 500/Pkg	42-11192	<b>23.49</b>
Wykle Wedges Assorted, 500/Pkg	42-11193	<b>23.49</b>
Wykle Wedges Refill Green 17mm, 500/Pkg	42-11194	<b>23.49</b>

# Matrix & Wedges



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## FenderWedge (Directa)



### Protects adjacent teeth during preparation!

In almost 2/3 of cases adjacent teeth are damaged during class II preparation. Avoid this with Fender Wedge™ as you establish and maintain correct positioning for your contact point.

As easy to insert as a wedge.

Jams in and Stays put! Avoids interference.

Placing a Fender Wedge™ prior to preparation facilitates easy placement of matrix and establishment of contact points.

Large Yellow	Size 1.0mm wide	100 pcs/Box	42-11174	<b>75.75</b>
X-Small Purple	Size 1.0mm wide	100 pcs/Box	42-11173	<b>75.75</b>
Small Orange	Size 1.4mm wide	100 pcs/Box	42-11172	<b>75.75</b>
Medium Green	Size 1.9mm wide	100 pcs/Box	42-11171	<b>75.75</b>
Assortment Kit	144 pcs/box (36pcs ea size)		42-11170	<b>127.95</b>

## Wizard Wedges® Matrix Wedges (Waterpik)



These proven wedges ensure better access and stable holding of matrix systems during your critical restorative procedures. Specially crafted from birch wood for greater tooth separation and faster, better tooth position recovery upon removal.

Wizard Wedges, Wide Base, 500/Pkg	42-11185	<b>37.95</b>
Wizard Wedges, Large, 500/Pkg	42-11187	<b>37.95</b>
Wizard Wedges, Slim Jim, 500/Pkg	42-11186	<b>37.95</b>
Wizard Wedges, Assorted, 500/Pkg	42-11189	<b>37.95</b>
Wizard Wedges, Small, 500/Pkg	42-11188	<b>37.95</b>

## Matrix Retainer No. 8 (Miltex)



42-20014 **54.95**

## Tofflemire™ Matrix (DSC)



Universal Retainer	42-20000	<b>10.99</b>
Contra Angle Retainer	42-20001	<b>10.99</b>
Junior Contra Angle Retainer	42-20002	<b>10.99</b>
#1 .0015 Bands 144/pkg	42-20003	<b>13.99</b>
#2 .0015 Bands 144/pkg	42-20004	<b>13.99</b>
#3 .0015 Bands 144/pkg	42-20005	<b>13.99</b>
#13 .0015 Bands 144/pkg	42-20006	<b>13.99</b>
#1 .002 Bands 144/pkg	42-20007	<b>13.99</b>
#2 .002 Bands 144/pkg	42-20008	<b>13.99</b>
#3 .002 Bands 144/pkg	42-20009	<b>13.99</b>
#13 .002 Bands 144/pkg	42-20010	<b>13.99</b>
#1 .001 Bands 144/pkg	42-20011	<b>13.99</b>
#2 .001 Bands 144/pkg	42-20012	<b>13.99</b>
#13 .001 Bands 144/pkg	42-20013	<b>13.99</b>

## Matrix Retainer No. 9 (Miltex)



42-20015 **40.95**





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Matrix & Wedges

### ElastoWedge (Danville)



Elastique Radiopaque Interdental Wedge. Provides complete gingival margin adaptation and seal. Simultaneously wedges from both buccal and lingual.

ElastoWedge Assortment Kit (255)	42-20020	<b>97.95</b>
Contains: 85 each of small, medium, and large.		
ElastoWedge Refill Small (85)	42-20120	<b>37.95</b>
ElastoWedge Refill Medium (85)	42-20121	<b>37.95</b>
ElastoWedge Refill Large (85)	42-20122	<b>37.95</b>

### Clear Matrix Bands

Standard, pk100

42-20021 **8.99**



### Contour Matrix Bands (Waterpik)

Pkg. of 50

.0015 Ultra Thin	42-20022	<b>36.95</b>
.002 Regular	42-20023	<b>36.95</b>



### Matrix Dispenser (Zirc)

- Contains one roll of 12yd. x 3/8" matrix tape, 0.002ga.
- Used with Matrix mylar rolls
- Allows matrix to be contoured (U-shaped as it dispenses)
- Stainless steel cutter (allows for desired length)
- Plastic dispenser: 2" x 2" x 1/8"
- Color: White

42-20024 **21.95**



### Strip Holder (Miltex)

<b>Regular</b>	42-20025	<b>21.95</b>
<b>Right Angle</b>	42-20026	<b>21.95</b>



### Wizard Wedges, Anatomical (Water Pik)

Box 400

Small, Natural	42-20027	<b>35.95</b>
Medium, Yellow	42-20028	<b>35.95</b>
Large, Purple	42-20029	<b>35.95</b>
X-Large, Red	42-20030	<b>35.95</b>
Assorted	42-20031	<b>35.95</b>



### Contact Matrix System (Danville)



- Makes tight contacts easy
- Contoured to restore the natural anatomy
- Proximal finishing is greatly reduced or eliminated
- Unique in resilience

Trial Kit

Contains: Mega Grip forceps (with serrated jaws to make the removal of matrices easy), 2 Contact rings (one inward and one outward), 50 small and 50 large Thin Flex matrices.

42-20032 **59.95**



Prestige Dental Products, Inc.

www.prestigedentalproducts.com

Order Form

Date:

[Date input box]

Ordered By

Deliver To

Same as Above

Business Name:

[Business Name input box]

Address:

[Address input box]

State/Province:

[State/Province input box]

Zip/Postal Code:

[Zip/Postal Code input box]

Phone:

[Phone input box]

Fax:

[Fax input box]

Contact Name:

[Contact Name input box]

Business Name:

[Business Name input box]

Address:

[Address input box]

State/Province:

[State/Province input box]

Zip/Postal Code:

[Zip/Postal Code input box]

Phone:

[Phone input box]

Fax:

[Fax input box]

Contact Name:

[Contact Name input box]

Item Number	Description	Quantity	Unit Price	Amount

Payment

Credit Card

Online Payment (Bank to Bank)

E-mail

[E-mail input box]

American Express

MasterCard

Visa

Check Payable To

[Check Payable To input box]

Card Number:

[Card Number input box]

Expiration Date:

[Expiration Date input box]

Cardholder Name:

[Cardholder Name input box]

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

Dental License #

[Dental License # input box]

Year Licensed

[Year Licensed input box]

Expiration Date

[Expiration Date input box]

DEA #

[DEA # input box]

Year Licensed

[Year Licensed input box]

Expiration Date

[Expiration Date input box]



**Prestige Dental Products, Inc**

3031 E. Coronado St., #D  
Anaheim, Ca 92806  
Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

