

# Surgical Products



**Prestige Dental Products, Inc**

**PrestigeDental.com**  
**877-772-3888**

# Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

**As an authorized Dealer of many brands in North America**, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

**Prestige Dental Products** offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

**Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and**

## Customer Information

### Order

**PHONE** To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

**FAX** You may FAX your orders to 714-632-8821.

**INTERNET** [www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



*Note:* New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



### Customize order for special customer

Email us to [mass@prestigedentalproducts.com](mailto:mass@prestigedentalproducts.com) a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

## Payment

**CREDIT CARD** We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

### Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

### MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

## Shipping

### GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

**EXPEDITED SHIPPING** UPS Next Day and Second Day Service for which additional shipping fees will be charged.

**SHIPPING HAZARDOUS ITEMS** All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

## State and federal laws

### DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog  
Are Subject To Change Without Prior Notice.**



# Surgical Products



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## COE-Pak (GC America)

Periodontal Dressing, Regular set is ready in 3 minutes, Hard & Fast set is ready in 1 minute.  
Regular, Standard Pkg.  
Each package contains: 90 gm each base and catalyst



<b>Regular Set:</b> 90 gm each base and catalyst	21-11127	<b>74.95</b>
<b>Hard/Fast:</b> 90 gm each base and catalyst	21-11128	<b>74.95</b>

## Perio-Care Periodontal Dressing (PulpDent)



For Post Surgical Tissue Management  
Does not support the growth of bacteria.  
Tissues are clean upon removal or the dressing.  
A highly elastic periodontal dressing that sets resiliently hard and does not chip or fall apart in the mouth. Provides durable post surgical protection to tissues.  
Two-paste, highly elastic periodontal dressing. Assists in tissue placement after surgery.  
7 Minute working time; sets in 15 minutes.  
Neutral odor and taste.  
Contains no eugenol.

<b>Standard Pack Contains:</b> 3 oz. tube paste, 3 oz. tube gel & mixing pad	21-11223	<b>38.84</b>
---	----------	--------------

## Dry Socket Paste (Sultan)




Used to treat alvolitis. Gives instant pain relief. Paste is introduced into socket with flat-bladed instrument, tamped down to cover exposed bone & allowed to remain in socket heal. A one-visit treatment with no packing to remove.

**Jar, 1 oz. 52-11114 79.95**

## FG Surgical Carbide Burs

FGSL 6, 5/Pkg W/O Label	19-11380	<b>9.95</b>
FGSL 4, 5/Pkg W/O Label	19-11381	<b>9.95</b>
FGSL 2, 5/Pkg W/O Label	19-11382	<b>9.95</b>
FGSL 8, 100/Pkg W/O Label	19-11384	<b>159.95</b>
FGSL 6, 100/Pkg W/O Label	19-11385	<b>159.95</b>
FGSL 2, 100/Pkg W/O Label	19-11386	<b>159.95</b>
FGSL 4, 100/Pkg W/O Label	19-11387	<b>159.95</b>
FGSL 8, 5/Pkg W/O Label	19-11388	<b>9.95</b>

## FG Surgical Straight Fissure Crosscut Carbide Burs



FGSL 556, 100/Pkg W/O Label	19-11389	<b>139.95</b>
FGSL 557, 100/Pkg W/O Label	19-11390	<b>139.95</b>
FGSL 558, 100/Pkg W/O Label	19-11391	<b>139.95</b>
FGSL 559, 100/Pkg W/O Label	19-11392	<b>139.95</b>
FGSL 556, 5/Pkg W/O Label	19-11393	<b>8.95</b>
FGSL 557, 5/Pkg W/O Label	19-11394	<b>8.95</b>
FGSL 558, 5/Pkg W/O Label	19-11396	<b>8.95</b>
FGSL 559, 5/Pkg W/O Label	19-11395	<b>8.95</b>
FGSL 702, 100/Pkg W/O Label	19-11397	<b>139.95</b>
FGSL 702, 5/Pkg W/O Label	19-11398	<b>8.95</b>

## Impact Air 45 Surgical Handpiece (Palisades)



IMPACT AIR 45 Autoclavable Lubricant	35-11770	<b>78.99</b>
Model 9 Push-Button Head Cap with Cap Wrench	35-11771	<b>48.99</b>
Model 9 Push-Button Chuck	35-11772	<b>99.95</b>
Model 9 Push-Button Turbine (replacement)	35-11773	<b>266.99</b>
Model 9 Push-Button Conversion Kit (includes turbine, head cap and cap wrench)	35-11774	<b>299.95</b>
2-Hole Surgical Handpiece	35-11775	<b>789.95</b>
Fiber Optic 5-Hole Surgical Handpiece	35-117761	<b>869.99</b>
4-Hole Surgical Handpiece	35-117771	<b>719.95</b>

## Nano Loupe Light (Dentlight)

The most power per pound, and the smallest overall light on the market with the most power at its size of only 6 grams on headlamp and 87 grams on battery pack.



Kit Includes slide-on composite filter, two battery packs, easy-to-clean battery pouch, universal charger, and loupe adapter with articulating arm that fits on all major loupes.	43-11157	<b>534.95</b>
Headlamp Only	43-11156	<b>168.95</b>
Battery Pack	43-11155	<b>119.95</b>
Power Adapter	43-11154	<b>36.95</b>



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Surgical Products

### LED Headlight System (LumaDent)



**Increase efficiency, lighted field wherever you look. without shadow  
Save your back, The better you see, the faster you can perform  
Complete LED Headlight**

**Weight: 0.19 oz (5.5 gr).** One of the lightest headlamp on the market, it is lighter than a quarter

**Diameter: 0.64" (16.3 mm), Length: 0.79" (20 mm).** Smallest in the market.

**Light Output: 40,000 Lux at 14" working distance.**

**Spot Size: 2.5" at 14", 3.5" at 18"**

**Battery Life: 8 hours continuously at maximum intensity, with Low-Battery indicator light.**

**Battery Charging Time: 3 hours**

**Battery size: 3" x 2" x 0.9" (8.8 cm x 5.5 cm x 2.3 cm)**

**Battery Weight: 5.8 oz (155 gr)**

***Package includes: LED headlamp, Battery Pack, Charger, Mounts for your loupes/frame, Wire accessory.***

If you are interested in Student Group Buy, please, email [sam@prestigedentalproducts.com](mailto:sam@prestigedentalproducts.com) for info.

Please download the Order Form from our website and select the maker name & the frame model, so we can send you the right one.

Complete System, Headlamp, Composite Filter, Mount, Battery Pack, Charger & carrying Case	43-11128	<b>489.95</b>
Battery Pack - Lenser H7, Includes: Battery Pack 500mA AC/DC Charger Carrying Case with Swivel Belt Clip	43-11131	<b>94.95</b>
Composite Resin Orange Filter.	43-11130	<b>34.95</b>
Battery Pack, 1x Lithium Ion Battery Pack 1x AC/DC Charger	43-11129	<b>94.95</b>



# Surgical Products



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## BenaCel (Unicare Biomedical)

BenaCel Resorbable, Dental Dressing, Plug for deep wounds biopsy sites and post extraction sockets  
Benacel absorbable hemostatic gauze acts instantly to control bleeding. Upon contact with blood, BenaCel hemostat will absorb the blood and transform into a gelatinous material.



50mm x 50mm Plug, 10/Pkg	52-11111	57.95
15mm x 15mm Plug, 12/Pkg	52-11113	39.95
5mm x 7mm, 10/Pkg	52-11112	29.95

## CytoFlex Membrane NonResorbable (Unicare Biomedical)



CytoFlex® is a non-resorbable membrane for use as a space-making barrier in guided tissue regeneration procedures.  
Composed of a proprietary micro porous polytetrafluoroethylene (PTFE) thin film, CytoFlex® has a surface texture and porosity suitable for adhesion of host

cells to the material while preventing passage and integration of bacteria within the interstices of the material.  
The supple membrane conforms easily to tissue contours, and yet offers sufficient stiffness to maintain a space over the bony defect.  
The smooth, micro-texture surface is designed to improve flap attachment but prevent tissue and bacteria ingrowth.  
[FDA has approved Cytoflex® for use as a barrier membrane in the treatment of periodontal defect using guided tissue regeneration procedures.](#)

25mm x 30mm, 5/Pkg	18-11112	194.95
--------------------	----------	--------

## Cytoflex Resorb (Unicare Biomedical)



Cytoflex® Resorb is an absorbable membrane for use as a space-making barrier in guided tissue regeneration procedure. Cytoflex® Resorb membranes are designed to enhance the adhesion of host cells, and at the same time provide a favorable environment for neo-vascularization and repopulation of bone cells to regenerate osseous

tissue. The flexible barrier membranes are easily adaptable to tissue contours, yet offer sufficient strength to protect the bony defect from the intrusion of bacteria and soft tissue. Composed of synthetic polylactide and polyglycolide copolymers, Cytoflex® Resorb barriers are resorbable and do not require a second retrieval procedure.

30x40mm Pkg of 2	18-11114	184.95
30x40mm, Each	18-11115	109.95

## UniGraft Bioactive Glass (Unicare Biomedical)



Unigraft is made of the highest quality bioactive glass by a proprietary process. The synthetic bioactive granules are optimally sized for ease of administration and enhancement of tissue growth. The product is conveniently packaged in a sterile vial, which is stored within a hermetically sealed foil pouch to ensure product sterility prior to implantation. Unigraft not only offers an osteo-genic framework for new bone growth, but also bolsters a higher level of osteoblastic activity to enhance bone regeneration. FDA has cleared this product for use in treating periodontal defects, extractions sites, cranial and facial augmentation and augmentation of the alveolar ridges.

### Easy to use:

Hemostatic -- quickly stabilizes the wound site  
Compatible with autologous and freeze-dried bone

### Bioactive :

Accelerates bone regeneration  
Transforms into natural bone with time

### Safe Alternative to Freeze-Dried Bone:

Consistent high quality  
Free of disease transmission  
No complicated release forms or restrictions for blood-donors

200-420 Micron, 1.0gm, 5/Pkg	18-11133	259.95
200-420 Micron, 0.37gm, 5/Pkg	18-11134	159.95

## Cytoflex Mesh (Unicare Biomedical)



Cytoflex® Mesh is made of a pure, implantable Titanium mesh with precision mesh profile. This product is designed for use to ensure three-dimensional reconstruction of alveolar bone defects and to facilitate bone augmentation through adequate fixation of the augmentation material.

The thin, supple mesh conforms easily to tissue contours, and yet presents sufficient stiffness to maintain a space over the bony defect.  
Each CytoFlex® mesh is double pouched and provided sterile.

25mm x 30mm, 5/Pkg	18-11113	174.95
--------------------	----------	--------




Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### NeoBur FG Surgical Carbide Burs (Microcopy)




MAX 1, 25/pkg	19-11517	54.95
MAX 2, 25/pkg	19-11518	54.95

### NeoBur FGSL Surgical Carbide Burs (Microcopy)



FGSL 1/2, 25/pkg	19-11549	54.95
FGSL 701, 25/pkg	19-11550	54.95
FGSL 558, 25/pkg	19-11551	54.95
FGSL 557, 25/pkg	19-11552	54.95
FGSL 8, 25/pkg	19-11553	54.95
FGSL 6, 25/pkg	19-11554	54.95
FGSL 4, 25/pkg	19-11555	54.95
FGSL 2, 25/pkg	19-11556	54.95
FGSL 1, 25/pkg	19-11557	54.95
FGSL 702, 25/pkg	19-11558	54.95

### NeoBur RASL Surgical Carbide Burs (Microcopy)



RASL 2, 25/pkg	19-11579	69.95
RASL 4, 25/pkg	19-11580	69.95
RASL 6, 25/pkg	19-11581	69.95
RASL 8, 25/pkg	19-11582	69.95

## Surgical Products

### Surgical Mallets (ASA Italy)



Surgical Mallets Partsch Cm.19	40-13897	39.95
Surgical Mallet Mead Cm.19,5	40-13898	39.95
Spare Parts For 1881	40-13899	5.95

### Surgical Scissors (ASA Italy)



Gum Scissors		
Metzenbaum-Delikat Straight Cm.14,5	40-13900	32.95
Ligature Gum Scissors Spencer Cm.13	40-13901	20.95
Ligature Gum Scissors Spencer Cm.9	40-13902	17.95
Ligature Gum		
Scissors Metzenbaum Curved Cm.14,5	40-13903	29.95
Gum Scissors		
Metzebaum-Delikat Curved Cm.14,5	40-13904	34.95

### Surgical Vacuum (Tech West)



Tech West's innovative Surgical Vacuum leads the industry in performance and value. Quiet, durable and reliable by design, the Surgical Vacuum features a thermally-protected, fan-cooled motor with rubber-mounted heads to minimize noisy vibrations and maximize energy efficiency. The fast recovery time

of the Tech West Surgical Vacuum enables smooth and continuous operation. Requiring virtually no maintenance, the unit is factory set to run at 17 to 23 inches of mercury for optimum performance. To enable considerable savings by the dentist who does not require NFPA certification, the pump is uncertified. However, Tech West can also provide NFPA-certified pumps for those offices requiring certification. Complete with a two year warranty, the Surgical Vacuum Pump is offered in both dual and triple configurations with plumbing, electrical and installation assistance available upon request.

3 Users, 1.5 HP, 2-Head, 208-230V	30-113961	Call for Price
4 Users, 2.25 HP, 3-Head, 208-230V	30-113971	Call for Price

# Surgical Products



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Suture Removal Kit

Kit Contains: Sharp scissor, forceps and sponge.



52-88961

6.95

## Gelfoam (Upjohn)



Sterile absorbable gelatin sponge, for use in oral surgery providing hemostasis. May be used either dry or moistened. Should not be used in closure of incisions because of interference with healing.

12 Single-use  
Sterile Dental Sponges  
Absorbable gelatin sponge, USP (Size 4)  
4sq cm (2 X 2 cm)  
1/2 sq in (3/4 X 3/4 in)

### Size 4 (12)

Contains Carton of 6 envelopes x 2 sponges each

52-11115 104.95

12-7mm Sterile, 12/Pkg

52-81461 149.95

## 0.9% Sodium Chloride Irrigation U.S.P (B Braun)



For Irrigation Only.  
Not for Injection.  
Each 100 mL contains: Sodium Chloride USP 0.9 g;  
Water for Injection USP qs  
pH adjusted with Hydrochloric Acid NF  
pH: 5.0 (4.5-7.0) Calculated Osmolarity: 310  
mOsmol/liter Concentration of Electrolytes  
(mEq/liter): Sodium 154; Chloride 154

1000 ml Each

52-12000

5.50

1000 ml 12/Case

52-12001

54.95

## Implant-Vance (ATW)

2 x 50 ml cartridge /  
3 Mixing tip / 3 Intra-oral tip / 1x bag



Light Body

36-11156

19.95

Heavy Body

36-11155

19.95

## Helimend Absorbable Collagen Membranes (Miltex)



Made from purified bovine tendon. Provides wound stabilization and creates space during guided tissue regeneration procedures where a longer absorption rate is desired.

30mm x 40mm, bx,

52-88946 140.95

20mm x 30mm,

52-88945 96.95

15mm x 20mm,

52-88944 79.95

30mm x 40mm,

52-88943 136.95

20mm x 30mm,

52-88942 84.95

15mm x20m,

52-88941 67.95

## Heliplug Control Bleeding (Miltex)

This cylindrical, soft, white, pliable collagen wound dressing is used to control bleeding and develop blood clots in extraction sites and biopsy sites.



Sterile Wound Dressing, 3/8"x3/4", 10/Pkg 52-88921 101.95

## Helitape Repair the Schniderian Membrane (Miltex)



Sterile Wound Dressing, 1" x 3", 10/Pkg 52-88911 157.95

## Surgifoam Gelatin Sponge (Ethicon)



Surgifoam Gelatin Sponge, USP, Oral Sponge is a sterile, water-soluble, malleable, porcine gelatin absorbable sponge intended for hemostatic use by applying to a bleeding surface.

They are individually packaged to prevent cross-contamination and avoid product waste.

Size 12-7, 12/Box

52-88971

104.95

1cm x 1cm x 1cm, 24/Box

52-88972

109.95





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Dynomite Scented Single Use Nasal Hoods (Midmark)

	Small	Medium	Large
Bubblegum	30-9001	30-9015	30-9016
Strawberry	30-9012	30-9002	30-9018
Vanilla	30-9011	30-9003	30-9017
Plain	30-9004	30-9005	30-9006
Orange	30-9013	30-9014	30-9019
Assorted	30-9007	30-9008	30-9009

**Assorted, Assorted, bag 24** 30-9010 **76.95**  
**Bag of 24** **76.95**

### Personal Inhaler Plus Nasal Hoods (Accutron)



Single patient use nasal hood gives your patients reliable, comfortable, high-quality N2O delivery. Reduces risk of cross-contamination. Thermoplastic elastomer

reduces concern for nasal hood collapse

Non-latex material

8 soothing scents/3 comfortable sizes

Most sizes available in the following scents: Groovy Grape, Birthday Bubble Gum, Sassy Strawberry, Outlaw Orange, Unscented Whisper Gray, Fresh Mint, Perfect Peach, French Vanilla, Piña Colada, and Variety Packs

	Small <b>78.95</b>	Medium <b>87.95</b>	Large <b>90.95</b>
	Small	Medium	Large
Bubblegum	30-88213	30-88222	
Strawberry	30-88212	30-88221	
Vanilla		30-88217	30-88231
Orange	30-88211	30-88220	
Peach		30-88218	30-88230
Mint		30-88219	30-88229
Pina Colada		30-88223	
Grape	30-88216		
Unscented Gray	30-88215	30-88224	30-88228
Variety pk 1 (Orange, Strawberry, and Bubble Gum)	30-88214	30-88225	
Variety pk 2 (Pina Colada, Vanilla, and Orange)		30-88226	
Variety pk 3 (Mint, Peach, and Vanilla)		30-88227	30-88233

### Reservoir Bag (Accutron)



**3 liter, Each** 30-89081 **39.95**

## Surgical Products

### Scavenging Circuit (Accutron)



Low profile scavenging assembly. Clear exhalation valve bulb. Easy-to-read vacuum control gauge. Autoclavable, X-ray compatible circuitry. Scavenging Circuit II Includes: Circuit, reservoir bag, corrugated tube and complimentary sampling of PIP+Nasal Hoods. Scavenging Circuit I Includes: Circuit and complimentary sampling of PIP+Nasal Hoods.

**Scavenging Circuit II, Each** 30-88240 **524.95**  
**Scavenging Hub, Each** 30-88241 **41.95**

### ART-E1 Electrosurgery / Cutting Unit (Bonart)



One of the best selling and most popular products. State of the art radio frequency electrosurgical unit proudly manufactured and offered by Bonart. Unit is comparable to major brands. Great quality & economical. Comes with: 7 piece variety electrodes set, 1 hand piece, 1 indifferent or grounding plate, 1 power cord, 1 foot switch. Operation mode of cut, coagulation, cut-coagulation. Operation mode is monopolar. Output power is up to 70 watts. Operation frequency is at 1.5 to 1.7 MHz with 10 different power settings. Unit is a monopolar mode but bipolar converter is optional and available upon request. Electrodes are autoclavable, 1/16" diameter shaft. Made universal to other major electrosurgical brands. Ideal for human (dental and dermatology use) and veterinary use.

**Complete Unit with 7 Pc. Electrodes** 27-11111 **429.99**  
**Electrodes Kit 7pc/Box** 27-11112 **69.95**



# Surgical Products



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## X-CUBE Brushless Surgical Implant System (Saeshin)



powerful, user-friendly surgical motor system used for general oral surgery, implantology as well as for all kinds of micro surgery procedures.

This **implant kit** is the best value in the industry and comes with everything needed to get started tomorrow.



- Push Button & Lever Type
- Reduction Gear Ratio 20:1 (16:1, 32:1, 64:1 also available)
- Torque: 20:1 (5.0 ~ 55 Ncm)
- RPM range with 20:1 Angle: 30 - 2500
- Autoclavable
- Motor speed control from 600 to 50,000 rpm with high power and high torque
- No noise and no vibration for precise surgical control
- With 32:1 geared angle for ideal speeds of 18-1562 rpm and up to 55 Ncm torque
- **Eye-appealing LCD display** with simple user interface
- **Powerful 230W output**
- **10 types of gear ratio setting** 1:5 / 1:4 / 1:2 / 1:1 / 16:1 / 20:1 / 27:1 / 32:1 / 64:1 / 256:1
- **10 program memory for setting** the speed, torque, operation direction, water pump
- **Quiet and strong water pump** ( 5 stage - control )
- **Automatic overload protection:** If the load on the bur is higher than torque. the motor is stopped automatically after 2 seconds. and to relas the overload function, push foot control pedal

**Ergonomic foot control pedal:** Control all the functions.

<b>X-Cube Implant Set V2.0 with 32:1 Handpiece Lever Type, Max torque: 65Ncm</b>	36-11301	<b>1954.00</b>
<b>X-Cube Implant Set V2.0 with 20:1 Handpiece Push Button Type , Max torque: 55Ncm</b>	36-11300	<b>1954.00</b>
<b>X-Cube Implant Set V2.0 with 20:1 Handpiece Lever Type, Max torque: 55Ncm</b>	36-11299	<b>1954.00</b>
<b>E-Type Handpiece</b>	36-11297	<b>329.95</b>
<b>Irrigation Tube, 10/Pkg</b>	36-11296	<b>74.95</b>
<b>X-Cube Implant Set V2.0 with 32:1 Handpiece Push Button Type, Max torque: 65Ncm</b>	36-11298	<b>1954.00</b>



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Surgical Cement (Prime Dental)

Consists of a powder and liquid that is a multi-purpose periodontal dressing and cement.

-Sets slowly out of the mouth but more quickly with contact of water or saliva.

-Promotes healing and has a soothing



**Kit** 4 oz. Powder & 2oz. Liquid      21-11236      **Call for price**

## Surgical Aspirator Tips (TC Dental)



Disposable. Made of sturdy poly material. Color coded for different tip opening sizes

**Surgical Tips White.** 6 1/4" x 1/8" tip 25pcs/Box      31-11172      **4.98**

**Surgical Tips Green.** 6 1/4" x 1/4" tips 25pcs/Box      31-11171      **4.98**

**Surgical Tips,** 6 1/4" x 1/16" Tip 25pcs/Box      31-11173      **4.98**

**Surgical Tips,** 7 1/2" x 5/16" Tip 25pcs/Box      31-11174      **4.98**

## Telescopic Loupes (ATW)



Extra wide field of view, large depth of focus. High light efficiency. Advanced mounting system with 3 adjustments for individualizing posture, working angles and achieve perfect eye alignment. Lightweight frames with soft silicone nose pad. Adjustable head strap. Removable side shields. you can

put the glass to fit your eye need.

**Magnification +250,** 17" Working Distance & 4.5" Field of View      43-11167      **249.95**

**Magnification +300,** 17" Working Distance & 3.5" Field of View      43-11168      **259.95**

**Magnification +350,** 17" Working Distance & 2.5" Field of View      43-11169      **269.95**

## Surgical Products

### Traus Implant Brushless Surgical Motor (Saeshin)



High performance brushless surgical motor 110V **with Fiper Optic hand-piece**

**Full Set** with Non-Optics - 20:1 PB Angle handpiece - 110V      36-17121      **3495.95**

**Full Set** with Optics 20:1 PB Angle handpiece - 110V      36-17122      **3795.95**

### Twin Power Surgical Handpiece (J-Morita)



- \*Safe, rear-facing exhaust vents
- \*High torque; up to 20 watts of power
- \*Quiet & long-lasting ceramic ball bearing

**Handpiece** 4H PAR-4HEX-O-45 (Morita Couplings)      35-98821      **979.95**

**Handpiece** 4H PAR-4HEX-O-KV-45 (Kavo MULTiflex LUX Coupling)      35-98822      **964.95**

**Turbine** 4HEX-45 Capsule Cartridge (Ceramic Bearing)      35-98823      **379.95**

### Surgical Push Button Handpiece- Pana-MAX

45 Degree Angle




**Swivel adapter included**      35-16521      **394.95**

## Surgical Products



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Piezo Bone Surgery Tips (Vista)



<b>BT1</b> (Acteon Style)	56-11160	<b>112.95</b>
<b>BT3</b> (Acteon Style)	56-11161	<b>112.95</b>
<b>ST1</b> (Acteon Style)	56-11159	<b>112.95</b>
<b>ST2</b> (Acteon Style)	56-11158	<b>112.95</b>
<b>ST3</b> (Acteon Style)	56-11157	<b>112.95</b>
<b>ST4</b> (Acteon Style)	56-11156	<b>112.95</b>
<b>LT1</b> (Acteon Style)	56-11155	<b>112.95</b>
<b>BT1</b> (Mectron Style)	56-11154	<b>112.95</b>
<b>BT3</b> (Mectron Style)	56-11153	<b>112.95</b>
<b>ST1</b> (Mectron Style)	56-11152	<b>112.95</b>
<b>ST2</b> (Mectron Style)	56-11151	<b>112.95</b>
<b>ST3</b> (Mectron Style)	56-11150	<b>112.95</b>
<b>ST4</b> (Mectron Style)	56-11149	<b>112.95</b>
<b>LT1</b> (Mectron Style)	56-11148	<b>112.95</b>

### Surgical push button Handpiece 45 Degree

6 Month Warranty



### 4-Hole Surgical Handpiece

with swivel

35-121131 **299.95**

### Delicate Post-Surgical Brush (Butler)



Delicate, with ultra gentle 0.004 bristles Ideal for post-surgical cleaning or in the presence of aphthous ulcers, mouth irritations, extractions, implants and grafts.

Brush w/31 Tufts, 12/Pkg 54-98201 **14.95**

### Myco Surgical Blades

### QuickSmart Blade Removal System (Myco)



- Removes all common blades from most size No3, No4, No5, No7 and No9 scalpels.
- Immediate containment of blade upon removal
- For single handed use with either left or right hand


- OSHA Standard Compliant
- Bench or Wall mounted

Each 52-88981 **30.85**



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Scalpel Blade Handles (Miltex)



No. 4, Stainless Steel 5.375" graduated in mm and cm, For Blades 20-25	51-87073	<b>8.95</b>
No. 4, Nickel 5.375" graduated in mm and cm, For Blades 20-25	51-87074	<b>17.95</b>
No. 7, Stainless Steel 6.5" For Blades 10-15C	51-87075	<b>22.95</b>
No. 9, Stainless Steel 4.5" For Blades 10-15C	51-87076	<b>22.95</b>
SIEGEL, Stainless Steel round handle, 6" For Blades 10-15C	51-87077	<b>35.95</b>
No. 5, plastic handle with stainless tip 5", green For Blades 10-15C	51-87078	<b>27.95</b>
No. 6, plastic handle with stainless tip 5.25", red For Blades 20-25	51-87079	<b>27.95</b>
No. 3, E-Z-Off Reusable Plastic Surgical Knife Handles, 5.5", gray For Blades 10-15C	51-870710	<b>15.95</b>
No. 4, E-Z-Off Reusable Plastic Surgical Knife Handles, 5.5", blue For Blades 20-25	51-870711	<b>15.95</b>
Blade-Safe Surgical Blade Remover	51-11122	<b>17.95</b>
No. 3, Stainless Steel 5" graduated in mm and cm, For Blades 10-15C	51-87071	<b>8.95</b>
No. 3, Nickel 5" graduated in mm and cm, For Blades 10-15C	51-87072	<b>17.95</b>

### Soft Tissue Trimmer (Axis)



For soft tissue surgery, the Soft Tissue Trimmer can be used as an alternative to electrode surgery and surgical blades. The biocompatible hard oxide point cuts tissue without burning. Recommended to be used without water/spray cooling. Autoclavable. Friction Grip only.

**FG No STT, Size 0.016mm,**  
Head Length 5.0mm, each **52-88991 82.95**

## Surgical Products

### Generic Surgical Blades

#### Surgical Blade Handles

**Blade Handle #3 #991** (Fits #10-15C)




**Blade Handle #4 #992** (Fits #20 - 25)



<b>Blade Handle #3 fits # 10-15C</b>	51-11123	<b>5.99</b>
<b>Blade Handle #4 fits # 20-25</b>	51-11124	<b>5.99</b>

#### Surgical Stainless Steel Blades



<b>Surgical Blades #23</b>	51-11134	<b>9.99</b>
<b>Surgical Blades #22</b>	51-11133	<b>9.99</b>
<b>Surgical Blades #10</b>	51-11128	<b>9.99</b>
<b>Surgical Blades #21</b>	51-11132	<b>9.99</b>
<b>Surgical Blades #20</b>	51-11131	<b>9.99</b>
<b>Surgical Blades #15C</b>	51-11130	<b>9.99</b>
<b>Surgical Blades #15</b>	51-11129	<b>9.99</b>
<b>Surgical Blades #12D</b>	51-11127	<b>9.99</b>
<b>Surgical Blades #12</b>	51-11126	<b>9.99</b>
<b>Surgical Blades #11</b>	51-11125	<b>9.99</b>
<b>Surgical Blades #25</b>	51-11135	<b>9.99</b>

### Miltex Surgical Blades

#### Blade-Safe Blade Remover (Miltex)

An instrument for the safe removal of scalpel blades from handles.



- Protects hand against cuts
- Economically priced
- Stainless steel
- Autoclavable
- Easy to use

**Each** **52-88951 18.95**



## Surgical Products



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Carbon Surgical Blades (Miltex)



Blades # 25 100/Pkg	51-87521	<b>54.95</b>
Blades # 22 100/Pkg	51-87519	<b>54.95</b>
Blades # 23 100/Pkg	51-87520	<b>74.95</b>
Blades # 21 100/Pkg	51-87518	<b>54.95</b>
Blades # 20 100/Pkg	51-87517	<b>54.95</b>
Blades # 15C 100/Pkg	51-87516	<b>99.95</b>
Blades # 15 100/Pkg	51-87515	<b>54.95</b>
Blades # 12 100/Pkg	51-87513	<b>74.95</b>
Blades # 12B 100/Pkg	51-87514	<b>99.95</b>
Blades # 11 100/Pkg	51-87512	<b>54.95</b>
Blades # 10 100/Pkg	51-88510	<b>54.95</b>

### Stainless Steel Surgical Blades (Miltex)



Blades #25 100/Pkg	51-87511	<b>49.95</b>
Blades #22 100/Pkg	51-87509	<b>49.95</b>
Blades #23 100/Pkg	51-87510	<b>59.95</b>
Blades #21 100/Pkg	51-87508	<b>49.95</b>
Blades #20 100/Pkg	51-87507	<b>49.95</b>
Blades #15C 100/Pkg	51-87506	<b>76.95</b>
Blades #15 100/Pkg	51-87505	<b>49.95</b>
Blades #12 100/Pkg	51-87503	<b>69.95</b>
Blades #12B 100/Pkg	51-87504	<b>76.95</b>
Blades #11 100/Pkg	51-87502	<b>49.95</b>
Blades #10 100/Pkg	51-87501	<b>49.95</b>

### Sterile Carbon Steel Blades (Hu-Friedy)



#12 100/pkg	51-2000	<b>84.95</b>
#12D 100/pkg	51-2001	<b>150.95</b>
#15 100/pkg	51-2002	<b>84.95</b>
#15C 100/pkg	51-2003	<b>137.95</b>

### Disposable Sterile Scalpels With Handles (DA)



#10	51-20010	<b>10.95</b>
#11	51-20011	<b>10.95</b>
#12	51-20012	<b>10.95</b>
#15	51-20013	<b>10.95</b>
#12B	51-20014	<b>11.95</b>

### Scalpels With Handles (Miltex)



Sterile  
Disposable  
Box of 10

#10	51-20015	<b>14.95</b>
#11	51-20016	<b>14.95</b>
#12	51-20017	<b>14.95</b>

### Disposable Scalpels (Myco)



Sterile Scalpel #10	51-20018	<b>7.49</b>
Sterile Scalpel #11	51-20019	<b>7.49</b>
Sterile Scalpel #12	51-20020	<b>7.49</b>
Sterile Scalpel #12B	51-20021	<b>7.49</b>
Sterile Scalpel #15	51-20021	<b>7.49</b>
Sterile Scalpel #20	51-20022	<b>7.49</b>



**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

