

Ultrasonic Cleaners



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**



Ultrasonic Cleaners



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Medflex™ Premium Ultrasonic Cleaners (First Medica)

In order to pursue safe infection control practices and prevent cross contamination the first step is effective cleaning. First Medica's MedFlex™ Premium Ultrasonic Cleaners -provide a crucial step in ridding your instruments of infectious bio-burden. Our high quality stainless steel construction will provide years of reliable service and with the -standard heating element effective cleaning is easy and cost effective. Not only will our Ultrasonic cleaners do an effective job of cleaning, they will look good in your -professional setting!



Stainless steel construction
Acoustic pulse generators produce powerful cavitations for thorough and effective cleaning
Quiet operation
Stainless steel drain for care-free removal of cleaning solution
Scratch and corrosion resistant
Adjustable temperature range for precision cleaning
Adjustable timer
Instrument basket included

3 Liter	55-11122	369.95
5 Liter	55-11121	469.95
7 Liter	55-11120	739.95
10 Liter	55-11119	929.95
Hanging Basket, 3 liter	55-11118	89.95
Round Basket	55-11117	69.95
Hanging Basket, 5 liter	55-11116	95.95
Standing Basket, 5 liter	55-11115	95.95
Wire Rack, 7 liter	55-11114	129.95
Hanging Basket, 10 liter	55-11113	199.95

Ultrasonic Cleaner (ATW)



3L, with Timer, Basket, Cover w/o Drainage, 120 Watt, Inside Dimension 240×140×90 mm	55-11129	279.95
5L, with Heater, Timer, Basket, Cover & Drainage, 120 Watt, Inside Dimension 240×140×150 mm	55-11130	449.95

10L, with Heater, Timer, Basket, Cover & Drainage, 120 Watt. Inside Dimension 305×245×150 mm	55-11131	679.95
--	----------	---------------

Ultrasonic Cleaner (Tuttnauer)

Tuttnauer's Clean & Simple Ultrasonic Cleaning Systems are manufactured with only the highest quality components. These machines provide 60 kHz of quiet cleaning power which is less aggressive to instruments while still removing all surface contamination. The contemporary design and quiet operation compliment any office environment. Clean & Simple Ultrasonic Systems are available in 2 sizes with a full selection of accessories to meet your specific needs. The 1 gallon system holds a standard selection of instruments or 2 half size cassettes. The 3 gallon holds 3 full size cassettes.



The Tuttnauer CSU1 Table Top Ultrasonic Cleaner

Features:

- High efficiency 60 kHz cavitations for a thorough and quieter cleaning
- Durable vinyl clad steel housing... easy to clean and resists scratching
- Full line of accessories available
- 2 1/2, 5 & 10 year Limited Warranty**

Unit size: 1 Gallon	55-11132	494.95
Unit size: 3 Gallon	55-11134	1,094.95
Unit size: 3 Gallon with Heater	55-11135	1,394.95
Unit size: 1 Gallon with Heater	55-11133	529.95
SS Basket for 1 Gallon size	55-11136	104.95
Accessory Kit for 1 Gallon size: includes:positioning cover and 2-600ml beakers with lids and bands	55-11137	179.95
SS Basket for 3 Gallons size	55-11138	214.95
Accessory Kit for 3 Gallons size: includes:1/2 size ss basket, positioning cover and 2-600ml beakers with lids and bands	55-11139	269.95



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Ultrasonic cleaner Unit (TPC)

- Powerful 40 KHz of cavitations or thorough, quieter cleaning
- Stainless steel cover and tank with welded drain
- Durable vinyl clad steel housing is scratch resistant and easy to clean
- Digital timer for quick and easy process
- Sweep technology ensures exceptional cleaning efficiency



4.0 qt, with basket, Stainless steel lid,
Timer without Drain Kit 55-45681 **374.95**

10.6 qt, with basket, Stainless steel lid,
Timer & Drain Kit 55-45682 **694.95**

Brite Shield (Premier)



Enzymatic action is up to 35% faster than other brands
Chemo-mechanically debrides clogged dental burs, removes alginate and stone from impression trays. Also removes set glass ionomer, carboxylate and zinc phosphate cements. Ideal for removing bio-burden from hand instruments.

A protective shield coats the instrument through a unique passivation process to restrict rust and corrosion while imparting a lustrous surface to stainless-steel instruments.

800gm/Jar, Container makes up to 34 gallons of finished solution. 55-4921 **49.95**

Pro One SonicTab (Cory Lab)



These easy-to-use enzymatic tablets dissolve in a gallon of water in under a minute to create a powerful and effective ultrasonic cleaning solution.
Made from a biodegradable, non-chlorine, non-corrosive formulation that contains fast acting protease enzymes.

The solution will greatly assist in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents.
Virtually eliminates handling and storage of heavy liquid concentrates.

64 Tablets/Box 55-24891 **32.95**

Ultrasonic Cleaners

Efferzyme Ultrasonic Cleaner Tablets (Crosstex)



The first effervescent tablet for all your general purpose instrument and endoscopic cleaning.
Fast acting prolease enzymes formulated specifically for cleaning soils from medical instruments.
Just add 2 tablets plus one gallon of water to produce an effective pre-soak or

ultrasonic cleaning solution.
Also used as an evacuation system cleaner.

52 Tablets/Box 55-24881 **39.95**

Enzymax Powder (Hu-Friedy)

A lighter-weight concentrated powder allows Enzymax to be cost effective while saving valuable storage space. Enzymax powder is specially formulated with a fresh scent.



800gm Makes 55 Gallons 55-24951 **46.95**

EZ-Zyme (Miltex)



- Removes mucous and tissue from hemostats.
- Concentrated, Each 3/4 ounce packet makes 1 gallon of solution.
- Highly effective instrument presoak, ultrasonic and evacuation cleaning solution.

Box of 32 Packets 55-24871 **32.95**

Ultrasonic Cleaner (Yates & Bird)

Cleans any object that can be immersed in liquid. This versatile unit is ideal for both dental office and laboratory use. Compact enough (5" cube overall size) for convenient chair side use, but large enough with its 31/2" diameter x 2" deep tank to accommodate full dentures and many other laboratory items. Model 42 generates 42,000 high-intensity cleaning pulses every second that instantly clean any dental diamond, bur, instrument or dental appliance.



Unit Model 42 41-11648 **144.95**

Ultrasonic Cleaners



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Ultrasonic Enzymatic Tablets (Defend)



Ultrasonic Enzymatic Tablets

Easy-to-use tablets dissolve in a gallon of water in under a minute to create a powerful, effective cleaning solution. Greatly assists in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents. Contains fast-acting protease enzymes and has a pleasant scent. Also use as an instrument presoak or evacuation system cleaner.

- Safe for amalgam separators
- Biodegradable, neutral pH
- Non-chlorine • Non-corrosive

64 Tablets/Box

55-11142

32.95

Citrizyme (Pascal)

Ultra Concentrated Enzyme Cleaner and Instrument Presoak

- Fast-acting dual enzymes remove protein and biofilms quickly
- Leaves no residue
- Premium non-foaming formula
- Powdered concentrate dissolves instantly in warm tap water-no mixing required



300gm, Jar

55-24861

29.95

900gm, Canister

55-24862

89.95

All in One (Enzyme Industries)

- Ecologically sound
- Evacuation system cleaner
- Instrument pre-soak
- Organic stain cleaner



Gallon, Makes 32 Gallons

55-24901

18.95

Clean and Simple Ultra Tablet (Tuttnauer)

Clean & Simple is an ultrasonic / enzymatic solution dispensed in a unique tablet form. Drop 2 tablets into your ultrasonic tank already filled with tap water and the tablets will dissolve in about a minute. This is a non-ionic solution which is non-corrosive, non-toxic, biodegradable and non-caustic.

FEATURES:

- Easy to use tablet form.
- Eliminates handling of heavy 1 gallon containers.
- No bulky storage.
- No pre-mixing required.
- Enzymatic solution effectively breaks down protein.



64 tablets/Box

55-11140

54.95

Cetyl-Zyme Pro-Am (Cetylite)



Cetyl-Zyme Pro-Am is a low foaming, bacteriostatic, protease amylase dual enzymatic detergent for cleaning delicate medical and dental instruments/equipment.

It dissolves both protein and starch-containing soils and is useful for cleaning instruments with hard to reach areas such as flexible endoscope channels, box locks

and serrated edges. It may be used as a presoak or in manual, ultrasonic and automated cleaning systems or as an

- One gallon concentrated formula makes 128 gallons of use solution
- Neutral pH Formula
- Non-Corrosive
- Includes no mess easy-to-use 1oz metered pump dispenser
- Pleasant scent
- Use as a presoak or in manual, ultrasonic and automated cleaning systems

1Gallon Concentrate, 1 oz makes

1gallon of ready to use solution

55-4931

44.95

Enzymax Liquid (Hu-Friedy)

Enzymax is the only specially formulated enzymatic line on the market with dual enzymes and steel guard agents for added corrosion protection & sparkling instruments restoration



Box of 40 1/3 oz Packets,

1 Packet makes 1 Gal

55-24941

77.95

1 Quart Pump Bottle, 1 Packet

makes 1 Gal, 4 Pumps makes 1 Gal

55-24942

39.95

1 Gallon, No Pump

55-24943

77.95

Enzymax PAX (Hu-Friedy)

The dissolvable Enzymax PAX are the ultimate in convenience, just drop them in the water. There's no measuring or scooping needed, just use one packet per gallon of warm water.



Box of 32 powder packets, makes

32 gallons

55-24961

48.95

Box of 96 powder packets, makes

96 gallons

55-24962

112.95



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

General Purpose Ultrasonic Solutions (DA)

Concentrated solution that dilutes
1:10 with water.

Formulated especially to remove
debris from a variety of operatory
and laboratory items.

1 Gallon	55-89601	19.95
4 Gallon	55-89602	78.95

Good Vibrations™ (Certol)

Effective Cleaning Action

- Supports ultrasonic cleaning action.
- Removes organic blood debris and inorganic dental materials such as cement.
- Non-Foaming.

Protects Valuable Instruments

- Most effective anti-corrosive agents available.
- Enhances instrument appearance and shine.
- Safe for delicate surgical items and jointed orthodontic pliers.

Economical

- Concentrated Formula.
- 1 ounce makes 1 gallon solution

1 oz. Unit Dose – 24/Box	55-11111	20.50
1Gallon with Pump	55-11112	22.99



IMS Daily Clean (Hu-Friedy)

Concentrated powder detergent specifically formulated for safe, efficient ultrasonic cleaning and soaking of all dental/surgical instruments and accessories.



1 lb. (Makes 35 gal.)	55-89161	24.95
5 lbs. (Makes 160 gal.)	55-89162	77.95

Biosonic Enzymatic (Coltene Whaledent)

Enzymatic Ultrasonic Cleaning
Concentrate.



8oz. Concentrate	55-24911	33.95
-------------------------	----------	--------------

Ultrasonic Cleaners

ProEZ 1™ (Certol)

Single Enzymatic Detergent and ultrasonic
cleaning concentrate, instrument cleaner



**Single Enzymatic Detergent. Protease enzymes
attack blood based soils.**

- Most economical choice.
- Concentrated: One ounce makes one to two gallons of solution.
- May be used as holding solution or in ultrasonic cleaners.
- Neutral pH and safe for delicate instruments.
- Powerful surfactants added into the formula aid the enzyme cleaning action.
- Water softening and chelating agents work to prevent the minerals in your office tap water from interfering with the cleaning process and prevent corrosion on your valuable instrument.

1Gallon with Pump	55-11123	29.50
--------------------------	----------	--------------

ProEZ foam™ (Certol)

Ready-to-Use Foaming Enzymatic Instrument Pre-
cleaning Spray

ProEZ foam is an easy to use instrument pre-
cleaning spray

Powerful quad enzyme formula quickly dissolves
sticky blood soils while unique anti-corrosive system
protects and brightens delicate metal surfaces.



ProEZ foam is especially valuable in the dental setting when
soiled instruments must be held for extended periods of time.
Simply spray cleaning foam on tray of instruments.

24 oz. Spray Bottle	55-11127	14.00
1Gallon Refill	55-11128	42.00

**It's easy to join our mailing
list!**

**Just send your email address
by text message:**

**Text
PRESTIGEDENTAL
to 22828 to get started.**



Message and data rates may apply.

Ultrasonic Cleaners

Sonic Cleaner (DA)

- Detergent Concentrate.
- Safe, efficient, non-corrosive.
- 5 lbs. makes 160 gallons.



5 Lb. Can

55-89651 **44.95**

Tartar and Stain Remover (DA)



Efficiently removes calculus, tartar, tobacco, and food stains from removable appliances. Use at full strength.

1 Gallon

55-89621 **19.95**

4Gallon

55-89622 **78.95**

Temporary Cement Remover (DA)

Ideal for effectively removing temporary cement from a variety of surfaces. Use at full strength.

1 Gallon

55-89631 **19.95**

4Gallon

55-89632 **78.95**

Ultrasonic Concentrated Solution (Defend)



Ultrasonic Solutions General Purpose Cleaner (#1)

Non-ammoniated, 10:1 concentrated formula effectively cleans dried blood, tissue, pumice, food, rouge, Tripoli & buffing compounds from the surface of any appliance while providing a rust inhibitor.

1 Gallon Bottle

55-11144 **19.99**

1 Case of 4 Gallons

55-11143 **69.99**



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

Ultrasonic Heavy-Duty Cleaner (DA)



Super-concentrate general purpose cleaner designed for most cleaning applications in the dental laboratory.

1 Gallon

55-89641 **19.95**

4Gallon

55-89642 **78.95**

Plaster and Stone Remover (DA)

Removes plaster and stone residue from metal trays and articulators.


1 Gallon

55-89611 **19.95**

4Gallon

55-89612 **78.95**

www.PrestigeDentalProducts.com



**Check our
website for more
items and
variation**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Ultrasonic Cleaners

Just to mention a few of our products. To see more visit us online at:
www.prestigedental.com

Best Dental
Supplies!!

Tartar & Stain
Remover

DA



Enzymax
Liquid

Hu-Friedy



Coltene Whaledent

Biosonic Enzymatic



Good Vibrations

Certol



*We Work Hard
So You Don't Have To.*



Call Us at 877-772-3888



www.prestigedental.com

Ultrasonic Cleaner Tuttnauer

Features:

- High efficiency 60 kHz cavitation for a thorough and quieter cleaning
- Durable vinyl clad steel housing... easy to clean and resists scratching
- Full line of accessories available
- 2 1/2, 5 & 10 year Limited Warranty



Follow Us On Social Media!!!

facebook

twitter

LinkedIn

Call us for Quarterly Promotion

Prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

