

# Ultrasonic Scalers



**Prestige Dental Products, Inc**

**PrestigeDental.com**  
**877-772-3888**

# Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

**As an authorized Dealer of many brands in North America**, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

**Prestige Dental Products** offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

**Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and**

## Customer Information

### Order

**PHONE** To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

**FAX** You may FAX your orders to 714-632-8821.

**INTERNET** [www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



*Note:* New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



### Customize order for special customer

Email us to [mass@prestigedentalproducts.com](mailto:mass@prestigedentalproducts.com) a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

## Payment

**CREDIT CARD** We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

### Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

### MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

## Shipping

### GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

**EXPEDITED SHIPPING** UPS Next Day and Second Day Service for which additional shipping fees will be charged.

**SHIPPING HAZARDOUS ITEMS** All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

## State and federal laws

### DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog  
Are Subject To Change Without Prior Notice.**



# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Advance Piezo Ultrasonic Scaler & Endodontics (TPC)



### Standard Scaler Features:

- Compact unit housing for easy placement anywhere in the operatory
- Extremely gentle on patient, causes little or no sensitivity
- No heat build up in handpiece tip with extended use
- Autoclavable handpiece
- Extensive power range to allow widest choice of clinical applications
- Fine water spray eliminates excessive overspray onto patient
- Thin sub gingival tips available
- 5 tips included with unit
- 1 year warranty

### Optional Endodontics Unit Features:

- Endodontic file holder available
- Easy and rapid operating for root canal treatment
- Achieve enlargement of canal quickly specifications
- Input rating: 115/230V, 50/60 Hz
- Dimensions: 160mm(W) x 155mm(D) x 70mm(H)
- Weight: 2.1kg
- Frequency: 32 KHz
- 1 year warranty (tips excluded)

<b>Piezo Ultrasonic Scaler &amp; Endodontic Unit</b>	56-17031	<b>394.95</b>
<b>Endo files and holder (21mm #25)</b>	56-17032	<b>69.95</b>
<b>Piezo scaler tip #1</b>	56-17033	<b>25.95</b>
<b>Piezo scaler tip #2</b>	56-17034	<b>25.95</b>
<b>Piezo scaler tip #3</b>	56-17035	<b>25.95</b>
<b>Piezo scaler tip #4</b>	56-17036	<b>25.95</b>
<b>Piezo scaler tip #5</b>	56-17037	<b>25.95</b>
<b>Piezo scaler handpiece</b>	56-17038	<b>195.95</b>
<b>Tip wrench</b>	56-17039	<b>20.95</b>
<b>Piezo scaler O-rings (5 pack)</b>	56-17040	<b>9.95</b>
<b>Replacement water line</b>	56-17041	<b>23.95</b>
<b>Replacement scaler foot control</b>	56-17042	<b>79.95</b>

### ART-MB3 Ultrasonic Scaler Built in type (Bonart)

Power Package includes: PC Board, LCD digital display screen, VR Switch Set, Terminal Block, Magnet Hand Piece w/ Cable Assembly, 2 Inserts. Optional parts includes: Solenoid Valve, Transformer, and Foot Control at additional costs. Available in either 25Khz or 30Khz system.

**Available now in all regions.**

### Warranty Terms:

- **5 year limited** warranty on the internal PC board and electrical components
- **1 year limited warranty** on the hand piece assembly, foot control, power cord, LCD display screen etc.
- 120 days on the expendable part like inserts



<b>Unit, 25 KHZ with 2 Insert</b>	56-11119	<b>499.00</b>
<b>Unit, 30 KHZ with 2 Insert</b>	56-11120	<b>499.00</b>





Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Ultrasonic Scalers

### ART-M1 Ultrasonic Scaler (Bonart)

Powerful easy to use scaling performance that outlasts its competitors. Excellent quality at an affordable price.

Specifications: Operation Frequency: AVAILABLE ONLY IN 25KHZ. Output Power: From 10 to 30 Watts of output power during normal scaling speed. Turbo Function: Increase power during scaling with one touch of a button. C25W to 35W output power. Water Supply: 25-60 PSIG (172-414 Kpa) Optional accessories: 25KHz inserts ISO Certified; CE Marked; 510K registered. Warranty Terms: 18 month limited warranty on the internal electronic module; 1 year warranty for hand piece, foot control; 120 days warranty for all insert tips.



<b>Complete Unit with 2 Insert 30KHZ</b>	56-11111	<b>344.95</b>
<b>Complete Unit with 2 Insert 25KHZ</b>	56-11112	<b>344.95</b>

### ART-M4 Ultrasonic Scaler (Bonart)



Blue is the only System is available in either 25KHz or 30KHz stand-alone system. Up to a maximum power setting of 50 with auto tune system that locks the frequency and stabilize power setting The well-equipped feature called "AGC" (Auto Gain Control) stabilizes the power setting and prevents stalling during scaling operation. System comes with a built in LCD digital screen that displays the mode and power setting, and allows the user to choose from Scaling, Perio and Endo mode with one turn of the knob, and an easy to access water control knob at the front end of the hand piece cable.

System includes 1 plastic grip internal flow insert and 1 grip handle external flow insert. ISO certified; CE marked; 510K registered. Warranty Terms: 5 year limited warranty on the internal PC board and electrical components 1 year limited warranty on the hand piece assembly, foot control, power cord etc. 120 days on the expendable part like insert

<b>Complete Unit with 2 Insert 30KHZ</b>	56-11117	<b>424.95</b>
<b>Complete Unit with 2 Insert 25KHZ</b>	56-11118	<b>424.95</b>

### ART-M3 II Ultrasonic Scaler (Bonart)

Has our unique "AGC" (Auto Gain Control) LCD digital display screen that shows scaling mode and power setting. System is available in either 25KHz or 30KHz stand-alone system. Up to a maximum power setting of 50 with auto tune system that locks in the frequency and stabilize power setting. The well-equipped feature called "AGC" (Auto Gain Control) stabilizes the power setting and prevents stalling during scaling operation. System

comes with a built in LCD digital screen that displays the mode and power setting, and allows the user to choose from Scaling, Perio and Endo mode with one turn of the

knob, and an easy to access water control knob at the front control panel. System includes 1 plastic grip internal flow insert and 1 metal grip external flow insert. ISO certified; CE marked; 510K registered. Warranty



#### Complete Unit with 2 Insert 30KHZ

Inside the box & Another Universal Insert Free 56-11115 **449.95**

#### Complete Unit with 2 Insert 25KHZ

& Another Universal Insert Free 56-11116 **449.95**

### Biosonic Ultrasonic Scaler System (Coltene Whaledent)

#### New Features:

- Water intensity and power control knobs now rotate in the same direction
- On demand turbo button
- Significantly less expensive than competitors units
- Operates with either 25KHz or 30KHz inserts
- Turbo button for an operator-controlled power boost
- Accepts most competitive magnetostrictive inserts
- Units small footprint saves space in the operatory
- No bulky external power supply
- 115V Model



**Complete Unit** 56-51247 **724.95**



# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## ART-P3 II Piezo Electric Scaler (Bonart)

Up to a maximum power setting of 50 with auto tune system that locks the frequency. The well-equipped feature called "AGC" (Auto Gain Control) stabilizes the power setting and prevents stalling during scaling operation. System comes with a built in LCD digital screen that displays the mode and power setting, allows the user to choose from Scaling, Perio and Endo mode with one turn of the knob, and has an easy to access water control knob at the front control panel. System includes 4 various tips (2 pointed tips & 2 flat tips), 4 individual plastic autoclavable wrenches, 1 aluminum sterilization box. ISO certified; CE marked; 510k registered. Warranty Terms: 5 year limited warranty on the internal PC board and electrical components 2 year limited warranty on detachable piezo hand piece 1 year limited warranty on the hand piece cable, foot control, power cord etc. 120



**Complete Unit with 4 Piezo Tips** 56-11121 **479.95**

## ART-P6 Piezo Electric Scaler (Bonart)

ART-P6 Piezo Electric Compact Sized Scaler

Unit is designed and manufactured to meet the most demanding scaling environment. System is portable and lightweight weighing only 3.52 lbs; it comes with a built in power setting knob (1-20) that allow users to set power output at their desired level of operation. Technical specifications are as follows:



- Power Supply: Available in either 115V or 230V
  - Working Freq: 28KHz to 32KHz
  - Water Supply: 25-60Psi (172-414kPA)
  - Dimension: 7.68"(L)x3.38"(W)x3.27"(H)
  - Weight: 3.52 lbs
  - Handpiece Cable Length: 8.52 feet
  - Footswitch Cable Length: 8.52 feet
  - Operation Environment Temperature: 0-140 degrees F
  - Relative Humidity: 10-90% (non-condensing) system
- includes 4 tips and 1 plastic wrench** in US; 3 tips and 1 plastic wrench in other regions. **5 year power unit warranty**

**Complete Unit with 4 Piezo Tips** 56-11125 **339.95**

## ART-PB3 Piezo Electric Built in Type Scaler (Bonart)



Power Package includes: PC Board, LCD display screen, VR Switch Set, Piezo Titanium Hand Piece w/ Cable Assembly, 4 Tips, 4 wrenches, 1 sterilization box. Optional parts includes: Solenoid Valve, Transformer, Foot Control at additional costs.

### Warranty Terms:

**5 year limited warranty** on the internal PC board and electrical components


**2 year limited warranty** on detachable Piezo hand piece

**1 year limited warranty** on the hand piece cable, foot control, power cord etc.

120 days on the expendable part like Piezo tips

**Complete Unit with 4 Piezo Tips** 56-11126 **395.95**

## NiTi Endo U-Files



Color Coded used by Ultrasonic scaler units for root canal cleaning

Tip with excellent tenacity is used for root canal cleaning  
File size: 32.5mm±0.5mm  
Handpiece size: Ø0.68mm-0.71mm  
K-file taper: 2%

<b>#10, 33mm, 6/Pkg</b>	29-12247	<b>12.95</b>
<b>#15, 33mm, 6/Pkg</b>	29-12246	<b>12.95</b>
<b>#20, 33mm, 6/Pkg</b>	29-12245	<b>12.95</b>
<b>#25, 33mm, 6/Pkg</b>	29-12244	<b>12.95</b>
<b>#30, 33mm, 6/Pkg</b>	29-12243	<b>12.95</b>
<b>#35, 33mm, 6/Pkg</b>	29-12242	<b>12.95</b>
<b>#40, 33, 6/Pkg</b>	29-12241	<b>12.95</b>
<b>Full Set, 42/Pkg 6 each size</b>	29-12248	<b>79.95</b>



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Cavitron Jet Combination System SPS (Dentsply)



#### Cavitron Jet Combination System SPS

56-11128 3399.95

### Piezo Advantage (Vista)



Little or no hand piece vibration virtually eliminates hand fatigue and patient sensitivity, while still powerful enough to remove heavy calculus.

- Compact, functional design with removable and autoclavable control covers
- It can easily be wall, pole or table-top mounted
- Extensive range of Hygiene/Perio tips and diamond coated Endodontic tips
- Compatible with Vista's Hygiene Irrigator™ System
- Table Top or Modular Unit Styles

#### Table Top Unit Includes:

- Piezo Advantage™ Ultrasonic Unit
- 3-Autoclavable hand piece sleeves
- 3-Tip Wrenches
- 6-Ultrasonic Tips of your choice (\*HST, UCT, Safety 20-4D and 40-4D Tips not included in this offer.)
- Choice of wall or pole mount
- Modular unit also includes insulation material

**Suggested Tips: SU1, HA, HC, HFL, HFR, PFL, PFR, GE20, GE40**

<b>Piezo Advantage Table Top Unit</b>	56-11145	<b>1069.95</b>
<b>Piezo Advantage Modular Unit</b>	56-11144	<b>999.95</b>
<b>Piezo Tip Sterilization Tray</b> (for P5 & EMS Style Tips)	56-11143	<b>89.95</b>
<b>Tip Wrench</b>	56-11142	<b>39.95</b>
<b>Replacement Sleeve</b>	56-11141	<b>36.95</b>
<b>Quick-a Tip Sleeve</b>	56-11140	<b>36.95</b>
<b>Piezo Tip Sterilization Tray</b> (for Piezo+/Advantage Tips)	56-11146	<b>89.95</b>

### Focused Spray Slimline 1000 Cavitron Insert (Dentsply)

Cavitron FSI SLI 1000 Ultrasonic Insert offers improved access and familiar adaptation so you can efficiently remove deposits while maintaining patient comfort

- Triple bend design
- Slimline tip design
- 30KHz
- Coined edge design
- Focused Spray technology
- Can be used at all power levels

**30k Insert each** 56-91953 **158.95**

### Jetpolisher 2000– Air Polisher (Deldent)



Advanced compact, non-clogging airpolishing system for about half the price of standard countertop polisher

*Recommended to clean tooth surfaces prior to bonding procedures, placing sealants, and bleaching.*

- Non-clogging nozzle due to exclusive H.S.T. technology.
- Harmless to titanium implants.
- No installation required, completely portable, and can be moved easily from room to room
- Measure: 3.5" x 3.5" x 7.5"
- Connects directly to existing hand piece tubing (30-35 psi) via 2-3 hole or 4 hole connectors (*kavo, W&H, Sirona connectors*)
- Using existing dental unit's foot pedal
- Lightweight and easy on the wrist
- Fully detachable autoclavable hand piece and water control knob to prevent cross-infection
- Large, easy-view powder reservoir.
- **FDA cleared for use around implant.**
- **2 years warranty.**

<b>Midwest Type with 4-Hole connector</b>	56-89691	<b>1,195.95</b>
<b>Kavo Type connector</b>	56-89692	<b>1,295.95</b>
<b>Sirona/Siemens Type connector</b>	56-89693	<b>1,295.95</b>
<b>W&amp;H Type connector</b>	56-89694	<b>1,295.95</b>
<b>JetStream Standard handpiece</b>	56-89695	<b>164.95</b>
<b>Booster Handpiece</b> (for heavy stain)	56-89696	<b>179.95</b>
<b>Spray Head Tubes</b> (50/Pkg)	56-89697	<b>28.95</b>

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## DelSonic 2000 Piezo Scaler (Deldent)



- Sleek, compact, quiet, powerful, and portable.
- Linear tip vibration, for more efficient and more comfortable scaling.
- High frequency 30 kHz design provides greater patient comfort than the competition (less vibration against teeth).
- Cool running, high torque hand piece (minimal tip warmth).
- Fully autoclavable hand piece sleeves, so purchasing additional hand pieces is not necessary.
- Range of tips available for all clinical needs.
- Tips are designed to introduce spray even into deep pockets.
- May be used with any independent fluid supply system dispensing purified water or medicaments.
- 2 year warranty (compare with competitive brands' one year warranty on handpiece).

<b>Full Unit</b>	56-25841	<b>995.95</b>
<b>#1 Universal Piezo Tip</b>	56-25842	<b>69.95</b>
<b>#2 Round Straight Piezo Tip</b>	56-25843	<b>89.95</b>
<b>#2 Round Right Piezo Tip</b>	56-25844	<b>89.95</b>
<b>#2 Round Left Piezo Tip</b>	56-25845	<b>89.95</b>
<b>#3 Flat Horizontal Piezo Tip</b>	56-25846	<b>69.95</b>
<b>#4 Perio Straight Piezo Tip</b>	56-25847	<b>89.95</b>
<b>Autoclavable Handpiece Sleeves (set of 3)</b>	56-25848	<b>38.95</b>
<b>Tip Changing Tool</b>	56-25849	<b>53.95</b>

## Jetsonic 2000M Combo - With Cavitron - Type Scaler (Deldent)

Advanced combination scaling and polishing unit with patent H.S.T. air polishing system and advanced scaling design features.

The air polishing portion of the JetSonic 200M has all the features of the JetPolisher 2000 **with the addition of warm water flow.**

The scaler side of this unique unit also offers advanced magneto restrictive technology, including

- Unique sensor technology recognizes and activates the selected hand piece (scaling or polishing)
- Advanced patented digital circuitry **automatically tunes any "stack-type" insert, whether 25K, or 30K.**
- Advanced microchip feedback system maintains optimum scaling power for all settings.
- Power-Surge Mode technology increases power by 15-20% with a "double-click" on the foot switch for removal of tenacious deposits. The pre-selected power setting returns automatically when the foot control is released
- Large, easy-view powder reservoir (no need to stop to add powder).
- Convenient rinse mode available
- Water heating feature



<b>Combo Unit</b>	56-88681	<b>3084.00</b>
<b>Stream Selector 2000 (Fluid supply system): Unit, 2 Bottles and 2 Quick Disconnect</b>	56-88682	<b>244.95</b>
<b>Bottle for Stream Selector 2000, each</b>	56-88683	<b>11.95</b>





Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Piezo Co-Pilot Scaler (Vista)



The piezo **Co-Pilot™** is a beautifully designed and powerful self-contained dual reservoir scaling system. The Co-Pilot™ accepts a wide variety of piezo tips from major manufacturers and can be ordered with either a standard (Acteon\*) or metric (EMS\*) style hand piece.

- Two large, easy-to-fill 350ml reservoir tanks with easy-to-see liquid levels
- Accepts a wide variety of ultrasonic scaling tips for General Scaling, Periodontal, Endodontic and Apical Surgery applications
- Advanced circuitry, touchpad controls, and push-button controlled variable irrigation flow settings
- Compact design for improved productivity and portability (needs NO operatory air/water connections)
- Soft-lit dome and reservoir. Autoclavable handpiece sleeves available in six contemporary colors

#### System Includes:

- 1 - Piezo Pilot" ultrasonic scaler unit,
- 4 - Ultrasonic tips and 4 torque wrenches (Either Hygiene or Endo tip packages)
- 1 - Detachable handpiece (either Universal Acteon\* or Metric EMS\* thread)
- 3 - Color-matched autoclavable handpiece sleeves
- 1 - Quick-disconnect water line and replaceable in-line water filter
- 1 - Omni-directional foot pedal

<b>Acteon Style Unit, Grey Color</b>	56-11162	<b>1999.95</b>
<b>Acteon Style Unit, Purple Color</b>	56-11163	<b>1999.95</b>
<b>Acteon Style Unit, Blue Color</b>	56-11164	<b>1999.95</b>
<b>Acteon Style Unit, Yellow Color</b>	56-11165	<b>1999.95</b>
<b>Acteon Style Unit, Blue/Green Color</b>	56-11166	<b>1999.95</b>
<b>Acteon Style Unit, Pink Color</b>	56-11167	<b>1999.95</b>
<b>EMS Style Unit, Grey Color</b>	56-11168	<b>1999.95</b>
<b>EMS Style Unit, Purple Color</b>	56-11169	<b>1999.95</b>
<b>EMS Style Unit, Blue Color</b>	56-11170	<b>1999.95</b>
<b>EMS Style Unit, Yellow Color</b>	56-11171	<b>1999.95</b>
<b>EMS Style Unit, Blue/Green Color</b>	56-11172	<b>1999.95</b>
<b>EMS Style Unit, Pink Color</b>	56-11173	<b>1999.95</b>

### Piezo Handpiece (Bonart)



Piezo handpiece truly defines the innovative modern piezo technology developed by Bonart. The outer casing of the hand piece is made of a specialized material that can withstand steam sterilization of 135 degrees Celcius. Its newly improved transducer design offers a more powerful scaling, endo and perio vibrating frequency with less pain and great ease. It comes with a protective cap to be installed at the bottom portion of the hand piece during autoclave session to prolong the life of the entire hand piece. ISO certified; CE Pending. 510k registered!. 2 year limited warranty.

**Each** 56-11184 **129.95**

### ProSonic Air Scaler (Pac-Dent)



The ProSonic™ air scaler is a pneumatically driven device that generates an elliptical oscillation of the working tips in the sound range (< 6.5 kHz) with longitudinal and transverse movements. It is as powerful as it gets without being an Ultra Sonic, and has been clinically proven to be efficient and gentle with a wide range of uses, including: Hygiene, Periodontics, Minimally invasive caries therapy, Cavity shaping and finishing, Endodontics and Cementation. It features a controlled water spray that prevents the "drowning effect" often experienced by patients during deep scaling. The scaler is entirely autoclavable / chemically, and the tips can also be dry heat sterilized. To further safeguard the operator and patient from cross contamination, it also comes with replaceable tactile silicone grips which also relieves finger fatigue. The ProSonic™ has a slim and lightweight design that permits complete scaling even in tight interproximal areas. Packaged with 5 scaling tips - Universal, Sickle, Perio, Perio Left and Perio Right

- Sonic technology effectively removes calculus deposits and dental stains.
- Connects to existing air driven tubings and foot control
- Lightweight - helps to relieve wrist and finger fatigue
- Can be safely used on patients with pacemakers
- Fully autoclavable
- Ergonomic & color coded silicone grips
- Selection of tips for different scaling needs

#### ProSonic™ Air Scaler Unit: 1x

ProSonic air scaler, 1 x universal scalping tip, 1 x sickle scaling tip, 1 x perio scaling tip

56-55401 **414.70**

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Piezo Pilot Scaler Unit (Vista)



Introducing the Piezo Pilot™ a New Generation of ultrasonic scalers. Featuring state-of-the-art circuitry and a compact modern design, the Pilot is a lightweight powerhouse. The soft-lit dome and removable / autoclavable hand piece sleeves are available in six contemporary colors to adapt to any operatory decor. The unit also features a magnetic handpiece holder that can be positioned on either side of the Pilot or onto any magnetic surface.

A complement to any office decor, the Piezo Pilot™ is a powerful, lightweight, and beautifully designed piezo scaler that allows the hygienist to work quickly and effortlessly. The Pilot's state-of-the-art circuitry offers a wide range of amplitude settings for regular hygiene procedures and all ultrasonic endodontic procedures, including certain bone surgery applications. It comes with a choice of either four ultrasonic hygiene or four ultrasonic endodontic tips, your choice of either a Universal (P5) or Metric (EMS) handpiece, and is compatible with tips from a variety of major tip manufacturers, including Vista's full line of piezo tips.

<b>Acteon Style Unit, Grey Color with Hygiene 4 Tips Package</b>	56-11186	<b>1099.95</b>
<b>Acteon Style Unit, Purple Color with Hygiene 4 Tips Package</b>	56-11187	<b>1099.95</b>
<b>Acteon Style Unit, Blue Color with Hygiene 4 Tips Package</b>	56-11188	<b>1099.95</b>
<b>Acteon Style Unit, Yellow Color with Hygiene 4 Tips Package</b>	56-11189	<b>1099.95</b>
<b>Acteon Style Unit, Blue/Green Color with Hygiene 4 Tips Package</b>	56-11190	<b>1099.95</b>
<b>Acteon Style Unit, Pink Color with Hygiene 4 Tips Package</b>	56-11191	<b>1099.95</b>
<b>Acteon Style Unit, Grey Color with Hygiene 4 Tips Package</b>	56-11192	<b>1099.95</b>
<b>Acteon Style Unit, Purple Color with Endodontic 4 Tips Package</b>	56-11193	<b>1099.95</b>
<b>Acteon Style Unit, Blue Color with Endodontic 4 Tips Package</b>	56-11194	<b>1099.95</b>
<b>Acteon Style Unit, Yellow Color with Endodontic 4 Tips Package</b>	56-11195	<b>1099.95</b>
<b>Acteon Style Unit, Blue/Green Color with Endodontic 4 Tips Package</b>	56-11196	<b>1099.95</b>
<b>Acteon Style Unit, Pink Color with Endodontic 4 Tips Package</b>	56-11197	<b>1099.95</b>
<b>EMS Style Unit, Grey Color with Hygiene 4 Tips Package</b>	56-11198	<b>1099.95</b>
<b>EMS Style Unit, Purple Color with Hygiene 4 Tips Package</b>	56-11199	<b>1099.95</b>
<b>EMS Style Unit, Blue Color with Hygiene 4 Tips Package</b>	56-11200	<b>1099.95</b>
<b>EMS Style Unit, Yellow Color with Hygiene 4 Tips Package</b>	56-11201	<b>1099.95</b>
<b>EMS Style Unit, Blue/Green Color with Hygiene 4 Tips Package</b>	56-11202	<b>1099.95</b>
<b>EMS Style Unit, Pink Color with Hygiene 4 Tips Package</b>	56-11203	<b>1099.95</b>
<b>EMS Style Unit, Grey Color with Endodontic 4 Tips Package</b>	56-11204	<b>1099.95</b>
<b>EMS Style Unit, Purple Color with Endodontic 4 Tips Package</b>	56-11205	<b>1099.95</b>
<b>EMS Style Unit, Blue Color with Endodontic 4 Tips Package</b>	56-11206	<b>1099.95</b>
<b>EMS Style Unit, Yellow Color with Endodontic 4 Tips Package</b>	56-11207	<b>1099.95</b>
<b>EMS Style Unit, Blue/Green Color with Endodontic 4 Tips Package</b>	56-11208	<b>1099.95</b>
<b>EMS Style Unit, Pink Color with Endodontic 4 Tips Package</b>	56-11209	<b>1099.95</b>



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### Scalex 800 Ultrasonic Scaler (Dentamerica)



The scalex 800 Ultrasonic Scaler is 25,000 scaling strokes per second, solid state power unit.

- Streamlined compact size and slim light weight handpiece
- Convenient foot control
- No power loss at maximum

pressure

- Unit is available in 100/110V 50-/60 Hz, or 220/240V 50 Hz operation
- Dimensions 22L x 19W x 7.5H cm, weight 3.6kg

Unit 25K, Tips not Included 56-11267 **354.95**

F4 Scaler Universal, External Water Tip 25K each 56-11268 **79.95**

F1 Scaler Universal, Internal Water Tip 25K each 56-11269 **49.95**


### Scalex 880 Piezo Scaler (Dentamerica)

The Scalex 880 Piezo Scaler has State of the art piezoelectric technology is self-tuning and load adjusting for maximum smooth operation

- Streamlined compact size and slim lightweight handpiece.
- Innovative water-cooled circuitry designed for greater stability during operation

- Extremely gentle on the patient
- The Scalex 880 Piezo Scaler has Available a wide variety of high quality tips to correspond with specific clinical applications

- Three universal tips included. **1 year warranty**



Unit with 3 Universal A Tips 56-11276 **444.95**

Piezo Handpiece only 56-11277 **249.95**

Piezo Scaler Tip (Universal A) each 56-11278 **24.95**

Piezo Scaler Tip (Universal B) each 56-11279 **24.95**

Piezo Scaler Tip (Chisel) each 56-11280 **24.95**

Piezo Scaler Tip (Perio) each 56-11281 **44.95**

Endo Chuck each 56-11282 **84.95**

## Ultrasonic Scalers

### Scalex 850 Self Contained Ultrasonic Scaler (Dentamerica)

Scalex™ 850 offers self-contained solution tank for total portability, adding therapeutic medicaments



The **Scalex 850 Self-Contained Ultrasonic Scaler** has all the cleaning power

you need with performance guaranteed for one full year. Self-contained solution tank allows for convenient portability and the addition of medicated solutions to scaler therapy. Automatic tuning activates inserts in a magnetostrictive cleaning pattern for feather-touch removal of calculus and stains under all loading conditions—no power loss under pressure. Each system includes base unit with graduated 1000ml self-contained solution tank, handpiece with extra-long 7 ft. cord and one 30kHz insert. 110V.

**1 year manufacturer's warranty.**

Unit, 30K with 1 Tip 56-11273 **849.95**


F6 Scaler Universal External Water Tip 30K each 56-11274 **94.95**

F3 Scaler Universal Internal Water Tip 30K each 56-11275 **59.95**

Water Tank 56-12275 **59.95**

### Scalex 830 Ultrasonic Scaler (Dentamerica)

- 30,000 scaling strokes per second, solid state power unit.
- Streamlined, compact size and slim lightweight hand piece
- Innovative water-cooled circuitry designed for greater stability during operation
- No power loss at maximum pressure
- Available selection of a wide variety of high quality inserts t correspond to specific clinical applications
- Unit for 25K scaling strokes is available
- Specifications : 110/220V 50/60Hz AC
- Dimensions : L140 x W140 x H50 mm
- Weight : 1 kg
- Accessories : Control Box, Hand piece, Foot Control, Water Connector, Tube, Power Cable



Unit 30K with 1 Tip 56-11272 **469.95**

F6 Scaler Universal External Water Tip 30K each 56-11271 **69.95**

F3 Scaler Universal Internal Water Tip 30K each 56-11270 **49.95**



# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Scalex 890 Self-Contained Piezo Scaler (Dentamerica)

Scalex 890 offers self-contained solution tank for total portability, adding therapeutic medicaments



The **Scalex 890 Self-Contained Piezo Scaler** has all the cleaning power you need with performance guaranteed for one full year. Self-contained solution tank allows for convenient portability and the addition of medicated solutions to scaler therapy. State of the art Piezo Electric Self-tuning And load adjusting technology for smooth and effortless operation.

Innovative water cooled circuitry designed for greater stability during operation.

**1 year manufacturer's warranty.**

<b>Unit with 3 Universal A Tips</b>	56-11289	<b>779.95</b>
<b>Piezo Handpiece only</b>	56-11288	<b>229.95</b>
<b>Piezo Scaler Tip (Universal A) each</b>	56-11287	<b>24.95</b>
<b>Piezo Scaler Tip (Universal B) each</b>	56-11286	<b>24.95</b>
<b>Piezo Scaler Tip (Chisel) each</b>	56-11285	<b>24.95</b>
<b>Piezo Scaler Tip (Perio) each</b>	56-11284	<b>49.95</b>
<b>Endo Chuck each</b>	56-11283	<b>79.95</b>

## Sonic Air Scaler (TPC)

### Features:

- Includes 3 tips AS1, AS2, AS3
- Tip wrench included
- Light weight
- Portable
- Saves time
- Autoclaveable
- Available in 2 or 4 hole

### Package includes the following:

- AS5000 Air sonic scaler handpiece
- ASI Universal Tip
- AS2 Sickle Tip
- AS3 Perio Tip
- A776 Tip Wrench
- **6Month warranty**



<b>4 Hole Kit</b>	56-17051	<b>198.95</b>
-------------------	----------	---------------

<b>2 Hole Kit</b>	56-17052	<b>198.95</b>
<b>ASI Universal Tip</b>	56-17053	<b>19.95</b>
<b>AS2 Sickle Tip</b>	56-17054	<b>19.95</b>
<b>AS3 Perio Tip</b>	56-17055	<b>19.95</b>


## Ultrasonic Inserts

### After Five Plus Ultrasonic Inserts (Hu-Friedy)

#### Sub-gingival

- External fixed water delivery
- Thin, probe-sized design
- Provides smooth, deep access into periodontal pockets
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use

AFTER FIVE PLUS - For Subgingival



<b>Left, 25 kHz, Teal Handles, each</b>	56-90975	<b>136.95</b>
<b>Straight, 25 kHz, Blue Handle, each</b>	56-90976	<b>136.95</b>
<b>Right, 25 kHz, Red Handle, each</b>	56-90977	<b>136.95</b>
<b>Set, 25 kHz, each</b>	56-90978	<b>374.95</b>
<b>Left, 30 kHz, Teal Handles, each</b>	56-90979	<b>136.95</b>
<b>Straight, 30 kHz, Dark Blue Handle, each</b>	56-90980	<b>136.95</b>
<b>Right, 30 kHz, Red Handle, each</b>	56-90981	<b>136.95</b>
<b>Set, 30 kHz, each</b>	56-90982	<b>374.95</b>

### After Five Ultrasonic Inserts (Hu-Friedy)

#### Sub-gingival

- External water delivery
- Provides smooth, deep access into periodontal pockets
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use
- Thinner Tip, Slender Taper with Metal

HU-FRIEDY

AFTER FIVE - For Subgingival



<b>Straight, 25 kHz, each</b>	56-90983	<b>127.95</b>
<b>Set of 3:</b>		
Left, Straight, Right, 25 kHz, each	56-90984	<b>344.95</b>
<b>Straight, 30 kHz, each</b>	56-90985	<b>127.95</b>
<b>Set of 3:</b>		
Left, Straight, Right, 30 kHz, each	56-90986	<b>344.95</b>
<b>Left, 3. kHz, each</b>	56-90987	<b>127.95</b>



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Biosonic Ultrasonic Inserts (Whaledent)

BioSonic Inserts - variety of most popular tip designs in both 25 kHz and 30 kHz

Better performance reliability and longevity through patented manufacturing process  
 Ergonomically designed grip for operator's comfort and reduced hand stress

Maximum comfort in popular tip designs



<b>Biosonic Ultrasonic Scaler Tip Slim #10 (25kHz)</b>	56-21110	<b>98.95</b>
<b>Biosonic Ultrasonic Scaler Tip Slim #10 (30kHz)</b>	56-21111	<b>98.95</b>
<b>Biosonic Ultrasonic Scaler Tri-Bend #1000 (25kHz)</b>	56-21112	<b>98.95</b>
<b>Biosonic Ultrasonic Scaler Tri-Bend #1000 (30kHz)</b>	56-21113	<b>98.95</b>
<b>Biosonic Ultrasonic Scaler Tip Opti-Mist (25kHz)</b>	56-21114	<b>114.95</b>
<b>Biosonic Ultrasonic Scaler Tip Opti-Mist (30kHz)</b>	56-21115	<b>114.95</b>
<b>Biosonic Ultrasonic Scaler Tip Slim OptiMist #10 (25kHz)</b>	56-21116	<b>119.95</b>
<b>Biosonic Ultrasonic Scaler Tip Slim OptiMist #10 (30kHz)</b>	56-21117	<b>119.95</b>
<b>Biosonic Scaler O-Rings Pkg of 12</b>	56-21118	<b>19.95</b>

### Satin Swivell Ultrasonic Inserts (Hu-Friedy)

- Exclusive 360° rotation
- Color-coded silicone grips
- Large handle diameter reduces finger pinching
- Internal water delivery
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use



<b>10-25K, Satin Swivel Insert, each, UI25SS10</b>	56-90950	<b>147.95</b>
<b>100-25K, Satin Swivel Insert, each, UI25SS100</b>	56-90951	<b>147.95</b>
<b>1000-25K, Satin Swivel Insert, each UI25SS1000</b>	56-90952	<b>147.95</b>
<b>10-30K, Satin Swivel Insert, each UI30SS10</b>	56-90947	<b>147.95</b>
<b>100-30K, Satin Swivel Insert, each UI30SS100</b>	56-90948	<b>147.95</b>
<b>1000-30K, Satin Swivel Insert, each, UI30SS1000</b>	56-90949	<b>147.95</b>

## Ultrasonic Scalers

### Sonix Air Scaler with 3 Tips (Dentamerica)



- 6000 Hz frequency
- Fits onto any highspeed handpiece tubing
- Lightweight, portable, and space-saving
- Autoclavable/Chemoclavable
- Available in 2 or 4 hole configuration

<b>Sonix Air Scaler, 4 Holes with 3 Tips</b>	56-11291	<b>194.95</b>
<b>Sonix Air Scaler, 2 Holes with 3 Tips</b>	56-11290	<b>194.95</b>
<b>Air Scaler Tips (Universal)</b>	56-63631	<b>24.95</b>
<b>Air Scaler Tips (Perio)</b>	56-63632	<b>24.95</b>
<b>Air Scaler Tips (Sickle)</b>	56-63633	<b>24.95</b>

### Streamline Ultrasonic Inserts (Hu-Friedy)

- Internal water delivery
- The #3 Beavertail design is excellent for removal of gross calculus and orthodontic cement



- Use #10 Universal for the removal of all types of debris
- Use #100 for removing light deposits
- Use #1000 for removing tenacious stains and calculus
- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use.

<b>Lavender, 25 kHz, #10 Universal</b>	56-90953	<b>123.95</b>
<b>Black, 25 kHz, #100 Thin</b>	56-90954	<b>123.95</b>
<b>Orange, 25 kHz, #1000 Triple Bend</b>	56-90955	<b>123.95</b>
<b>Yellow, 25 kHz, #3 Beavertail</b>	56-90956	<b>123.95</b>
<b>Lavender, 30 kHz, #10 Universal</b>	56-90957	<b>123.95</b>
<b>Black, 30 kHz, #100 Thin</b>	56-90958	<b>123.95</b>
<b>Orange, 30 kHz, #1000 Triple Bend</b>	56-90959	<b>123.95</b>
<b>Yellow, 30 kHz, #3 Beavertail</b>	56-90960	<b>123.95</b>



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

## Ultrasonic Scalers

### Swivel Direct Flow Ultrasonic Inserts (Hu-Friedy)

- Patented 360° rotation for improved access and adaptability

- Exceptional water delivery focuses the water flow directly to the tip for improved visibility and increased patient comfort

- Lifetime warranty against defects in materials and workmanship, and guaranteed against breakage, joint failure, and corrosion under normal use



Univ #10, Swivel, 25K, each	56-90931	147.95
Univ #10, Swivel, 30K, each	56-90932	147.95
Thin #100, Swivel, 25K, each	56-90933	147.95
Thin #100, Swivel, 30K, each	56-90934	147.95
After Five, Left, 25K, each	56-90935	147.95
After Five, Straight, 25K, each	56-90936	147.95
After Five Set, 25K, 3/Pkg	56-90937	374.95
After Five, Left, 30K, each	56-90938	147.95
After Five, Right, 30K, each	56-90939	147.95
After Five, Straight, 30K, each	56-90940	147.95
After Five 30K, 3/Pkg	56-90941	374.95
#1000 Triple Bend, 25K, each	56-90942	147.95
#1000 Triple Bend, 30K, each	56-90943	147.95

### Endo Chuck Peizo Tips (ATW)



E1 120 Degree Angle, each (Endo File Not Included)	29-11294	39.99
E2 90 Degree Angle, each (Endo File Not Included)	29-11295	39.99

### Ultrasonic Scaler Inserts (Bonart)



SPA Series Insert  
Internal water flow



P-10 Series Insert  
Universal external water flow




P-100L Insert  
Left angle slim external water flow



P-100R Insert  
Right angle slim external water flow



P-100 Series Insert  
Straight angle slim external water flow



IF-50 Streamline Insert  
Universal internal water flow



IF-5 Beavertail Insert  
Beavertail internal water flow



IF-100 Series Insert  
Slim internal water flow

### SPA Series– Internal water flow insert:

Bonart Medical Technology, Inc. holds the patent to the SPA series of internal water flow inserts. The SPA inserts spray water close to the tip utilizing a silicone sleeve that pushes the water from the back of the tip to the front. This feature allows for a more focused water spray. **Autoclavable, Made in USA**

### P-10 Series– Universal external water flow insert

The Bonart universal external water flow insert series (P-10) is a straight angled universal insert with an external water tubing that sits above the tip. The P-10 series is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**

### P-100L - Left angle slim external water flow insert

The Bonart P-100L insert is a left-angled insert with a slim tip. Water flows from an external water tube that sits above the tip. The P-100L insert is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**





Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

### P-100R– Right angle slim external water flow inserts

The Bonart P-100R insert is a right-angled insert with a slim tip. Water flows from an external water tube that sits above the tip. The P-100L insert is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**

### P-100 Series– Straight angle slim external water flow insert

The P-100 series is a straight-angled insert with a slim tip. The slim tip is best used for perio scaling of calculus on sensitive areas of the mouth. Water flows from an external water tube that sits above the tip. The P-100 insert is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**

### IF-50 Streamline Insert

The IF-50 Streamline #10 ultrasonic insert has a universal sized tip. Water flow comes from the back of the tip. This tip is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**

### IF-5 Beavertail Insert

The IF-5 Beavertail insert has a spatula-style tip for removing hidden calculus more easily. Water flows from the back of the tip. This tip is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**

### IF-100 Slim Series Insert

The FS slim series insert has a slim tip for better perio scaling of calculus on sensitive areas of the mouth. Water flow comes from the back of the tip. This tip is designed to be compatible with most major brand magnet scalars. **Autoclavable, Made in USA**

Internal water flow insert 30KH, each	56-11306	49.95
Internal water flow insert 25KH, each	56-11305	49.95
Universal External water flow insert, 30KH each	56-11304	43.95
Universal External water flow insert, 25KH each	56-11303	43.95
Left angle slim external water flow insert 30KH, each	56-11302	49.95
Left angle slim external water flow insert 25KH, each	56-11301	49.95
Right angle slim external water flow insert 30KH, each	56-11300	49.95
Right angle slim external water flow insert 25KH, each	56-11299	49.95
Straight angle slim external water flow insert 30KH each	56-11298	43.95
Straight angle slim external water flow insert 25KH each	56-11295	43.95
Universal internal water flow insert 30KH, each	56-11297	43.95
Universal internal water flow insert 25KH, each	56-11294	43.95

## Ultrasonic Scalers

### Slim internal water flow insert

30KH, each 56-11293 **43.95**

### Slim internal water flow insert

25KH, each 56-11292 **43.95**

### Beavertail internal water flow insert

30KH, each 56-11296 **43.95**

### Beavertail internal water flow insert

25KH, each 56-11307 **43.95**

### Ultrasonic Scaler Inserts (Dentamerica)

#### ※ULTRASONIC INSERTS

Reorder No.80005

Universal Insert 25K  
(External Water)



Reorder No.83005

Universal Insert 30K  
(External Water)



Reorder No.80006

Universal Insert 25K  
(Internal Water)



Reorder No.83006

Universal Insert 30K  
(Internal Water)



Reorder No.80007

Chisel Type Insert 25K  
(External Water)



Reorder No.83007

Chisel Type Insert 30K  
(External Water)



Reorder No.80008

Chisel Type Insert 25K  
(Internal Water)



Reorder No.83008

Chisel Type Insert 30K  
(Internal Water)



Reorder No.80011

Endodontic Insert 25K  
(Canalsonic Tip)



#### SPECIFICATIONS :

- \*INPUT : 115V 60Hz or 230V 50Hz 3A
- \*OUTPUT : AC 24V 2.4A
- \*FREQUENCY : 25,000 or 30,000 Hz/sec.
- \*DIMENSION : W140xD140xH45mm
- \*WEIGHT : 1000g

F4 External Water, 25K Universal Insert 56-11313 **69.90**

F1 Internal Water, 25K Universal Insert 56-11312 **49.95**

External Water, 25K Chisel Type Insert 56-11311 **69.95**


Internal Water, 25K Chisel Type Insert 56-11310 **49.95**

F6 External Water, 30K Universal Insert 56-11309 **69.95**

F3 Internal Water, 30K Universal Insert 56-11308 **49.95**

### Endo Chuck Ultrasonic Insert (Bonart)

The Bonart endo chuck insert is designed to support endodontic procedures. The end of the tip has a slot to hold a



U-file. U-files are not included with the insert and must be purchased separately. Compatible with the ART-M4 and ART-M3II. **Autoclavable Materials Made in USA**

Operating Frequency: 25KHZ only

90 degree Endo Chuck Insert,

25KH each 56-11129 **59.95**



Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888


## Ultrasonic Scalers

### Piezo Ultrasonic Periodontic Tips (ATW)



Each 56-11133 **39.99**

### Piezo Ultrasonic Scaler Tips (ATW)



G1 Scailing tips, each 56-11135 **39.99**

G2 Scailing tips, each 56-11134 **39.99**

G4 Scailing tips, each 56-11136 **39.99**

### Piezo Bone Surgery Tips (Vista)

- Acteon\* Style
- Mectron\* Style



BT1 (Acteon Style) 56-11160 **112.95**

BT3 (Acteon Style) 56-11161 **112.95**

ST1 (Acteon Style) 56-11159 **112.95**

ST2 (Acteon Style) 56-11158 **112.95**

ST3 (Acteon Style) 56-11157 **112.95**

ST4 (Acteon Style) 56-11156 **112.95**

LT1 (Acteon Style) 56-11155 **112.95**

BT1 (Mectron Style) 56-11154 **112.95**

BT3 (Mectron Style) 56-11153 **112.95**

ST1 (Mectron Style) 56-11152 **112.95**

ST2 (Mectron Style) 56-11151 **112.95**

ST3 (Mectron Style) 56-11150 **112.95**

ST4 (Mectron Style) 56-11149 **112.95**

LT1 (Mectron Style) 56-11148 **112.95**

### Piezo Endodontic Tips (Vista)

**GE20 20mm Tip from Bend Angle** Same geometry as the GE20D without the diamond coating. This tip reacts less aggressively and is not recommended for cements.

**GE20D Diamond Coated 20mm Tip from Bend Angle** Used in a variety of endodontic procedures, including the

removal of dental pulp, re-treatment, retrieval of separated instruments, identifying canal openings and the removal of cements. Vista Dental recommends usage of irrigation as a cooling agent when removing crowns and posts to avoid injury.



#### 4 Acteon Style Tips

(GE20, GE20D, GEBD-2 & GEFD-2) 56-11174 **269.95**

**Piezo (Acteon, Vista P5, Tulsa, Spartan, NSK, Hu - Friedy) Style Tips (GE20) each**

56-11175 **84.95**

**Piezo (Acteon, Vista P5, Tulsa, Spartan, NSK, Hu - Friedy) Style Tips (GE20D) each**

56-11176 **84.95**

**Piezo Acteon Style Tips (GEBD-2) each**

56-11177 **84.95**

**Piezo Acteon Style Tips (GEFD-2) each**

56-11178 **84.95**

#### 4 EMS Style Tips

(GE20, GE20D, GEBD-2 & GEFD-2) 56-11179 **269.95**

**Piezo EMS Style Tips (GE20) each**

56-11180 **84.95**

**Piezo EMS Style Tips (GE20D) each**

56-11181 **84.95**

**Piezo EMS Style Tips (GEBD-2) each**

56-11182 **84.95**

**Piezo EMS Style Tips (GEFD-2) each**

56-11183 **84.95**

**Piezo (Amdent\* (US30), Piezo Advantage) Style Tips (GE20) each**

56-12183 **84.95**


**Piezo (Amdent (US30), Piezo Advantage) Style Tips (GE20D) each**

56-12184 **84.95**

### Piezo Hygiene Tips (Vista)

**SU1 Light to Medium - Heavy Cases of Supra-Gingival Calculus,**

Utilizing the side radius of the initial 2-3mm of the top of the tip in a sweeping motion on the tooth enamel provides a smooth scaled surface with minimal need for secondary hand instrumentation.



**SU6 Designed for Subgingival Scaling** Long and thin design allows for a less evasive entry into subgingival pockets resulting in minimal patient sensitivity and reduced bleeding



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**HC Universal Curette Design for Sub and Supragingival Scaling** Designed for deep pocket access and powerful enough for medium to heavy calculus removal. Use the edge of the tip to guide the motion for smooth end results.

**PTU Subgingival Pocket Scaling** With a design to simulate the natural curvature of a tooth, this universal tip provides for easy access into pockets with minimal patient discomfort. Extremely thin design helps address patient sensitivity and allows tip to create a smooth scaled surface.

**4 Tips (SU1, SU6, HC & PTU) 56-11185 269.95**


### Piezo Scaler Tips (Dentamerica)

※PIEZO SCALER TIPS



<b>Endo Chuck, each</b>	56-11214	<b>69.95</b>
<b>Piezo Scaler Tip (Universal A) each</b>	56-11213	<b>24.95</b>
<b>Piezo Scaler Tip (Universal B) each</b>	56-11212	<b>24.95</b>
<b>Piezo Scaler Tip (Chisel) each</b>	56-11211	<b>24.95</b>
<b>Piezo Scaler Tip (Perio) each</b>	56-11210	<b>49.95</b>

### Piezo Ultrasonic Tip Kits (Vista)



Tip kits offer great value and provide the clinician with the tips most likely used for a specific procedure in one complete package in both Acteon\* Style and EMS\* Style. Vista Dental's piezo tips kits include:

## Ultrasonic Scalers

- General Scaling (SU1, SU2, SU3, SU4, SU5, SU6)
- Apical Surgery (S70D, S90D, SPL, SPR, SA, SAD)
- Gracey Mini Curette (HC-1, HC-2, HC-3, HC-4, HC-7, HC-8)
- Periodontal (DST, HC, PTU, PTL, PTR, SU6)
- Retreatment (GE20, GE20D, GE40D, PRT, CT4, VPro™)
- Restorative (GEBD-1, GEBD-2, GEFD-1, GEFD-2)
- Bone Surgery (BT1, BT3, ST1, ST2, ST3, ST4, LT1)
- Hygienist Kit: (DST, HC, PTU, PTL, PTR, SU6)

**EMS\* Style:** Mini Piezon\*, Scalex 880\*, Bonart\*, Super-Son EIE\*, Master 400\*, Parkell (Piezo)\*, Dent America\*

**Amdent\* Style:** Piezo +

**Satelec\* Style:** Vista P5\*, NSK\*, ENAC\*, Sybron Endo\*, Obtura Spartan\*, ADEC\*, J.Morita\*, Dentsply/Tulsa\*, Amdent\*

**Cavity Prep Kit (Acteon Style) 56-11264 399.95**

**General Scaling Tip Kit:** SU1, SU2, SU3, SU4, SU5, SU6 one each Plus 6 Torque Wrenches (Acteon Style) 56-11263 **369.95**

**Perio Scaling Kit (Acteon Style) 56-11262 399.95**

**Hygienist Tip Kits:** DST, HC, PTU, PTL, PTR, SU6 one each Plus 6 Torque Wrenches (Acteon Style) 56-11261 **429.95**

**Retreatment Tip Kit:** GE20, GE20D, GE40D, PRT, CT4, VPro™ one each Plus 2 Wrenches (Acteon Style) 56-11260 **399.95**

**Restorative Tip Kit:** GEBD-1, GEBD-2, GEFD-1, GEFD-2 one each Plus 2 Wrenches (Acteon Style) 56-11259 **239.95**

**Apical Surgery Tip Kit:** S70D, S90D, SPL, SPR, SA, SAD one each (Acteon Style) 56-11258 **399.95**

**Bone Surgery Tip Kit:** BT1, BT3, ST1, ST2, ST3, ST4, LT1 one each (Acteon Style) 56-11257 **759.95**

**Gracey Mini Curette Tip Kit:** HC-1, HC-2, HC-3, HC-4, HC-7, HC-8 one each (Acteon Style) 56-11256 **469.95**

**General Scaling Tip Kit:** SU1, SU2, SU3, SU4, SU5, SU6 one each (Amdent Style) 56-11255 **399.95**

**Perio Scaling Kit (Amdent Style) 56-11254 399.95**

**Cavity Prep Kit (EMS Style) 56-11253 399.95**

**Bone Surgery Tip Kit:** BT1, BT3, ST1, ST2, ST3, ST4, LT1 one each (EMS Style) 56-11252 **759.95**

**Gracey Mini Curette Tip Kit:** HC-1, HC-2, HC-3, HC-4, HC-7, HC-8 one each (EMS Style) 56-11251 **469.95**



# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Bone Surgery Tip Kit:

BT1, BT3, ST1, ST2, ST3, ST4, LT1  
one each (Mectron Style) 56-11250 **759.95**

**General Scaling Tip Kit:** SU1, SU2,  
SU3, SU4, SU5, SU6 one each Plus 6  
Torque Wrenches (EMS Style) 56-11249 **399.95**

**Perio Scaling Kit (EMS Style)** 56-11248 **399.95**

**Hygienist Tip Kit:** DST, HC, PTU,  
PTL, PTR, SU6 one each Plus 6  
Torque Wrenches (EMS Style) 56-11247 **429.95**

**Retreatment Tip Kit:** GE20, GE20D,  
GE40D, PRT, CT4, VPro™ one each  
Plus 2 Wrenches (EMS Style) 56-11246 **399.95**

**Restorative Tip Kit:** GEBD-1, GEBD-  
2, GEFD-1, GEFD-2 one each Plus 2  
Wrenches (EMS Style) 56-11245 **224.95**

**Apical Surgery Tip Kit:** S70D, S90D,  
SPL, SPR, SA, SAD one each (EMS  
Style) 56-11244 **399.95**


**Bone Surgery Tip Kit:** BT1, BT3,  
ST1, ST2, ST3, ST4, LT1 one each  
(Piezotome Style) 56-11265 **679.95**

## Ultrasonic Piezo Tips (Vista)

Hygiene/Perio

### SU1 Light to Me- dium - Heavy Cases of Supra- Gingival Calculus

Utilizing the side ra-  
dius of the initial 2-  
3mm of the top of the  
tip in a sweeping mo-  
tion on the tooth  
enamel provides a  
smooth scaled surface  
with minimal need for  
secondary hand instru-  
mentation.



### SU2 Scaling of High Density Supra-Gingival Calculus Deposits

Utilizing the flat edge of the top of the tip and the side of the tip  
perpendicular to the tooth surface, use an up-and-down motion  
on the tooth surface to shear off calculus.

### SU3 Removal of Stains (Tobacco, Coffee & Other Simi- lar Stains)

Utilizing both the round edge of the top of the tip and the side of  
the tip on the tooth surface, use an up and down motion on the  
tooth surface to shear off stains.

## SU4 General Scaling of Sensitive Areas

Less abrasive than the SU1 - SU3 tips, this tip is ideal for use in  
the periodontal pockets and with other sensitive locations.

## SU5 Designed for Easy Interproximal Scaling

The thin design of the SU5 tip allows for entry into most inter-  
proximal spaces.

## SU6 Designed for Subgingival Scaling

Long and thin design allows for a less evasive entry into subgingi-  
val pockets resulting in minimal patient sensitivity and reduced  
bleeding.

**EMS\* Style:** Mini Piezon\*, Scalex 880\*, Bonart\*, Super-  
Son EIE\*, Master 400\*, Parkell (Piezo)\*, Dent America\*

**Amdent\* Style:** Piezo +

**Satelec\* Style:** Vista P5\*, NSK\*, ENAC\*, Sybron Endo\*,  
Obtura Spartan\*, ADEC\*, J.Morita\*, Dentsply/Tulsa\*, Am-  
dent\*

## SU1 - Advantage / Amdent™

Thread Style, each 56-11328 **86.95**

SU1 - P5 Thread Style, each 56-11327 **86.95**

SU1 - EMS™ Thread Style, each 56-11325 **86.95**

SU2 - P5 Thread Style, each 56-11323 **86.95**

## SU2 - Advantage / Amdent™

Thread Style, each 56-11324 **86.95**

SU2 - EMS™ Thread Style, each 56-11322 **86.95**

SU3 - P5 Thread Style, each 56-11321 **86.95**

## SU3 - Advantage / Amdent™

Thread Style, each 56-11320 **86.95**

SU3 - EMS™ Thread Style, each 56-11319 **86.95**

SU4 - P5 Thread Style, each 56-11318 **86.95**

## SU4 - Advantage / Amdent™

Thread Style, each 56-11317 **86.95**

SU4 - EMS™ Thread Style, each 56-11316 **86.95**

SU5 - P5 Thread Style, each 56-11315 **86.95**

## SU5 - Advantage / Amdent™

Thread Style, each 56-11314 **86.95**

SU5 - EMS™ Thread Style, each 56-11329 **86.95**

SU6 - P5 Thread Style, each 56-11330 **86.95**

## SU6 - Advantage / Amdent™

Thread Style, each 56-11331 **86.95**

SU6 - EMS™ Thread Style, each 56-11332 **86.95**




Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

## Ultrasonic Scalers


### Piezo Tips (Bonart)

**ART Piezo Tip: BE-1 / BS-1:** is a pointed universal tip best used for general deposit removal in all quadrants.

**ART Piezo Tip: BE-2 / BS-2:** is a flat crescent tip best used for removing resin-free ortho cements. This tip may also be used for supragingival calculus removal.



ART Piezo Tip: BE-1 / BS-1  
Universal pointed tip



ART Piezo Tip: BE-2 / BS-2  
Flat chisel tip



ART Piezo Tip: BE-3 / BS-3  
Rounded flat tip

**ART Piezo Tip: BE-3 / BS-3:** is a rounded flat tip. The rounded edge is designed for gentle treatment.

**ART Piezo Tip: BS-PFR:** is a right-angled tip designed for perio procedures and reaching difficult-to-reach spots near the posterior teeth.



ART Piezo Tip: BS-PFR  
Right-angled perio tip



ART Piezo Tip: BS-P  
Slim perio tip



ART Piezo Tip: BS-PFL  
Left-angled perio tip

**ART Piezo Tip: BS-P:** is a slim perio tip designed for supra and subgingival scaling.



ART Piezo Tip: BS-PFU  
Straight-angled perio tip



ART Piezo Tip: BS-10P  
45° angled perio tip



ART Piezo Tip: BS-H2L  
Left-angled curved perio tip

**ART Piezo Tip: BS-PFL:** is a left-angled tip designed for perio procedures and reaching difficult-to-reach spots near the posterior teeth.

**ART Piezo Tip: BS-PFU:** is a straight-angled universal tip specially designed for perio procedures and reaching difficult-to-reach spots near the posterior teeth.



ART Piezo Tip: BS-H2R  
Right-angled curved perio tip



ART Piezo Tip: BS-10Z  
Flat headed perio tip

**ART Piezo Tip: BS-10P:** is a 45° angled tip recommended for procedures on shallow pockets.

**ART Piezo Tip: BS-H2L:** is a left-angled curved tip specially designed for perio functions and reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for treatment on premolars, molars, furcation, and narrow spaces.

# Ultrasonic Scalers



Prestige Dental Products, Inc.  
Prestigedentalproducts.com  
Tel: 877-772-3888

**ART Piezo Tip: BS-H2R:** is a right-angled curved tip specially designed for perio functions and reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for treatment on premolars, molars, furcation, and narrow spaces.

**ART Piezo Tip: BS-10Z:** is a wide straight-angled tip specially designed for perio functions and reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for work on medium pockets less than 4mm.

**BE tips** are designed to be **EMS** compatible like **ART-P3, ART-SP2**

**BS tips** are designed to be compatible with **Satelec Acteon** like **ART-P3II, ART-P4, ART-P6**

<b>BE-1 Universal Pointed Piezo Tip for EMS type, each</b>	56-11230	<b>39.95</b>
<b>BS-1 Universal Pointed Piezo Tip for Satelec Acteon type, each</b>	56-11231	<b>39.95</b>
<b>BE-2 Flat chisel Piezo Tip for EMS type, each</b>	56-11223	<b>39.95</b>
<b>BS-2 Flat chisel Piezo Tip for Satelec Acteon type, each</b>	56-11229	<b>39.95</b>
<b>BE-3 Rounded Flat Piezo Tip for EMS type, each</b>	56-11228	<b>39.95</b>
<b>BS-3 Rounded Flat Piezo Tip for Satelec Acteon type, each</b>	56-11227	<b>39.95</b>
<b>BS-PFR Right-Angled perio Tip for Satelec Acteon type, each</b>	56-11226	<b>39.95</b>
<b>BS-P Slim Perio Tip for Satelec Acteon type, each</b>	56-11225	<b>39.95</b>
<b>BS-PFL Left-Angled Perio Tip for Satelec Acteon type, each</b>	56-11224	<b>39.95</b>
<b>BS-PFU Straight-Angled Perio Tip for Satelec Acteon type, each</b>	56-11222	<b>39.95</b>
<b>BS-10P 45° Angled Perio Tip for Satelec Acteon type, each</b>	56-11221	<b>39.95</b>
<b>BS-H2L Left-Angled Curved Perio Tip for Satelec Acteon type, each</b>	56-11220	<b>39.95</b>
<b>BS-H2R Right-Angled Curved Perio Tip for Satelec Acteon type, each</b>	56-11219	<b>39.95</b>
<b>BS-10Z Flat Headed Perio Tip for Satelec Acteon type, each</b>	56-11218	<b>39.95</b>

## Piezo Tip Wrenches (Vista)



## Tip Wrenches - 3 Styles

- Torque - Acteon\* Style (Frosted Beige)
- Torque - EMS\* Style (Blue)
- Universal - Acteon\* EMS\*, Amdent\* Style

<b>Universal - Acteon, EMS, Amdent</b>	56-11216	<b>9.95</b>
<b>Torque - Acteon (Frosted Beige)</b>	56-11215	<b>18.95</b>
<b>Torque - EMS (Blue)</b>	56-11217	<b>18.95</b>

## Sonic Air Scalers (MTI)



### ScaleMaster

- The perfect choice for air-driven scalers.
- Gentle, yet powerful to remove stains and calculus.
- Attach to any 2-or 4-hole connections.

### Black color With:

4 Holes & Three differently shaped scaling tips and wrench	56-89811	<b>279.95</b>
--	----------	---------------

### Blue color With:

4 Holes & Three differently shaped scaling tips and wrench	56-89812	<b>279.95</b>
--	----------	---------------

### Black color With:

2 Holes & Three differently shaped scaling tips and wrench	56-89815	<b>279.95</b>
--	----------	---------------

### Blue color With:

2 Holes & Three differently shaped scaling tips and wrench	56-89816	<b>279.95</b>
--	----------	---------------




Prestige Dental Products, Inc.  
 Prestigedentalproducts.com  
 Tel: 877-772-3888

### Sonic Scalers Tips (MTI)



<b>C-Perio Tip</b>	56-89821	<b>14.95</b>
<b>B-Sickle Tip</b>	56-89822	<b>14.95</b>
<b>A-Universal Tip</b>	56-89823	<b>14.95</b>
<b>Sonic Scaler Tip Wrench</b>	56-89824	<b>12.95</b>

### Star Titan Scalers Type Tips (MTI)



<b>Perio Scaling Tips</b>	56-89831	<b>49.95</b>
<b>Sickle Scaling Tips</b>	56-89832	<b>49.95</b>
<b>Universal Scaling Tips</b>	56-89833	<b>49.95</b>
<b>Star Type Tip Wrench</b>	56-89834	<b>12.95</b>

### Foot Control Assembly Cavitron Select SPS BOOST Type (Chapman-Huffman)



- A footswitch assembly heavy enough to stay where you put it, no chasing it around with your foot
- Same light foot pressure to operate as original equipment switch
- Round pin male plug as illustrated with 8 foot cable
- Two position switch
- Fits DENTSPLY CAVITRON Select SPS BOOST
- DENTSPLY cross ref. # 565162003

<b>Fits Dentsply Cavitron Select SPS Boost</b>	56-88091	<b>195.95</b>
--	----------	---------------

## Ultrasonic Scalers

### Perio Flow Ultrasonic Adapter (Vista)



- Transforms Piezo ultrasonic handpiece into a versatile subgingival pocket irrigator!
- Attaches any style of luer-lock irrigation tip to your piezo ultrasonic scaler
- Autoclavable
- Available in three styles to fit all Vista Piezo Ultrasonic units, Satelec and Tulsa P5 units, EMS and Amdent™ Piezo ultrasonic systems
- Great for subgingival irrigation after scaling

<b>Perio-Flow™ - P5 / Satelec Style Adapter</b>	56-11138	<b>99.95</b>
---	----------	--------------

<b>Perio-Flow™ - Piezo Advantage™ / Amdent Adapter</b>	56-11137	<b>99.95</b>
--	----------	--------------

<b>Perio-Flow™ - EMS Style Adapter</b>	56-11139	<b>99.95</b>
--	----------	--------------

### O-Ring for Cavitron Insert Replacement (Plasdent)



<b>Black For Cavitron Insert, 12/Pkg</b>	56-11132	<b>11.23</b>
--	----------	--------------

<b>Green For Cavitron Insert, 12/Pkg</b>	56-11131	<b>11.23</b>
--	----------	--------------

### O-Ring Ultrasonic Insert (Hu-Friedy)

The O-Ring acts as a water seal when the insert is placed into the open end of the handpiece.



<b>30K, 12/Pkg</b>	56-90961	<b>11.95</b>
--------------------	----------	--------------





**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

