

Vacuum Forming



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.

Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.


**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**



Vacuum Forming

Vacuum Forming Machines

Air Vac XQ Vacuum Former (Keystone)



The Keystone Air Vac XQ Dental Vacuum Forming Machine is a work horse Dental vacuum forming machine. It's price, features and reliability make it one of the best dental vacuum formers on the market. The Air Vac XQ has greater accuracy and definition when processing bleaching trays, sports mouth guards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.

- The first and only patented 'Xtra Quiet', American made vacuum forming machine in the industry
- Utilizes compressed air to create all types of dental appliances in less time, with precise adaptation and superior accuracy
- **5 year warranty**
- UL, CUL, and CE rated

Air Vac XQ Vacuum Former 110V 57-11111 **459.95**

Air Vac XQ Vacuum Former 220V 57-11112 **479.95**

Sta-Vac Premium Electric Vacuum Forming Machine (Buffalo)

Features/Benefits:

Versatile New ConvertTable Combination Vacuum Table and Bead Basket unique new table design can be used in quick and easy flat table mode or flipped over and used in the versatile bead basket mode for infinite model positioning and improved vacuum flow

High Performance Electric Vacuum Forming System Dual-Stage electric vacuum motor provides powerful suction and highly accurate adaptations—no compressed air connections required.

Accurate Frame and Alignment System spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table

Heavy Duty Construction Solid Construction (7"W x 16" D x 11" H, 32 lbs) for maximum durability and long life

Convenient and Safe Design Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation



Sta-Vac Vacuum Former, 120V
AC, Low Profile, High Powered 57-11469 **539.95**

Sta-Vac Vacuum Former, 220V
AC, Low Profile, High Powered (CE) 57-11470 **549.95**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Sta-Vac II Premium Electric Vacuum Forming Machine (Buffalo)

Features/Benefits:

Versatile New ConvertTable Combination Vacuum Table and Bead Basket unique new table design can be used in quick and easy flat table mode or flipped over and used in the versatile bead basket mode for infinite model positioning and improved vacuum flow

High Performance Electric Vacuum Forming System Dual-Stage electric vacuum motor provides powerful suction and highly accurate adaptations—no compressed air connections required.

Accurate Frame and Alignment System spring loaded, self-leveling, gasket-sealed frame handles a variety of 5" x 5" sheet materials, and the square center post assures smooth operation and guarantees alignment of frame and table

Heavy Duty Construction Solid Construction (7"W x 16" D x 11" H, 32 lbs) for maximum durability and long life

Convenient and Safe Design Unique housing design prevents motor contamination, shields heater, and promotes ergonomic operation



Sta-Vac II Vacuum Former, 120V

AC, w/2 ConvertTable Tables 57-11467 **589.95**

Sta-Vac II Vacuum Former, 220V

AC, w/2 ConvertTable Tables 57-11466 **599.95**

Replacement ConvertTable Table for

Sta-Vac II Vacuum Former 57-11465 **37.95**

Replacement Steel Beads for Sta-Vac

II (and Accu-Vac), 3 lb. Pkg. 57-11468 **79.95**

Econo-Vac Vacuum Forming Machine Heavy Duty (Buffalo)

Dual Stage vacuum motor
Square center post
Self-leveling frame
Heavy-duty construction
Frame (7 1/2" wide x 10 1/2" deep x 15" high, 32 lbs)



Econo-Vac Vacuum Former,
120V AC, High Powered 57-11332 **398.95**

Econo-Vac Vacuum Former,
220V AC, High Powered (CE) 57-11333 **449.95**



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

The Machine III Vacuum Former (Keystone)

This is the most popular Dental Vacuum Forming machine in the dental industry. Its price, features and reliability make it the first choice in dental vacuum formers. More power means greater accuracy and definition when processing bleaching trays, sports mouth guards, implant stints, temporary splints, custom impression trays, denture bases, copings, temporary crowns and bridges and orthodontic appliances.



- Uses 5" x 5" resin sheets
- 22 lbs

The Machine III 110V 57-11471 **369.95**

The Machine III 220V CE Euro Cord 57-11472 **389.95**

Tray-Vac Vacuum Forming Machine (Buffalo)



Excellent Performance Heavy duty electric vacuum motor provides powerful and immediate suction, creating accurate adaptations, no compressed air connections required.

Gasket-Free Easy Frame System spring loaded, self-leveling, gasket-free frame with clever quick-release "fingernail-saving" latch securely handles a variety of 5" x 5" sheet materials and maintains excellent alignment with the flat vacuum table

Solid Construction Solid construction (7 1/2" W x 9 1/2" D x 15" H, 17lbs) for improved durability, long life and virtually no maintenance in a lighter-weight design

Convenient and Safe Design heater hood rotates behind the center post for safety and convenience

Wide Range of Available Accessories Please see the complete line of Buffalo Dental vacuum forming equipment and supplies including a wide range of Vacuum Forming Sheets, Block-out Compound, Silicone Spray, Tray Handles, Knives and Burs, the ThermaKnife Thermal Cutting Knife, and PolyShoK Sports Mouth guard Materials.

Tray-Vac Vacuum Former, 120V AC, Light Weight, Small, Powerful 57-11479 **349.95**

Tray-Vac Vacuum Former, 220V AC, Light Weight, Small Powerful (CE) 57-11472 **389.95**

Vacuum Forming

Vacuum Forming Materials

Base Plate Material (Keystone)

- Three thicknesses (.060, .080 & .100)
- Pure, high impact styrene 5" x 5" sheets
 - Low melt temperature used for denture base or wax bite plates
 - Packages of 100 or 625 sheets



Also available in round

| | | |
|--|----------|---------------|
| .060 Base Plate 5" x 5" 100/Pkg | 57-11117 | 34.95 |
| .080 Base Plate 5" x 5" 100/Pkg | 57-11116 | 48.95 |
| .100 Base Plate 5" x 5" 100/Pkg | 57-11115 | 52.95 |
| .060 Base Plate 5" x 5" 625/Pkg | 57-11114 | 209.95 |
| .080 Base Plate 5" x 5" 625/Pkg | 57-11113 | 279.95 |
| .100 Base Plate 5" x 5" 625/Pkg | 57-11118 | 299.95 |

Coping Material (Keystone)

- Excellent for temporary bridges
- Get good Definition
 - 5" x 5" sheets
 - Three thicknesses (.020, .030 & .040)



Also available in round

| | | |
|------------------------------------|----------|---------------|
| .020 Coping 5" x 5" 50/Pkg | 57-11123 | 17.95 |
| .030 Coping 5" x 5" 50/Pkg | 57-11122 | 19.95 |
| .040 Coping 5" x 5" 50/Pkg | 57-11121 | 22.95 |
| .020 Coping 5" x 5" 625/Pkg | 57-11120 | 119.95 |
| .030 Coping 5" x 5" 625/Pkg | 57-11119 | 139.95 |
| .040 Coping 5" x 5" 625/Pkg | 57-11127 | 179.95 |

Crown & Bridge Material (Keystone)

- Thin, clear rigid material
- 5" x 5" sheets
- Acrylic does not bond to material
- For temporary bridges
- Packages of 50 or 625 sheets



Also available in round

| | | |
|--|----------|---------------|
| .020 Crown & Bridge 5" x 5" 50/Pkg | 57-11128 | 24.95 |
| .020 Crown & Bridge 5" x 5" 625/Pkg | 57-11129 | 259.95 |



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Vacuum Forming

Duraclear Sheets Round (Dentsply)

Best material for combination of strength, durability, and clarity. Unsurpassed clarity when formed. Can be cleaned with almost any method without staining or wear. Extremely durable, will not crack, bend out of shape, or break. Will bond to acrylic, making it easy to modify for other appliance uses. Specific uses include invisible retainers, minor tooth movement, temporary bridges, aligners, and bite planes.



| | | |
|---------------------------|----------|-------|
| .030" 125mm 20 sheets/box | 57-11139 | 34.95 |
| .030" 120mm 20 sheets/box | 57-11138 | 34.95 |
| .040" 125mm 20 sheets/box | 57-11137 | 34.95 |
| .040" 120mm 20 sheets/box | 57-11136 | 34.95 |
| .060" 125mm 10 sheets/box | 57-11135 | 34.95 |
| .060" 120mm 10 sheets/box | 57-11134 | 34.95 |
| .080" 125mm 10 sheets/box | 57-11133 | 34.95 |
| .080" 120mm 10 sheets/box | 57-11132 | 34.95 |

Duraclear Sheets Square (Dentsply)

Best material for combination of strength, durability, and clarity. Unsurpassed clarity when formed. Can be cleaned with almost any method without staining or wear. Extremely durable, will not crack, bend out of shape, or break. Will bond to acrylic, making it easy to modify for other appliance uses. Specific uses include invisible retainers, minor tooth movement, temporary bridges, aligners, and bite planes.



| | | |
|---------------------------|----------|-------|
| .030" 5" 20 sheets/box | 57-11147 | 34.95 |
| .030" 125mm 20 sheets/box | 57-11146 | 34.95 |
| .040" 5" 20 sheets/box | 57-11145 | 34.95 |
| .040" 125mm 20 sheets/box | 57-11144 | 34.95 |
| .060" 5" 10 sheets/box | 57-11143 | 34.95 |
| .060" 125mm 10 sheets/box | 57-11142 | 34.95 |
| .080" 5" 10 sheets/box | 57-11141 | 34.95 |
| .080" 125mm 10 sheets/box | 57-11140 | 34.95 |

Duraforce Sheets (Dentsply)

Superior conformity to the model when heated, little or no shrink. Cooling spray is not required.



Great long term wear. Strong and durable, will not crack or break. Strong enough to create divots and initiate tooth movement if needed. Can be easily cleaned without heavily staining the material. Good clarity. Not bondable to acrylic. Specific uses include final retainers, bruxism splints, temporary bridges, and minor tooth movement.

| | | |
|-------------------------------|----------|-------|
| Square .030 Box of 100 sheets | 57-11155 | 56.95 |
| Round .030 Box of 100 sheets | 57-11154 | 56.95 |
| Square .040 Box of 100 sheets | 57-11153 | 56.95 |
| Round .040 Box of 100 sheets | 57-11152 | 56.95 |
| Square .060 Box of 50 sheets | 57-11151 | 41.95 |
| Round .060 Box of 50 sheets | 57-11150 | 41.95 |
| Square .030 Bag of 10 sheets | 57-11149 | 8.95 |
| Round .030 Bag of 10 sheets | 57-11159 | 8.95 |
| Square .040 Bag of 10 sheets | 57-11156 | 8.95 |
| Round .040 Bag of 10 sheets | 57-11157 | 8.95 |
| Square .060 Bag of 5 sheets | 57-11158 | 10.95 |
| Round .060 Bag of 5 sheets | 57-11148 | 10.95 |

Duraguard Sheets Round (Dentsply)



Specially formulated EVA heats and forms easily, with great conformity to the model. Specific uses include customer mouth guards, soft splints, and whitening trays. Material forms soft and flexible, with good durability. 25 sheets/box

| | | |
|----------------------|----------|-------|
| Aqua .040" 5" | 57-11206 | 28.95 |
| Black .040" 5" | 57-11223 | 28.95 |
| Blue .040" 5" | 57-11224 | 28.95 |
| Blue Light .040" 5" | 57-11177 | 28.95 |
| Clear .040" 5" | 57-11176 | 28.95 |
| Green .040" 5" | 57-11178 | 28.95 |
| Orange .040" 5" | 57-11179 | 28.95 |
| Purple .040" 5" | 57-11180 | 28.95 |
| Red .040" 5" | 57-11182 | 28.95 |
| Silver .040" 5" | 57-11181 | 28.95 |
| White .040" 5" | 57-11183 | 28.95 |
| Yellow .040" 5" | 57-11184 | 28.95 |
| Multi Color .060" 5" | 57-11165 | 28.95 |
| Aqua .080" 5" | 57-11167 | 28.95 |
| Black .080" 5" | 57-11168 | 28.95 |
| Blue .080" 5" | 57-11169 | 28.95 |



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Vacuum Forming

| | | |
|----------------------|----------|-------|
| Blue Light .080" 5" | 57-11170 | 28.95 |
| Clear .080" 5" | 57-11172 | 28.95 |
| Green .080" 5" | 57-11171 | 28.95 |
| Orange .080" 5" | 57-11173 | 28.95 |
| Purple .080" 5" | 57-11190 | 28.95 |
| Red .080" 5" | 57-11191 | 28.95 |
| Silver .080" 5" | 57-11192 | 28.95 |
| White .080" 5" | 57-11209 | 28.95 |
| Yellow .080" 5" | 57-11211 | 28.95 |
| Multi Color .080" 5" | 57-11210 | 28.95 |
| Aqua .120" 5" | 57-11212 | 37.95 |
| Black .120" 5" | 57-11213 | 37.95 |
| Blue .120" 5" | 57-11214 | 37.95 |
| Blue Light .120" 5" | 57-11215 | 37.95 |
| Clear .120" 5" | 57-11216 | 37.95 |
| Green .120" 5" | 57-11217 | 37.95 |
| Orange .120" 5" | 57-11218 | 37.95 |
| Purple .120" 5" | 57-11219 | 37.95 |
| Red .120" 5" | 57-11220 | 37.95 |
| Silver .120" 5" | 57-11221 | 37.95 |
| White .120" 5" | 57-11222 | 37.95 |
| Yellow .120" 5" | 57-11208 | 37.95 |
| Multi Color .120" 5" | 57-11207 | 37.95 |
| Aqua .160" 5" | 57-11193 | 43.95 |
| Black .160" 5" | 57-11194 | 43.95 |
| Blue .160" 5" | 57-11195 | 43.95 |
| Blue Light .160" 5" | 57-11196 | 43.95 |
| Clear .160" 5" | 57-11197 | 43.95 |
| Green .160" 5" | 57-11199 | 43.95 |
| Orange .160" 5" | 57-11198 | 43.95 |
| Purple .160" 5" | 57-11200 | 43.95 |
| Red .160" 5" | 57-11201 | 43.95 |
| Silver .160" 5" | 57-11202 | 43.95 |
| White .160" 5" | 57-11204 | 43.95 |
| Yellow .160" 5" | 57-11203 | 43.95 |
| Multi Color .160" 5" | 57-11205 | 43.95 |


Duraguard Sheets Scented Square (Dentsply)



Specially formulated EVA heats and forms easily, with great conformity to the model. Specific uses include customer mouthguards, soft splints, and whitening trays. Material forms soft and flexible, with good durability. 25 sheets/box

| | | |
|--------------------|----------|-------|
| .040" Clear/Mint | 57-11264 | 40.95 |
| .040" Green/Mint | 57-11253 | 40.95 |
| .040" Purple/Grape | 57-11252 | 40.95 |
| .040" Red/Cherry | 57-11251 | 40.95 |
| .060" Clear/Mint | 57-11250 | 40.95 |
| .060" Green/Mint | 57-11249 | 40.95 |
| .060" Purple/Grape | 57-11248 | 40.95 |
| .060" Red/Cherry | 57-11247 | 40.95 |
| .080" Clear/ Mint | 57-11246 | 40.95 |
| .080" Green/Mint | 57-11254 | 40.95 |
| .080" Purple/Grape | 57-11255 | 40.95 |
| .080" Red/Cherry | 57-11263 | 40.95 |
| .120" Clear/Mint | 57-11262 | 45.95 |
| .120" Green/Mint | 57-11261 | 45.95 |
| .120" Purple/Grape | 57-11260 | 45.95 |
| .120" Red/Cherry | 57-11259 | 45.95 |
| .160" Clear/Mint | 57-11258 | 53.95 |
| .160" Green/Mint | 57-11257 | 53.95 |
| .160" Purple/Grape | 57-11256 | 53.95 |
| .160" Red/Cherry | 57-11245 | 53.95 |

Duraguard Sheets Scented Round (Dentsply)



Specially formulated EVA heats and forms easily, with great conformity to the model. Specific uses include customer mouthguards, soft splints, and whitening trays. Material forms soft and flexible, with good durability. 25 sheets/box

| | | |
|--------------------|----------|-------|
| .040" Clear/Mint | 57-11244 | 40.95 |
| .040" Green/ Mint | 57-11233 | 40.95 |
| .040" Purple/Grape | 57-11232 | 40.95 |
| .040" Red/Cherry | 57-11231 | 40.95 |
| .060" Clear/Mint | 57-11230 | 40.95 |
| .060" Green/Mint | 57-11229 | 40.95 |
| .060" Purple/Grape | 57-11228 | 40.95 |
| .060" Red/Cherry | 57-11227 | 40.95 |
| .080" Clear/Mint | 57-11226 | 40.95 |
| .080" Green/Mint | 57-11234 | 40.95 |
| .080" Purple/Grape | 57-11235 | 40.95 |

Vacuum Forming

| | | |
|--------------------|----------|-------|
| .080" Red/Cherry | 57-11243 | 40.95 |
| .120" Clear/Mint | 57-11242 | 45.95 |
| .120" Green/Mint | 57-11241 | 45.95 |
| .120" Purple/Grape | 57-11240 | 45.95 |
| .120" Red/Cherry | 57-11239 | 45.95 |
| .160" Clear/Mint | 57-11238 | 53.95 |
| .160" Green/Mint | 57-11237 | 53.95 |
| .160" Purple/Grape | 57-11236 | 53.95 |
| .160" Red/Cherry | 57-11225 | 53.95 |

Mouthguard Resin Sheets (Keystone)

Soft, for energy absorption
 •5"x 5" sheets, .160 thick
 •Mouthguards & Niteguards
 •Packages of 25 or 300 sheets



| | | |
|-------------------------------|----------|--------|
| .160 Clear 5" x 5" 25/Pkg | 57-11351 | 27.95 |
| .160 Red 5" x 5" 25/Pkg | 57-11349 | 27.95 |
| .160 Blue 5" x 5" 25/Pkg | 57-11350 | 27.95 |
| .160 Green 5" x 5" 25/Pkg | 57-11348 | 27.95 |
| .160 Yellow 5" x 5" 25/Pkg | 57-11347 | 27.95 |
| .160 Black 5" x 5" 25/Pkg | 57-11346 | 27.95 |
| .160 Orange 5" x 5" 25/Pkg | 57-11345 | 27.95 |
| .160 Eva White 5" x 5" 25/Pkg | 57-11344 | 27.95 |
| .160 Purple 5" x 5" 25/Pkg | 57-11343 | 27.95 |
| .160 Clear 5" x 5" 300/Pkg | 57-11342 | 229.95 |
| .160 Red 5" x 5" 300/Pkg | 57-11341 | 229.95 |
| .160 Blue 5" x 5" 300/Pkg | 57-11340 | 229.95 |
| .160 Green 5" x 5" 300/Pkg | 57-11339 | 229.95 |
| .160 Yellow 5" x 5" 300/Pkg | 57-11338 | 229.95 |
| .160 Black 5" x 5" 300/Pkg | 57-11337 | 229.95 |
| .160 Orange 5" x 5" 300/Pkg | 57-11336 | 229.95 |
| .160 White 5" x 5" 300/Pkg | 57-11335 | 229.95 |
| .160 Purple 5" x 5" 300/Pkg | 57-11334 | 229.95 |

Mouthguard Straps (Keystone)

Keystone has developed an easy-to-attach helmet strap used in the production of custom-made mouthguards. These straps have a safety release feature that allows the mouthguard to



Prestige Dental Products, Inc.

Prestigedentalproducts.com

Tel: 877-772-3888

break away under 12 pounds of pressure. Available in black and clear. Sold in packages of 12 and 25.

| | | |
|--------------|----------|-------|
| Black 12/Pkg | 57-11354 | 11.69 |
| Black 25/Pkg | 57-11353 | 17.95 |
| Clear 12/Pkg | 57-11352 | 11.69 |
| Clear 25/Pkg | 57-11355 | 17.95 |

Duraguard Sheets Square (Dentsply)

Specially formulated EVA heats and forms easily, with great conformity to the model. Specific uses include customer mouthguards, soft splints, and whitening trays. Material forms soft and flexible, with good durability.



| | | |
|-------------------------------------|----------|-------|
| Aqua .040" 5", 25 sheets/box | 57-11329 | 44.95 |
| Black .040" 5", 25 sheets/box | 57-11277 | 44.95 |
| Blue .040" 5", 25 sheets/box | 57-11279 | 44.95 |
| Blue Light .040" 5", 25 sheets/box | 57-11280 | 44.95 |
| Clear .040" 5", 25 sheets/box | 57-11282 | 44.95 |
| Green .040" 5", 25 sheets/box | 57-11281 | 44.95 |
| Orange .040" 5", 25 sheets/box | 57-11283 | 44.95 |
| Purple .040" 5", 25 sheets/box | 57-11284 | 44.95 |
| Red .040" 5", 25 sheets/box | 57-11285 | 44.95 |
| Silver .040" 5", 25 sheets/box | 57-11286 | 44.95 |
| White .040" 5", 25 sheets/box | 57-11287 | 44.95 |
| Yellow .040" 5", 25 sheets/box | 57-11288 | 44.95 |
| Multi Color .040" 5", 25 sheets/box | 57-11289 | 44.95 |
| Aqua .060" 5", 25 sheets/box | 57-11291 | 44.95 |
| Black .060" 5", 25 sheets/box | 57-11290 | 44.95 |
| Blue .060" 5", 25 sheets/box | 57-11278 | 44.95 |
| Blue Light .060" 5", 25 sheets/box | 57-11276 | 44.95 |
| Clear .060" 5", 25 sheets/box | 57-11328 | 44.95 |
| Green .060" 5", 25 sheets/box | 57-11327 | 44.95 |
| Orange .060" 5", 25 sheets/box | 57-11318 | 44.95 |
| Purple .060" 5", 25 sheets/box | 57-11266 | 44.95 |
| Red .060" 5", 25 sheets/box | 57-11265 | 44.95 |
| Silver .060" 5", 25 sheets/box | 57-11268 | 44.95 |
| White .060" 5", 25 sheets/box | 57-11267 | 44.95 |
| Yellow .060" 5", 25 sheets/box | 57-11269 | 44.95 |
| Multi Color .060" 5", 25 sheets/box | 57-11270 | 44.95 |
| Aqua .080" 5", 25 sheets/box | 57-11271 | 44.95 |
| Black .080" 5", 25 sheets/box | 57-11272 | 44.95 |
| Blue .080" 5", 25 sheets/box | 57-11273 | 44.95 |
| Blue Light .080" 5", 25 sheets/box | 57-11274 | 44.95 |
| Clear .080" 5", 25 sheets/box | 57-11275 | 44.95 |



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

| | | |
|---|----------|--------------|
| Green .080" 5" , 25 sheets/box | 57-11292 | 44.95 |
| Orange .080" 5" , 25 sheets/box | 57-11293 | 44.95 |
| Purple .080" 5" , 25 sheets/box | 57-11394 | 44.95 |
| Red .080" 5" , 25 sheets/box | 57-11311 | 44.95 |
| Silver .080" 5" , 25 sheets/box | 57-11312 | 44.95 |
| White .080" 5" , 25 sheets/box | 57-11313 | 44.95 |
| Yellow .080" 5" , 25 sheets/box | 57-11314 | 44.95 |
| Multi Color .080" 5" , 25 sheets/box | 57-11315 | 44.95 |
| Aqua .120" 5" , 25 sheets/box | 57-11316 | 49.95 |
| Black .120" 5" , 25 sheets/box | 57-11317 | 49.95 |
| Blue .120" 5" , 25 sheets/box | 57-11319 | 49.95 |
| Blue Light .120" 5" , 25 sheets/box | 57-11320 | 49.95 |
| Clear .120" 5" , 25 sheets/box | 57-11321 | 49.95 |
| Green .120" 5" , 25 sheets/box | 57-11322 | 49.95 |
| Orange .120" 5" , 25 sheets/box | 57-11323 | 49.95 |
| Purple .120" 5" , 25 sheets/box | 57-11324 | 49.95 |
| Red .120" 5" , 25 sheets/box | 57-11325 | 49.95 |
| Silver .120" 5" , 25 sheets/box | 57-11309 | 49.95 |
| White .120" 5" , 25 sheets/box | 57-11310 | 49.95 |
| Yellow .120" 5" , 25 sheets/box | 57-11295 | 49.95 |
| Multi Color .120" 5" , 25 sheets/box | 57-11296 | 49.95 |
| Aqua .160" 5" , 25 sheets/box | 57-11297 | 54.95 |
| Black .160" 5" , 25 sheets/box | 57-11298 | 54.95 |
| Blue .160" 5" , 25 sheets/box | 57-11299 | 54.95 |
| Blue Light .160" 5" , 25 sheets/box | 57-11300 | 54.95 |
| Clear .160" 5" , 25 sheets/box | 57-11301 | 54.95 |
| Green .160" 5" , 25 sheets/box | 57-11302 | 54.95 |
| Orange .160" 5" , 25 sheets/box | 57-11303 | 54.95 |
| Purple .160" 5" , 25 sheets/box | 57-11304 | 54.95 |
| Red .160" 5" , 25 sheets/box | 57-11305 | 54.95 |
| Silver .160" 5" , 25 sheets/box | 57-11306 | 54.95 |
| White .160" 5" , 25 sheets/box | 57-11307 | 54.95 |
| Yellow .160" 5" , 25 sheets/box | 57-11308 | 54.95 |
| Multi Color .160" 5" , 25 sheets/box | 57-11326 | 54.95 |

Pro-Flex & Dual Combo (Keystone)

An excellent way to make niteguards. First, vacuum form a dual laminate onto the patient's cast. This gives a hard bite plate with a built-in soft liner. Then, modify the occlusal surface by adding the semi-flexible Pro-flex NFC acrylic resin. The final result is a hard bite plate modified to the patient's needs that is slightly flexible with a soft inner liner for patient comfort. Two great products combined into one terrific final appliance.



| | | |
|--|----------|---------------|
| Dual Laminate Kit | 57-11364 | 59.95 |
| Clear NFC 160gm Powder w 4oz Liquid Kit | 57-11365 | 139.95 |

Vacuum Forming

Pro-Form Fun Mouthguards (Keystone)

Pro-form has added some new designs to its extensive line of laminated mouthguards. They are all made using the same laminating procedure to insure the same Pro-form quality.




Also available in round

| | | |
|---------------------------------------|----------|--------------|
| Blue Polkadot 5"x 5" , 12/Pk | 57-11391 | 59.95 |
| White Polkadot 5" x 5" , 12/Pk | 57-11390 | 59.95 |
| Red Polkadot 5" x 5" , 12/Pk | 57-11389 | 59.95 |
| Fang 5" x 5" , 12/Pk | 57-11388 | 59.95 |
| Blood 5" x 5" , 12/ Pk | 57-11392 | 59.95 |

Pro-Form Dual Laminates - E- Gasket (Keystone)

The Pro-form Dual Laminate is a unique combination of two materials formed into a single laminated sheet. It was made especially for use in bruxing cases. The soft liner assures patient comfort and the hard surface assures long-lasting wear. Also, acrylic will bond to the hard surface to allow the dental technician to modify the occlusal surface. Available in 3 different thicknesses (2, 3, & 5 mm). 12 laminated sheets per package.



Also available in round

| | | |
|---|----------|--------------|
| 2mm E-Gasket 5"x 5" , 12/Pk | 57-11130 | 57.95 |
| 3mm Dual Laminates 5"x 5" , | 57-11131 | 62.95 |
| 5mm Anti Snore Device 5"x 5" , 12/Pk | 57-12132 | 68.95 |

Pro-Form Niteguards (Keystone)

This laminate or layered soft/soft material was made specifically for nighttime bruxers. The layering gives the thermoplastic a more abrasion-resistant characteristic that will withstand grinding better than traditional soft sheet materials. Available in packages of 12 laminates.



Also available in round

| | | |
|-----------------------------------|----------|--------------|
| Niteguards 5" x 5" , 12/Pk | 57-11363 | 54.95 |
|-----------------------------------|----------|--------------|

Vacuum Forming



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Pro-Form Tie-Dye Mouthguards (Keystone)



Also available in round

Pro-form has added some new Tie Dye designs to its extensive line of laminated mouthguards. They are all made using the same laminating procedure to insure the same Pro-form quality.

| | | |
|----------------------------------|----------|--------------|
| Wilderness 5" x 5", 6/Pkg | 57-11360 | 32.95 |
| Midnight 5" x 5", 6/Pkg | 57-11359 | 32.95 |
| Shamrock 5" x 5", 6/Pkg | 57-11358 | 32.95 |
| Flamingo 5" x 5", 6/Pkg | 57-11357 | 32.95 |
| Rainbow 5" x 5", 6/Pkg | 57-11356 | 32.95 |
| All In 5" x 5", 6/Pkg | 57-11362 | 32.95 |
| Patriot 5" x 5", 6/Pkg | 57-11361 | 32.95 |

Pro-Form Mouthguard Laminates (Keystone)

Pro-form Mouth-guard Laminate is a strong two layer material made of high-energy absorbing EVA vinyl's. 0.150" thick in dual-color sheets. Dual-color Mouth guard Laminate in assorted two-color combos. 5" x 5", 12/Pk.



Also available in round

| | | |
|----------------------|----------|--------------|
| Blue/Green | 57-11386 | 56.95 |
| Black/White | 57-11385 | 56.95 |
| Blue/Red | 57-11384 | 56.95 |
| Black/Yellow | 57-11382 | 56.95 |
| Black/Red | 57-11383 | 56.95 |
| Purple/Yellow | 57-11381 | 56.95 |
| Maroon/Yellow | 57-11380 | 56.95 |
| Blue/Yellow | 57-11379 | 56.95 |
| Red/White | 57-11378 | 56.95 |
| Orange/White | 57-11377 | 56.95 |

| | | |
|--------------------------------------|----------|--------------|
| Blue/White | 57-11376 | 56.95 |
| Orange/Blue | 57-11375 | 56.95 |
| Green/Yellow | 57-11374 | 56.95 |
| Black/Orange | 57-11367 | 56.95 |
| Purple/White | 57-11368 | 56.95 |
| Green/White | 57-11369 | 56.95 |
| Red/Yellow | 57-11370 | 56.95 |
| Maroon/White | 57-11371 | 56.95 |
| Red/Clear | 57-11372 | 56.95 |
| Black/Blue | 57-11373 | 56.95 |
| Orange/Clear | 57-11366 | 56.95 |
| Assorted 2 Color | 57-11387 | 56.95 |
| Clear 5" x 5", 12/Pk | 57-11403 | 56.95 |
| Orange 5" x 5", 12/Pk | 57-11402 | 56.95 |
| Yellow 5" x 5", 12/Pk | 57-11401 | 56.95 |
| Maroon 5" x 5", 12/Pk | 57-11400 | 56.95 |
| White 5" x 5", 12/Pk | 57-11399 | 56.95 |
| Green 5" x 5", 12/Pk | 57-11398 | 56.95 |
| Black 5" x 5", 12/Pk | 57-11397 | 56.95 |
| Red 5" x 5", 12/Pk | 57-11396 | 56.95 |
| Blue 5" x 5", 12/Pk | 57-11395 | 56.95 |
| Purple 5" x 5", 12/Pk | 57-11394 | 56.95 |
| Pink 5" x 5", 12/Pk | 57-11394 | 56.95 |
| Assorted Solid 5" x 5", 12/Pk | 57-11404 | 56.95 |
| Red/White/Blue | 57-12388 | 56.95 |
| Canada Red/White/Red | 57-12389 | 56.95 |

Pro-Form Mouthguard Ti-color Laminates (Keystone)

Pro-form Mouth-guard Laminate is a strong two layer material made of high-energy absorbing EVA vinyls. 0.150" thick in tri-color sheets.



Also available in round

| | | |
|--|----------|--------------|
| Red/White/Black 5" x 5", 12/Pk | 57-11432 | 59.95 |
| Red/White/Blue 5" x 5", 12/Pk | 57-11431 | 59.95 |
| White/Green/Black 5" x 5", 12/Pk | 57-11430 | 59.95 |
| Red/White/Red Canada 5" x 5", 12/Pk | 57-11429 | 59.95 |
| Blue/Yellow/Blue 5" x 5", 12/Pk | 57-11428 | 59.95 |
| Red/Black/Red 5" x 5", 12/Pk | 57-11427 | 59.95 |



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

| | | |
|------------------------------------|----------|-------|
| Black/White/Green 5" x 5", 12/Pk | 57-11426 | 59.95 |
| Green/Clear/Yellow 5" x 5", 12/Pk | 57-11425 | 59.95 |
| Orange/White/Orange 5" x 5", 12/Pk | 57-11421 | 59.95 |
| Orange/Blue/Orange 5" x 5", 12/Pk | 57-11422 | 59.95 |
| Purple/Blue/Orange 5" x 5", 12/Pk | 57-11423 | 59.95 |
| Red/White/Red 5" x 5", 12/Pk | 57-11424 | 59.95 |
| Red/Yellow/Blue 5" x 5", 12/Pk | 57-11420 | 59.95 |
| Blue/White/Yellow 5" x 5", 12/Pk | 57-11419 | 59.95 |
| Black/Blue/Red 5" x 5", 12/Pk | 57-11418 | 59.95 |
| Yellow/Blue/Yellow 5" x 5", 12/Pk | 57-11417 | 59.95 |
| Orange/Black/Orange 5" x 5", 12/Pk | 57-11416 | 59.95 |
| Purple/White/Purple 5" x 5", 12/Pk | 57-11415 | 59.95 |
| Yellow/Red/Yellow 5" x 5", 12/Pk | 57-11414 | 59.95 |
| Green/Black/Green 5" x 5", 12/Pk | 57-11413 | 59.95 |
| Blue/White/Blue 5" x 5", 12/Pk | 57-11412 | 59.95 |
| Black/White/Yellow 5" x 5", 12/Pk | 57-11411 | 59.95 |
| Purple/Black/Purple 5" x 5", 12/Pk | 57-11410 | 59.95 |
| Black/Yellow/Green 5" x 5", 12/Pk | 57-11409 | 59.95 |
| Blue/Black/Yellow 5" x 5", 12/Pk | 57-11408 | 59.95 |
| Clear/Green/Clear 5" x 5", 12/Pk | 57-11407 | 59.95 |
| Clear/Red/Clear 5" x 5", 12/Pk | 57-11406 | 59.95 |
| Black/Green/Yellow 5" x 5", 12/Pk | 57-11405 | 59.95 |
| Assorted 3 Color 5" x 5", 12/Pk | 57-11433 | 59.95 |

Retainer Material (Keystone)

- Stiff, clear material
- 5" x 5", sheets
- Acrylic sticks to material
- Invisible orthodontic retainers
- Packages of 50 or 625 sheets



Also available in round

| | | |
|---|-----------|--------|
| .030 Super Clear Retainer 5" x 5", 25/Pk | 57-11435 | 44.95 |
| .030 Super Clear Retainer 5" x 5", 50/Pk | 57-11434 | 74.95 |
| .030 Super Clear Retainer 5" x 5", 625/Pk | 57-11436 | 795.95 |
| Clear Retainer 5" x 5", 25/Pk | 57-114371 | 23.95 |
| Clear Retainer 5" x 5", 625/Pk | 57-114381 | 549.95 |

Vacuum Forming

Soft EVA Sheets (Keystone)

- Soft, clear material
- 5" x 5" sheets
- Easily formed and trimmed
- Bleaching trays, soft brushing
- Custom fluoride trays, model duplication
- Packages of 25 or 300 sheets



Also available in round

| | | |
|-------------------------------|----------|--------|
| .040 Soft EVA 5" x 5", 25/Pk | 57-11443 | 18.95 |
| .060 Soft EVA 5" x 5", 25/Pk | 57-11442 | 19.95 |
| .080 Soft EVA 5" x 5", 25/Pk | 57-11441 | 21.95 |
| .120 Soft EVA 5" x 5", 25/Pk | 57-11440 | 26.95 |
| .040 Soft EVA 5" x 5", 300/Pk | 57-11439 | 129.95 |
| .060 Soft EVA 5" x 5", 300/Pk | 57-11438 | 179.95 |
| .080 Soft EVA 5" x 5", 300/Pk | 57-11437 | 199.95 |
| .120 Soft EVA 5" x 5", 300/Pk | 57-11444 | 209.95 |

Splint Material (Keystone)

- Stiff, clear material
- 5" x 5", sheets
- Acrylic sticks to material
- Packages of 25, 50 or 625 sheets
- For hard splints, thin splints, surgical trays, orthodontic retainers and bleaching trays



Also available in round


| | | |
|------------------------------|----------|--------|
| 1 .020 Splint 5" x 5", 50/Pk | 57-11464 | 16.95 |
| 3 .040 Splint 5" x 5", 50/Pk | 57-11463 | 23.95 |
| 4 .060 Splint 5" x 5", 50/Pk | 57-11461 | 28.95 |
| 5 .080 Splint 5" x 5", 50/Pk | 57-11460 | 34.95 |
| 6 .100 Splint 5" x 5", 50/Pk | 57-11459 | 39.95 |
| .020 Splint 5" x 5", 625/Pk | 57-11452 | 194.95 |
| .040 Splint 5" x 5", 625/Pk | 57-11450 | 269.95 |
| .060 Splint 5" x 5", 625/Pk | 57-11449 | 339.95 |
| .100 Splint 5" x 5", 625/Pk | 57-11447 | 439.95 |

Vacuum Forming



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Tuff Form Sheets (Dentsply)



Patented arch form indentation. Provides additional material surface area, which eliminates thinning along the incisal edge, and provides a retainer with uniform thickness throughout the appliance. Will bond to acrylic, making it useful for creating modified appliances. Superior conformity to the model when heated. Material has little or no shrink. Material forms to the contours of the model, and cooling spray is not required. A very secure fit to the model is assured. Great clarity when formed, creating a true invisible retainer for patients. Will not get brittle and crack. Material forms thin, making it comfortable for the patient, and increasing patient compliance. Specific uses include invisible retainers, minor tooth movement, temporary bridges, aligners, and bite planes.

| | | |
|--|----------|--------------|
| Clear Square .040" , Box of 50 sheets | 57-11481 | 79.95 |
| Clear Round .040" , Box of 50 sheets | 57-11482 | 79.95 |

Tray Material (Keystone)



Also available in round

- Two thickness
- Pure, high impact styrene
- 5" x 5", sheets
- Custom impression trays
- Packages of 100 or 625 sheets
- Spacer foam included

| | | |
|--|----------|---------------|
| .125 Tray Material 5" x 5" , 100/Pk | 57-11477 | 69.95 |
| .150 Tray Material 5" x 5" , 100/Pk | 57-11476 | 74.95 |
| .125 Tray Material 5" x 5" , 625/Pk | 57-11475 | 389.95 |
| .150 Tray Material 5" x 5" , 625/Pk | 57-11478 | 469.95 |

Specialty Breast Cancer Pink Mouthguard (Keystone)



Specialty breast cancer pink mouthguard

Breast Cancer Pink Mouthguard

| | | |
|------------------------|----------|--------------|
| 5" x 5" , 12/Pk | 57-11446 | 59.95 |
|------------------------|----------|--------------|

Spacer Foam (Keystone)



- 5" x 5" square sheets
- Use as a spacer material
- Makes cast easier to remove
- Packages of 25
- One size works for all materials
- Reusable

| | | |
|------------------------------------|----------|--------------|
| Foam Spacer 5" x 5" , 25/Pk | 57-11445 | 19.95 |
|------------------------------------|----------|--------------|

Thermoform Sheet Intro Kit (Dentsply-Glenroe)



Special introductory kit featuring 10 sheets each of Duraclear, Duraforce, Tuff Form, and Duralight. 40 sheets/box

| | | |
|---------------|----------|--------------|
| Square | 57-11474 | 57.95 |
| Round | 57-11473 | 57.95 |



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Vacuum Forming

Youth Mouthguard (Keystone)

- 5" x 5", square sheets, .120 thick
- Soft EVA, for energy absorption
- Great for little kids in sports
- Packages of 25 or 300 sheets



Also available in round

| | | |
|--------------------------------------|----------|---------------|
| .120 Soft EVA 5" x 5", 25/Pk | 57-11494 | 27.95 |
| .120 Red 5" x 5", 25/Pk | 57-11493 | 27.95 |
| .120 Blue 5" x 5", 25/Pk | 57-11492 | 27.95 |
| .120 Green 5" x 5", 25/Pk | 57-11491 | 27.95 |
| .120 Yellow 5" x 5", 25/Pk | 57-11490 | 27.95 |
| .120 Black 5" x 5", 25/Pk | 57-11489 | 27.95 |
| .120 White 5" x 5", 25/Pk | 57-11488 | 27.95 |
| .120 Soft EVA 5" x 5", 300/Pk | 57-11487 | 209.95 |
| .120 Soft EVA 5" x 5", 300/Pk | 57-11487 | 209.95 |
| .120 Red 5" x 5", 300/Pk | 57-11486 | 209.95 |
| .120 Blue 5" x 5", 300/Pk | 57-11485 | 209.95 |
| .120 Green 5" x 5", 300/Pk | 57-11484 | 209.95 |
| .120 Yellow 5" x 5", 300/Pk | 57-11483 | 209.95 |
| .120 Black 5" x 5", 300/Pk | 57-11495 | 209.95 |

The Machine III Vacuum Former (Keystone)



**Duraguard Sheets
 Scented Square
 (Dentsply)**

**Pro-Form Dual
 Laminates-E Gasket
 (Keystone)**



**To See More Products Visit
www.prestigedental.com**



Place Your Order NOW



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

| Item Number | Description | Quantity | Unit Price | Amount |
|-------------|-------------|----------|------------|--------|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

