

Bite Registration



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**



Bite Registration

BiteVance Bite Registration Material (ATW)

BiteVance is a vinyl Polysiloxane Bite Registration Material:

- Optimal hardness, 80+ durometer
- Thixotropic for easy application
- Non-slumping consistency
- Fluffy, mousse-like
- Excellent dimensional stability
- Improves fit & alignment of crowns & bridges
- Used for bite registration, quick die for temporary inlays & onlays, copying transfers, and primary impression with bite tray
- wont compress or flex
- 30 seconds working time, and 50 seconds setting time



Fast Set 2x 50ml Cartridges & 6 mixing tips	16-11129	19.95
Fast Set Bulk Pack 60x 50ml, no tips	16-11128	389.95
Fast Set 4x 50ml cartridges, No tips	16-11127	32.95

Blu-Mousse Bite Registration (Parkell)



Superfast 30-Second Set Split Cartridge Contains: 2 cartridges (100ml.) total, mixing tips	37-89901	34.95
Regular Classic 2-Minute Set Split Cartridge Contains: 2 HP Style cartridges, mixing tips	37-89902	34.95
Superfast 30-Second Set Cartridge for Standard Cartridge Gun Contains: 2 cartridges (100ml.) total, mixing	37-89903	34.95
Classic 2-Minute Set Hand mix Tubes Contains: 2 tubes	37-89904	34.95
Dispensing Gun-Split Cartridge	37-89905	99.95
Split Cartridge Mixing Tips w/Pink Base Pkg. of 30	37-89906	22.95



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888


Opti-Bite Bite Registration (Pac-Dent)



- Fast procedure time
- Generous work time
- Excellent handling and accuracy
- Included ribbon tips ensure easy dispensing

Fast set, 2x 50ml cartridges + 6 Tips	16-98601	22.95
Fast set Mint Flavor, 4 50ml cartridges	16-98602	34.95
Regular set, 2 x 50ml cartridges +6 Tips	16-98603	22.95
Super fast, 2 x 50ml cartridges + 6 Tips	16-98604	22.95
Teal Mixing Tips, 1:1, 6.5mm, 48/Pkg	16-98605	24.95

Thermobite Bite Registration Wafers (Dentsply)



- Specially-formulated resin provides highly accurate bite impressions of teeth and horizontal molar area
- Heat according to instructions, bite down, and trim
- Minor perforations allowed, for a correct occlusion with touching teeth surfaces
- Easy to store and transport to labs; resistant to temperature changes during transit

30 Pcs/Pkg	45-14069	27.95
-------------------	----------	--------------



Prestige Dental Products, Inc.
 Prestigedentalproducts.com
 Tel: 877-772-3888

Bite Registration

Bite Tray

Disposable Bite Tray (Winner)




Anterior Adult 35/Pkg	37-20017	6.50
Anterior Jr. 35/Pkg	37-20018	6.50
Posterior 35/Pkg	37-20019	6.50
Quadrant 35/Pkg	37-20020	6.50
Sideless 35/Pkg	37-20021	6.50
Assorted 35/Pkg	37-20022	6.50

Triple Trays (Premier)

Proven Accurate In Over 80 Million Impressions

All Premier Triple Tray designs have very thin, loose webbing which will not distort when the patient closes. Look for Premier on the handle; your assurance of a genuine Triple Tray. The superior performance of Triple Tray, combined with today's very accurate materials, provides you with an ideal impression in an easy technique. Triple Tray has been consistently rated an outstanding product by a leading independent research organization due to "least water absorption and thinnest inner occlusal wafer"



Anterior 35/Box	16-11142	42.95
Quadrant 40/Box	16-11141	42.95
Posterior 48/Box	16-11140	42.95
Sideless Posterior 48/Box	16-11143	42.95
Assorted Trays 40/Box	16-11144	42.95
3/4 Arch Partial 35/Box	16-11145	42.95
Ext Quadrant 40/Box	16-11146	42.95
Wide Body Posterior 48/Box	16-11147	42.95


Bite Registration Trays (TC Dental)

The right selection of durable, solid plastic bite trays for different impressions



Bite Trays-Anterior 35 Pcs/Box	16-11120	9.99
Bite Trays-Posterior 50 Pcs/Box	16-11119	9.99
Bite Trays-Quadrant 35 Pcs/Box	16-11116	9.99
Bite Trays-Sideless 50 Pcs/Box	16-11117	9.99
Bite Trays-Full Arch 30 Pcs/Box	16-11118	9.99

Excellent Bite Trays (Plasdent)



Excellent Anterior Trays 40 Pcs/Box	16-11126	10.95
Excellent Posterior Trays 50 Pcs/Box	16-11125	10.95
Excellent Sideless Posterior Trays 50 Pcs/Box	16-11124	10.95
Excellent Quadrant Trays 40 Pcs/Box	16-11123	10.95
Excellent Full Arch Trays 40 Pcs/Box	16-11122	10.95
Excellent Anterior, Posterior, Sideless Posterior, Quadrant 40pcs/box 10 pcs/each	16-11121	10.95

PolyBite Bite Trays

Impression, counter impression, bite registration 3 in one disposable tray



Anterior, 35/Pkg	16-11135	17.95
Posterior, 50/Pkg	16-11136	17.95
Quadrant, 35/Pkg	16-11137	17.95
Full Arch, 30/Pkg	16-11138	17.95
Sideless Posterior, 50/Pkg	16-11139	17.95

Bite Registration



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Colorbite D (Zhermack)

Accurate and stable impressions: zero thickness.

• **Changes color from green to yellow**

- Thixotropic, does not slump
- Extremely hard after setting (45 Shore-D)
- Very precise

• **Extra Fast**

• **Time in mouth less than 60 seconds**



Colorbite D, changes color from green to yellow, 2x 50ml cartridges/ 12x green mixing tips

16-11131 **33.95**

Colorbite Rock (Zhermack)

Thermochromic VPS for bite registration formulated with pigments that permit the visible control of the various working steps to minimize the margin of error



• **Changes color from purple to pink**

- Thixotropic, does not slump, mousse consistency

• **Rigid set (95 Shore-A)**

- Very precise
- Extra Fast

• **Time in mouth less than 45 seconds**

Purple to Pink To Purple

(Purple at room temperature)

12x green mixing tips 16-11132 **31.95**

Bite Registration Accessories (Zhermack)

D2 Dispenser (1:1) 16-11115 **78.95**

Replacement Plunger for Universal Dispenser 50ml 16-11114 **13.95**

D3 Multi-Dispenser

(1:1, 2:1, 10:1, 4:1) 1x Applicator / 1x 4:1 plunger/ 1x 1:1/2:1 plunger/ 1x 10:1 plunger/ 1x 1: 1/2:1 breech/ 1x 4:1/10:1 breech

16-11113 **73.95**

Large mixing tips

48x green mixing tips 16-11112 **39.95**

Intra Oral Bite tips

(for extruding a flat bead of material)

48x Clear Bite tips 16-11111 **16.95**

Occlufast CAD (Zhermack)

A-Silicone for bite registration developed to be used with all CAD CAM systems



Advantages:

- Compatible with CAD-CAM systems using optical, laser and contact reading
- Fast: 60 secs permanence in the oral cavity
- Imperceptible consistency for the patient who bites naturally
- After setting occlufast CAD can be finished with a cutter and/or bur
- Product easy to apply in the mouth, both with the bite technique and vestibular side applications

Characteristics:

Thixotropic, does not drip, hard (95 Shore A),Precise, dimensionally stable over time and to disinfection, can be re-positioned in the mouth and on the models after setting, light blue color

• **Rigid set (95 Shore-A)**

- Very precise

• Extra fast

• **Time in mouth less than 60 seconds**

Blue 2x 50ml cartridges/12x green mixing tips

16-11133 **33.95**

Occlufast Rock (Zhermack)

A-Silicone for bite registration. Maximum intercuspitation bite registration, centric, protusion or lateral registrations, gnathological evaluations, diagnostic evaluation of cranialcervicalmandibular disorders, orthodontic centric occlusion



Advantages:

- Imperceptible consistency for the patient Easy handling and bubble-free Positioning control
- Accurate and stable impressions: zero thickness
- Easy to trim for excess material removal
- Flow able: it avoids pressure and slipping during positioning
- Accurate: maximum dimensional variation of just -0.05% after 24 hours allows long-term stability
- Non-deformable: it is not affected by any thermal shocks during transport

Characteristics:

Final hardness 95 Shore-A, Fast: 60sec permanence in the oral cavity, Thixotropic purple color

Violet 2x 50ml cartridges/ 12x green mixing tips

16-11134 **33.95**



Prestige Dental Products, Inc.
Prestigedentalproducts.com
Tel: 877-772-3888

Bite Registration

You'll find Brands like:

Bite Vance ATW



And More!!!



Call us for Quarterly Promotion

Prestigedentalproducts.com



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

