

# Cameras Photography



**Prestige Dental Products, Inc**

**PrestigeDental.com**  
**877-772-3888**

# Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

**As an authorized Dealer of many brands in North America**, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

**Prestige Dental Products** offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

**Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and**

## Customer Information

### Order

**PHONE** To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

**FAX** You may FAX your orders to 714-632-8821.

**INTERNET** [www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



*Note:* New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



### Customize order for special customer

Email us to [mass@prestigedentalproducts.com](mailto:mass@prestigedentalproducts.com) a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

## Payment

**CREDIT CARD** We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

### Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

### MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

## Shipping

### GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

**EXPEDITED SHIPPING** UPS Next Day and Second Day Service for which additional shipping fees will be charged.

**SHIPPING HAZARDOUS ITEMS** All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

## State and federal laws

### DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog  
Are Subject To Change Without Prior Notice.**



# Cameras - Photography



Prestige Dental Products, Inc.  
Tel: 877-772-3888

## CU2 Digital Video USB intra oral camera (MyRay)



Digital USB 2.0 Camera with Progressive Scan Technology

The MyRay C-U2 Intraoral Camera generates clear, explicit images and features the highest degree of illumination.

The digital video system is fully integrated into the handpiece, and progressive scanning allows perfect image-freeze and guarantees top-quality definition.

The viewing angle is partially inclined away from the extremity of the handpiece for easy access to distal and rear tooth areas.

The C-U2's 12 micro-light provide the highest degree of illumination in its category.

Communication made simple: Efficient communication thanks to sharp images. Encourage patients to undergo treatment by involving them at each step of the plan. Visual communication is a simple and hands-on experience. C-U2 brings this type of communication to life: straightforward and simple to use, totally automatic without the need for any manual adjustments, it generates clear and explicit images. Technologically advanced beneath the surface, it is easy and intuitive to use.


C-U2 is the communication tool designed for professionals who have chosen to increase their efficiency through the use of a computer in the operating area, right beside where they work. The use of a computer makes acquisition and filing of dental images a simple task, organising them in nominal clinical folders, and ensures long-term rapid access for consultation purposes.

C-U2, digital by nature.

12 micro-lights: The highest degree of illumination in its category, enabling an efficient spread of light examining the oral cavity, thus producing bright images w/sharp defined colors. Automatic brightness setting, Progressive scan, Communication made simple, Focus Free, Touch Sensitive Control, Retro Flex Vision, Easy to use and Quick connect

<b>ItraOral Camera USB2,</b>	20-11122	<b>1694.95</b>
<b>Active holster</b> for usage in additional room	20-11123	<b>49.95</b>
<b>Hygienic Sheaths</b> for C-U2 camera (50 pcs)	20-11124	<b>29.95</b>

## MultyCam Intraoral Camera USB Type (Dentazon)



### Multy Cam USB Type Full Set,

include: MultyCam, USB Cable, Camera Holder and 300 Pcs. Camera Sheath

20-11125

**669.95**

## Cammy USB Intraoral Camera (Dentamerica)



Directly connects to a computer via the USB port.

Designed with a freeze button to capture images and a joystick to control digital zoom.

Self-sufficient, CCD technology, digital image processing.

Four white LED lamp illumination to provide maximum realistic color reproduction.

built-in lens protector designed to prevent dust and scratches to optical components.

Mini-head design allows for easy access to third molar.

Lightweight (only 92 grams). ergonomically designed handpiece.

User friendly software. Easy to use and integrate into patient management software.

Digital images can be easily stored, retrieved and transmitted electronically.

Disposable sleeves included to prevent cross contamination.

### Cammy USB Intraoral Camera

With Software, Digital Zoom

20-11112

**594.95**

**Disposable Sleeves, 500/Pkg**

20-11113

**39.95**





Prestige Dental Products, Inc.  
Tel: 877-772-3888

### Camrex 191A Intraoral Camera (Dentamerica)




#### Intraoral Camera with Multiple Image Memory

The CAMREX handpiece is small and lightweight, weighing in at only 89 grams. It contains 4 white LED lamps for illumination to provide maximum color reproduction. The CAMREX has a 4 freeze frame memory stored in its video process unit. Single images can be viewed full screen or 4 images can be viewed simultaneously in a split screen format. The toggle switch is conveniently located on the handpiece so that the user can rotate through the captured images with ease. A built-in lens protector is designed to prevent dust and scratches to the optical components. Also included are disposable sleeves designed to fit snugly on the handpiece to prevent patient cross-contamination.

The Camrex is truly universal, being offered in the NTSC and PAL video formats as well as coming with a universal AC power source (AC 100V~240V). Live images are captured by using either the handpiece freeze button or the included foot pedal. The foot pedal has a special circular design so that it can be comfortably activated from all 360°. The CAMREX 191A is incredibly easy to setup. Just four connections are needed: the AC power supply, the CAMREX handpiece, the foot pedal, and the television video cable. You will be up and running in no time with the simplicity of the CAMREX 191A.

Camera 191A with America NTSC System	20-11117	739.95
Camera 191AP with Europe PAL System	20-11116	739.95
Adapter & softwear	20-11118	139.95
Zoom Lens each	20-11114	48.95
Disposable Sleeves, 500/Pkg	20-11115	42.95

### Camrex 192 Wireless Intraoral Camera (Dentamerica)



Built-in wireless receiver with 4 separate frequencies for optimal signal reception. Videoutput connects directly to TV or computer monitor through AV or VGA input. Self-sufficient, CCD technology, digital image processing. Four lamp white LED illumination to provide maximum realistic color reproduction.

Built-in lens protector designed to prevent dust and scratches to optical components. Mini-head design allows for easy [access](#) to the third molar region. Disposable sleeves included to prevent cross contamination. NTSC or PAL available.

Camera with PAL Europe System	20-11120	984.95
Camera with NTSC System	20-11121	984.95
Zoom Lens each	20-11119	48.95


## Cameras - Photography

### Any Multipurpose Intraoral Camera (GoodDrs)




Anything Anypurpose  
'ANY'

With Any you can check your oral health at home  
You can manage your oral health and your family's oral health easily if you only invest one minute a day with ANY  
Save time and money with early detection  
Easy to install and use. Produces quality images similar to professional microscopes used in dentistry.  
Also can enhance your dental office's image by increasing the chance of early detection of cline's oral problems.  
Disposable plastic covers so you can share ANY with your whole family.



Any Multipurpose Intraoral Camera	20-11111	99.99
-----------------------------------	----------	-------

### Photographic Mirrors (Ortho Technology)



double-ended mirrors are perfect imaging solution. High reflectivity and are virtually indestructible. Autoclavable up to 273F/134C.

<b>F-</b> Adult Buccal/Child Buccal, 2 7/16" (62.0mm) x 7 3/8" (189.0mm) x 1 13/16" (47.0mm)	41-11404	59.95
<b>E-</b> Lingual Upper/Lingual Lower, 1 13/16" (47.0mm) x 7 11/16" (195.0mm) x 1 13/16" (47.0mm)	41-11399	59.95
<b>D-</b> Child Occlusal/Child Buccal, 2 9/16" (65.0mm) x 6 15/16" (176.0mm) x 1 13/16" (47.0mm)	41-11400	59.95
<b>C-</b> Adult Occlusal/Child Occlusal, 2 7/8" (73.0mm) x 7 1/4" (185.0mm) x 2 9/16" (66.0mm)	41-11401	59.95
<b>B-</b> Adult Occlusal/Adult Buccal, 2 13/16" (72.0mm) x 7 1/2" (190.0mm) x 2 1/8" (34.0mm)	41-11403	59.95
<b>G-</b> Buccal Lateral, 1 5/8" (42.0mm) x 6 11/16" (170.0mm) x 1 5/8" (42.0mm)	41-11405	59.95
<b>A-</b> Adult Occlusal/Child Occlusal, 2 3/4" (70.0mm) x 7" (178.0mm) x 2 3/16" (55.0mm)	41-11402	59.95





**Prestige Dental Products, Inc.**  
www.prestigedentalproducts.com

## Order Form

Date:

### Ordered By

### Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

### Payment

### ☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

**Total**

Dental License #	
DEA #	

Year Licensed	
Year Licensed	

Expiration Date	
Expiration Date	



## Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

[www.prestigedentalproducts.com](http://www.prestigedentalproducts.com)

