

Curing Light Units



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**




Curing Light Units

Ledex Violet LED Curing Light Unit (Dentmate)

- Radiant intensity up to 1800 mw/cm²
- Color: Black only
- LEDs type: 10 watt violet & blue LEDs
- LEDs wavelength: 390 ~ 480 nm, peak: 405, 460 nm
- Efficient optics: Up to cure resin 2mm in 3 seconds (High Mode)
- Handpiece size Ø 38 mm (max.) x L190 mm
- Handpiece weight 180g
- Cradle size Ø 115 mm (max.) x H68 mm
- Battery type Lithium-ion, 3.7V, 2500mAh
- Battery charge time Apprx. 2.5Hours
- Operating time 720 times in 10 seconds standard mode
- Power supply Input: 1~ 240 V, 50/60 Hz, output: DC 5 V / 1.2 A

Combo of Curing Light & Caries Detector



Full Unit: Hand piece, Battery, Cradle, Power supply (AC100~240V/DC5V), Optical fiber light guide rod (Ø 11>8mm), Optical fiber light guide rod (Ø 11>5mm), Filter, FINDTM Optical Filter, Anti-glare shield, Light Guide Sleeve, Curing Light Sleeve and User*

25-99971 **599.00**

Battery Fit

WL-90+WL, 3.7V/2500mAh 25-99972 **89.95**

Optical Fiber Light guide Ø 11>5mm, black, sterilize-able 25-99973 **89.95**

Optical fiber light guide Ø 11>3mm, black, sterilize-able 25-99974 **89.95**

Optical fiber light guide Ø 11>8mm, black, sterilize-able 25-99975 **79.95**

Curing Light Sleeve Fit

LEDX WL-90+, 100Pcs/Box 25-99976 **19.95**

ART-L3 LED Cordless Curing Light (Bonart)



3 WATTS CORDLESS CURING LIGHT SYSTEM

At long last, introducing the latest innovative dental technology from Bonart, the portable and light weight ART-L3(PS) pen type cordless LED curing light. Great quality, economical, and compact. Built to outlast its competitors.

Blue Color 25-11112 **224.95**



Prestige Dental Products, Inc.
Tel: 877-772-3888

ART-L5 LED Cordless Curing Light (Bonart)



System includes following accessories:
Light Guide Fiber and Protector
Rechargeable Base
AC Power Adapter Cord
User Manual

ART-L5 Blue Color 25-11113 **324.95**
Curing Light Tip 25-11114 **79.95**
Charge Base 25-11115 **54.95**

Radii Plus LED Curing Light Unit (SDI)



Cordless LED Curing Light.
High Intensity:
1500m W/cm²., 6mm depth of cure
360 Degree turning tip.
Lightweight 6.2 oz/177g.
1200 ten second cures before needing to be recharged.
Built-in radiometer, No noisy fan

Unit 25-11178 **868.95**

Radii-Cal LED Curing Light Unit (SDI)




Cordless lightweight
Radically powerful
Intensity of 1200mW/cm²,
Stays charged for 720 ten second cures
No noisy fan
Built in radiometer

Unit 25-11179 **499.95**

Ledex Blue LED Curing Light Unit (Dentmate)

- Radiant intensity Up to 2400 mW/cm²
- Color:Blue
- LEDs type: 10 watt blue LEDs
- LEDs wavelength: 440 ~ 480 nm, peak: 460 nm
- Effient optics: up to cure resin 2mm in 2 seconds (Turbo mode)
- Handpiece size Ø 38 mm (max.) x L190 mm



Full Unit: Hand piece, Battery, Cradle, Power supply (AC100~240V/DC5V), Optical fiber light guide rod (Ø 11>8mm), Filter, Anti-glare shield, Light Guide Sleeve, Curing Light Sleeve and User's Manual 25-99961 **449.00**

Battery Fit WL-90/WL, 3.7V/2500mAh 25-99962 **79.95**

Optical Fiber Light Guide Ø 11>8mm, black, sterilize-able 25-99963 **79.95**

Curing Light Units

Curing Light Sleeve Fit LEDEX WL-90/ 25-99964 **19.95**
WL 100Pcs/Box

Fusion Curing Light Unit (Dentlight)



Dentlight has developed patient pending focus beam technology- FBT. FBT allows fusion to unleash high power intensity, light beam and more importantly maintaining high power intensity over distance. This high intensity beam can be seen when testing the unit. Fusion Cures 2mm A3 and lighter shades in less than 5 seconds and can cure bulk composites as deep as 10mm in 20-30 seconds.

Silver Color	25-11137	694.95
Pink Color	25-11136	694.95
Blue Color	25-11135	694.95
Gold Color	25-11134	694.95
Battery Replacement	25-11133	109.95
Safety Glasses	25-11132	15.95

Hilux 250 Curing Light (First Medica)



HILUX™ HALOGEN LIGHTS have been a market leader for over 10 years. Hilux™ has earned a reputation for reliability and exceptional value. Hilux™ Lights will cure all light activated materials. All Hilux™ Curing Lights are pistol shaped for good ergonomics. Hilux™ Lights are corded so the user never needs to worry about low battery charges.

Unit 25-11139 **584.95**

Clever Curing Light Unit (GoodDrs)



System Includes following accessories:

- * Wave Length: 440~490 (nm)
- * Power Source: DC 6~12V
- * Light Source: High Power LED
- * Light Weight (handpiece Type): Approx. 104 g & much more

Clever 25-11118 **449.95**



Prestige Dental Products, Inc.
Tel: 877-772-3888

Dr's Light Wireless Curing Light (GoodDrs)



High-Tech, High-Performance Wireless, LED Curing Light, Full Color Display and 6 Curing Modes
Specifications: 5 Seconds are enough for composite cure

- * High power 1600mW/cm² (+10%)
- * Various curing program modes
- * Ergonomic design (gun type & pen type)
- * Color LCD screen
- * Unique thermal control system
- * Long time use
- * Warranty Terms: 18 Month limited warranty.

Curing Light Unit	25-11131	499.95
Guide tip 8mm	25-11129	129.95
Guide tip 11mm	25-11128	109.95
Battery for Dr's Light	25-11127	99.95
Anti-glare Shield	25-11126	19.95

Ledex LED Curing Light Unit (Dentmate)



White Color	25-11153	199.95
Blue Color	25-11152	199.95
Orange Color	25-11151	199.95
Apple Green Color	25-11150	199.95
Black Color	25-11149	199.95
Yellow Color	25-11148	199.95
Pink Color	25-11147	199.95
8 mm Optical fiber light guide rod, sterilizable.	25-11145	89.95
Light Cure Probe	25-11144	99.95
Battery 3.7 V / 800 m Ah	25-11146	79.95



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

