

Equipment



Prestige Dental Products, Inc

PrestigeDental.com
877-772-3888

Welcome to Prestige Dental Products, Inc.

located in Anaheim, California USA is recognized for its exceptional customer service, high quality products, and its competitive low prices.

As an authorized Dealer of many brands in North America, we have provided dentists with advanced dental supplies and equipment, along with excellent customer service. Our bulk purchases allow us to provide products to our clients with highly competitive low prices, thus ensuring that our customers' needs are met with full satisfaction.

Prestige Dental Products offers all Dental Specialties a comprehensive choice of the products and equipment that are necessary to operate a successful dental practice in the fields of General Dentistry, Cosmetic Dentistry, Orthodontics, Oral Surgery, Endodontic, Periodontics, Pedodontics and Implant Dentistry.

Our Mission is to provide innovative, integrated Dental products and excellent services; and to be trusted advisors and consultants to our customers, enabling them to deliver the best quality patient care and enhance their practice efficiency and

Customer Information

Order

PHONE To place an order by phone, please call toll free Number 877-772-3888 Monday through Friday 9.00 a.m. -5:00 p.m. (PST).

FAX You may FAX your orders to 714-632-8821.

INTERNET www.prestigedentalproducts.com

You Can Place Orders Online any time in our secure website.

Our Online Catalog offers access to our Product Line, Monthly Specials, Manufacturer Specials, and Promotions.



Note: New Customers/First time buyers do not need User ID and Password to shop from us.

Simply add items to the shopping cart. At CHECKOUT, you will have the option to pay with a Credit Card through our Secure Servers.



Customize order for special customer

Email us to mass@prestigedentalproducts.com a list of the items you order on a regular basis, or contact our sales representative to the Toll Free (877) 772-3888 and we will prepare a special Price list just for you.

Payment

CREDIT CARD We accept Visa, MasterCard, American Express and PayPal. Credit card order will ship only to credit card holders name and address. If you want to ship your order to different address than the billing address, please pay through PayPal or through Bank to Bank method.

Bank to Bank payment through High Secure Intuit Payment-Network

Pay electronically through high secure payment network, after you place your order, Prestige Dental Products will email you the invoice which has a link to guide you through a one-time set-up process to send your payment, and it is free for you. The payment transaction will be processed from your bank account to Prestige dental Products bank account securely. Nobody can see or know your bank information through this high secure Payment-Network.

MONTHLY BILLING "Net 30 days payment"

To get a net 30 days payment terms, please fill out our credit application form. These accounts are subject to our open account terms pending credit approval. Satisfactory credit references are required. Your account payment will be due upon receipt of your monthly statement.

Shipping

GROUND DELIVERY

Free shipping for all Orders \$400.00 or more "Tax not included \$400.00" in the state of the west coast (Excluding: Equipment, Solutions, Large Heavy Cases of disposable Goods, Casting Materials, Hazardous and Refrigerated Items).

EXPEDITED SHIPPING UPS Next Day and Second Day Service for which additional shipping fees will be charged.

SHIPPING HAZARDOUS ITEMS All Hazardous products will have an Additional \$20.00 charge due to Dept. of Transportation Regulations. Hazardous materials may not be shipped by air (for example: Alcohol, Methanol, Ethanol, Laboratory Reagents, or anything flammable). No returns will be accepted.

All Refrigerated Items are non-returnable and are shipped Next Day Air only and charges are added to Invoice.

State and federal laws

DEA CERTIFICATE & LICENSE

To enable us to ship to you in compliance with the law, we mandatory request that you supply us with a copy of your DEA Registration Certificate and Physician License prior to your first order of any purchase items such as needles, syringes, I.V. Products, anesthetics and pharmaceuticals.

Prestige Dental Products does not take responsibility for any typographical errors.

**All Prices And Specifications In This Catalog
Are Subject To Change Without Prior Notice.**



Equipment



Prestige Dental Products, Inc.
Tel: 877-772-3888



Electromechanical Patient Operatory Chair (DCI)

Standard Features:

One operator-programmable position, one preset exit position
Electronic foot switch preposition and auto return switch
Built-in fuse/surge protection
Toe guard protection to prevent scuffing
All-metal construction
Safety shut-off base plate
Double-articulating headrest, adjustable for height and wheelchair positioning
Locking swivel, 60°
Right-hand or left-hand chair mount adapter
Pivoting arms for easy entry to and exit from chair
Seamless upholstery
Coordinated leg rest tilt
Input voltage 110V or 220V

Electric Patient Operatory Chair 30-11191 **3199.95**



Hydraulic Patient Operatory Chair (DCI)

Standard Features:

Five operator-programmable positions
Electronic foot switch preposition and auto return switch
Built-in fuse/surge protection
Toe guard protection to prevent scuffing
All-metal construction
Safety shut-off base plate
Double-articulating headrest, adjustable for height and wheelchair positioning
Locking swivel, 60°
Right-hand or left-hand chair mount adapter
Pivoting arms for easy entry to and exit from chair
Seamless upholstery
Coordinated leg rest tilt
Input voltage 110V or 220V

Hydraulic Patient Operatory Chair 30-11233 **3699.95**

Orthodontic Chair (Boyd)

STANDARD FEATURES: All-Steel Frame Construction,
Low Voltage
DC Motors, Cantilevered Style Fixed Height Pedestal Base,
Motorized Back, Three (3) Programmable Positions,
Auto Return to Home Switch, Snap-In/Snap-Out Upholstery,
Winged Back with Contoured, One-Piece Seat, Convenient
Integrated Foot Control.



Hydraulic Patient Operatory Chair 30-11297 **3299.00**



Prestige Dental Products, Inc.
Tel: 877-772-3888

Equipment



Ceiling Mount Operatory Light (DCI)

On off switch. Input voltage 110V or 220V

Ceiling Mount Operatory Light 30-11150 **Call for Price**

Top Mount Chair Operatory Light (DCI)

On off switch. Input voltage 110V or 220V

Top Mount Chair Operatory Light 30-114021 **Call for Price**

Curverd Post Operatory Light (DCI)

On off switch. Input voltage 110V or 220V

Curverd Post Operatory Light 30-11158 **Call for Price**

Track Mount Operatory Light (DCI)

On off switch. Input voltage 110V or 220V

Track Mount Operatory Light 30-114071 **Call for Price**



Celux Operatory Light (Dentamerica)

Post Mount, Ivory color 30-11152 **629.95**

Replacement Bulb H3, 12V, 55W 30-11151 **15.95**

Clear Safety Shield 30-11153 **24.95**

Operatory Light Bulb (Parts)



14.95

30-11289

E-1552 for

* Belmont
(040, 046, 047, 048, 049, Excaliber)
* Adec Performer
(OEM #041-513-00)
* Chayes Virginia
* Healthco
(Aristocrat & Celebrity)
* Celux
12V, 55 Watt



30-11288

9.95

E-1550 for

Adec Cascade and older
(OEM #041-179-01)
* J. Morita
* LSM
* Eurodent
* Midmark
* Planmecca Delight
* Faro Sunlight 70
(Original Units)
* Siemens A1/E1 Sirona
* SDS 1350
* Knight Model L
24 V, 150 Watt



E-1551 for

22.95

Pelton & Crane LF Plus,
LF-1, LF-II
24 Volt, 150
Watt



E-1554 for


17.95

Faro
(Kool Beam, Sunlight S-90, S2000)
Faro Sunlight 70
(Upgraded units)
Ritter K Super Starlight-
Schein (Crusader, Star-brite)
Porter Daystar S-90B
Marus

30-11286

Prestigedentalproducts.com

www.PrestigeDentalProducts.com



**Check our
website for
more items
and
variation**

Call us for Quarterly Promotion

Equipment



Prestige Dental Products, Inc.
Tel: 877-772-3888

Alliance Ensemble Dental Operatory Package (DCI)

DC1235 hydraulic chair

Standard Chair Features:

Cast bracket withswing-out arms

Five operator-programmable positions

One operator-programmable position, one preset exit position

Electronic foot switch preposition and auto return switch

Built-in fuse/surge protection

Toe guard protection to prevent scuffing

All-metal construction

Safety shut-off base plate

Double-articulating headrest, adjustable for height and

wheelchair positioning

Locking swivel, 60°

Right-hand or left-hand chair mount adapter

Input voltage 110V or 220V

Coordinated leg rest tilt

BioFree tubing standard on all units

Straight aseptic tubing for three handpieces

Convenient brake release for flex arm

Master on-off control

Individual non-retracting water coolant adjustment

Purge system for handpiece tubing

Plastic B-size tray

Pre-wired with 8 conductor cable for accessory addition

Wet-dry foot control, Floor utility center with regulators and shut-off valves

Post-mount utility center, Post-mount utility center with solids collector

Self-contained water system with city water/bottled water selector, Patented quick-switch system for water bottle

QD water connection with flow control, QD air connection



Alliance Ensemble Dental Operatory Package

30-11140 **Call for Price**

Alliance Swing Mount Ensemble Operatory Package (DCI)

Rotating system design quickly and easily converts from left- to right-hand position

Individual drive air adjustment

Straight aseptic tubing for three handpieces

Convenient brake release for flex arm

Master on-off control

Individual non-retracting water coolant adjustment

Purge system for handpiece tubing

Plastic B-size tray

Pre-wired with 8 conductor cable for accessory addition

Wet-dry foot control

Floor utility center with regulators and shut-off valves

Individual drive air adjustment

Individual non-retracting water coolant adjustment

Purge system for handpiece tubing

Convenient brake release for flex arm

Self-contained water system with city water/bottled water selector

Patented quick-switch system for water bottle

DC1235 hydraulic chair

Standard Chair Features:

Cast bracket withswing-out arms, Five operator-programmable positions

One operator-programmable position, one preset exit position

Electronic foot switch preposition and auto return switch

Built-in fuse/surge protection, Toe guard protection to prevent scuffing, All-metal construction, Safety shut-off base plate Double-articulating headrest, adjustable for height and wheelchair positioning

Locking swivel, 60°, Right-hand or left-hand chair mount adapter

Input voltage 110V or 220V, Coordinated leg rest tilt



Alliance Swing Mount Ensemble Operatory Package

30-11144 **Call for Price**



Prestige Dental Products, Inc.
Tel: 877-772-3888

Equipment

Pro Ensemble Operatory Package (DCI)



30-11311

Call for Price

Package:

DC1335 Electric Operatory Chair

PRO30 Delivery System

1244 Curved Operating Light

Traditional Doctor's Stool

Traditional assistant's stool

EC-DPROE

C1335 Chair, PRO30 Delivery System, 1244 Curved Operating Light, Traditional Doctor's Stool & Traditional assistant's stool

Pro Hygiene Operatory Package (DCI)

30-11317

Call for Price

Individual non-retracting water coolant control
Hinged control head

Standard three-way syringe with autoclavable tip

Standard junction box, includes preset regulators
Master on/off switch

Purge system

Wet/dry foot control

Plastic B-size tray

High volume evacuator and hose

Saliva ejector and hose

DC1335 Electromechanical Chair

Standard Chair Features:

One operator-programmable position, one preset exit position

Electronic foot switch preposition and auto return switch

Built-in fuse/surge protection

Toe guard protection to prevent scuffing

All-metal construction

Safety shut-off base plate

Double-articulating headrest, adjustable for height and wheelchair positioning

Locking swivel, 60°

Right-hand or left-hand chair mount adapter

Input voltage 110V or 220V

Coordinated leg rest tilt



A10 Operatory System (Flight Dental)

30-11117 Call for Price

The A10 is one of our premium packages featuring a large doctor work surface and the elegant one-piece cuspidor. The generous size of the A10 will accommodate your largest patients comfortably while providing ideal working access for the doctor and assistant.



- Combine with your choice of A6 or A10 patient chair
- 3 Handpiece automatic control with air/water syringe
- Rotating utility center with one piece cuspidor
- Adjustable timer setting for bowl rinse and cup fill
- Unit mounted treatment light with two intensity settings
- Touchpad control for doctor and assistant controls chair, light, cuspidor and bowl rinse
- Choose from large or compact doctor's handpiece control
- City/Closed water system standard on all units
- Foot controls handpieces and offers water on/off toggle and chip air blower

Options:

Built in ultrasonic or peizo scaler, Handpiece fibre-optics system , Multi-media monitor mount, Integrated intra-oral camera, Continental style delivery, Glass Bowl Cuspidor Optional

A12 Operatory System (Flight Dental)

30-11123 Call for Price

A12H Hydraulic Lift Chair is integrated with cutting edge technology improving patient and doctor comfort. This new ergonomically designed chair provides unparalleled comfort and functionality for the dentist, making it the perfect chair for any dental operatory.



The exceptional range of motion coupled with its durability and sensible value, is the result of the perfect union between man and space.

Hydraulic motor with 380mm as the lowest position capability of the chair.

Dual 180 rotatable arms allows for convenient placement of cuspidor and tray units for doctor and assistant.

Left-hand and right-hand capability optimizes efficiency and space for doctor and patient.

New designed dental light providing superb luminosity.

Thin backrest provides the optimal space and access for the dentist and close positioning to the patient.

12 seamless or plush upholstery colors.

4 programmable chair position settings to make the work more efficient for the dentist.

Synchronized backrest and seat movements to provide the most comfortable patient seating

Optional LCD monitor mount, Optional G.Comm LED Light.



Prestige Dental Products, Inc.
Tel: 877-772-3888

Equipment

A2 Operator System (Flight Dental)

30-11130 Call for Price

The A2 chair is the most economical chair in our line of products. The backbone to a reliable chair is contingent on both the engineering and quality of the materials used. The chair is produced using some of the highest standards and as such has received praise as one of the best chairs on the market. Furthermore, the entire unit is finished using a Powder-coated baked finish to enable the highest level of asepsis control and with a cast iron frame the chair provides excellent durability. The company continues to exceed expectations, satisfying customer's desire for an economically priced unit built using the highest standards.



- Combine with the A2 patient chair
- 3 way air/water syringe with fully autoclavable tips
- Rotating one piece cuspidor
- Automatic bowl rinse and timed-cup fill functions accessible by patient and assistant side
- Unit mounted treatment light with two intensity settings
- Touchpad control of patient chair for doctor and assistant
- Choose from large or compact doctor's handpiece control
- City/Closed water system standard on all units
- Foot controls handpieces and offers water on/off toggle and chip air blower
- 12 seamless and plush upholstery colors

Options:

Built in ultrasonic or piezo scaler, Handpiece fibre-optics system, Multi-media monitor mount, Integrated intra-o

A3 Operator System (Flight Dental)

30-11132 Call for Price

The A3 is Flight Dental's economy chair with a European inspired design. The workmanship seen on this chair is incomparable in the marketplace today, and has been one of Flight Dental's most successful endeavours. It combines versatility, aesthetics and superior quality to produce an operator package that is suitable for hygienists and dentists alike. The chair is built with the patient and dentist in mind, with its thin backrest and ergonomically correct seat positioning, it produces the ultimate synergy between man and environment.



- 3 Handpiece automatic control with air/water syringe
- Utility center with one piece rotating unique coloured glass bowl cuspidor
- Adjustable timer setting for bowl rinse and cup fill
- Unit mounted treatment light with two intensity settings
- Touchpad control for doctor and assistant controls chair, light, cuspidor and bowl rinse
- Choose from large or compact doctor's handpiece control
- City/Closed water system standard on all units

Foot controls handpieces and offers water on/off toggle and chip air blower

Equipment



Prestige Dental Products, Inc.
Tel: 877-772-3888

Autoclave, Sterilizer Unit (Tuttnauer 2340M)

Manual autoclave is designed to compliment any healthcare facility. The bright new easy to use panel was designed with the operator in mind. Our 2340M autoclave will satisfy all of your sterilization needs without compromising on quality, safety or reliability.



2340M

30-11145

Call for Price

Elara 11 Class B Automatic Autoclave Sterilizer (Tuttnauer)

The Elara 11 Sterilizer is designed for sterilization of wrapped, porous and hollow loads. The chamber volume is large, external dimensions are small, making it especially suitable for smaller clinics requiring a high turnover.



11 Class B

30-11190

Call for Price

EZ Automatic Autoclave/Sterilizer (Tuttnauer)



With the simplicity of one touch design, all your sterilization and drying needs are fulfilled.

EZ9	EZ10	EZ10K
30-11205	30-11204	30-11206

Call for Price

Large Capacity Manual Autoclave/Sterilizer (Tuttnauer 3870M)

Our largest manual table top autoclave. Perfect for multi-doctor practices, clinics or surgery



3870M

30-11257

Call for Price

Tuttnauer Autoclave/Sterilizer (Tuttnauer 2540M)



The 2540M Tuttnauer manual autoclave is designed to compliment any healthcare facility.

2540M

30-114111

Call for Price

Valueklave 1730 Sterilizer (Tuttnauer)



The Valueklave 1730 Autoclave will satisfy all of your basic sterilization needs at an economical price.

1730

30-114141

1949.95

Manual Kwiklave Sterilizer (Tuttnauer 2540MK)



The Tuttnauer 2540MK Series is also available in a quick cycle model, the Kwiklave (MK).

2540MK

30-11274

3749.95



Call us for Quarterly Promotion

Prestigedentalproducts.com



Prestige Dental Products, Inc.
Tel: 877-772-3888

Equipment

Dental Equipment

RXDC HyperSphere X-Ray Unit (MyRay)



HyperSphere Plus incorporates an ergonomically advanced touch-activated positioning technology and the freedom of wireless control.

HyperSphere technology

Locked and unlocked by 2 touch-sensitive areas, the tube head revolves freely around the spherical ball-joint, enabling freedom and precision in terms of positioning.

Easy controls

The wireless control device, which can be attached to any available surface in the surgery, offers a full range of intuitive exposure scenarios designed to make correct X-ray acquisition a straightforward process.

X-Pod Wireless Dental Digital X-Ray Sensor (MyRay)



Stand-Alone Wireless Digital X-Ray sensor with Hand-held Touchscreen

Self-sufficient diagnostics

The sensors have cut corners, a two fold filtering system to give crisp high detailed image boasting 25 lp/mm.

X-POD is smaller than an IPAD and contains its own imaging software. **Call for best deal!**

X-Pod with sensore size #1	30-11469	Call for Price
X-Pod with sensore size #2	30-11468	Call for Price
X-pod Kit with sensors size 1 and 2, touchscreen, Bluetooth dongle, charger, rayMage imaging software, sensor positioning kits, holster	30-11470	Call for Price
Bluetooth dongle	30-11471	99.00

ZENX Digital X-Ray Sensor (MyRay)



The Zen sensors have cut corners, a two fold filtering system to give crisp high detailed image boasting 25 lp/mm.

Zen-X has been developed to simplify the whole process of taking radiographs thanks to the outstanding ergonomic shape, which ensures maximum patient comfort.

Versatility and speed

Powerful USB 2.0 connectivity makes the system extraordinarily convenient and mobile and minimizes the time between X-ray exposure and on-screen viewing.

The sensors also have an inexpensive option to extend their 2 year warranty to 5 years at time of purchase.

Zen-X Kit with sensor size 1 (active area 20x30mm), control unit, rayMage imaging software, sensor positioning kit	30-11482	4994.00
Zen-X Kit with sensor size 2 (active area 26x34mm), control unit, rayMage imaging software, sensor positioning kit	30-11481	6194.00
Zen-X Kit Duo both #1 & #2 Sensors with control unit, rayMage imaging software, sensor positioning kit	30-11483	9994.00



Eva Digital X-Ray Sensor (Dent-X)

EVA Dental is a filmless, digital dental imaging radiology system that produces highest resolution, diagnostic quality images. Utilizing patented "Active Column" CMOS technology, EVA produces high efficiency, low noise images which are digitized and transmitted to your computer via a standard USB (Universal Serial Bus) port and is compatible with laptop or desktop systems.

Digital imaging delivers significantly increased revenue to today's high-tech practice. Requiring less labor, it increases patient throughput and is safer for the patient with reduced X-Ray exposure.

Easy-to-use, state-of-the-art digital imaging eliminates the inconvenience of film processing and the 3 to 6 minute wait for a single X-Ray. Patient archiving is fast and simple for patient data base archives. EVA interfaces with Dent-X ProImage or other image and practice management software suites. The operational simplicity of digital radiography consolidates workflow as it builds productivity.

Size #2 Sensor W/ Docking Station and Basic Software	30-97792	Call for price
Size #1 Sensor W/ Docking Station and Basic Software	30-97791	Call for price



Prestige Dental Products, Inc.
Tel: 877-772-3888

Equipment

REX Dental Digital X-Ray Sensor



2 Sizes Available

Specifications	Size 1 (600 mm ²)	Size 2 (900 mm ²)
Pixel Size	19 µm x 19 µm	
Pixel Pitch	19 µm in both directions	
Area Array Size	19.95 mm x 30.02 mm	36.00 mm x 25.99 mm
Number of Lines	1580	1896
Number of Columns	1050	1368
Cut Corners (front Side)	130 pixels	235 pixels
Cut Corners (back side)	53 pixels	137 pixels
Cable Outer Material	polyurethane	
Cable Length	3 m +/- 6 cm	

This sensor is designed with the most innovative features including unique grooves, and easy-to-use sensor holder system with high-resolution

Connection eliminates the need for a bulky control box

Size 1	30-97131	3995.00
Size 2	30-97132	4395.00

RXDC eXTend X-Ray Unit (MyRay)

We pleased to announce the eXTend DC x-ray unit is now available for sale. The intra-oral DC market is extremely competitive when it comes to pricing, and the eXTend allows MyRay to be much more competitive.

Simple wireless control


One handheld wireless digital control device allows you to program the unit from wherever you are in the surgery. There are no other wall-mounted or wired control panels to worry about. Easy-to-use and easy-to-handle, the control device offers a full range of intuitive exposure scenarios designed to make correct X-ray acquisition a straightforward process.

No more complex programming or button-crowded panels. RXDC defines automatically the correct exposure by selecting the region of interest.



RXDC eXTend with remote control, 0.7mm focal spot, 16" (40cm) extension arm	30-11346	Call for Price
RXDC eXTend with remote control, 0.7mm focal spot, 24" (60cm) extension arm	30-11347	Call for Price
RXDC eXTend with remote control, 0.7mm focal spot, 35" (90cm) extension arm	30-11344	Call for Price
Wall mounting plate for single stud mounting	30-11345	99.00

RXDC HyperSphere X-Ray Unit (MyRay)



Call for best deal!

HyperSphere technology

Locked and unlocked by 2 touch-sensitive areas, the tube head revolves freely around the spherical ball-joint, enabling freedom and precision in terms of positioning.

Easy controls

The wireless control device, which can be attached to any available surface in the surgery, offers a full range of intuitive exposure scenarios designed to make correct X-ray acquisition a straightforward process.

MyRay Multi-Mode The most appropriate combination of technical factors (kV, mA and exposure time) is automatically determined according to patient size and anatomic region of interest

0.4mm focal spot, remote control, 16" (40cm) extension arm	30-11350	Call for Price
0.4mm focal spot, remote control, 24" (60cm) extension arm	30-11351	Call for Price
0.4mm focal spot, remote control, 35" (90cm) extension arm	30-11348	Call for Price
Wall mounting plate for single stud mounting	30-11349	99.00

Hyperion Panoramic Imager Morphology Recognition Technology (MyRay)



Call for best deal!

Panoramic Imager with Morphology Recognition Technology (MRT)

Hyperion is a highly-advanced panoramic imager, which employs the latest technologies to enhance radiography acquisition and diagnostics. Via laser-guide multi-motor kinematics, Hyperion positions itself around patients, reads subjects and scans accordingly, Morphology Recognition Technology (MRT) automatically identifies patient size and all parameters required to ensure correct X-Ray exposure.

No need to program exposure times, KV or mA technical factors.

15 different diagnostic programs are available to satisfy any possible diagnostic need.

Hyperion Panoramic Imager With RayMage Image Software	59-11215
Hyperion X7 Panoramic Imager, Var. Collimator Ceph Upgradable	59-11214
Hyperion X7-C1 Pan/Ceph Imager, Movable Pan/Ceph Sensor	59-11213
Hyperion X7-C2 Pan/Ceph Imager, two Dedicated Sensors	59-11212
X7 Upgrade to X7-C2: Ceph Support Arm, right side, and Pan/Ceph sensor	59-11216





Prestige Dental Products, Inc.
Tel: 877-772-3888

Equipment

Portable Digital X-Ray Unit D3

It has the following features:

- *Convenience of being cordless using rechargeable Li-polymer battery pack
- *Benefit of having a great reduction of X ray exposure to the patient and dentist i-e, about 1/7 radiation of the conventional X ray machines
- *Has a high frequency circuit of inverter type incorporated in it
- *Has a most suitable and diverse sensor system
- *Very easy to use with a comfortable handle for holding it



Full unit plus 1. Shield, Two batteries, Charger, Bag, Manual CD, Software

30-97161

3,695.00

Manufacturer extended warranty

30-97162

500.00

Portable Digital All In One X-Ray Systems D4



Cordless Portable Digital All In One X-Ray Systems D4

helping out your professional skills in dentistry. It is now ready to replace the previous versions of mounted X-Ray machines.

Portable X-ray D4 provides you with

- *Convenience of being cordless
- *Charging battery pack of loading/unloading type
- *Minimum chances of radiation exposure as it emits 1/7th radiations as compared to previous systems
- *More compact in size and shape
- *has inverter type high frequency circuit incorporated in it
- *Applicable to both digital and film type sensor

Unit/Sensor Size 1, Build in Wifi module

30-97221

8594.00

Unit/Sensor Size 2, Build in Wifi module

30-97222

9994.00

Manufacturer extended warranty

30-97223

500.00

Portable X-Ray Unit REXTAR-X (Meta)

Compact & Light

- *One-hand operation to make exposure
- *Only 1.6Kg with battery providing up to 300 exposures with one time charge

Automatic and Manual exposure control system

- *With smart exposure control system, proper exposure vales and time can be memorized to get the optimized darkness and contrast for future to same time to reset every single case and patient



Unit with Carrying Bag

30-97231

3,494.00



Prestige Dental Products, Inc.
www.prestigedentalproducts.com

Order Form

Date:

Ordered By

Deliver To

☐ Same as Above

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Business Name:

Address:

State/Province:

Zip/Postal Code:

Phone:

Fax:

Contact Name:

Item Number	Description	Quantity	Unit Price	Amount

Payment

☐ Credit Card

☐ Online Payment (Bank to Bank)

E-mail

☐ American Express

☐ MasterCard

☐ Visa

☐ Check Payable To

Card Number:

Expiration Date:

Cardholder Name:

Sub Total

C.O.D. Charges

Freight Charges

Sales Tax

Total

Dental License #

Year Licensed

Expiration Date

DEA #

Year Licensed

Expiration Date



Prestige Dental Products, Inc

3031 E. Coronado St., #D

Anaheim, Ca 92806

Tel:877-772-3888 Fax:714-632-8821

www.prestigedentalproducts.com

