

Contents

2-9 Elastomerics

- 2-6 Ligature Chain and Ties, Storage
- 7 Thread
- 8-9 Modules, Wedges, Separators

10-12 Adhesives & Accessories

13-15 Wire & Accessories

16-19 Brackets

20-21 Instruments

22-23 Elastics

24-27 Procedure Aids

- 24 Abrasive Strips, Arch Gauge, Ortho Zip Ties
- 25 Tongue Blocks
- 26-27 Expanders, Retractors, Propy Aids

28 Lab Products

29-31 Impressions

32-34 Plastic

- 32 Plastic Comparison Chart
- 33 Duraclear™ Plastic, Duraforce™ Plastic
- 34 Duraguard™ Plastic

35 Sports Mouthguards

36-39 Patient Aids

- 36-37 Wax Boxes, Silicone Boxes
- 38-39 Pontic Wax, Chewies™ Aligner Tray Seaters, Retainer Retrievers, Kits


40-41 Appliance Cases





Elastomeric Ligature Ties

- Low profile fits under bracket tie wings to avoid irritating the gingiva
- Individual Ties – Configured to help eliminate cross contamination
- Mini Ties – Allow the use of a unique color on each arch, without contamination issues or waste
- Choose from 61 colors



Available in **BULK PACKS!**
Call for details.

		"O" Ties	"S" Ties	Mini Ties	Individual Ties		Long Ties
	Diameter	.120"	.120"	.120"	.110"	.120"	.120"
Single Color	Standard, Glow, Jelly, Fluorescent	TT-3(color)	TT-1(color)	TT-4(color)	TTS-2(color)	TT-2(color)	TT-5(color)
	Metallic	TTM-3(color)	TTM-1(color)	TTM-4(color)	TTSM-2(color)	TTM-2(color)	TTM-5(color)
	Jewel	JT-3(color)	JT-1(color)	JT-4(color)	JTS-2(color)	JT-2(color)	JT-5(color)
Multi-Color	Standard (B, G, HP, OR, R, SIL)	—	—	TT-4MC	—	TT-2MC	TT-5MC
	Metallic (B, BLK, C, GRN, P, R)	—	—	TTM-4MC	—	TTM-2MC	—
	Fluorescent	—	—	TTF-4MC	—	TTF-2MC	—
	Glow	—	—	TT-4GLOWMC	—	TT-2GLOWMC	—
	Jelly	—	—	TT-4JMC	—	TT-2JMC	—
	Masculine (B, BLK, BRN, BY, TG)	—	—	TT-4MMC	—	TT-2MMC	—
		100 per dispenser; 10 dispensers per bag		10 per dispenser; 100 dispensers per bag	22 per dispenser; 45 dispensers per bag		52 per dispenser; 20 dispensers per bag



Elastomeric Ligature Chain

Tuff Chain™ Elastomeric Chain

- Chain available in 62 colors
- Spool color can match chain color or be clear

Ultra Tuff™ Elastomeric Chain

- .020" thickness provides clearance for tying chain on top of ligature ties and piggybacking with other force systems
- Spool color matches chain color


Thin Chain Elastomeric Chain

- Perfect for lingual treatment, self ligating brackets and overall patient comfort
- .015" thickness for easy placement and piggybacking over other appliances; Ligates easily over tie wings and self ligating brackets
- Spool color matches chain color

Thin Chain not available in Open Medium style.

Available in **BULK PACKS!**
Call for details.

	Tuff Chain™ Elastomeric Chain		Ultra Tuff™ Chain	Thin Chain
15ft Spool	.024" Thickness		.020" Thickness	.015" Thickness
Chain Color	Colored Spool	Clear Spool	Clear Spool	Clear Spool
Standard, Glow, Jelly, Fluorescent	GCC-(color)(style)	GC-(color)(style)	—	—
Metallic	GCCM-(color)(style)	GCM-(color)(style)	—	—
Jewel	GCCJ-(color)(style)	GCJ-(color)(style)	—	—
Clear	—	GC-C(style)	UT-C(style)	GCT-C(style)
Dark Gray	GCC-DG(style)	GC-DG(style)	UT-DG(style)	—
Silver	GCC-SIL(style)	GC-SIL(style)	UT-SIL(style)	GCT-SIL(style)
Pearl	GCC-PEARL(style)	GC-PEARL(style)	—	GCT-PEARL(style)



Elastomeric Ligature Ties & Chain Color Options

Standard

Aqua A	Fire Engine Red FIRE	Light Pink LP	Red R
Azure AZ	Fuchsia F	Maroon M	Spring Green SG
Blue B	Gold GLD	Navy Blue NB	Silver SIL
Black BLK	Green GRN	Orange OR ties, O chain	Sky Blue SKY
Champagne CHPN	Hot Pink HP	Orchid ORCH	Teal T
Clear C	Lavendar L	Pearl PEARL	Tooth TC
Cobalt Blue CB	Light Blue LB	Pearl Blue PEARLB	Violet V
Dark Green DG ties, DGRN chain	Light Green LG ties, LGRN chain	Pearl Green PEARLG	White W
Dark Grey D ties, DG chain	Light Grey LG chain only	Purple P	Yellow Y
Emerald Green E			

Metallic

Black M-BLK
Blue M-B
Copper M-C
Green M-GRN
Ocean M-OG
Purple M-P
Red M-R

Jewel

Jade J-J
Sapphire J-S

Fluorescent

Red FR
Green FG
Yellow FY
Orange FO
Pink FP

Jelly

Jelly Blue JB
Jelly Green JG
Jelly Purple JP
Jelly Red JR
Jelly Yellow JY ties only

Glow

Glow GLOW
Glow Blue GLOWB
Glow Green GLOWG
Glow Orange GLOWO
Purple GLOWP
Red GLOWR

Certain colors require specific minimum quantities when ordered. Call for more information.
Colors are for illustration purposes only. Actual color of product may vary.

Loose & Caned Elastomeric Ligature Ties



Loose Ties


Diameter	.110"		.120"	
Standard	PET(color) 1000SB	PET(color) 5000SB	PET(color) 1000B	PET(color) 5000B
Multi	PETMC 1000SB	PETMC 5000SB	PETMC 1000B	PETMC 5000B
Ultra Slide™ Dark Gray	—	—	PETUS 1000B	—
	1000 ties per bag	5000 ties per bag	1000 ties per bag	5000 ties per bag

Caned Ties

Diameter	.110"	.120"
Standard	PET(color) 1000S	PET(color) 1000
Multi	PETMC 1000S	PETMC 1000
Ultra Slide™ Dark Gray	—	PETUS 1000
	50 ties per cane; 20 canes per pack	

Standard

Azure AZ	Glow GLOW	Light Pink LP	Silver SIL
Black BLK	Glow Green GLOWG	Light Yellow LY	Smoke SMK
Blue B	Gold GLD	Maroon M	Spring Green SG
Clear C	Green G	Orange OR	Tooth Color TC
Dark Grey D	Hot Pink HP	Purple P	White W
Emerald Green E	Lavendar L	Red R	Yellow Y
Fuchsia F	Light Blue LB		



Certain colors require specific minimum quantities when ordered. Call for more information. Colors are for illustration purposes only. Actual color of product may vary.













ChromaTies™

ChromaTies™ Ligature Ties

- Two colors in a single ligature tie
- Popular sports team and fashion colors are fun for patients

		.120" Diameter
Loose Ties	EDT-(color) 1000B	1000 ties per bag
Loose Ties	EDT-(color) 5000B	5000 ties per bag
Caned Ties	EDT-(color)	50 ties per cane; 20 canes per pack

Colors are for illustration purposes only.
Actual color of product may vary.



Available in
BULK PACKS!
Call for details.

Ultra Slide™ Elastomeric Chain and Ties


- Slick surface makes placement on ceramic brackets easy, and provides reduced friction with wires.
- Available in Clear (C), Dark Gray (DG), Silver (SIL) or Pearl (PEARL)

Ultra Slide™ Chain

- 15ft Spool

GCUS-(color)(style)

Style



Mini Patient Ties

- .120" diameter
- Bag of 100 dispensers, 10 ties per dispenser

Not available in Pearl.

TTUS-4(color)

Individual Patient Ties

- .120" diameter
- Bag of 45 dispensers, 22 ties per dispenser


TTUS-2(color)

"S" Ties

- .120" diameter
- Bag of 10 dispensers, 100 ties per dispenser

Not available in Pearl.

TTUS-1(color)



Storage

Acrylic Chain Dispenser

- Dispenses spools of chain easily and quickly
- Heavy duty acrylic with slide top
- Holds up to 10 spools of chain or thread

Clear	CD-1AC
Smoke	CD-1AS

Ligature Tie Dispenser

- Heavy duty clear acrylic
- Stackable, compact design
- Each drawer holds up to 50 dispenser sticks or canes

Clear	TD-1A
-------	--------------



Elastomeric Thread

Lip Bumper Tubing

- Slides over the wire without snagging
- Fits up to .046" wire
- 10ft Spool

Clear	LBT-C
Gray	LBT
Opaque Tooth	LBT-TC

Archwire Sleevings

- Great patient relief aid; helps compliance
- Fits to archwire to help prevent food from entering tube
- 10ft Spool

	.018" Diameter	.027" Diameter
Clear	AWS-18C	AWS-27
Gray	AWS-18G	AWS-27G
Opaque Tooth	AWS-18TC	AWS-27TC



Tuff™ and Ultra Tuff™ Elastomeric Thread

- Tuff™ thread has a smooth surface
- Ultra Tuff™ thread has a rough surface so knots will not slide or come unraveled
- Different diameters provide varying degrees of force

	Diameter	Tuff™ Square	Tuff™ Round	Tuff™ Round - Hollow	Ultra Tuff™ Round	Spool
Gray	.020"	—	ETT-020	ETH-020	UTT-020	25ft spool
	.025"	ETT-2525	ETT-025	ETH-025	UTT-025	
		—	ETT-025x50FT	ETH-025x50FT	—	50ft spool
	.030"	ETT-3030	ETT-030	ETH-030	UTT-030	25ft spool
		—	ETT-030x50FT	ETH-030x50FT	—	50ft spool
	.020"	—	ETTC-020	ETHC-020	UTTC-020	25ft spool
Clear	.025"	ETTC-2525	ETTC-025	ETHC-025	UTTC-025	
		—	ETTC-025x50FT	ETHC-25x50FT	—	50ft spool
	.030"	ETTC-3030	ETTC-030	ETHC-30	UTTC-030	25ft spool
		—	ETTC-030x50FT	ETHC-30x50FT	—	50ft spool



Shrink Tubing

- Creates added protection and comfort when used over archwire and lip bumpers
- Slide over the wire and heat to create a smooth crack-resistant surface

6 per tube **ST-045**

Modules & Wedges

Tuff Link™ Elastomeric Modules

- Ligature ties with connecting bridges, for use with closing intra arch space
- Excellent for rotating individual teeth and connecting bracket hooks
- Molded from urethane that provides continuous gentle force




	4.3mm	6.0mm	8.5mm	11.5mm	Assorted
Dark Gray	TM-1	TM-2	TM-3	TM-4	TM-5A
Clear	TM-1C	TM-2C	TM-3C	TM-4C	TM-5AC
	100 per pack				

Rotation Wedges

- Unique wedge shape provides increased contact, allowing more force to be transmitted to the tooth

	Intro Kit	Standard	Increased
Gray	RW-KIT	RW-1	RW-2
Clear	RW-KITC	RW-1C	RW-2C
	100 per pack; 50 each style	100 per pack	




Standard
Wedge



Increased
Wedge




Separators



Clear Anterior Separators

- Thin for easy insertion; .054"-.055" thickness; .080"-.082" inner diameter; .190"-.192" outer diameter


10 per dispenser;
50 dispensers per pack **OS-2GCA**



Jumbo Separators

- Use for correcting partial molar impaction
- 2mm thickness and 6mm diameter allows maximum extension and memory

4 per dispenser;
32 dispensers per pack **OS-2XL**




Resep™ Extra Large Posterior Separators

- Great for hard-to-separate cases, or for when appliances such as a palatal expander are used
- Flat sides and round corners for easy insertion



10 per dispenser;
50 dispensers per pack **RS-2**

Orthoseps™ Elastomeric Separators

- Strong with a high resistance to tearing
- Rounded corners and a flat top help make insertion into tight contacts easy
- Blue for easy identification inside the mouth
- **Individual** configuration helps prevent cross contamination
- **Individual Radiopaque** is identifiable by x-ray, to help avoid loss underneath gum tissue




Individual	Individual Radiopaque	Molded Posterior
Blue	Dark Blue	Light Blue
OS-2	SS2	OS-1
10 per dispenser; 100 dispensers per bag		70 per dispenser; 15 dispensers per bag



Duraseps™ Elastomeric Separators

- Vertical slot makes it easy to insert with separator pliers
- Bars on both ends allow maximum stretch without breakage and help prevent slipping under gum tissue

		Standard Size
70 Shore A	Purple	DS8-2
85 Shore A	Orange	DS-2
10 per dispenser; 50 dispensers per pack		



Loose Separators

- Radiopaque style is identifiable by x-ray, to help avoid loss underneath gum tissue

	Anterior	Posterior	Radiopaque
Color	Green	Blue	Dark Blue
Thickness	.047"	.047"	.050"
Interior Diameter	.067"	.090"	.090"
Outer Diameter	.147"	.170"	.200"
100 per bag	ESS100	ESL100	ESR100
1000 per bag	—	ESL1000	ESR1000

Etching



Bonding Sponges

- Great method for applying etch and primers
- Conserves material and provides controlled application

200 per Jar

GBS-200

Mini Etch Sticks

- Use with etchant, bonding agents and sealants
- Bends easily into any angle for controlled placement in areas with limited access
- Reduces excess solution

100 per Tube

PES-100

Gel Etch Syringe Kits

- Syringes with disposable tips allow controlled etchant placement, without cross contamination
- Contains 37% phosphoric acid
- Thick viscosity gel for no slump or running

12g syringe,
(40) disposable tips

GE-1

(4) 3ml syringe,
(20) disposable tips

GE-1KIT

60g syringe,
(40) disposable tips

GE-JKIT

Individual Patient Etchant

- Helps eliminate cross contamination
- Dispenser allows accurate placement on each tooth, for less waste
- Each dispenser has enough etchant for approximately 28 brackets
- Contains 36% phosphoric acid

Available for US distribution only.

(25) 1.5g dispensers **ME-25**

Band Cement

Fluoride Glass Ionomer Band Cement

- Formula chemically bonds to non-precious metals and enamel
- Setting time formulated for orthodontic use

Intro Kit	15g powder, 15ml liquid	GI-15
Large Kit	50g powder, 55ml liquid	GI-50
Econo Kit	100g powder, 55ml liquid	GI-100

Fluoride Zinc Phosphate Band Cement

- Contains fluoride
- Smooth, creamy mix for easy application
- Light gold color helps distinguish it from the tooth
- High resistance to washout

Intro Kit	100g jar of powder, 55ml liquid	BCFI
Econo Kit	1lb (453.6g) jar of powder, (4) 55ml liquid	BCFI

Bonding Systems

Lite Bond Orthodontic Bonding System

- Unlimited working time enables adjustment of brackets as much as desired prior to final placement
- Approximately 30 sec. set time when used with most visible lights
- Color stable, no bracket float and resistance to abrasion

(4) 5g syringes, 6cc sealant, 9g etchant	LB-1
--	-------------

Paste & Paste Orthodontic Bonding System

- Adhesive system provides controlled application, color stability and bond strength
- Micro-filled to withstand abrasion
- Paste mixture can be adjusted to vary working time

18g base, 18g catalyst, 7.5ml etchant, 7g base resin, 7g catalyst resin	P-PP
--	-------------

Adhesives

Orthodontic Adhesives

- Cohesively stronger than the bond to the tooth, allowing easy one-piece removal without enamel scarring
- Tacky viscosity and good flow helps eliminate bracket float
- **One-Step Ortho Adhesive Kit** - sets on contact
- **VLC Orthodontic Adhesive Kit** - sets with a light cure source

One-Step Ortho Adhesive Kit	PNMG-1S
(4) 5g syringes, 12ml primer, 7ml etchant Available for US distribution only.	

VLC Orthodontic Adhesive Kit	LBG-1
(4) 5g syringes, 7ml resin, 7ml etchant	




Mixing Accessories

Mini Mixing Sticks

- 2.2" long
- Small design for easy adhesive mixing control and application

100 per pack **MMS-100**



Mixing Sticks

- Perfect for mixing banding and bonding cements
- 4.4" length makes it comfortable to hold and control

100 per pack **GMS-100**



Tuff Heat™ Spatula

- Excellent for glass ionomer cements
- Flexible for excellent control when mixing
- Rigid, handle grip helps prevent slippage
- Autoclavable; 7" long

1 per pack **TS-2**



Ultra Cement Spatula

- Tapered end is sturdy, yet flexible for easy mixing
- Long wide handle is easy to grip and hold
- 5.75" long

10 per pack **TS-3B**



Stay Cold Mixing Disc

- Allows convenient chairside mixing
- Ideal for loose band appointments
- Blue ice inner core stays cold up to one hour
- 3.125" diameter is large enough for one or two mixes
- Includes a plastic Tuff Heat™ Spatula

MD-1

Multi-Strand Wire

Super Three Strand Twist Wire


- Ideal initial archwire provides gentle force
- Will not unravel when cut
- Loops can be formed



Size	Natural Arch	14" Straight Length
.0155"	TSPA015	TSS015
.0175"	TSPA175	TSS175
.0195"	TSPA195	TSS195
	10 per pack	10 per tube

Super Smooth Cable Wire

- Five strands are wrapped around a core wire
- Excellent first wire – gentle and easy to tie
- Excellent for unraveling and leveling teeth
- Great resiliency and smooth surface reduces friction
- Will not fray when cut




Size	Natural Arch	14" Straight Length
.0155"	PCWA015	UCWS015
.0175"	PCWA175	UCWS175
.0195"	PCWA195	—
	10 per pack	10 per tube

Wire Accessories

Archwire Markers

- Disposable markers help prevent cross contamination
- Thick wax tips of sturdy dowels provide precise marking, high visibility, and superior adhesion
- Dispenser includes Velcro mounting

100 Archwire Markers	AWM-1
Archwire Marker Dispenser	AWMD-1




Archwire Holder





- Convenient chairside dispenser
- Holds 6 sizes of wire in both upper and lower forms; 12 slots total
- Ample room for labels next to each slot
- Stainless steel
- 7.875" long

AWH-N

Stainless Steel Wire

- Smooth surface for reduced friction
- Heat treated to help prevent corrosion
- Laser etched midline mark



	Size	Natural Arch		Standard Arch		14" Straight Length
		Upper	Lower	Upper	Lower	
Round 	.012"	P012U	P012L	—	—	—
	.014"	P014U	P014L	RWS14U	RWS14L	RWL14
	.016"	P016U	P016L	RWS16U	RWS16L	RWL16
	.018"	P018U	P018L	RWS18U	RWS18L	RWL18
	.020"	P020U	P020L	RWS20U	RWS20L	RWL20
		100 per pack		100 per pack		20 per tube
Retainer Round 	.025"	—	—	—	—	RET25
	.028"	—	—	—	—	RET28
	.030"	—	—	—	—	RET30
	.032"	—	—	—	—	RET32
	.036"	—	—	—	—	RET36
		—		—		10 per tube
Rectangular 	.016" x .022"	P622U	P622L	S622U	S622L	RWL622
	.017" x .022"	—	—	—	—	RWL722
	.017" x .025"	P725U	P725L	S725U	S725L	RWL725
	.018" x .025"	P825U	P825L	S825U	S825L	RWL825
	.019" x .025"	P925U	P925L	S925U	S925L	RWL925
	.021" x .025"	P215U	P215L	—	—	RWL215
Square 		10 per pack		10 per pack		10 per tube
	.016" x .016"	P616U	P616L	S616U	S616L	RWL616
	.018" x .018"	P818U	P818L	—	—	—
		10 per pack		10 per pack		10 per tube



Weldable Lingual Sheaths

- Designed for use with quad helixes and palatal bars
- .036" slotted and horizontal with latching indent
- Fits .072" x .036" appliances
- Stainless steel

4 per pack **LS36**

Nickel Titanium Wire

Marsenol Archwire


Nickel titanium alloy uses shape memory effect and is more elastic than stainless steel

- Delivers consistent force over a long period of time
- Exhibits a high load deflection characteristic when displaced over a short distance
- Smooth polished surface with laser-etched midline mark

Marsenol Archwire - Thermal

Heat-activated nickel titanium wire accepts a bend, then slowly returns to its original shape when placed intraorally and elevated to mouth temperature

- Smooth polished surface with laser-etched midline mark



	Natural Arch		"E" Arch		Natural Arch - Thermal		"E" Arch - Thermal	
Size	Upper	Lower	Upper	Lower	Upper	Lower	Upper	Lower
.012"	M12U	M12L	ME12U	ME12L	TM12U	TM12L	TME12U	TME12L
.014"	M14U	M14L	ME14U	ME14L	TM14U	TM14L	TME14U	TME14L
.016"	M16U	M16L	ME16U	ME16L	TM16U	TM16L	TME16U	TME16L
.018"	M18U	M18L	ME18U	ME18L	TM18U	TM18L	TME18U	TME18L
.020"	M20U	M20L	ME20U	ME20L	—	—	—	—
.016" x .016"	M616U	M616L	ME616U	ME616L	—	—	—	—
.016" x .022"	M622U	M622L	ME622U	ME622L	TM622U	TM622L	TME622U	TME622L
.017" x .025"	M725U	M725L	ME725U	ME725L	TM725U	TM725L	TME725U	TME725L
.018" x .025"	M825U	M825L	ME825U	ME825L	TM825U	TM825L	TME825U	TME825L
.019" x .025"	M925U	M925L	ME925U	ME925L	—	—	—	—
.021" x .025"	M215U	M215L	ME215U	ME215L	—	—	—	—

10 per pack — Also available in bulk 100 pack: add "B" after item code



Nickel Titanium Reverse Curve Archwire

- Excellent for opening and closing a bite
- Also used for cases with minor crowding

Size	Shape	Upper	Lower
.016"	Round	RCA16U	RCA16L
.016" x .022"	Square	RCA622U	RCA622L

10 per pack



Closed



Open

Nickel Titanium Coil Springs

- Excellent for opening or closing spaces with constant, gentle force

Size	Length	Style	Force	
.009" x .030"	7"	Open	Light	CSO-9NT
.010" x .030"			Medium	CSO-10NT
.012" x .030"			Heavy	CSO-12NT
.010" x .030"	9mm	Closed	—	CSC-9NT
.010" x .030"	12mm		—	CSC-12NT

3 springs per pack

Clear Composite Brackets

- Smaller size gives a better appearance, without loss of wings and rotational control
- Rounded everywhere for patient comfort and less chance of breakage
- Continuous slot helps eliminate trapped food and adds strength
- Non-polycarbonate plastic resists staining and water absorption for longer life
- Bottom locking ridges allow use of color-stable adhesive without a plastic primer



Roth Prescription		Torque	Angulation	L .018" Slot R		L .022" Slot R	
Maxillary Arch	Central	12	5	FRU118L	FRU118R	FRU122L	FRU122R
	Lateral	8	9	FRU218L	FRU218R	FRU222L	FRU222R
	Cuspid - Hook	0	11	FRU318LHK	FRU318RHK	FRU322LHK	FRU322RHK
	Cuspid	0	11	FRU318L	FRU318R	FRU322L	FRU322R
	Bicuspid - Hook	-7	0	FRU18BIHK		FRU22BIHK	
	Bicuspid	-7	0	FRU18BI		FRU22BI	
Mandibular Arch	Lower Anterior	0	0	FRL118		FRL122	
	Lower Cuspid - Hook	-11	5	FRL318LHK	FRL318RHK	FRL322LHK	FRL322RHK
	Lower Cuspid	-11	5	FRL318L	FRL318R	FRL322L	FRL322R
	Lower 1st Bicuspid - Hook	-17	0	FRL418HK		FRL422HK	
	Lower 1st Bicuspid	-17	0	FRL418		FRL422	
	Lower 2nd Bicuspid - Hook	-22	0	FRL518HK		FRL522HK	
	Lower 2nd Bicuspid	-22	0	FRL518		FRL522	
Standard Edgewise		Torque	Angulation	.018" Slot		.022" Slot	
Upper Central		0	0	FSU118		FSU122	
Upper Lateral		0	0	FSU218		FSU222	
Lower Anterior		0	0	FSL118		FSL122	
Cuspid		0	0	FS318		FS322	
Bicuspid		0	0	FS18BI		FS22BI	


Ceramic Aesthetic Brackets

- Small size, smooth surface and rounded corners are designed for patient comfort and aesthetics
- A true mechanical bonding base allows easy removal
- Dovetail base design
- Ample tie wing space makes ligation easy



Roth Prescription		Torque	Angulation	L	.018" Slot	R	L	.022" Slot	R
Maxillary Arch	Central	+12	+5	CRU118L		CRU118R	CRU122L		CRU122R
	Lateral	+8	+9	CRU218L		CRU218R	CRU222L		CRU222R
	Cuspid with Hook	0	+11	CRU318LHK		CRU318RHK	CRU322LHK		CRU322RHK
	Cuspid	0	+11	CRU318L		CRU318R	CRU322L		CRU322R
	Bicuspid	-7	0	CR18BI			CR22BI		
Mandibular Arch	Lower Anterior	0	0	CRL118			CRL122		
	Lower Twin Anterior	0	0	CRL118TWN			CRL122TWN		
	Cuspid with Hook	-11	5	CRL318LHK		CRL318RHK	CRL322LHK		CRL322RHK
	Cuspid	-11	5	CRL318L		CRL318R	CRL322L		CRL322R
	1st Bicuspid	-17	0	CRL418			CRL422		
	2nd Bicuspid	-22	0	CRL518			CRL522		

Standard Edgewise	Torque	Angulation	.018" Slot	.022" Slot
Upper Central	0	0	CSU118	CSU122
Upper Lateral	0	0	CSU218	CSU222
Lower Anterior	0	0	CSL118	CSL122
Cuspid	0	0	CS318	CS322
Bicuspid	0	0	CS18BI	CS22BI



Bracket Height Gauge

- Used to measure from incisal edge to bracket slot for precise bracket placement
- Instrument slot guides at 3.5mm, 4mm, 4.5mm, and 5mm
- Designed without right angles for user comfort
- Autoclavable

	.018" Slot	.022" Slot
Blue Gauge	BHG-18B	BHG-22B
White Gauge	BHG-18	BHG-22

Stainless Steel Brackets



Econo Line Stainless Steel Bracket

- High quality stainless steel bracket
- Torque and angulation are built into the base, not cut in like other metal brackets
- Laser welded hooks for extra strength
- One piece bonding base for better adaption to tooth
- Micro particle blasting in bonding base
- Micro size for excellent aesthetic

Roth Prescription		Torque	Angulation	L .018" Slot	R	L .022" Slot	R
Maxillary Arch	Central	12	5	ERU118L	ERU118R	ERU122L	ERU122R
	Lateral	8	9	ERU218L	ERU218R	ERU222L	ERU222R
	Cuspid - Hook	0	11	ERU318LHK	ERU318RHK	ERU322LHK	ERU322RHK
	Bicuspid - Hook	-7	0	ERU18BILHK	ERU18BIRHK	ERU22BILHK	ERU22BIRHK
Mandibular Arch	Lower Anterior	0	0	ERL118		ERL122	
	Cuspid - Hook	-11	7	ERL318LHK	ERL318RHK	ERL322LHK	ERL322RHK
	1st Bicuspid - Hook	-17	0	ERL418LHK	ERL418RHK	ERL422LHK	ERL422RHK
	2nd Bicuspid - Hook	-22	0	ERL518LHK	ERL518RHK	ERL522LHK	ERL522RHK

Standard Edgewise		Torque	Angulation	.018" Slot	.022" Slot
Upper Central		0	0	ESU118	ESU122
Upper Lateral		0	0	ESU218	ESU222
Lower Anterior		0	0	ESL118	ESL122
Cuspid		0	0	ES318	ES322
Bicuspid		0	0	ES18BI	ES22BI



Ultra Line Stainless Steel Bracket

- Colored arch ID
- Compound curvature to better fit tooth morphology
- Stainless steel foil mesh for bonding strength
- Precision-milled arch slot
- Rounded edges for patient comfort

Roth Prescription		Torque	Angulation	L .018" Slot	R	L .022" Slot	R
Maxillary Arch	Central	+12	+5	MRU118L	MRU118R	MRU122L	MRU122R
	Lateral	+8	+9	MRU218L	MRU218R	MRU222L	MRU222R
	Cuspid - Hook	0	+10	MRU318LHK	MRU318RHK	MRU322LHK	MRU322RHK
	Bicuspid - Hook	-7	0	MRU18BI		MRU22BI	
Mandibular Arch	Lower Anterior	-1	0	MRL118		MRL122	
	Cuspid - Hook	-11	7	MRL318LHK	MRL318RHK	MRL322LHK	MRL322RHK
	1st Bicuspid - Hook	-17	0	MRL418LHK	MRL418RHK	MRL422LHK	MRL422RHK
	2nd Bicuspid - Hook	-22	0	MRL518LHK	MRL518RHK	MRL522LHK	MRL522RHK

Standard Edgewise		Torque	Angulation	.018" Slot	.022" Slot
Upper Central		0	0	MSU118	MSU122
Upper Lateral		0	0	MSU218	MSU222
Lower Anterior		0	0	MSL118	MSL122
Cuspid		0	0	MS318	MS322
Bicuspid		0	0	MS18BI	MS22BI

Bracket Trays



No Light Bonding Bracket Tray

- Excellent for storing one complete set of brackets, second bicuspid to second bicuspid, upper and lower
- Smooth sliding removable top for easy access
- No light trays allow pre-loading of light-activated adhesives for storage until needed

Clear	—	BK-T
Transparent	No light	BK-NLT
Blue	No light	BK-NL
Black	No light	BK-BLKNL

Disposable Bracket Tray

- Individual bracket trays with adhesive backing for temporary storage and setup
- Deep wells help prevent bracket spillage
- Space for patients name, for easy identification
- Available with or without a mixing well
- Disposable; Single use only

	Tray	Tray with Mixing Well
Blue	BK-B2	BK-B2W
Teal	BK-T2	BK-T2W
White	BK-W2	BK-W2W

Bracket Stacker Trays

- 20 individual case compartments provide convenient, stackable storage
- Clear plastic lid is designed for smooth sliding and dispensing
- 4 colors for easy bracket prescription identification

Blue	BTS-B
Green	BTS-G
Red	BTS-R
Yellow	BTS-Y

Bracket Tray Holder

- Made from strong acrylic
- Holds up to 10 bracket trays
- Keeps setups away from contaminants and light
- Excellent for transporting bracket cases
- Attractive enough for any office décor

Smoke	BK-THSM
Clear	BK-TH




Bite Sticks

Bite Sticks

- Handle is made of high-impact plastic; will not flex or bend
- Sure grip handle is easy to hold, even when wet
- Stainless steel tip is swirled to grip bands easily and efficiently
- Autoclavable
- Short-tipped styles are easier to maneuver inside the mouth




		Standard-Tipped	Short-Tipped
Star		TS-3S Cobalt Blue	TS-3ST Cobalt Blue
Regular		TS-1S Cobalt Blue	TS-1(color)
Reverse		TS-2S Cobalt Blue	TS-2ST Cobalt Blue
Round		TS-1R Cobalt Blue	TS-RST Cobalt Blue



Ergo Bitestick

- Unique wave design and smooth surface for ergonomic comfort
- Short-tipped for easy maneuvering inside the mouth
- Stainless steel tip is swirled to grip bands easily and efficiently
- Autoclavable

	Short-Tipped
Regular	 TS-1ECB Cobalt Blue



Instruments

Lingual Contouring Instrument

- Ideal for placing lingual 3 to 3 wire retainers
- Angled head for easy access to lingual of anterior teeth
- Flat plastic tip makes cement contouring easy; will not adhere to bonding cement
- Hold wire in place while placing bonding cement
- Autoclavable

LCCI



Instrument Code Rings

- For use with pliers, scalars, band seaters, hemostats, and other orthodontic instruments
- High temperature silicone can withstand autoclaving

100 Code Rings CR-100(color)

Instrument Code Ring Kit
50 each of 12 colors, dispenser CR-KIT

 Blue B	 Green G	 Red R
 Black BLK	 Orange OR	 Tan T
 Brown BRN	 Pink PNK	 White W
 Dark Gray DG	 Purple P	 Yellow Y



Intra Oral Elastics

Latex Elastics

- Force in ounces is measured at 3x passive diameter
- Simple to dispense and easy for patient to use
- Packaged in color-coded tamper-proof patient bags

Available in
BULK PACKS!
Call for details.




Diameter	Light 2.5oz		Medium 3.5oz		Heavy 4.5oz		Extra Heavy 6oz	
1/8" / 3.2mm	Red A	G731	Green A	G733	Blue A	G735	Brown A	G736
3/16" / 4.8mm	Red B	G741	Green B	G743	Blue B	G745	Brown B	G746
1/4" / 6.4mm	Red C	G751	Green C	G753	Blue C	G755	Brown C	G756
5/16" / 7.9mm	Red D	G761	Green D	G-763	Blue D	G765	Brown D	G766
3/8" / 9.5mm	Red E	G771	Green E	G773	Blue E	G775	Brown E	G776

100 elastics per patient bag; 50 bags per dispenser box

Non Latex Elastics

- Eliminates latex allergic reaction concerns
- Use with aligner cases or ceramic and composite brackets
- Force in ounces is measured at 3x passive diameter
- Simple to dispense and easy for patient to use
- Packaged in color-coded tamper-proof patient bags




Diameter	Light 2.5oz		Medium 3.5oz		Heavy 4.5oz		Extra Heavy 6oz	
1/8" / 3.2mm	Red A	UNLN731	Green A	UNLN733	Blue A	UNLN735	Brown A	UNLN736
3/16" / 4.8mm	Red B	UNLN741	Green B	UNLN743	Blue B	UNLN745	Brown B	UNLN746
1/4" / 6.4mm	Red C	UNLN751	Green C	UNLN753	Blue C	UNLN755	Brown C	UNLN756
5/16" / 7.9mm	Red D	UNLN761	Green D	UNLN763	Blue D	UNLN765	Brown D	UNLN766
3/8" / 9.5mm	Red E	UNLN771	Green E	UNLN773	Blue E	UNLN775	Brown E	UNLN776

100 elastics per patient bag; 50 bags per dispenser box

Neon Non-Latex Elastics

- Eliminates latex allergic reaction concerns
- Six bright colors in each patient bag
- Force in ounces is measured at 3x passive diameter
- Simple to dispense and easy for patient to use
- Packaged in color-coded tamper-proof patient bags



Diameter	Light 2.5oz		Medium 3.5oz		Heavy 4.5oz		Extra Heavy 6oz	
1/8" / 3.2mm	Red A	UNLN731	Green A	UNLN733	Blue A	UNLN735	Brown A	UNLN736
3/16" / 4.8mm	Red B	UNLN741	Green B	UNLN743	Blue B	UNLN745	Brown B	UNLN746
1/4" / 6.4mm	Red C	UNLN751	Green C	UNLN753	Blue C	UNLN755	Brown C	UNLN756
5/16" / 7.9mm	Red D	UNLN761	Green D	UNLN763	Blue D	UNLN765	Brown D	UNLN766
3/8" / 9.5mm	Red E	UNLN771	Green E	UNLN773	Blue E	UNLN775	Brown E	UNLN776

100 elastics per patient bag; 50 bags per dispenser box


PERSONALIZE YOUR BAGS!

Keep your name in front of customers every day and generate referrals.

Minimum order sizes vary. Bags printed with "teeth" graphic are also available. Please call for details.



Elastics Accessories



Extra Oral Latex Elastics

- Activated force is reached when elastic is stretched to three times its original size

Diameter	8 oz	14oz	16oz
1/4"	G457	G458	G456
5/16"	G467	G465	G466
3/8"	G474	G475	G476
1/2"	G484	G485	G486

1000 elastics per bag



Elasti™ Placers

- Unique hook and pusher design
- Excellent tool to encourage wearing and regularly changing elastics

Vivid Multi Pack	100 per pack	EP-100
Individual Colors	100 per pack	EP-(color)

Individual		Vivid Multi Pack	
	Dark Blue DB		Glow-in-the-Dark GLOW
	Dark Green DG		Red/Cherry Scented RWC
	Dark Purple DP		Violet
			Magenta
			Yellow
			Azure

Colors are for illustration purposes only. Actual color of product may vary.

Elastics Dispenser

- Heavy duty acrylic
- Holes on the back for wall mounting
- Holds up to 5 boxes of elastics

Clear	ED-1AC
-------	--------



Procedure Aids



Abrasive Strip Holder

- Manual operation for safer, more accurate control
- Incisal guide helps prevent tissue laceration
- Holds stainless steel abrasive strips up to 8mm thick
- Also accepts many linen, cloth or plastic abrasive strips
- Autoclavable


Abrasive Strip Holder	ES-1000
Abrasive Strip Holder with Sample Strips	ESK



Abrasive Strips

- 4.75" long with extra-strong bonded grit
- Autoclavable


Width	Medium Grit	Medium Grit
	Single-Sided	Double-Sided
4mm	AS-S4M	AS-D4M
6mm	AS-S6M	AS-D6M
8mm	AS-S8M	AS-D8M
12 per pack		



Ortho Zip Ties

- Adjustable elastomeric ring can be used for a variety of applications that require tying or cinching
- Anchoring orthodontic appliances
- Ligating impacted cuspids
- Tying exposed cuspids
- Ligating overlay wires onto main wires
- Latex Free

Pearl	ZT-2PEARL
Silver	ZT-2SIL
8 per dispenser; 5 dispensers per bag	



Arch Gauge

- Use during IPR treatment planning to measure the width of teeth and the arch
- Marked with graduated 5mm lines
- Disposable, flexible plastic can be cut into sections for many orthodontic uses

Orange	AMDOR-25
White	AMD-25
25 per tube	

Tongue Blocks

Disposable Tongue Blocks

- Allows patient to relax their jaw while keeping their tongue away from their teeth during bonding and banding procedures
- Soft plastic for easy trimming with scissors
- Saliva ejectors may be attached to holes for maintaining a dry field
- Single use only




	Soft Unscented	Soft Cherry Scented	Super Soft Unscented
	Red	Red	Blue
Small	TBU-DS	TB-DS	TB-SSDS
Medium	TBU-D	TB-D	TB-SSD
50 per pack			

Saliva Ejector

- Saliva is removed through holes at end of vinyl tubing
- Metal insert allows tubing to bend in any direction
- Quick disconnect fitting for easy placement and removal
- Adaptor tubes for connection to slow speed ejector systems

Saliva Ejectors	3 per pack	SE-3
Adapter Tubes	3 per pack	AT-3




Cheek Expanders and Retractors

Spand Ezz™ Cheek Expander

- Flexible plastic cheek expanders allow access to the whole mouth by holding it open and depressing the tongue

	Disposable Single use only	Autoclavable
Small	PCKE-0D Green	PCKE-0 Clear
Medium	PCKE-1D Blue	PCKE-1 Clear
Large	PCKE-2D Red	PCKE-2 Clear
	10 per pack	Sold individually



Saliva Ejector System

- Non-intrusive design easily connects to Spand Ezz™ Cheek Expander to create a dry working field
- Disposable; Single use only


5 Saliva Ejectors, 5 Sets of Tubing	SE-PCKE
3 Adapter Tubes	AT-3



Double-Ended Retractors

- Take intra-oral photos without an assistant
- Elongated ends allow unrestricted access to anterior and posterior regions in the mouth
- Autoclavable

2 Adult	BPR-2A
2 Pedo	BPR-1A
Combo Pack 2 each size	BPR-KITA



Two Piece Photo Retractor Set

- Buccal extension holds lips and cheeks away from mouth
- Elongated handle for easy gripping and holding

	Autoclavable
Small	PHOTO-1
Large	PHOTO-2



Prophy Aids

Extra Oral Cheek Expanders

- Designed to prevent breakage, with extra long buccal extension that allows greater access to posterior region
- Optional tabs allow cheeks to be pulled further out

	Disposable Single use only	Autoclavable
Small	EXO-1DNT	EXO-1NT
Small with Tabs	EXO-1D	EXO-1
Regular	EXO-2DNT	EXO-2NT
Regular with Tabs	EXO-2D	EXO-2
	10 per pack	Sold individually



Wyred™ Cheek Expander

- Stainless steel spring wire support for lip and cheek retraction provides consistent and adjustable force for improved access to posterior teeth
- Autoclavable

Small Retractor	CKE-1A
Medium Retractor	CKE-2A
Wire	CKE-KITWIRE
Endcaps	CKE-ENDCAPS
Small Retractor Kit Retractors, wire and endcaps	CKE-KIT1A
Medium Retractor Kit Retractors, wire and endcaps	CKE-KIT2A



Prophy Cup Holder

- Great for fast individual tooth rebonding
- Easy single tooth cleaning

Disposable Holder and Prophy Cups 25 per pack; Single use only, Red	PS-1D
Autoclavable Prophy Cup Holder 1 Holder, Blue; 20 Disposable prophy cups; Single use only	PS-1
Refill Disposable Prophy Cups 144 per pack; Single use only	PROPHYCUPS

Aerosol Reduction Device

- Fits on CAVITRON® PROPHY-JET® Prophylaxis System to capture aerosol spray
- Use during tooth surface prep prior to bonding orthodontic brackets and other appliances
- Use to clean teeth after de-bonding brackets and after in-office bleaching procedures



Adapter Kit Reusable suction and assembly parts, 10 disposable cups	AR-K
10 Disposable Refill Cups	AR-R

Lab Products



Disposable Model Base Formers

- Create consistent model bases while saving cleanup time
- Durable plastic will not deform during model creation
- Flexible plastic makes separation easy
- Single use only

60 tray pack; 20 each small, medium, large	EBD-60
120 tray pack; 40 each small, medium, large	EBD-120



Shipping Cases

- Durable plastic will not deform during shipment
- Interior foam material provides protection and moisture retention; can be adjusted to fit

Model Shipping Case	MSB-52
Impression Tray Shipping Case	MSBIT-52
52 per box	

Custom Positioner Blanks

- Made of standard (78 Shore A) or soft (73 Shore A) EVA plastic
- Available in scented and unscented

Standard	78 Shore A	PB-(color)
Soft	73 Shore A	PBS-(color)

*Certain colors require a minimum order of 100.
Colors are for illustration purposes only.
Actual color of product may vary.

Soft

	Clear/Unscented C
	Clear/Bubble Gum CBG
	Clear/Mint CM
	Transparent Blue/ Mint* BM

Standard Solid

	White/Unscented* W
	Black/Licorice BL
	Pink/Bubble Gum PBG
	White/Licorice* WL

Standard Transparent


	Blue/Unscented* B
	Black/Unscented* BLK
	Blue/Mint BM
	Dark Blue/Mint DBM
	Green/Apple* GA
	Purple/Grape PG
	Red/Strawberry RS
	Yellow/Pineapple* YPIN

Standard Clear

	Unscented C
	Bubble Gum CBG
	Cherry CWC
	Chocolate* CCHOC
	Cinnamon* CCIN
	Grape* CG
	Licorice* CL
	Mint CM
	Pineapple* CPIN
	Raspberry* CR
	Strawberry* CS
	Vanilla* CV
	Watermelon* CW



Impression Trays



Disposable Impression Trays

- Raised ends to capture terminal molar and provide all-around detail
- Rounded edges provide patient comfort
- Four easy-to-use numbered sizes
- Single-use only

	Upper	Lower
Adult Medium/Large (#4)	GITN-U4	GITN-L4
Adult Small (#3)	GITN-U3	GITN-L3
Pedo Medium (#2)	GITN-U2	GITN-L2
Pedo Small (#1)	GITN-U1	GITN-L1
50 per bag		
Intro Kit	GITN-KIT	
Mixed upper & lower trays, 48 per bag, 6 of each size		



Impression Trays

Tuff Lock™ Disposable Impression Trays

- A unique mechanical locking system eliminates the need for spray adhesives
- High sides provide excellent margin and gingival detail
- Single use only



		Upper	Lower
Adult Large	Yellow	GIT-U1	GIT-L1
Adult Medium	Green	GIT-U2	GIT-L2
Adult Small	Blue	GIT-U3	GIT-L3
Pedo Medium	Orange	GIT-U4	GIT-L4
Pedo Small	Red	GIT-U5	GIT-L5
50 per bag			
Intro Kit	5 upper and 5 lower of each size (50 total)		GIT-KIT



Accessories

Alginate Flavoring

- Create great tasting alginate
- Concentrated flavor
- Sugar free
- 2oz bottle

Banana	ALGF-BAN
Bubble Gum	ALGF-BG
Cherry	ALGF-C
Cinnamon	ALGF-CIN
Cola	ALGF-COLA
Fruit Punch	ALGF-FP
Grape	ALGF-G
Lemon	ALGF-L
Orange	ALGF-O
Pina Colada	ALGF-PC
Peppermint	ALGF-PEP
Raspberry	ALGF-R
Strawberry	ALGF-S
Chocolate	ALGF-CHOC
Watermelon	ALGF-WM



Alginate Flavoring Holder

- Heavy duty acrylic
- Holds up to six bottles
- Smoke color

ALGF-HOLDER

Alginate Mixing Bowl

- Extremely smooth, pliable plastic for easy mixing
- Smooth interior makes removing mixed alginate easy
- Does not contain PVC
- 4" diameter

Super Soft Vinyl Alginate Mixing Bowl

- Molded from extremely flexible, smooth and pliable vinyl for controlled mixing and cleaning

Alginate Spatula

- Made of plastic; won't scratch bowls surface
- Flexibility and stiffness for good control
- Easy grip handle reduces slipping

Bowl	AMB-(color)
Super Soft Vinyl Bowl - Yellow	AMB-Y
Spatula	TS-(color)
Autoclavable Spatula - Blue	TS-ABA
Autoclavable up to 260°	

Colors are for illustration purposes only.
Actual color of product may vary.

	Dark Blue (no code)
	Red R
	Pearl PEARL
	Pearl Blue PEARLB
	Ocean OCEAN
	Kiwi KS
	Splash SPLASH
	Tropical Blue TB
	Tropical Green TG
	Tropical Orange TO
	Tropical Pink TP
	Tropical Yellow TY



Thermobite Bite Registration Wafers

- Specially-formulated resin provides highly accurate bite impressions of teeth and horizontal molar area
- Heat according to instructions, bite down, and trim
- Minor perforations allowed, for a correct occlusion with touching teeth surfaces
- Easy to store and transport to labs; resistant to temperature changes during transit

30 per pack **TB-30**

Impression Guards


- Take impressions without removing archwire
- Thin elastomer strip fits over brackets, keeping impression material from getting behind bracket tie wings
- Impression can be removed easily with minimal damage
- Each strip fits one full arch, and can be trimmed for patients who are not fully bonded
- Individually wrapped to prevent patient cross contamination

50 per box **IMP-100**



The Right Plastic for Each Application

Clarity, thickness, strength, aesthetics and specific application are all important things to consider when selecting thermoforming materials. This side-by-side comparison chart can help you select the right plastic for the right application.



Duraclear™ Plastic

Duraforce™ Plastic

Duraguard™ Plastic

Bondable to acrylic	X		
Creating clear retainers	X	X	
Creating thin retainers for better patient comfort and compliance	X	X	
Creating temporary retainers	X	X	
Creating tooth movement aligners	X	X	
Creating bruxism splints		X	X
Creating customized sports mouthguards			X
Creating temporary bridges	X	X	
Creating bite planes	X		
Creating whitening trays			X
Available in round and square sizes	X	X	X
Available in colors		X	X
Latex and BPA free	X	X	X

Plastic

Duraclear™ Plastic

- Use for retainers, tooth movement aligners and temporary anterior bridges
- Good combination of strength, durability and clarity
- Durable; resistant to cracking, bending out of shape or breaking
- Bondable to acrylic; store in a cool, dry place



Thickness	5" Square	125mm Circle	Quantity
.030" / ¾mm	DCS030X20S	DCS030X20R	20 per box
.040" / 1mm	DCS040X20S	DCS040X20R	20 per box
.060" / 1½mm	DCS060X10S	DCS060X10R	10 per bag
.080" / 2mm	DCS080X10S	DCS080X10R	10 per bag

Duraforce™ Plastic

- Average life: 24 months
- Use for bruxism splints, retainers, temporary bridges and tooth movement aligners
- Conforms to model when heated with minimal shrink
- Extremely durable; resistant to cracking or breaking
- Not for use with acrylics/bonding agents
- Store in a cool, dry place



Thickness	5" Square	125mm Circle	Quantity
.030" / ¾mm	DFS030X100S	DFS030X100R	100 per box
	DFS030X10S	DFS030X10R	10 per bag
.040" / 1mm	DFS040X100S	DFS040X100R	100 per box
	DFS040X10S	DFS040X10R	10 per bag
.060" / 1½mm	DFS060X50S	DFS060X50R	50 per box



Composite Aligner Button

- Made from strong plastic resin
- Bonds with acrylic resin or other bonding agent to most plastic appliances to attach elastics
- Grooves in base provide bonding retention
- Smooth low-profile design minimizes patient discomfort

10 per bag **PAB**

Sports Mouthguard Plastic

Duraguard™ Plastic

- Use for custom sports mouthguards and soft splints
- .040" thickness is good for creating whitening trays
- Durable EVA material forms soft and flexible, with great conformity to model
- Offer patients their favorite team colors
- Store in cool, dry place

	Yellow Y		Clear C		Green G
	White W		Orange O		Black BLK
	Silver SIL		Purple P		Blue B
	Red R		Light Blue LB		Aqua A
Multi Color MC					

Thickness	5" Square	125mm Circle	Quantity
.040" / 1mm	DGS(color)040X25S	DGS(color)040X25R	25 per box
.060" / 1½mm	DGS(color)060X25S	DGS(color)060X25R	
.080" / 2mm	DGS(color)080X25S	DGS(color)080X25R	
.120" / 3mm	DGS(color)120X25S	DGS(color)120X25R	
.160" / 4mm	DGS(color)160X25S	DGS(color)160X25R	

NOTE: Certain colors require specific minimum quantities when ordered. Please call for additional information.
 Colors are for illustration purposes only. Actual color of product may vary.



Sports Mouthguards and Cases

Tuff Guard™ Breathable Sports Mouthguard

- Innovative design provides upper and lower arch protection from braces while playing non-contact sports
- Three breathing vents in front allow optimum airflow
- Fits comfortably between bite pads with clearance for the tongue
- Choose colors from chart below

	Tuff Guard™	Tuff Guard™ with Strap
Single Mouthguard with Case	TG-(color)	—
5 per dispenser box	TGD-(color)	—
10 per bag	TGP-(color)	TGPS-(color)

Super SportsGuard

- Designed for upper and lower arch protection
- Soft enough to protect the gingival, but tough enough to stand up to hard use
- Choose colors from chart below

	Super SportsGuard	Super SportsGuard with Strap
Small	SMG-(color)	SMGS-(color)
Regular	MG-(color)	MGS-(color)
	10 per pack	

 Aqua A	 Clear C	 Purple P	 Violet V
 Azure AZ	 Green GRN	 Red R	 White W
 Blue B	 Maroon M	 Silver SIL	 Yellow Y
 Black BLK	 Orange OR	 Teal T	 Multi Color MC



Multi Color pack includes A, B, BLK, GRN, OR, P, R, SIL, W, and Y.
Colors are for illustration purposes only. Actual color of product may vary.



Sports Mouthguard Cases

- 1.5" deep x 3" wide x 2.5" long
- Strong hinge and secure lock
- 0.625" x 0.4" opening for mouthguard strap

Sparkle	RC-5MS(color)
Tropical Sparkle	RC-5MTS(color)
	12 per pack

Sparkle	 Aqua A	 Blue B	 Green GRN	 Purple P
	 Black BLK	 Cool Mint JCM	 Purple P	 Plum JPLM
Tropical Sparkle	 Blue B	 Orange O	 Green G	 Pink P

Wax and Silicone Boxes

Orthosil™ Silicone Patient Boxes

- Silicone used for temporary relief from appliance irritation
- Transparent and waterproof
- Resistant to hot liquid

50 per pack

SBMC



Orthosil™ Bracket Protector

- Thicker strips of silicone for covering an entire arch
- Ideal for patients who play musical instruments
- Provides added protection when used with a sports mouthguard
- One piece covers cuspid to cuspid

50 per pack

SB-P



Scented Relief Wax

- Clear box showcases bright colors
- Scented wax appeals to kids and adults

White	Unscented	WB-C
White	Mint	WBC-MINT
Yellow	Lemon	WBC-YL
Purple	Grape	WBC-PG
Green	Lime	WBC-GL
Orange	Orange	WBC-OR
Red	Cherry	WBC-RC
Multi Color	Assorted	WBC-MC
50 per pack		



Wax and Silicone Box Holder

- Made from strong, high quality clear acrylic
- Removable cover slides on and off for easy dispensing and re-stocking
- Holds 100 boxes

WB-HOLDER

Patient Relief Wax

- Ideal for fast, temporary relief from fixed brackets and appliances
- 5 wax strips per case
- 1.375" long x 2.125" wide case

Colors are for illustration purposes only.
Actual color of product may vary.

50 per pack **WB-(color)**

PERSONALIZE YOUR CASES!

Keep your name in front of customers every day and generate referrals.

Imprint ink color is Silver Metallic.
Minimum order sizes vary. Please call for details.



UNSCENTED

Silver
SIL

Kiwi
KS

Ocean
OCEAN

Glow-in-the-Dark
GLOW

Sparkle Cool Mint
SJCM

Sparkle Plum
SJPLM

Splash
SPLASH

Red
R

USA
USA

White
W

Tropical
Sparkle Blue
TSB

Tropical
Sparkle Green
TSG

Tropical
Sparkle Pink
TSP

Tropical
Sparkle Orange
TSO

Tropical Sparkle
Multi-Color
TSMC

Cherry
RWC

Blueberry
BB

Chocolate
CHOC

Apple
GA

Mint
GM

Orange
OR

Bubble Gum
PBG

Cotton Candy
PCC

Grape
PG

Cinnamon
RCIN

Coconut
WC

Vanilla
VV

Pineapple
YPIN

SCENTED

Take-Home Aids



Bruxism "S" Splint

- Protects teeth from abrasion due to bruxism
- Use over braces without affecting tooth movement
- Easily adjusted to fit each patient
- Includes appliance box and detailed instructions

Clear	Unscented	S2-BS
Pink	Bubble Gum Scented	S2-BSBG



Chewies™ Aligner Tray Seaters

- Same size (1.375" long x 0.4375" wide) as standard dental cotton roll – but better! Made of soft plastic, they are hygienic and easier to use.

White	Unscented	AC-25WPP
Green	Mint Scented	AC-25GMPP
Pink	Bubble Gum Scented	AC-25PBGPP
Yellow	Pineapple Scented	AC-25YPINPP
Purple	Grape Scented	AC-25PGPP


2 seaters per patient pack; 10 patient packs per bag



Patient Relief Aid Kit

- Great kit filled with essential items for patient relief
- Give to patients at start of orthodontic treatment

Patient Relief Wax, Elasti™ Bite Therapy Wafer, Ortho-LipGuard™ Strips	PRK-10
10 kits per pack	



Retainer Retrievers

- Remove appliances cleanly and efficiently
- Store in appliance case
- Packaged in individual tamper-proof patient bags

White	RR-1W
Tropical Blue	RR-1TSB
Tropical Green	RR-1TSG
Tropical Orange	RR-1TSG
Tropical Pink	RR-1TSP

10 patient packs per bag



Hook Cushions

- Helps prevent patient discomfort
- Use over bracket hooks to provide a buffer between patient's cheek and hook
- Soft plastic does not swell or absorb water
- Tight fit helps prevent slipping

14 per dispenser, 10 dispensers per pack	HC-2
---	-------------



APK-1

APK-2

APK-3

Aligner Patient Kits

Key aligner/retainer products used during treatment, packaged in these handy kits.

Sparkle Black Econo Lock™ Retainer Case, Sparkle Green Mini Case, 2 White Retainer Retrievers, 2 Green Mint Chewies™ Aligner Tray Seaters	APK-1
Tropical Sparkle Pink Econo Lock™ Retainer Case, Tropical Sparkle Blue Mini Case, 2 Tropical Green Retainer Retrievers, 2 Green Mint Chewies™ Aligner Tray Seaters	APK-2
Tropical Blue Retainer Retriever, Pontic Wax, 3 Green Mint Chewies™ Aligner Tray Seaters	APK-3



Retainer Case Labels

- Helps eliminate lost retainers
- Attach to the bottom of retainer cases

Lips, Round, 2.625" diameter	RC-L
Toothbrush, Round, 2.625" diameter	RC-L2
Toothbrush, Square, 2" wide x 1.5" tall	RC-L3
Text, Half Circle, 2.25" wide x 1.75" tall	RC-L4
100 per roll	



Ortho-LipGuard™ Strips

- Provides temporary relief of direct pressure to the lips caused by orthodontic appliances
- Applications include playing instruments and non-contact sports
- Thermoplastic material can be heated and formed to contours of the patient's mouth
- 2 guards per sealed patient pack; 5 packs per bag

Patent #6,080,923
Ortho-LipGuard is not a trademark of DENTSPLY.

LG-1



Elasti™ Bite Therapy Wafer

- Helps reduce discomfort associated with wearing orthodontic braces
- Chewing on wafer may increase blood circulation around recently moved teeth by stimulating and massaging gum tissue

Clear	Cherry Scented	EB25
Green	Mint Scented	EB-25M
Purple	Grape Scented	EB-25G
Pink	Bubble Gum Scented	EB-25PBG
25 per pack		



Retainer Case Clips

- Easily clips to Econo Lock™, Euro, and Square retainer cases
- Multicolored neon pink, green, blue, and yellow plastic

10 per pack **RC-KR**

Appliance Cases

Appliance Case Options

One-piece cases are designed with a locking mechanism that stands up to heavy use. Two-piece cases have a positive locking system that prevents accidental opening. Cases can be personalized with your information.

Germguard™ Retainer Cases

An antimicrobial agent is molded into the plastic to provide protection in a wet, moist environment that's susceptible to germ and bacteria growth.

Available in select styles and colors. Look for the green GermGuard™ seal.

PERSONALIZE YOUR CASES!

Keep your name in front of customers every day and generate referrals.

Imprint ink color is Silver Metallic.
Minimum order sizes vary. Please call for details.



Sparkle Aqua A Black BLK Blue B Burgundy BRG Cool Mint JCM Green GRN Plum JPLM Purple P Red R Multi Color MC			Tropical Sparkle Blue B Green G Orange O Pink P Multi Color MC Neon Green G Orange O Red R			Splash Blue BSPLASH Green GSPLASH Kiwi KS Ocean OCEAN Purple PSPLASH Red RSPLASH Splash SPLASH Tangerine TGNS USA USA White WSPLASH		
Unscented Black BLK Blue B Glow GLOW Purple P Pearl Green PRL Red R White W Yellow Y			Scented Apple GA Banana YB Blueberry BB Bubble Gum PBG Cherry RWC Cinnamon RCIN Coconut WC Cotton Candy PCC Licorice BL Mango MNGO Mint GM Orange OR Purple/Grape PG Raspberry PR Strawberry RS Tutti Frutti TF Vanilla WV Wintergreen WG			Two Color Candy Cane CANDY Lemon Lime LL		

Colors are for illustration purposes only.
Actual color of product may vary.





Econo Lock™ Retainer Cases

- 1.2" deep x 3.4" wide
- Big enough for most retainers, small enough to fit in a pocket or purse


Sparkle	RC-6S(color)
Sparkle Multi (SA, SB, SBLK, SGRN, SJCM, SJPLM, SP, SR)	RC-6SMC
Tropical Sparkle	RC-6TS(color)
Neon	RC-6N(color)
Unscented Multi (BLK, R, W, B)	RC-6MC
Unscented	
Scented	RC-6(color)
Splash	
20 per pack	



Euro Cases

- 1.5" deep x 3" wide x 2.5" long
- Deeper to allow for slightly larger appliances, including TMJ splints and nightguards

Not available in Pearl Green, Glow, Tropical Sparkle Multi, Germguard Sparkle Cool Mint or Plum. Also available: case with opening for mouthguard strap, see page 39.

Sparkle	RC-5S(color)
Sparkle Multi (SA, SB, SBLK, SGRN, SJCM, SJPLM, SP, SR)	RC-5SMC
Tropical Sparkle	RC-5TS(color)
Neon	RC-5N(color)
Unscented	RC-5(color)
GermGuard 	GG-5S(color)
12 per pack	



Square Cases

- 1" deep x 2.25" wide x 2.25" long
- Compact and easy to carry

Not available in Yellow or Germguard Sparkle Plum.

Splash	RC-2SPASH
USA	RC-2USA
Sparkle	RC-2S(color)
Sparkle Multi (SA, SB, SBLK, SGRN, SJCM, SJPLM, SP, SR)	RC-2SMC
Tropical Sparkle	RC-2TS(color)
Neon	RC-2N(color)
Unscented	RC-2(color)
Scented	
GermGuard 	GG-2S(color)
24 per pack	



Mini Cases

- 1" deep x 2" wide x 3.375" long
- Excellent for smaller appliances such as aligners & whitening trays
- Great fit in pockets and backpacks

Not available in Burgundy, Sparkle Aqua, Cool Mint, Plum or Purple.


Sparkle	RC-7S(color)
Sparkle Multi (SA, SB, SBLK, SGRN, SJCM, SJPLM, SP, SR)	RC-7SMC
Tropical Sparkle	RC-7TS(color)
Splash Blue	RC-7SPASHB
Splash Purple	RC-7SPASHP
25 per pack	



Two-Piece Cases

- 1.0625" deep x 3.25" wide
- Slim design for clear aligners and bleach trays
- Positive locking system prevents accidental opening

Not available in Germguard Sparkle Black, GermGuard Sparkle Red, Coconut or Pineapple.

Sparkle	RC-1S(color)
Sparkle Multi (SA, SB, SBLK, SGRN, SJCM, SJPLM, SP, SR)	RC-1SMC
Tropical Sparkle	RC-1TS(color)
Neon	RC-1N(color)
Unscented	
Scented	RC-1(color)
Two Color	
GermGuard 	GG-1S(color)
18 per pack	




Lab Cases

- Two sizes: 1.5" deep x 3.25" wide or 2" deep x 3.25" wide
- Fits most lab appliances including bionators, positioners and nightguards
- Positive locking system prevents accidental opening
- Top vents help reduce odor

Not available in Tropical Sparkle Blue.

	1.5" deep	2" deep
Sparkle Blue	LC-1B	LC-3SB
Tropical Sparkle Multi	LC-1TSMC	LC-3TSMC
Tropical Sparkle	LC-1TS(color)	LC-3TS(color)
12 per pack		



**Elastomeric
Ligature Chain**
Page 3

Distributed by:

Prestige Dental Products
3031 E. Coronado St. #D
Anaheim, CA 92806
877-772-3888

Catalog expires December 31, 2014. Not responsible for typographical errors. Colors are for illustration purposes only. Actual color of product may vary. All descriptions and items listed herein are subject to change.

© 2014 DENTSPLY RTE-105-14

DENTSPLY Raintree Essix, 7290 26th Court East, Sarasota, FL 34243, USA