



Prestige Dental Products

www.prestigedental.com

Toll Free: 877-772-3888 Fax: 714-632-8821

Exp: 06/30/2016



Brackets

MIM Brackets 20/Set \$29.95
Buy 5 Get 2 Free - Net \$21.39
Buy 10 Get 5 Free - Net \$19.96
Buy 25 Get 25 Free - Net \$14.97

Econo Brackets 20/Set \$19.95
Buy 5 Get 2 Free - Net \$14.25

Ortho-x Brackets 20/Set \$11.95
Buy 5 Get 1 Free - Net \$9.95

Ceramic Brackets 20/Set \$94.95
Buy 5 Get 2 Free - Net \$67.82

Ultra Clear Brackets 20/Set \$139.95
Buy 5 Get 2 Free - Net \$99.96

Composite Brackets 20/Set \$49.95
Buy 5 Get 1 Free - Net \$41.62



Tubes

1st Molar or 2nd Molar Tubes
 UR/UL, LR/LL 5/Pkg \$14.95

Molar Band

Buccal Tube



+




+

CLEAT



Assembly Bands + Tube + Cleat
 1st or 2nd Molar \$4.74
Buy 5 Get 1 Free - Net \$3.95b



ArchWire

Nickel Titanium 10/Pk \$8.95
Buy 5 Get 2 Free - Net \$6.39

Coated Nickel Titanium 10/Pk \$21.95
 Tooth Tone, Red, Blue, Gold or Green
Buy 5 Get 1 Free - Net \$18.29

Thermal Nickel Titanium 10/Pk \$8.25
Buy 5 Get 2 Free - Net \$6.39

RCS — Reverse Curve of Spee
 Nickel Titanium 10/Set \$19.95
Buy 5 Get 1 Free - Net \$16.62

Twist 3 Strand NiTi 10/Pk \$6.95
Buy 5 Get 1 Free - Net \$5.79

Stainless Steel 10/Pk \$3.95
Buy 5 Get 2 Free - Net \$2.82

NiTi Open Coil Spring 3x7"/Pk \$19.95
Buy 5 Get 1 Free - Net \$16.62



Cheek Retractor
 Adult or Child
 2/Pk \$7.95



Multi Adjustable
 Facemask
 Red, Teal & Blue
 \$37.99
 Face Mask
 (Reverse Pull)
 \$39.99

Pliers and Instruments

Distal End Cutter



Safety hold is ideal for cutting wires within the patient's mouth. Plier shear cuts hard wires close to the buccal tube and holds loose wire. Maximum cutting capacity: Archwires up to .022" x .028".

• 800-205

Small Distal End Cutter



Compact tip for easier access. Safety hold is ideal for cutting wires within the patient's mouth. Plier shear cuts hard wires close to the buccal tube and holds loose wire. Maximum cutting capacity: Archwires up to .021" x .025".

Regular Handle • 800-215

Long Handle • 800-215XL

Lab Hard Wire Cutter



The extra strong edges make it one of the most popular hard wire cutters. Designed to cut all sizes and types of archwires, it features a 15° cutting angle for easy intraoral access. Cuts archwires up to .022" x .028".

• 800-149

Pin & Ligature Cutter



This cutter's super-sharp diamond honed edges are effective for cutting all soft wires, pins and elastics up to .016" diameter. 15° cutting angle.

• 800-107

Mini Pin & Ligature Cutter



Fine tips to easily access hard to reach areas. Cuts wire ligatures, pins and elastics. Maximum cutting capacity: Soft wire up to .012". 7° cutting angle.

• 800-146

Angled Direct Bond Bracket Remover



Designed to quickly and safely remove direct bond brackets. The 60° angle allows access to anterior and posterior brackets.

• 800-312

Straight Direct Bond Bracket Remover



Safely removes all direct bond brackets and attachments. Nylon tip rests on incisal/occlusal edges of teeth, while metal edge fits under bracket flange. May remove remaining adhesive.

• 800-313

3/16" Replacement Tips (3/pk)

• 800-313T

1/4" Replacement Tips (3/pk)

• 800-313TT

Replacement Blade (1/pk)

• 800-313B

Posterior Band Remover



Safely and easily removes posterior bands.

• 800-105

3/16" Replacement Tips (3/pk)

• 800-105T

1/4" Replacement Tips (3/pk)

• 800-105TT

How Plier



Matching serrated pads assure a firm, positive grip. Rounded tips for patient comfort and safety. 3/32" diameter pads.

• 800-106

Angle How Plier



Efficient utility plier with serrated beaks for a firm grip. 45° angled tips to access hard to reach areas.

• 800-301

Weingardt Plier



Serrated tips hold wires firmly at a convenient working angle. Tapered beaks fit easily between brackets and hard to reach areas. Our most popular utility plier.

• 800-303

Tweed Loop Forming Plier



Convex/concave beak makes precise Omega loops. Tip has 3 forming sections: .045", .060", and .075". Etched tip surfaces prevent wire slippage.

• 800-406

This Page Only

Buy Any 5 of the Same Pliers and Get 1 Free!

\$59.95

Bird Beak Plier



Versatile plier for working with round wire up to .030". Round and pyramid beaks are smooth to prevent wire scoring. Can also be used to form springs.

• 800-426

\$59.95

Wire Contouring Plier



Smooth working area is ideal for forming and contouring archwires. Useful in placing accentuated curves in wires. Excellent plier for forming nickel titanium archwires.

• 800-401

\$59.95

Light Wire Plier



Designed for working on small wires and forming wire springs. Tapered beaks make it easy to bend small diameter loops. Bends wires up to .020".

• 800-425

\$59.95

3-Prong Clasp Adjusting Plier



Precision tips for accurate wire and clasp adjusting and contouring. Also an excellent plier for bending of orthopedic appliances. Superior strength for wires up to .030".

• 800-315

\$59.95

Jarabak Plier



Excellent for precise bending and forming loops in light wires up to .020". Three sets of precision grooves assure a firm grip.

• 800-302

\$59.95

Elastic Separating Plier



Angulated beaks for accessibility, ease and accuracy. Grooved beaks assure positive hold for positioning elastics.

• 800-408

\$29.95

Arch Torquing & Bending Plier



Ideal for forming square or rectangular archwires. Beveled edges prevent scoring of archwire. Bends wires up to .022" x .028".

• 800-209

\$59.95

Mathieu Pliers



\$29.95



- Wide Tip • 800-101
- Narrow Tip • 800-121
- 1/2 Fine Tip • 800-123
- Hook Tip • 800-110






Amber Elastics



New Stylish Packaging!

- **Manufactured from Surgical Grade Latex**
- **Precision Cut for Reliable Forces**
- **Color-Coded Bags for Easy Identification**
- **Available in Amber or Neon Colors**
- **100 Bands per Bag, 50 Bags per Box**

\$29.95

					
DIAMETER	1/8"	3/16"	1/4"	5/16"	3/8"
mm Size	3.2	4.8	6.4	7.9	9.4
Bag Color	Purple	Green	Red	Blue	Orange
Light 2.5 oz. A	350-102AS	350-202AS	350-302AS	350-402AS	350-502AS
Medium 4.5 oz. B	350-104AS	350-204AS	350-304AS	350-404AS	350-504AS
Heavy 6.5 oz. C	350-106AS	350-206AS	350-306AS	350-406AS	350-506AS

This Page Only
Buy Any 5 of the Same Pliers and Get 1 Free!



Bracket Placer/Slot Aligner **\$8.95**
Design features a bracket holder on one end and a slot aligner on the other. This end fits into the slot for accurate bracket placement. US-332



Direct Bond Bracket Holder **\$8.95**
Unique reverse action tweezers with finely serrated beaks for a nonslip grip. Easy release of both metal and plastic direct bond brackets. (fits narrow to wide twin) US-309



Deluxe Bracket Holder **\$8.95**
Stainless steel holder is designed for accurate placement of anterior and posterior direct bond brackets and buccal tubes. Excellent field of vision. USA-44-302



Bracket Placing Tweezers **\$8.95**
Slim, precision-curved tips allow for an excellent field of vision while maintaining a firm grip to assure accurate bracket placement. US-44-321




Posterior Direct Bond Bracket Holder **\$19.98**
Designed for direct bonding of posterior teeth. Small head for excellent field of vision. Serrated tips for a positive grip. USA-44-311



Bracket Height Gauge
Similar to the Boone Gauge but in a straight, double-ended design. Accurately measures height of bracket on tooth from the incisal edge. Gauge has 3.5mm, 4.0mm, 4.5mm and 5.0mm.
For .018 Slots • For .022 Slots
Also available in Plastic:
For .018 Slots • For .022 Slots



Boone Bracket Position Gauge
This unique gauge is placed on the incisal/occlusal surface of the tooth to simplify bracket positioning. Measures 3.5mm, 4mm, 4.5mm and 5mm.



Ligature Director (Double-Sided) **\$14.95**
Design features ligature "tie-spin" tip on one end; ligature tucking on the other. All stainless steel construction. USA-44-201



Schure Band Seater **\$14.95**
This design features a finely serrated band pusher on one end, all purpose scaler on the opposite end. USA-44-208



Band Biter **\$9.95**
AUTOCLAVABLE! Multi-purpose seater with hardened stainless steel tip. Molded handle also serves as bite stick.

Ketac Cem Radiopaque



Triple Pack: 36ml liquid and 99gm powder

\$169.95

Intro Kit: 12ml liquid, 33gm powder & mix pad.

\$69.95

Permanent Glass Ionomer Luting Cement for Crowns, Bridges, Posts, Inlays, Onlays and Orthodontic bands. Adheres to both tooth structure and metal oxides with no conditioning for an easier technique than hybrid cements.

Fuji Ortho Band Paste Pak



Starter Pak : 1 refill, 1 mixing pad, 1 paste pak dispenser & mixing spatula. **\$199.95**

Refill Contains: 2 - 13.3g paste pak cartridges & mixing pad **\$129.95**

Fluoride Releasing Glass Ionomer Band Cement
A blue-colored glass ionomer for bonding orthodontic bands. It can also be used for the cementation of appliances (e.g. Herbst or expansion appliances). The Paste Pak delivers consistent mixes every time, with a 3½ minute working time, and a 3 minute final set time.

Prestige Dental Products

Exp: 06/30/2016

\$7.95/Pk

\$8.95/Pk

Order Form














































Power Chain Ties. 15"/Spool

Patient Stick Ties 1000/Pk

Fax: 714-632-8821

Cont. Qty Short Qty Long Qty

Regular Qty Mini Qty


	Hot Pink	771-001		771-101		771-201		771-501		771-701	
	Light Pink	771-002		771-102		771-202		771-502		771-702	
	Green	771-003		771-103		771-203		771-503		771-703	
	Light Green	771-004		771-104		771-204		771-504		771-704	
	Dark Green	771-005		771-105		771-205		771-505		771-705	
	Blue	771-006		771-106		771-206		771-506		771-706	
	Sky Blue	771-007		771-107		771-207		771-507		771-707	
	Light Blue	771-008		771-108		771-208		771-508		771-708	
	Cobalt Blue	771-009		771-109		771-209		771-509		771-709	
	Orange	771-010		771-110		771-210		771-510		771-710	
	Yellow	771-011		771-111		771-211		771-511		771-711	
	Purple	771-012		771-112		771-212		771-512		771-712	
	Fuschia	771-013		771-113		771-213		771-513		771-713	
	Teal	771-014		771-114		771-214		771-514		771-714	
	Aqua	771-015		771-115		771-215		771-515		771-715	
	Met. Green	771-016		771-116		771-216		771-516		771-716	
	Bright Red	771-017		771-117		771-217		771-517		771-717	
	Fluor Green	771-018		771-118		771-218		771-518		771-718	
	Fluor. Pink	771-019		771-119		771-219		771-519		771-719	
	Emerald	771-020		771-120		771-220		771-520		771-720	
	Tooth Tone	771-021		771-121		771-221		771-521		771-721	
	Metallic Gold	771-022		771-122		771-222		771-522		771-722	
	Burgundy	771-023		771-123		771-223		771-523		771-723	
	Pearl	771-024		771-124		771-224		771-524		771-724	
	Aquamarine	771-025		771-125		771-225		771-525		771-725	
	Lime	771-026		771-126		771-226		771-526		771-726	
	Ocean Blue	771-027		771-127		771-227		771-527		771-727	
	Watermelon	771-028		771-128		771-228		771-528		771-728	
	Tangerine	771-029		771-129		771-229		771-529		771-729	
	Violet	771-030		771-130		771-230		771-530		771-730	
	Pearl Waterm.	771-031		771-131		771-231		771-531		771-731	
	Grape	771-032		771-132		771-232		771-532		771-732	
	Met. Ocean Blue	771-033		771-133		771-233		771-533		771-733	
	Lilac	771-034		771-134		771-234		771-534		771-734	
	Navy	771-035		771-135		771-235		771-535		771-735	
	Dark Gray	771-036		771-136		771-236		771-536		771-736	
	Sparkle Blue	771-037		771-137		771-237		771-537		771-737	
	Sparkle Green	771-038		771-138		771-238		771-538			
	Sparkle Red	771-039		771-139		771-239		771-539			
	Periwinkle	771-040		771-140		771-240		771-540			
	Clear	771-041		771-141		771-241		771-541			
	Silver	771-042		771-142		771-242		771-542			
	White	771-001		771-101		771-201		771-543			
	Black	771-044		771-144		771-244		771-544			
	Red	771-045		771-145		771-245		771-545			

Intraoral Auxiliaries

This Page Only

Buy 5 of the same item, get 1 free!

Preformed Ligature Wire




(1,000 per pack)
 .009 • 270-009
 .010 • 270-010
 .011 • 270-011
 .012 • 270-012

Pretwist (500 per pack)
 .010 • 270-016

Toothtone (100 per pack)
 .012 • 270-017

Kobayashi Tie Hooks



(100 per pack)
 .012 • 270-002
 .014 • 270-004

Pretwist (100 per pack)
 .012 • 270-007
 .014 • 270-006

Stainless Steel Open & Closed Coil Spring



3 ft. Super Resilient

Open .010 x .030 • 270-040
 Open .010 x .036 • 270-050
 Closed .010 x .030 • 270-060

Pletcher Springs



Provides Continuous Force for Space Closure

.012 x .030 • 270-0281
 (50 per pack)


Australian H.T. Wire



Smooth Finish,
 Highly Resilient,
 25 Foot Spool

.014 • 270-081 .018 • 270-083
 .016 • 270-082 .020 • 270-084


Ball Clasps



1/2 Hard Retainer
 Wire with a
 Smooth Ball End

.028 • 270-0241 (100 per pack)
 .028 • 270-024 (1,000 per pack)
 .032 • 270-0261 (100 per pack)
 .032 • 270-026 (1,000 per pack)

Crimpable Stops & Hooks



a. Crowder Stops (100 per pack)

Up to .016 - 2mm • 270-090
 Up to .016 - 3mm • 270-091
 Up to .018 x .025 - 2mm • 270-092

b. Ball Hooks .018 x .025 (20 per pack) • 270-200

c. Crimp Stop .021 x .025 (10 per pack) • 270-201

d. Crimp Tube with Hook .017 x .025 (10 per pack) • 270-207

d. Crimp Tube with Hook .019 x .025 (10 per pack) • 270-209

Bondable Attachments



Button (20 per pack)
 • 270-206



Micro-Button (10 per pack)
 • 270-205



Eyelet (10 per pack)
 • 270-213



Cleat (20 per pack)
 • 270-212

3x3 Lingual Retainers



(5 per pack)
 22mm • 270-122
 24mm • 270-124
 26mm • 270-126
 28mm • 270-128
 30mm • 270-130
 Kit (20 pieces) • 270-150

Rotation Wedge




(100 per pack)
 Gray • 270-021
 Clear • 270-022

Niti Separators



(25 per pack)
 Small (bicusps) • 225-144
 Large (molars) • 225-146

Brass Separating Rings



(100 per pack)
 .020 • 270-216
 .025 • 270-217

Radio Opaque Single Separators



(1,000 per pack)
 Standard (1,000 per pack) • 325-302
 Large (1,000 per pack) • 325-303

Photographic Mirrors

\$59.95

Adult Occlusal/Child Occlusal



• 801-020

Buccal Lateral



• 801-021

Adult Occlusal/Child Occlusal



• 801-022

Adult Buccal/Adult Occlusal



• 801-023

Lingual Upper/Lingual Lower



• 801-024

Child Buccal/Adult Buccal



• 801-025

Child Occlusal/Child Buccal



• 801-026

Elastic Placers



Assists for easy placement of intra-oral elastics. Mixed colors.
Package of 100

\$9.95

Prophy Paste



Plain paste designed specifically for preparation of orthodontic bonding surfaces. 200 cups per box.
Prophy Paste

\$24.95

Relief Wax



For fast temporary relief of irritated tissue.
5 strips per patient pack. 50 packs per box.

Comes with 1 Flavor/Box

\$23.95

Utility Wax





3/16" diameter strips x 11"
90 strips per box.

-Clear
-Red

\$9.95

Rapid Expanders (Leone)

	 arch mm	 body mm	 expansion mm	activation turns for maximum expansion level			μg
	1.5	12	8 mm	35	0.8 mm	A0620-06	1
	1.5	14	9 mm	40	0.8 mm	A0620-09	1
	1.5	15	11 mm	50	0.8 mm	A0620-11	1
	1.5	18	13 mm	60	0.8 mm	A0620-13	1

8mm, Body 12 mm 1/pkg

9mm, Body 14 mm 1/pkg

11mm, Body 15 mm 1/pkg

13mm, Body 18 mm 1/pkg

\$19.99

Arch Markers



Wax coated wooden sticks for easy arch marking.

ARCH WIRE MARKING STICK
DISPENSER **\$16.00**

WAX ARCH WIRE MARKER
ORANGE 10 **\$17.00**


Retainer Cases



Protect your appliances with a quality case. This attractive case offers incredible strength and a variety of fun colors!
12 per box.

\$7.95

Excellent II Disposable Impression Trays



#1 Large-Upper, #2 Large-Lower,
#3 Medium-Upper, #4 Medium-Lower,
#5 Small-Upper, #6 Small-Lower,
#7 Upper-Right/Lower-Left,
#8 Upper-Left/Lower-Right,
#9 Anterior

Baby Blue
12/Bag

\$4.95

Micro Brush



Regular, Assortment
Package, Pkg. of 400

Fine, Pkg. of 400

Super Fine, Pkg. of 400

\$24.95

Package of 400

Ortho-Choice LC Bracket Adhesive



Kit: 2 X 4 gm syringes
adhesive, 5 mL
bonding resin, 6 mL
etchant + accessories

\$52.99

Light Cure - Fluoride Releasing Provides secure bonds, easy handling, and rapid curing on demand, saving valuable chairside time. Viscosity is designed to prevent bracket float. Fluoride release protects enamel against recurring caries and decalcification under brackets.

Ligature Ties Dispenser



Ligature Dispenser Sold with 24
assorted colors Patent Stick Ties
\$189.95

Ligature Dispenser Sold Empty
\$69.99



Digital Phosphor Plates Scanner

- INSTANT & INTELLIGENT
- LASER FAST IMAGE PROCESSING
- MORE FLEXIBILITY
- POWERFUL SOFTWARE INCLUDED



\$6299

FREE! 12 Imaging Plates
FREE! Imaging Software
FREE! 500 Barrier Envelope
FREE Shipping in USA

Ledex WL-90+ (Dentmate)

2400mW/cm²



\$424.95

ART-L5 (Bonart)

1200mW/cm²
LED Cordless
Curing Light



\$279.95

Cybird XD Cordless LED Curing Light (Dentazon)

2700mW/cm²
High power &
plasma emula-
tion modes



\$639.95

Radii Plus (SDI)

1500mW/cm²
Lightweight
Cordless
Curing Light
360° Degree
Rotate Tip for
Full Access



\$849.95

Dr's Light Wireless Curing Light (GoodDr's)

5 Seconds



\$749.95

Ledex LED Curing Light (Dentmate)



**10 Second
Cure 3mm Composite**
Available in 7 Colors:
BUY 2+ for \$199 ea.

\$219.95

Mini-Blast Cordless Curing Light (First Medica)



\$385.95

Blast-Lite LED Cordless Curing Light (First Medica)

5 Seconds



\$589.95

Litex 696 Cordless LED Curing Light (Dentamerica)

1200mW/cm²



Autoclavable 8.5mm Fiber optic Tip
360° Degree Tip for Full Access

\$449.95

LED Cordless Curing Light (Wi)



\$139.95

Litex 680A (Dentamerica)



600mW/cm²
\$289.95

iCure Cordless LED Curing Light (Pac-Dent)

1400mW/cm²
Light guides
compatible with
L.E. Demetron &
Demi Curing
Lights



\$409.95

Coltolux LED (Coltene Whaledent)



High Power
LED Curing Light

\$419.95

Fusion 4 Curing Light (DentLight)

2700mW/cm²
5 second cure on all com-
posite
360° Swivel head
90 minute recharge time
Color:
Silver, Pink, Blue, Gold



\$895.95

SPEC 3 LED Curing Light (Coltene & Whaledent)

High Intensity
-5 sec. Cure 2mm in Standard Mode
-1 sec. Cure 2mm in 3K Mode

\$989.95



Endomax Cordless Endo Handpiece (Dentamerica)

Compare to
Endo-Mate TC(NSK)

\$695.95



Endo A Class Cordless Endodontic Handpiece (Saeyang)

\$549.95

